

**Site Preparation Specification****Purpose of Procedure**

To assure that the installation of Agilent instruments and systems can be completed successfully by careful preparation and evaluation of the installation site and by ensuring the availability of appropriate utilities, consumables and supplies.

Customer Responsibilities

Customers should ensure that all necessary operating supplies, consumables and usage dependent items such as columns, vials, syringes and solvents required for the successful installation of instruments and systems are available. In addition, customers must provide the Agilent Support Provider with appropriate administrator-level access to the system during the installation. Installation sites should be prepared in accordance with the following specifications.

Important Information

If you have problems in providing any of the following, please contact your local Agilent office for assistance. Assistance with user specific applications may be provided but should be contracted separately. Users of the instrument should be present throughout the installation and familiarization otherwise important operational, maintenance and safety information may be missed.

Procedure Checklist**A) Method Validation Pack requirements****PC Hardware Requirements**

- ☐ Pentium II processor operating at 400 MHz or higher
- ☐ At least 96 MB of RAM for a single instrument configuration under Windows NT and 128 MB of RAM for Windows 2000 configurations.
- ☐ At least 128 MB of RAM for a multi instrument configuration (256 MB recommended)
- ☐ Hard-disk with at least 4 GB capacity
- ☐ CD-ROM drive
- ☐ Display: min. resolution 1024 x 768, more than 256 colors
- ☐ Microsoft Windows compatible pointing device
- ☐ A backup device or network-based backup strategy

Software Requirements Operating System

- ☐ Windows NT 4.0 with service pack 6a or higher and Internet Explorer 4.01 or higher
or
- ☐ Windows 2000 with Service Pack 2 or higher.
- ☐ User logon with Administrative capabilities

Software Requirements Chromatographic Data system

- ☐ ChemStation A.08.04 or higher for
 - A/D (G2172AA)
 - CE (G1601A)
 - CE/MSD (G2201AA)
 - GC (G2070AA)
 - LC (G2170AA)
 - LC/MSD (G2710AA)
- ☐ ChemStore B.02.01 (G2181BA)

B) Additional Requirements for Client/Server Systems**The following requirements only apply if Method Validation Pack needs to be installed in a client/server environment**

- ☐ Existing ChemStore Client/Server configuration
- ☐ Oracle Client software 8.1.7.0.0
- ☐ At least 8 GB free disk space on the server for the creation of a small database, 34 GB free disk space for a medium database and 96 GB MB free disk space for a large database, respectively.

Tick Boxes