

# Lower Operating Costs, Safer Operation, More Productive. Agilent 4100 MP-AES



**Steve Wall, Agilent Technologies**

# Common Challenges facing Laboratories doing Elemental Analysis Today

- Increased need for multi-element determination over a wide dynamic range
- Desire to reduce the overall cost of analysis due to rising costs (instrument supplies and consumables; power; labour etc)
- Difficulty in sourcing some gases – esp. in remote areas and emerging geographies
- Availability of suitably trained personnel to develop methods, perform sample measurement and interpret results
- Some laboratories under pressure to improve safety by removing flammable gases

# Agilent 4100 Microwave Plasma-Atomic Emission Spectrometer (MP-AES)

New technique for elemental determination using atomic emission

- Microwave excited plasma source
- Nitrogen based plasma - runs on air (using a N<sub>2</sub> generator)

## Improved performance compared with flame AA:

- Higher sample throughput with fast sequential measurement
  - More than 2x faster than conventional flame AA
- Superior detection limits and improved dynamic range

## Easy to use:

- New generation software featuring automated optimization and software applets that load a preset method
- One piece torch with easy torch removal and replacement – no alignment

## Reduced operating costs:


- Runs on air – eliminates need for Acetylene, Argon, etc.
- Eliminates need for source/hollow cathode lamps
- Simple installation – no chiller, 10 A supply

## Improved Safety:

- Eliminates need for flammable gases and cylinder handling
- Safe, reliable unattended multi-element overnight operation




# Schematic Diagram



1. Instrument exhaust
2. Pre-optics window
3. Torch loader
4. Peristaltic pump
5. Plasma enable button
6. High voltage power supply
7. Electronics (control PWB)
8. Pre-optics
9. Plasma (*Magnetron*)
10. Monochromator with CCD detector
11. External gas control module
12. Cooling air inlet
13. Inlet gas connections



# Microwave Plasma Emission Overview



# How Does MP-AES Work?

## Agilent's patented microwave waveguide technology

- Using nitrogen as the plasma gas gives a robust plasma with a conventional torch. Nitrogen can be supplied either via bottled gas or via nitrogen generator
- Magnetic excitation gives a toroidal plasma and an effective central zone for sample injection
- The microwave magnetically excited nitrogen plasma
  - Provides a robust, high temperature source in conventional torches (approx. 5000 K)
  - A cooler central channel suitable for sample atomization



# Performance – Typical Detection Limits

*DL's in ppb, clean water samples*


Element	Flame AA	MP-AES	Element	Flame AA	MP-AES
K	0.8	0.65	As*	60	45
Ca	0.4	0.05	Cd	1.5	1.4
Mg	0.3	0.12	Cr	5	0.5
Na	0.3	0.12	Mn	1.0	0.25
Au	5	1.8	Pb	14	4.4
Pt	76	4.5	Sb	37	12
Pd	15	3.8	Se*	500	70
Ag	1.7	0.5	V	20	0.2
Rh	4	0.5			

*\* 30 second integration time used for As and Se*

***3 sigma DLs using a 10 second integration time***

# Performance – Extended Working Range


Au (267.595 nm) Calibration



MP - normal AA range

Intensity = 15775.550 \* Concentration - 0.205  
Correlation coefficient: 1.00000

Au (267.595 nm) Calibration




MP - Extended Range

Intensity = (16408.841 \* Concentration + 0.057) / (1 + 0.001 \* Concentration)  
Correlation coefficient: 0.99998



# Reduced Cost of Ownership

## MP-AES cost savings versus flame AA



Cost savings based on reductions in gas expenses ONLY!.

Assumes measuring 9 elements (4 nitrous-oxide elements) in batches of 100 samples, 3 days/week using the N<sub>2</sub> generator. For a detailed cost comparison, consult your local Agilent representative

# Simple Torch Installation – No Alignment

## Torch installation in three easy steps

**1** Open the torch loader



**2** Insert the torch



**3** Close the torch loader



# New Generation MP Expert Software

- Windows 7, worksheet based software with significant enhancements
  - Fresh, clean look
- Provides capability for:
  - Applet style operation using preset methods, or
  - Access to full capabilities
- Automated optimization tools
- Innovative and simple to use background correction
  - No correction points to select
- Extensive help with hover tips



# New Generation MP Expert Software

- Retains familiar "worksheet" based approach

The screenshot displays the MP Expert software interface with the following components:

- Worksheet Summary:** A table showing results for various samples and standards. The 'Standard 1' row is highlighted.
- Results matrix:** A detailed table of results for Sb (217.582 nm) across 10 replicates.
- Spectra:** A plot of Intensity vs. Wavelength (nm) for Sb (217.582 nm), Standard 1, Rack S1 Tube 2, showing a sharp peak at 217.582 nm.
- Calibration:** A plot of Intensity vs. Concentration (ppm) for Sb (217.582 nm), showing a linear relationship with a correlation coefficient of 1.00000.
- Replicate data:** A table showing the concentration and intensity for 10 replicates of Standard 1.
- Summary Panel:** A small panel on the left showing the element name (Sb), wavelength (217.582 nm), concentration (5.00000000 x ppm), and RSD (%): N/A.

Rack: Tube	Sample Label	Sb 217.582 nm ppm	As 228.812 nm ppm	Pb 368.346 nm ppm	Cd 228.802 nm ppm	Sb 231.146 nm ppm	Sb 259.804 nm ppm	Hg 253.652 nm ppm
51:1	Blank	0.00000000 x	0.00000000	0.00000000	0.00000000	0.00000000	0.00000000	0.00000000
51:2	Standard 1	5.00000000 x	5.00000000 x	5.00000000 x	5.00000000 x	5.00000000 x	5.00000000 x	5.00000000 x
1:1	Sample 1 Blank							
1:2	Sample 2 Blank							
1:3	Sample 3 Blank							
1:4	Sample 4 Blank							
1:5	Sample 5 Blank							
1:6	Sample 6 Blank							
1:7	Sample 7 Blank							
1:8	Sample 8 Blank							
1:9	Sample 9 Blank							
1:10	Sample 10 Blank							

Replicate	Concentration	Intensity
1	5.00000000 x	147638.8081990
2	5.00000000 x	150627.4222504
3	5.00000000 x	149264.2280479
4	5.00000000 x	149789.0037702
5	5.00000000 x	149478.1343917
6	5.00000000 x	145676.7305612
7	5.00000000 x	146140.8856437
8	5.00000000 x	137623.2013332
9	5.00000000 x	137583.1623369
10	5.00000000 x	136929.3558030



# Accessory Options for the 4100 MP-AES


- Automate and simplify analysis with the SPS3 autosampler (required for unattended overnight operation)
- For organic applications, use the EGCM to bleed air into the plasma minimizing C build-up and reducing background
  - also requires the OneNeb inert nebulizer (incl. with the Organics kit)
- To enable low ppb level detection of As, Se or Hg, use the Multimode Sample Introduction System (MSIS)
  - also requires the 5 channel peristaltic pump option



# Applications for the Agilent 4100 MP-AES



**Geochemistry, Mining and Minerals**



**Chemical and Petrochemical**



**Food and Agriculture**



**Environmental**

# Applications MP-AES Can Do – by Market

## GEOCHEMICAL

Geochem samples in aqua regia digests  
Trace elements in geological samples  
Trace level gold in cyanide leach  
Analysis of high purity gold  
Platinum group elements in ore grade material  
Various elements in plating solutions

## CHEMICAL & PETROCHEMICAL

Additives in lubricating oils  
Wear metal contaminants in used oils  
Analysis of coolant  
Analysis of petroleum and diesel fuel  
Major elements in polymers

## FOOD & AGRICULTURE

Major elements in foods, beverages and agricultural samples  
Cations in soils  
Nutrients in soils  
Metals in soil extracts  
Metals in agricultural soil samples

## ENVIRONMENTAL

Hg, Pb, Cd and Cr in electronics and plastics (for WEEE/RoHS compliance)  
Heavy metals in soils  
As, Sb and Se in sediments and waste  
Analysis of waste waters, sediments and soils



# Applications for the Agilent 4100 MP-AES

Geochemical, Metals and  
Mining Applications



# Mining (Geochemistry) Applications

Aqua regia digests, fire assay, cyanide leaches

- Au, Ag, Pt, Pd, Cu, Ni, Zn


## Analytical Challenges

- Low detection limits
- High levels of dissolved solids
- Linearity
- Sample throughput



# Gold and Precious Metals by Fire Assay



Au, Pt and Pb extracted into a silver sphere



# Instrument Parameters

Parameter	Setting
Nebulizer	Glass concentric
Spray chamber	Single-pass glass cyclonic
Read Time	3 s
Replicates	3
Sample uptake delay	10 s
Stabilization time	5 s
Fast Pump	On (80 rpm)
Background correction	Auto
Nebulizer Pressure	140 to 240 kPa (auto optimized)

# Extended Linear Range for Platinum & Palladium



**Both Pt & Pd  
linear to  
120 ppm**

Analyte	Correlation
Pt 265.495	0.9999
Pd 363.470	0.9995

# Accuracy –Platinum & Palladium in Custom CRMs

## Platinum

CRM	Certified value (mg/L)	MP-AES result (mg/L)
6	0.74	0.75
7	35.6	35.9
8	9.0	9.5

## Palladium

CRM	Certified value (mg/L)	MP-AES result (mg/L)
6	3.21	3.4
7	44.4	44.0
8	35.0	36.5

# Superior Performance – Detection Limits

Analyte	Wavelength (nm)	Method DL (ug/L)	Flame AA IDL (ug/L)
Au	267.595	4	11
Pt	265.945	13	76
Pd	363.469	0.7	15

**MDL =  $3\sigma$ , 3s read time, 20 blanks**



# MP-AES Improves Throughput over Flame AA

Ore analysis can be volume limited

Typical sample volume available ~4mL

Automate analysis & reduce sample consumption using the autosampler

- Analysis time only 55s sample to sample
- Sample consumption ~1.8 mL
- Throughput 65 samples/hour, or 520 samples in 8 hours
- Plus...
- Safe unattended, multielement operation




**Agilent SPS 3 Autosampler**



# Applications for the Agilent 4100 MP-AES

Chemical and Petrochemical  
Applications



# Applications MP-AES Can Do – by Market

## GEOCHEMICAL

Geochem samples in aqua regia digests  
Trace elements in geological samples  
Trace level gold in cyanide leach  
Analysis of high purity gold  
Platinum group elements in ore grade material  
Various elements in plating solutions

## CHEMICAL & PETROCHEMICAL

Additives in lubricating oils  
Wear metal contaminants in used oils  
Analysis of coolant  
Analysis of petroleum and diesel fuel  
Major elements in polymers

## FOOD & AGRICULTURE

Major elements in foods, beverages and agricultural samples  
Cations in soils  
Nutrients in soils  
Metals in soil extracts  
Metals in agricultural soil samples

## ENVIRONMENTAL

Hg, Pb, Cd and Cr in electronics and plastics (for WEEE/RoHS compliance)  
Heavy metals in soils  
As, Sb and Se in sediments and waste  
Analysis of waste waters, sediments and soils

# Organic Applications – The Challenge

High vapour pressure from solvents

- Plasma instability
- May extinguish the plasma

Carbon build up in injector

- Poor precision and drift
- Down time - injector may require regular cleaning

Nebulizer blockage

- Poor precision
- Down time - nebulizer needs cleaning

# Organics Kit - External Gas Control Module (EGCM)



1. Prevents carbon build up on the injector when running organics
2. Reduces the background emissions from the plasma
3. A controlled flow of air is bled into the Auxiliary gas flow through the torch

# Determination of Impurities in Gasoline and Petro-Diesel Fuels



# Limits on Target Elements in Diesel Fuel - EN14538



Property	ASTM Method	Limits	Units
Calcium & Magnesium, combined	EN 14538	5 max	ppm (ug/g)
<b>Flash Point (closed cup)</b>	<b>D 93</b>	<b>93 min.</b>	<b>Degrees C</b>
Alcohol Control (One of the following must be met)			
1. Methanol Content	EN14110	0.2 Max	% volume
2. Flash Point	D93	130 Min	Degrees C
<b>Water &amp; Sediment</b>	<b>D 2709</b>	<b>0.05 max.</b>	<b>% vol.</b>
Kinematic Viscosity, 40 C	D 445	1.9 - 6.0	mm <sup>2</sup> /sec.
Sulfated Ash	D 874	0.02 max.	% mass
<b>Sulfur</b>			
<b>S 15 Grade</b>	<b>D 5453</b>	<b>0.0015 max. (15)</b>	<b>% mass (ppm)</b>
<b>S 500 Grade</b>	<b>D 5453</b>	<b>0.05 max. (500)</b>	<b>% mass (ppm)</b>
Copper Strip Corrosion	D 130	No. 3 max.	
Cetane	D 613	47 min.	
<b>Cloud Point</b>	<b>D 2500</b>	<b>Report</b>	<b>Degrees C</b>
Carbon Residue 100% sample	D 4530*	0.05 max.	% mass
<b>Acid Number</b>	<b>D 664</b>	<b>0.50 max.</b>	<b>mg KOH/g</b>
<b>Free Glycerin</b>	<b>D 6584</b>	<b>0.020 max.</b>	<b>% mass</b>
<b>Total Glycerin</b>	<b>D 6584</b>	<b>0.240 max.</b>	<b>% mass</b>
Phosphorus Content	D 4951	0.001 max.	% mass
Distillation, T90 AET	D 1160	300 max.	Degrees C
Sodium/Potassium, combined	EN 14538	5 max	ppm
<b>Oxidation Stability</b>	<b>EN 14112</b>	<b>3 min</b>	<b>hours</b>



**Workmanship** Free of undissolved water, sediment, & suspended matter  
**BOLD = BQ-9000 Critical Specification Testing Once Production Process Under Control**

\* The carbon residue shall be run on the 100% sample.

# Sample Preparation – Method EN14538

A commercial diesel sample was analysed

- 1:10 dilution with Shellsol

Standards made from Conostan S21+K

- 0.5 ppm, 1 ppm, 5 ppm, 10 ppm

All samples and standards matrix matched  
with blank oil



# Instrument Parameters

Parameter	Setting
Nebulizer	Inert OneNeb
Spray chamber	Double-pass glass cyclonic
Sample tubing	Orange/green solvent resistant
Waste tubing	Blue/blue solvent resistant
Read Time	3 s
Replicates	3
Stabilization time	15 s
Fast Pump (80 rpm)	On
Background correction	Auto
Pump speed	5 rpm



# Method Detection Limits

3 $\sigma$ , 3s read time, 10 blanks

Element	Wavelength (nm)	MDL (ppb)
Mg	285.213	2.7
Ca	422.673	8.2
Na	588.995	18.7
K	766.491	2.7

# Spike Recoveries on Diesel Fuel

Diesel fuel sample spiked with 0.55 ppm S21+K

Element	Sample (ppm)	Spiked Sample (ppm)	Recovery (%)
Mg 285.213 nm	< MDL	0.53	97
Ca 422.673 nm	< MDL	0.51	93
Na 588.995 nm	< MDL	0.51	93
K 766.491 nm	< MDL	0.51	93

# Applications for the Agilent 4100 MP-AES

Food and Agricultural  
Applications



# Applications the MP-AES can do – by Market

## GEOCHEMICAL

Geochem samples in aqua regia digests  
Trace elements in geological samples  
Trace level gold in cyanide leach  
Analysis of high purity gold  
Platinum group elements in ore grade material  
Various elements in plating solutions

## CHEMICAL & PETROCHEMICAL

Additives in lubricating oils  
Wear metal contaminants in used oils  
Analysis of coolant  
Analysis of petroleum and diesel fuel  
Major elements in polymers

## FOOD & AGRICULTURE

Major elements in foods, beverages and agricultural samples  
Cations in soils  
Nutrients in soils  
Metals in soil extracts  
Metals in agricultural soil samples

## ENVIRONMENTAL

Hg, Pb, Cd and Cr in electronics and plastics (for WEEE/RoHS compliance)  
Heavy metals in soils  
As, Sb and Se in sediments and waste  
Analysis of waste waters, sediments and soils

# Food Applications – Types of Measurement



## Product Labeling

- Nutrient Levels



## Food Quality

- Contamination
- Pollutants



## Food Production

- QA/QC



## Research

- Developing Better Products



# Which Elements are Monitored in Foods?

## Toxic trace elements

- Low level analysis – typically ppb levels
- As, Cd, Hg, Pb, Tl, Cr

## Common elements that are toxic in excess

- Typically range from ppb to ppm level
- Al, Ni, Cu, Zn, Se, Mo, Sn

## Essential minerals

- Can cover a wide range from ppm to % levels
- Na, Mg, P, S, K, Ca, Fe

# Sample Types

## Wide Range of Matrices

- Dairy, Meat, Fish, Grain, Vegetables, Fruit, etc
- Processed goods
  - Canned, frozen, dried and preserved food
  - Drinks (soda, tea, coffee, alcoholic beverages etc.)
  - Preservatives, “fillers” and additives such as calcium carbonate
- “Nutraceuticals”

## Other Important Samples

- Water (used for irrigation)
- Soil
- Fertilizer
- Food additives / flavors / colors
- Packaging material



# Food and Agriculture Applications – The Challenge

## Robust Sample Introduction

- Large number of samples
- Wide range of matrices
- High total dissolved solids (TDS)

## Many Elements Per Sample

- Potential for spectral interferences
- Labs. need high sample throughput
- Labs. also want low cost of analysis

## Wide Dynamic Range

- Trace (ppb) level to % level
- Avoid sample dilution (if possible)



# Major, Minor and Trace Analysis in Foodstuffs



A range of samples

1. NIES CRM No.7 Tea Leaves
2. NIES CRM No.10c Rice Flour
3. NIST SRM 1577 Bovine Liver:

Plus pre-prepared sample solutions in 4 % HNO<sub>3</sub> (High Purity Standards)

4. CRM-Wheat Flour
5. CRM-Milk Powder
6. CRM-Oyster Tissue

# Instrument Parameters

Parameter	Setting
Nebulizer	Glass concentric
Spray chamber	Single pass glass cyclonic
Sample tubing	White/white
Waste tubing	Blue/blue
Read Time	10 s
Replicates	3
Stabilization time	10 s
Fast Pump	On (80 rpm)
Pump speed	12 rpm
Nebulizer pressure	160 – 180 kpa (auto optimized)
Background Correction	Auto

# Measured Results for CRM Food Samples



## NIES No. 7 Tea Leaves



## High Purity Stds CRM Milk Powder

Element	Measured values	Certified values
	wt%	wt%
Ca	0.314 ± 0.013	0.320 ± 0.012
Mg	0.150 ± 0.004	0.153 ± 0.006
K	1.861 ± 0.074	1.86 ± 0.07
	mg/kg	mg/kg
Ba	5.76 ± 0.57	5.7*
Cd	nd	0.03 ± 0.03
Co	nd	0.12*
Cr	nd	0.15*
Cu	7.13 ± 0.81	7 ± 0.3
Pb	nd	0.8 ± 0.03
Ni	6.03 ± 0.63	6.5 ± 0.3
Sr	3.63 ± 0.43	3.7*
Zn	34 ± 3	33 ± 3

Element	Measured values (mg/kg)	Certified values (mg/kg)
Al	nd	0.020 ± 0.002
Ca	131 ± 9	130 ± 1
Co	nd	0.0004*
Cr	nd	0.0003*
Cu	0.006 ± 0.001	0.007 ± 0.001
Fe	0.018 ± 0.002	0.020 ± 0.001
K	178 ± 6	170 ± 2
P	98.7 ± 1.3	100 ± 1
Pb	nd	0.002*
Mg	11.9 ± 0.2	12 ± 0.1
Mn	0.003 ± 0.002	0.003*
Na	48.7 ± 2.6	50 ± 1
Zn	0.48 ± 0.05	0.50 ± 0.01

\* Reference values only

# Measured Results for CRM Food Samples



## NIES No. 10c Rice Flour

Element	Measured values	Certified values
	wt%	wt%
Mg	0.127 ± 0.006	0.125 ± 0.008
K	0.279 ± 0.012	0.275 ± 0.010
P	0.300 ± 0.010	0.335 ± 0.008
	mg/kg	mg/kg
Al	1.49 ± 0.13	1.5*
Ca	95.4 ± 7.0	95 ± 2
Cd	1.83 ± 0.14	1.82 ± 0.06
Co	nd	0.007*
Cr	nd	0.08*
Cu	4.03 ± 0.32	4.1 ± 0.3
Fe	10.6 ± 0.15	11.4 ± 0.8
Mo	nd	1.6 ± 0.1
Ni	nd	0.30 ± 0.03
Sr	0.2	0.2*
Zn	21.8 ± 1.0	23.1 ± 0.8

## High Purity Stds CRM Wheat Flour

Element	Measured values (mg/kg)	Certified values (mg/kg)
Al	0.83 ± 0.02	0.85 ± 0.01
Ca	9.64 ± 0.97	9.5 ± 0.1
Cd	nd	0.0015*
Co	nd	0.001*
Cr	0.013 ± 0.001	0.014*
Cu	0.09 ± 0.008	0.1 ± 0.002
Fe	0.81 ± 0.04	0.90 ± 0.01
K	62.5 ± 0.5	65 ± 0.7
P	61.1 ± 1.7	65 ± 0.7
Pb	0.05 ± 0.001	0.050 ± 0.003
Mg	20.8 ± 0.1	20.0 ± 0.2
Mn	0.36 ± 0.02	0.4 ± 0.008
Ni	nd	0.009 ± 0.001
Zn	0.47 ± 0.05	0.50 ± 0.01

\* Reference values only

# Measured Results for CRM Food Samples



## NIST 1577 Bovine Liver

## High Purity Stds CRM Oyster Tissue

Element	Measured values	Certified values
	wt%	wt%
Na	0.247 ± 0.006	0.243 ± 0.013
K	1.00 ± 0.08	0.97 ± 0.06
	mg/kg	mg/kg
Ca	131	123*
Cd	nd	0.27 ± 0.04
Co	nd	0.18*
Cu	185 ± 6	193 ± 10
Fe	266 ± 5	270 ± 20
Pb	nd	0.34 ± 0.08
Mg	625 ± 45	605*
Mn	10.4 ± 1.41	10.3 ± 1
Mo	nd	3.2*
Sr	0.15 ± 0.07	0.14*
Zn	125 ± 4	130 ± 10

Element	Measured values (mg/kg)	Certified values (mg/kg)
Al	2.92 ± 0.07	3*
Ca	15.0 ± 0.49	15*
Cd	nd	0.03*
Co	nd	0.004*
Cr	nd	0.007*
Cu	0.56 ± 0.05	0.6*
K	100 ± 0.96	100*
P	79.1 ± 0.9	80*
Pb	nd	0.005*
Mg	12.1 ± 0.2	12*
Mn	0.18 ± 0.01	0.2*
Na	48.9 ± 0.8	50*
Ni	nd	0.01*
Zn	8.3 ± 0.4	9*

\* Reference values only

# Applications the MP-AES can do – by Market

## GEOCHEMICAL

Geochem samples in aqua regia digests  
Trace elements in geological samples  
Trace level gold in cyanide leach  
Analysis of high purity gold  
Platinum group elements in ore grade material  
Various elements in plating solutions

## CHEMICAL & PETROCHEMICAL

Additives in lubricating oils  
Wear metal contaminants in used oils  
Analysis of coolant  
Analysis of petroleum and diesel fuel  
Major elements in polymers

## FOOD & AGRICULTURE

Major elements in foods, beverages and agricultural samples  
Cations in soils  
Nutrients in soils  
Metals in soil extracts  
Metals in agricultural soil samples

## ENVIRONMENTAL

Hg, Pb, Cd and Cr in electronics and plastics (for WEEE/RoHs compliance)  
Heavy metals in soils  
As, Sb and Se in sediments and waste  
Analysis of waste waters, sediments and soils

# Instrument Parameters

Parameter	Setting
Nebulizer	Glass concentric
Spray chamber	Single pass glass cyclonic
Sample tubing	Black/black
Waste tubing	Blue/blue
Read Time	3 s
Replicates	3
Stabilization time	10 s
Fast Pump (80 rpm)	On
Pump speed	15 rpm
Viewing position	-20 to 20 steps (auto optimized)
Nebulizer pressure	120 – 240 kpa (auto optimized)

# Measured Results for NIST SRM 2710 Montana Soil


Analyte	MP-AES Results Average (mg/kg)	NIST Reference Range (mg/kg)
Al	24300 ± 400	12000 - 26000
As	550 ± 20	490 - 600
Cr	21 ± 1	15 - 23
Cu	2800 ± 20	2400 - 3400
Fe	28000 ± 300	22000 - 32000
Mn	8500 ± 200	6200 - 9000
Ni	9.5 ± 0.6	8.8 - 15
Pb	5600 ± 300	4300 - 7000
Zn	6100 ± 200	5200 - 6900



# Measured Results for NIST SRM 2711 Montana Soil

Analyte	MP-AES Results Average (mg/kg)	NIST Reference Range (mg/kg)
Al	20000 ± 200	12000 - 23000
As	90 ± 15	88 - 110
Cr	21 ± 3	15 - 25
Cu	90 ± 1	91 - 110
Fe	23000 ± 2000	17000 - 26000
Mn	600 ± 10	400 - 620
Ni	17 ± 3	14 - 20
Pb	1400 ± 30	930 - 1500
Zn	300 ± 10	290 - 340

# Improving Detection Limits for As, Se, Hg with MSIS




Analyte	MSIS in Simultaneous Mode ug/L
As	1.0
Se	2.0
Hg	0.5

# Application Notes Available



## Geochemical

- New methodology for analysis of gold and precious metals
- Analysis of base metals in geochemical samples



## Chemical & Petrochemical

- Analysis of wear metals and contaminants in engine oils
- Analysis of lubricating oil additives
- Analysis of trace elements in petroleum and diesel fuels



## Food & Agriculture

- Analysis of soil extracts
- Cost-effective analysis of major, minor, and trace elements in foodstuffs



## Environmental

- Determination of metals in industrial wastewaters

# Agilent 4100 MP-AES



## *Runs on air – the most significant advance in atomic spectroscopy*

- **Lowest running cost** of any atomic spectroscopy technique due to capability to run on air – ideal for remote and at site operation
- **Improved safety** – capability to run on air means no flammable gases and no manual handling of cylinders
- **Easy-to-use** software with MP applets and plug and play torch which simplify operation and maximize uptime
- **Superior performance to flame AA**, with capability to run unattended overnight

# Agilent's Atomic Spectroscopy Portfolio

**ICP-OES**



**ICP-MS**




**ICP-QQQ**



**Flame AAS**



**Graphite Furnace AAS**



**4100 MP-AES**





# QUESTIONS?