



Update on Extractable/Leachable Analysis
Challenges and Solutions from Separation, Selectivity
to Identification in Pharma Industry



David A. Weil, Ph.D.
 Senior Applications Scientist
 Agilent Technologies Inc.
 Schaumburg, IL USA

Acknowledgements:

Diana Wong	GC/MS Applications Chemist
Roger Firor	Senior GC/MS Applications Chemist
Anthony Macherone	Senior GC/MS Applications Scientist
Andreas Tei	Pharma Segment Manager (sm mol)
Amir Libra	US SPSD AE Manager
Srividya Kailasam	Pharma Application Chemist (India)

Alan Hendricker	Becton Dickenson, Durham, NC
Douglas E. Kiehl	Eli Lilly, Indianapolis, IN
Dana Reed	3M CRL, St. Paul, MN

Outline

- **Background:**
 - Recalls, Definitions, Regulations, Risk Assessment
- **Challenges of Extractable/Leachable Analysis**
- **Agilent GC/MS Solutions:**
 - GC/MSD and GC/QTOF Example
- **Agilent LCMS Solution: Accurate Mass, MS/MS**
 - Sensitivity with Jet Stream
 - Streamline MS/MS Databases Creation
- **Questions**

Why Worry about Extractable/Leachable

Merck recalls cholesterol drug Liptruzet due to packaging issues: Food and Drug recalls

Some of the outer laminate foil pouches may allow in air and moisture, according to the release, which "could potentially decrease the effectiveness or change the characteristics of the product." The likelihood of either of these outcomes is "remote," the company

Ranbaxy recalling allergy drug in US over defective packaging

(Ref: Business Standard, The Wall Street Journal, livemint.com)
May 2nd, 2014

According to the FDA, Ranbaxy is recalling 29 790 packs of an allergy drug in the US, manufactured by its Ohms Laboratories unit, due to defective packaging. The agency said the voluntary Class II recall, which began in February, follows the discovery of "an unacceptable level of blister defects" affecting the loratadine and pseudoephedrine sulphate extended release tablets.

Johnson and Johnson recalls over 43 Over-the-counter (OTC) children's Medicines manufactured by McNeil Consumers Healthcare

NeoProfen (ibuprofen lysine) Injection: Recall and Shortage - Risk of Particulate Matter

[Posted 08/02/2010]

Nutraloid Labs ejaculoid XXTREME and stimuloid II: Undeclared Drug Ingredient

Coumadin 1 mg Tablet Blister Packs: Recall

FDA Drug Recalls Surges over 836 in 2014!

As FDA data show, the last two years have seen almost as many recalls (2,061) as the previous nine years combined (2,217)—and that's only counting the first seven months of 2014. *ref raps.org August 2014*

Teva Parenteral Medicines Initiates Voluntary Nationwide Recall of Select Lots of Adrucil® (fluorouracil Injection, USP) 5 g/100 mL (50 mg/mL) Due to Particulate Matter

Mylan Initiates Voluntary Nationwide Recall of Select Lots of Injectable Products Due to the Presence of Particulate Matter

Heritage Pharmaceuticals Initiates a Nationwide Voluntary Recall of Colistimethate for Injection USP, 150 mg and Rifampin for Injection USP, 600 mg/vial Due to a Lack of Sterility Assurance

Hospira Issues A Voluntary Nationwide Recall Of One Lot Of Bupivacaine HCl Injection Due To Potential Iron Oxide Particulate In Glass Vials

Aventis Pharma Limited March 2015
Defect with Epilim/Na valproate tablets: some batches have **an odor attributed to Foil Packaging**

Reckitt Benckiser Healthcare June 2014
Fybogel Hi-Fibre and Fybogel Orange potential risk of contamination **with metal particles**

ViroPharma SPRL July 2014
Buccolam Oromucosal Solutions
Recall **Chemical Contamination**

 Agilent Technologies

ASTS Hawaii Meeting April 2015

5

Growing concern about impact of E/L impurities on final API and biologic products

Particulates, Mold, Moisture Stability

- 06/30/2014 Bristol-Myers Squibb Issues Voluntary Nationwide Recall of COUMADIN (Warfarin Sodium) for Injection Due to Presence of Particulate Matter
- 06/17/2014 Hospira Announces Voluntary Nationwide Recall of One Lot of 0.5% Marcaine™ (Bupivacaine HCl Injection, USP), 30 mL, Single-Dose, Preservative-Free Vial Due to Visible Particulates
- 06/02/2014 Alexion Initiates Voluntary Nationwide Recall of Certain Lots of Soliris® (eculizumab) Concentrated Solution for Intravenous Infusion Due to the Presence of Visible Particulate Matter in a Single Lot
- 03/05/2014 Baxter Initiates U.S. Voluntary Recall of One Lot of Peritoneal Dialysis Solution Due to Container-Closure Non-Integrity
- November 01, 2013 Perrigo Recalls Acetaminophen Infant Suspension Liquid, 160 mg/5 mL, Due to a Potential Defect with the Co-packaged Oral Syringe
- 05/16/2014 Labetalol Hydrochloride Injection 100 MG/20 ML (5MG/ML), 20 ML, Multidose Vial by Hospira: Recall - Visible Particulates
- 05/15/2014 Dobutamine Injection (250mg/20mL)/Hospira: Recall - Visible Particulates
- 05/14/2014 Hospira Announces Voluntary Nationwide Recall Of One Lot Of Dobutamine Injection, USP, 250 MG, 20 ML, Single-Dose Flip-top Vial, Due To Visible Particulates
- 04/21/2014 Hospira Announces Voluntary Nationwide Recall of One Lot of 0.25% Marcaine™ (Bupivacaine HCl Injection, USP), 10 ml, Single-Dose, Preservative-Free Vial Due to Visible Particulates
- 04/18/2014 Hospira Announces Voluntary Nationwide Recall of One Lot of 1% Lidocaine HCl Injection, USP, Due To Visible Particulates
- 04/18/2014 Cubist Pharmaceuticals Issues Voluntary Nationwide Recall of One Lot of CUBICIN (daptomycin for injection) 500 mg in 10 mL single use vials Due to Presence of Particulate Matter
- 04/17/2014 Hospira Announces Nationwide Voluntary Recall of Seven Lots of Propofol Injectable Emulsion, USP, Due to Visible Particulates

 Agilent Technologies

BioPharma Road Show July 2014

6



Agilent Technologies

Federal Regulations on Extractable/Leachable



U.S. Food and Drug Administration
Protecting and Promoting *Your* Health

FDA Cosmetic Act

Section 501 (a)(3) of the ACT states t a drug is deemed to be adulterated “if it container is composed, in whole or in part, of any poisonous or deleterious substances that may render the contents injurious to health. ..



Guidance for Industry

Container Closure Systems for Packaging
Human Drugs and Biologics

May 1999

“Drug product containers and closures shall **not be reactive, additive, or absorptive** so as to **later the safety, identity, strength, quality or purity of the drug** beyond the official or establishment requirements”



Agilent Technologies

ASTS Hawaii Meeting April 2015

7



European Medicines Agency
Inspections

Dec 2005

Guideline on Immediate Plastic Packaging Materials

- **Excludes** elastomers, natural and synthetic rubbers
- No Data needed if materials previously approved by EP or member state
- Sec 4: Extraction solvent used should have the same propensity to extract substances as the active substance dosage.
- Sec 5.1: Migration Studies Required: with specific objectives
- Sec 6: Toxicological Information and Documents



Agilent Technologies

ASTS Hawaii Meeting April 2015

8



Agilent Technologies

Pharma Industry Expert Working Groups



PQRI (Product Quality Research Institute) is a working group established to develop **regulatory guidance** for Extractable/Leachable analysis, which is also recognized by the FDA




Guidelines and Assessments for E/L in Construct

USP (United States Pharmacopeia)


- USP<1663> Assessment of Extractables Associated with Pharmaceutical Packaging/Delivery Systems **Revised**
- USP<1664> Assessment of Drug Product Leachables Associated with Pharmaceutical Packaging/Delivery Systems **Revised**
- **EP (European Pharmacopeia)**
 - EP 3.1 Materials used for the manufacture of containers
 - EP 3.2 Containers for pharmaceutical use
- **ICH (International Conference of Harmonization)**
 - ICH Q1A, Q3A, Q3B, Q6A, Q7B Specifications: Test procedures and acceptance criteria for new drug substances and new drug products: chemical substance
- **ISO (International Organization of Standardization)**
 - ISO10993 Biological Evaluation of Medical Devices
- **PQRI (Product Quality Research Institute)**

New USP Chapters just released



Ref: Denise R. Jenke, Daniel L. Norwood, and Desmond G Hunt

Definitions: Extractable, Leachable, Migrants



Extractables

Chemical compounds that can be extracted out of packaging component

- Analyze **packaging component** at
 - High-temperatures:** to obtain the worst case leachable profile
 - Solvent extraction:** polar and non-polar solvent to mimic similar properties as drug product

Leachable

- Chemical compounds from packaging component that leach into the drug product
- Analyze **drug product** at
 - Normal conditions
 - Stimulate extended storage conditions

Migrants

- Crossed a physical barrier primary packaging from secondary and tertiary packaging, accumulating in the drug product

Extractable and Leachable Applications

- **Pharma and BioPharma (Small Molecule/Biologics)**


Drug Delivery Technologies



- Single Dose Prefilled Syringes
- Medical Device Products
- Drug Packaging (Constructs, Films, Closures)
- Orally inhaled and nasal drug products (OINDPs)
- **Food Packaging**
 - Similar Materials and Concerns (BPA, Phthalates, Colorants, Dyes)

Plastics Industry – Sources of Extractables

Sources of extractables are plastic and elastomeric components (plasticizers, etc), ink and adhesives (label), and degradation products (processing, storage, and sterilization)



Cindy Zweiben, Pfizer, Inc., Characterization of Extractables and Leachable in Parenteral Drug Products

Common Sources of Extractable and Leachable

Plastic Components:	Polymers Additives, Phthalates, Lubricants, Fatty Acids, Nitrosamines
Elastomers/Rubber:	Vulcanizing Agent, PAHs, Accelerators Antioxidants, Carbon Black,
Inks/Labels:	Azo Dyes, Aromatic Amines
Adhesives:	Antioxidants (AO), Catalyst Residues, Heavy Metals, Silicones, Surfactants
Pigments:	Inorganic (TiO ₂ , FeO's), Organic Pigments
Cyclic Oligomers:	Polybutylene terephthalate (PBT) polyester Silicone Cyclic Oligomers

Risk categories of pharmaceutical packaging

Degree of concern associated with Route of Administration	Likelihood of interaction between packaging component and dosage form		
	High	Medium	Low
Highest	Inhalation aerosols and solution Injections and injectable suspensions	Sterile powders Injection powders Inhalation powders	
High	Ophthalmic solutions and suspensions Transdermal ointments and patches Nasal aerosols and sprays		
Low	Topical solutions and suspensions Topical and lingual aerosols Oral solutions and suspensions	Topical powders Oral powders	Oral tablets Oral hard capsules Oral soft gelatin capsules

Adapted from Guidance for Industry; *Container Closure Systems for Packaging Human Drug and Biologics*, US Department of Health and Human Services, Food and Drug Administration, Rockville, MD, May 1999

- Interested in **high-risk categories** packaging components
- What contributes to the high-risk in pharmaceutical packaging?
 - Prefilled syringe containing an injectable drug suspension
 - Interacts with multiple components in the packaging material (plastic barrel, rubber plunger, metal needle) with direct delivery to the bloodstream

Biologics Typically Differ from Traditional Small Molecule Drugs for the following Reasons

- Large MW and complex structure (glycosylated, multiple domains)
- Abundance of both hydrophilic and hydrophobic sites may render biologic drugs more efficient extraction media
- Could be heterogeneous mixtures
- Produced by a complex manufacturing process employing living organisms

Ref: Ingrid Markovic, CBER Presentation USP/PQRI E/L Workshop April 2014


Special Concerns about E/L Effect on Biologics

- **Leachable may affect protein products in the following ways**
 - Oxidation
 - Unfolding
 - Aggregation
 - Increase in particulates
 - Formation of clipped variants
 - Formation of Protein Adducts
 - Post translational events during fermentation (glycosylation)
 - Altered protein translation

Examples from Ref: Ingrid Markovic, CBER April 2014

Extractable Leachable Analytical Workflow

Objective: To detect low level quantities organic or inorganic impurities that result from interaction of API or Biologic with Drug Delivery



Simulated Extraction Chart

- Vigorous/exaggerated/-exhaustive conditions
- no material deformation
- „Hot“ extraction techniques but no sample dissolving solvents
- Solvents should be attributed to the expected universe of substances (cover wide range of polarity)
- Solvents should mimic drug product formulation

	Thermal	n-Hexane	Iso-propanol	Isopropanol/Water	Aqueous pH 2.5	Aqueous pH 9.5
Headspace	X					
Reflux	---	X	X	PC/PVC only	---	---
Soxhlet	---	X	X	---	---	---
Sealed Vessel	---	---	---	55°C/3d	(121°C/1hr) ¹	(121°C/1hr) ¹
Sonication	---	---	---	---	X	X

¹: autoclave conditions: (121°C/1hr)

Ref: T.Egert, A. Hendrick, C.Houston, 2011 Talk at Best Practices fro PODP



Analytical Techniques for Impurity Analysis

ICH Identification Limit	Genotoxic Impurities Typical Limits		
	0.01% 100 ppm	0.001% 10 ppm	0.0001% 1 ppm
NMR			
HPLC-UV			
GC-FID			
LC-MS			
GC-MS			
ICP-MS			

Hyphenated technique for better separation →

New Heavy Metal ICH and USP Methods/Regulations



Proposed new ICH and USP methods for elemental impurities: The application of ICP-MS and ICP-OES for pharmaceutical analysis

White paper

Authors
Amir Liba, Ed McCurdy and
Ross Ashdown
Agilent Technologies



Abstract

The United States Pharmacopoeial Convention (USP), in parallel with the International Conference on Harmonization (ICH), is developing new methods for inorganic impurities in pharmaceuticals and their ingredients. The current USP method, <231> "heavy metals limit test", is acknowledged to be inadequate and is due to be replaced with new General Chapters USP-232 (Limits) and <232> (Procedures) in December 2015. The new methods will address the limitations of the current method, extending the list of analytes, reducing maximum permitted exposure limits and taking account of the route of exposure. The new methods will also introduce the use of closed vessel sample digestion and modern instrumental techniques to ensure the accurate recovery and determination of individual analyte concentrations. This White Paper discusses the development of the new USP General Chapters and the ICH Guideline for Elemental Impurities (Q30) and how Agilent's 7700x ICP-MS and 5100x ICP-OES address the requirements of the proposed new methods.



pub 5990-9382EN 2014



Validating the Agilent 7700x ICP-MS for the determination of elemental impurities in pharmaceutical ingredients according to draft USP general chapters <232> / <233>

Application note

Pharmaceutical

Authors
Samira Hussain
Estiva
USA
Amir Liba and Ed McCurdy
Agilent Technologies
USA



Abstract

The United States Pharmacopoeia (USP) is developing new General Chapters relating to the determination of elemental impurities in pharmaceutical products and ingredients. USP-232 defines the analyte limits, while USP-233 defines sample preparation options including closed vessel microwave digestion, and recommends the use of modern instrumentation, such as multi-element ICP-MS and ICP-OES techniques. Analytical equipment qualification under USP-233 is based on performance testing, and includes requirements to demonstrate accuracy, repeatability, and the unequivocal identification of analytes. In this paper we present data to illustrate the successful validation of the Agilent 7700x ICP-MS for the measurement of elemental impurities in gelatine capsule samples, according to USP-232 / <233>.

pub 5990-9365EN, 2011



Inorganic Extractable Impurities

Inorganic Impurities can result from the manufacturing process and are generally known, identified and include:

- Reagents, Ligands, Catalysts
- Heavy Metals
- Inorganic Salts
- Manufacturing Aids (charcoal, filters)



Agilent 7700 Series ICP-MS

Determination of Toxic Elements in Traditional Chinese Medicine Using Inductively Coupled Plasma Mass Spectrometry

This Application Note describes a method for the analysis of Be, Cr, Mn, Ni, Cu, Zn, As, Ag, Cd, Ba, Hg, Tl and Pb in Traditional Chinese Medicines using Inductively Coupled Plasma Mass Spectrometry (ICP-MS). The samples were dissolved by microwave digestion. To validate the method, two certified reference materials were digested and measured. The method detection limits for all target elements were between 0.1-7.2 ng·g⁻¹.

Instrumentation for Extractable Leachable Analysis



**7010
GC/QQQ**



**5975E
GC/MS**



**7200
GC/QTOF**



SFC/MSD/QTOF

**Hi-DEF Q-TOF
6500 series**



**QQQ
6400 Series**



**TOF
6200 series**

Agilent GC and GC/MS



5977A with 7890B GC



7890B GC



5975T LTM SQ



Purge and Trap



7010 QQQ



7200 Q-TOF

 Agilent Technologies

GCMSD New Application Notes!

5991-5616EN

Analysis of Extractable/Leachable Compounds From Plastic Intravenous Bag Sets Using GC/MSD Systems

Application Note

Pharmaceuticals

Abstract

Two Agilent 5977A Series GC/MSD Systems were used for the analysis of extractable and leachable compounds in plastic IV bag sets. Two types of IV bags were investigated: 150 mL electronic bag (topring) and 1 L sodium chloride bag (w/armed). Potentially toxic additives, such as phthalate plasticizers, were shown to have migrated from the IV bag to its infusion solution using the complementation of headspace sampling and liquid injection techniques. High temperature analysis was accomplished using the 7897A Headspace and a 7890A GC coupled with a 5977A MSD. Solvent extracts were analyzed using the MMI 7890A GC coupled with a 5977A MSD. Single ion monitoring (SIM) was used to confirm compound migration.

 Agilent Technologies

5991-5605EN

Analysis of Extractable/Leachable Compounds From Transdermal Patches Using GC/MSD Systems

Application Note

Pharmaceutical

Abstract

A silicone adhesive patch and film release liner were used to investigate extractable and leachable compounds in transdermal drug delivery systems using two Agilent 5977A Series GC/MSD Systems. Plastic and adhesive additives were identified in acetone, dichloromethane, and hexane extracts using the large volume liquid injection technique. Pharmaceutical ingredients were also identified using high temperature headspace and liquid sampling techniques.

Introduction

Particular interest has been given to extraction techniques in container closure systems (CCS) used in the pharmaceutical industry. Regulators have become increasingly aware of the need to understand whether chemical species can be extracted from the primary packaging material (package with direct contact to the drug product), as well as whether the extracted species (from the package) will appear as leachable species in the drug product. Extractables analysis involves extracting compound from the packaging material using elevated temperatures and solvents related to the packaging composition. Leachables analysis involves identifying compounds in the drug formulation that may have leached from the primary packaging material.

The major source of extractables and leachables are additives that provide physical and protective properties to packaging material, such as flexibility, rigidity, stability, and barrier. Extractables include plastic and elastomeric components, inks and adhesives from coating, and degradation products during processing, storage, and sterilization. Leachables are usually a subset of extractables, however new compounds can form from the interaction between drugs and packaging material.

 Agilent Technologies

5991-5632EN

Analysis of Extractable/Leachable Compounds from Generic Liquid Drug Formulations Using GC/MSD Systems

Application Note

Abstract

Pharmaceutical liquid formulations are commonly stored in plastic containers at all risk categories. A pharmaceutical suspension was used as model for investigating compound migration from packaging material. Two Agilent 5977A Series GC/MSD Systems were used. Fatty acid plasticizers were identified using the 7897A Headspace Sampler and a 7890A GC coupled with a 5977A MSD. Phthalate plasticizers were found using the MMI 7890A GC coupled with a 5977A MSD. Single ion monitoring (SIM) confirmed the identification of these plasticizers.

 Agilent Technologies


 Agilent Technologies

Confidentiality Label
May 18, 2015
26

Extractable/Leachable from Plastic IV Bag Sets

Samples Investigated

150-mL dextrose IV bag (expired 8 years)
 1-L NaCl IV bag (2 months at 100°F)
 Kolliphor + NaCl solution
 IV tubing (expired 3 years)



Background

- PVC IV bags are a source of DEHP, which can leach into drug solutions
- DEHP extracts faster in non-polar solvents or when liquid is warmed or agitated
- Kolliphor EL is a nonionic polyethoxylated castor oil solubilizer/emulsifier that can enhance leaching of DEHP from PVC bags


Method

- Headspace temperatures: 85-275°C
- DCM and hexane extraction by sonication


Summary

- DEHP and benzophenone was identified in IV tubing using HS
- DEHP, BHT (antioxidant), benzothiazole (rubber), acetophenone (ink, coating, adhesives) were identified in 150-mL IV bag by HS
- Phthalates, BHT, acetophenone, and Metilox (plasticizer) were identified in 1-L IV bag using HS
- DEHA and DEHP migration was observed in DCM and hexane extracts of dextrose IV bags using MMI and SIM analysis
- Phthalates and DEHA were observed in DCM and hexane extracts of NaCl and Kolliphor-NaCl IV bags

IV tubing using the headspace GC/MS at 250 °C



Compound Migration due to Extended Storage Conditions




SIM analysis of DCM extracts from IV Bag and dextrose solution to confirm the identification and migration of DEHA and DEHP

- Hexadecane (1)
- Butylated Hydroxytoluene (BHT) (2) –antioxidant
- Kodaflex TXIB (3) –plasticizer
- Benzophenone (4) – UV stabilizer
- Isobutyl nonyl phthalate (5)
- 7,9-Di-tert-butyl-1-oxaspiro(4,5)deca-6,9-diene-2,8-dione (6)
- 2-Mercaptobenzothiazole (7) - Rubber**
- Palmitic acid (8) – plasticizer
- Isopropyl palmitate (9)
- Palmitic acid, butyl ester (10)
- 2-Ethylhexyl trans-4-methoxycinnamate (11)
- Benzyl butyl phthalate (12)
- DEHA (13) –plasticizer**
- di(oct-3-yl) phthalate (14)**
- DEHP (15) –plasticizer**
- Irgafos 168 (16) –antioxidant **
- Tris(2,4-di-tert-butylphenyl) phosphate (17)

5991-5616EN

Confidentiality Label
May 18, 2015
29


Compound migration due to extended heating & Kolliphor EL



MMI 5977 GC/MS of DCM extracts of NaCl and Kolliphor-NaCl IV bags

Kolliphor enhanced leaching of Kodaflex TXIB plasticizer, which may suggest that alternative plasticizers was utilized

At 27 min, migration can be bis(2-methylpropyl) phthalate or 6-ethyl-3-octyl butyl phthalate. At 28.5 min, migration can be dibutyl phthalate or butyl cyclobutyl phthalate at 28.5min




Kodaflex txib: 2,2,4-Trimethyl-1,3-pentanediol diisobutyrate
Metilox: benzenepropanoic acid, 3,5-bis(1,1-dimethylethyl)-4-hydroxy-, methyl ester

Confidentiality Label
May 18, 2015
30

Extractable/Leachable From Transdermal Patch

Schematic Diagram of Transdermal Drug Delivery System



Samples Investigated

- Adhesive material contains 5% lidocaine applied to a nonwoven polyester (PE) felt backing
- Film-release liner is made of polyethylene terephthalate (PET)

Method

- Headspace temperatures: 100-275°C
- Acetone, DCM, and hexane extraction by sonication

Summary

- **Identified** (active and inactive) **pharmaceutical** ingredients
- **Identified potentially toxic plasticizers**
- **Acetone** extraction: terephthalates, phthalates, fatty acids
- **Hexane** extraction: phthalates and fatty acid
- **DCM** extraction: Diethylhexylphthalate and benzophenone
- **Headspace:** Fatty acid plasticizers


Reference

- 5991-5605EN Application Note consisted of:
 - TIC and table listing extractables identified by HS and MMI
 - Summary Table of extractables organized by sampling techniques and origin
 - Detailed description of system parameters

5991-5605EN

Confidentiality Label
May 18, 2015
31

Analysis of Lidocaine Patch using MMI 5977 GC/MSD by DCM extraction




TIC is zoomed-in to emphasize low-level leachables

Active/Inactive ingredients produce distorted peaks due to high concentration

Phthalates and plasticizers may originate from the PE material in patch or from the PET in film

5991-5605EN

Lidocaine Patch using Headspace 5977 GC/MSD 250°C



Propylene glycol and glycerin produces distorted peaks due to high concentration of the inactive ingredient

Lidocaine produce broad peaks due to high concentration of an active ingredient

Agilent Technologies

New 7200B GC/Quadrupole Time-of-Flight (QTOF)

- Both EI and CI sources are Standard
- Backflush ready GC
- Resolving Power: **>12.5K at m/z 272**
13K to 15K typical
- Mass Accuracy: **<3 ppm** <2 ppm typical
- MS Sensitivity **<250 fg of OFN derived statistically**
- Dynamic Range: **> 3 orders of magnitude**
 - 4 orders with Dual Gain
- Extended Mass range: **m/z 20-3000**
- MS/MS Mass range: **m/z 20-1050**
- Spectral Rate: **1-50 Spectra/sec**



Agilent Technologies

34

Example: Fluorotelomer Alcohols in Biosolids

Source unknown

- Possible intermediate degradation products of fluorinated polymers
- Oxidize to form fluorinated carboxylic acid, some of which are **toxic**
- Structure: $F_3C(CF_2)_{N-1}(CH_2)_M OH$ N:M FTOH is the shorthand notation

Positive Chemical Ionization (PCI) 20% Methane

Acronym	Formula	Exact Mass + H	Observed Mass	Δ ppm
4:2 FTOH	C ₆ H ₅ F ₉ O	265.0269	265.0270	-0.38
6:2 FTOH	C ₈ H ₅ F ₁₃ O	365.0206	365.0206	0
8:2 FTOH	C ₁₀ H ₅ F ₁₇ O	465.0142	465.0140	0.43
10:2 FTOH	C ₁₂ H ₅ F ₂₁ O	565.0078	565.0078	0
7:2 sFTOH	C ₉ H ₅ F ₁₅ O	415.0174	415.0190	-3.85
5:1 FTOH	C ₆ H ₃ F ₁₁ O	301.0081	301.0079	0.66
6:1 FTOH	C ₇ H ₃ F ₁₃ O	351.0049	351.0050	-0.28
7:1 FTOH	C ₈ H ₃ F ₁₅ O	401.0017	401.0016	0.25
8:1 FTOH	C ₉ H ₃ F ₁₇ O	450.9985	450.9985	0
9:1 FTOH	C ₁₀ H ₃ F ₁₉ O	500.9953	500.9956	-0.60
10:1 FTOH	C ₁₁ H ₃ F ₂₁ O	550.9921	550.9922	-0.18
11:1 FTOH	C ₁₂ H ₃ F ₂₃ O	600.9889	600.9896	-1.16
MeFOSE	C ₁₁ H ₆ F ₁₇ NO ₃ S	558.0026	558.0042	-2.87
EtFOSE	C ₁₂ H ₁₀ F ₁₇ NO ₃ S	572.0183	572.0167	2.80

Targeted Quantitation vs Non-Target Screening

Targeted Quantitation

LC/QQQ

Benefits:

- Most Sensitive Technique
- Good Quantitative Results

Drawbacks:

- Limited number of analytes
- Non-target compounds undetected

Non-Target Screening

LC/TOF or LC/Q-TOF

Analysis Options:


- MS only screening
- MS/MS with accurate mass library for confident identification

Benefits:


- Detect, Quantify, and Identify hundreds or thousands of compounds in a single injection
- Reanalyze old sample data any time *without reinjection* upon discovery of new compounds

MS/MS using ALL Ions Approach

Fluorotelomer Neutral Loss Mechanism Inferred PCI, methane




11:1 FTOH Spiked in Biosolid Extract PCI, methane



The LCMS Match Game....

Choosing The Mass Spectrometer for the Task
Sensitivity, Selectivity, Identification..




Agilent Technologies

Impurity Profiling Dec 2013

Page 39


The Relationship of Mass Accuracy to the Number of Possible Molecular Formulas by Mass



Agilent Technologies

Polymer Characterization April 2013


Q(TOF) Innovations from Agilent



 Agilent Technologies

Polymer Characterization April 2013

How Does Agilent Jet Stream Technology Enhance ESI?




- Higher temperatures increase evaporation rate of mobile phase so formation of aerosol and evaporation of droplets will occur more quickly. Fewer large droplets are left in the spray.
- Concentric orientation of sheath gas with respect to spray helps reduce ion dispersion at normal LC flow rate.
- Introduces many more ions into the MS.
- Reduces number of neutral solvent clusters.

Result: stronger signals with lower relative standard deviation (RSDs) at the limit of detection.

 Agilent Technologies


Polymer Characterization April 2013

Irganox 1010: Jet Stream versus Dual ESI > 10X Enhancement in Ionization



Agilent Technologies

Irganox 1024 Accurate Mass Measurement





Agilent Technologies

Polymer Characterization April 2013

Molecular Feature Extraction (MFE)

Automated Data Reduction Software

Finds Features in TOF/QTOF Data



Data Reduced sum intensities of isotopes, adducts, clusters and multiply charges ions together.




species	RT	m/z	mass	abund.
M	8.162		342.1467	130843
M+H	8.165	343.1547	342.1474	11989
M+H+1	8.162	344.1581		2290
M+H+2	8.157	345.1748		363
M+H4N	8.164	360.1807	342.1469	8420
M+H4N+	8.156	361.1893		1227
M+Na	8.162	365.1359	342.1466	75678
M+Na+1	8.163	366.1394		15324
M+Na+2	8.162	367.1429		1901
2M+Na	8.164	707.2810	342.1459	4629
2M+Na+1	8.162	708.2850		1808
2M+Na+2	8.173	709.2895		336

Agilent Technologies

Polymer Characterization April 2013

Impurity ID Workflow: Untargeted Approach



Agilent Technologies

Impurity Profiling Dec 2013
Page 46



Agilent Technologies

Personal Compound Databases (PCDL)


Polymer Additive Database contains over 1550 Compounds Building MS/MS Searchable Library

Compound Name	Formula	Mass	Anion	Cation	RT (min)	CAS	ChemSpider	IUPAC Name	Num Spectra
Bis 2-ethyl hexyl adipate	C22H42O4	370.30831	<input type="checkbox"/>	<input type="checkbox"/>					0
Bis ethyl hexyl phthalate	C24H38O4	390.27701	<input type="checkbox"/>	<input type="checkbox"/>					0
Bis Phenol A	C15H16O2	228.11503	<input type="checkbox"/>	<input type="checkbox"/>					0
BIS(2-ETHYLHEXYL) 2-ETHYLHEXYLPHOSPH...	C24H51O...	418.35758	<input type="checkbox"/>	<input type="checkbox"/>					0
BIS(3,5-DI-T-BUTYL-4-HYDROXYBENZYL) SUL...	C30H48O...	470.32185	<input type="checkbox"/>	<input type="checkbox"/>	96.69.5			BIS(3,5-DI-T-BUTYL-4-HYDROXYBENZYL) SUL...	0
Bis(5-benzoyl-4-hydroxy-2-methoxyphenyl)methane	C29H26O6	466.19329	<input type="checkbox"/>	<input type="checkbox"/>				Bis(5-benzoyl-4-hydroxy-2-methoxyphenyl)methane	0
BIS(TRIDECYL) PHTHALATE	C34H68O4	530.43351	<input type="checkbox"/>	<input type="checkbox"/>					0
Bis 2-ethylbutyl Phthalate	C20H30O4	334.21441	<input type="checkbox"/>	<input type="checkbox"/>					0
BUTOXYETHYL OLEATE	C24H46O3	382.34470	<input type="checkbox"/>	<input type="checkbox"/>					0
Butyl Benzene Phthalate	C19H20O4	312.13616	<input type="checkbox"/>	<input type="checkbox"/>					0
Butyl hydroxyanisole	C11H16O2	180.11503	<input type="checkbox"/>	<input type="checkbox"/>	121.00.6	8145		Butyl hydroxyanisole	0
BUTYL OLEATE	C22H42O2	338.31848	<input type="checkbox"/>	<input type="checkbox"/>					0

Accurate Mass and MS/MS Library Searchable Databases

Analysis of IPA Reflux of Al Foil Closure


IPA Extract from Foil – MFE APCI (+)



Compounds Sorted by Abundance

Cpd	RT	Name	Formula	Score	Mass	Score (MFG)	DfH (MFG, ppm)	Height	Vol	ID Source	ID Techniques Applied
9	8.38		C20H32O8	39.16	400.20937	78.31	0.86	15545	69788	MFG	DBSearch-MFG
6	7	Cyclic TEDDEG	C24H24O10	87.98	472.12697	76.88	2.06	9111	53389	DBSearch-MFG	DBSearch-MFG
16	9.87		C23H38O8	97.71	442.25688	97.71	-0.47	13733	52129	MFG	MFG
21	10.21		C24H40O8	99.66	456.27236	99.66	-0.09	12219	44321	MFG	MFG
7	7.25	Cyclic TEDDEG	C24H24O10	49.63	472.13691			7245	43422	DBSearch	DBSearch-MFG
22	10.23	Proanthocyanidin A2	C30H24O12	98.05	576.12632	97.07	0.79	6870	36652	DBSearch-MFG	DBSearch-MFG
11	9.14		C29H36N2O8	99.41	540.24707	99.41	0.17	8423	32658	MFG	MFG
28	12.97	nDOP; Chissocizer nDOP; N...	C24H38O4	49.8	390.27698			9417	29114	DBSearch	DBSearch-MFG
20	10.13		C30H48O12	98.24	600.31369	98.24	1.48	7319	26720	MFG	MFG
23	11.47	Funghichromin	C35H58O12	96.21	670.39198	95.22	1.89	6324	26571	DBSearch-MFG	DBSearch-MFG
14	8.49		C22H38O8	49.75	428.24003	99.49	0.44	4146	13799	MFG	MFG
10	8.45	Dimethyl sebacate	C12H22O4	81.92	230.15187	77.93	-0.27	3612	12281	DBSearch-MFG	DBSearch-MFG
2	1.24	1-Heptene	C7H14	85.24	98.10912			2059	11617	DBSearch	DBSearch
30	13.08	1β-Hydroxylamin D3 / 1β-H...	C27H44O2	99.25	400.3343	99.45	-0.42	4010	11452	DBSearch-MFG	DBSearch-MFG
35	18.25	DGJP-14(0/161(S2))	C35H66O4	97.74	550.49534	98.11	1.4	2473	10693	DBSearch-MFG	DBSearch-MFG
8	7.61	N,N'-Dicyclohexylurea	C13H24N2O	86.26	224.18964	86.26	0.98	1629	10464	DBSearch-MFG	DBSearch-MFG
13	9.46	2H-1-Benzopyran-3,6-dicarb...	C15H22O2	84.54	234.16187	85.17	0.49	2645	9443	DBSearch-MFG	DBSearch-MFG
33	15.14	3-Oxoxylamin D3	C27H44	90.14	368.34432	90.14	-0.05	1716	8464	DBSearch-MFG	DBSearch-MFG
29	13.01	Diethyl adipate	C22H42O4	80.66	370.30809	82.87	0.6	2403	8168	DBSearch-MFG	DBSearch-MFG
18	10.1	Oxidized Laticin	C13H22O	87.27	194.16708	87.22	-0.06	2185	7919	DBSearch-MFG	DBSearch-MFG
31	13.24	Diamide L-200; Neutron-22...	C22H43N O	83.97	337.33437	84.5	0.28	2366	7886	DBSearch-MFG	DBSearch-MFG
26	11.94	1,2-Benzenedicarboxylic aci...	C20H30O4	84.23	334.21429	84.47	0.34	1774	7061	DBSearch-MFG	DBSearch-MFG
37	19.17		C37 H70 O4	93.8	578.50508	93.8	2.3	1334	6970	MFG	MFG
36	18.95	Squalene	C30H50	86.12	410.39047			1385	6790	DBSearch	DBSearch
15	9.6	Phenol, 2,4-dimethyl-6(2-m...	C15H22O	84.01	218.16725	85.57	-0.83	1675	6173	DBSearch-MFG	DBSearch-MFG
34	17.18				522.46369			1076	5632		
17	10.06	MEHP	C16H22O4	85.19	278.15206	85.88	-0.91	1922	5394	DBSearch-MFG	DBSearch-MFG
4	4.67				158.10935			1113	5241		
27	12.43	Sansocizer DHP; Chissocize...	C22H34O4	81.17	362.24515	81.67	1.53	1233	5143	DBSearch-MFG	DBSearch-MFG
19	10.1		C13H20	84.42	176.15638	84.42	0.67	1538	4839	MFG	MFG
12	9.46	Irganox degradate	C17H24O3	83.24	276.17256	83.68	-0.05	1793	4362	DBSearch-MFG	DBSearch-MFG
5	5.79				114.10354			1256	4248		
32	13.54	previtamin D3 / precholecal...	C27H44O	86.93	384.33629	95.17	2.43	1172	3635	DBSearch-MFG	DBSearch-MFG
24	11.84	4-oxo-9C,11Z,13E,15E-octad...	C18H26O3	84.86	290.18797	85.53	0.76	1574	3951	DBSearch-MFG	DBSearch-MFG
25	11.84	Photoinitiator 32; Phenylgly...	C10H10O3	84.74	178.06303	85.46	-0.19	1445	3024	DBSearch-MFG	DBSearch-MFG


Major Compound C₂₄H₂₄O₁₀ 0.19 ppm Error



Agilent Technologies

Extractable and Leachable Examples
Page 51


Isotope Pattern Matching



Agilent Technologies

Extractable and Leachable Examples
Page 52


APCI(+) Conc. Foil Extract AutoSMS Results



Agilent Technologies

Extractable and Leachable Examples
Page 53


MS/MS of [M+H]⁺ C₂₄H₂₄O₁₀ RT 6.9 minutes




Agilent Technologies

Extractable and Leachable Examples
Page 54

MS/MS Data $[M+NH_4]^+$ peak for $C_{24}H_{24}O_{10}$ at RT 7.0



MS/MS of $[M+H]^+$ $C_{24}H_{24}O_{10}$ at RT 7.27 mins



Molecular Structure Correlator of $[M+H]^+ C_{24}H_{24}O_{10}$

The screenshot displays the Agilent Molecular Structure Correlator software interface. The main window shows search results for the compound formula $C_{24}H_{24}O_{10}$. The interface is divided into several sections:


- Compound formula:** $C_{24}H_{24}O_{10}$
- Structure Search:** Parameters: 1/1, Add Structure, Set by Score, Set.
- Standard IChEM:** BKPYRILGRVZCR, UHFFFDYTSAN, Score: 14.85, MSC, Save, Delete.
- Fragments of structure #1 -- elucidated: 87.5% ions, 18 of 18 ions:**

Mass	Intensity	Weight%	No. of cand.	Best score	
1	73.0265	74.53	0.0	4	80.6
2	89.0587	334.04	0.0	2	92.9
3	148.0227	767.04	0.4	2	96.2
4	193.0470	776.90	2.3	4	79.3
5	237.0751	306.26	2.3	5	95.6
6	238.0837	113.00	1.1	0	98.0
7	281.1049	492.34	11.8	5	76.2
8	473.1906	137.72	81.8	0	93.0
- Fragment formulas for $C_{24}H_{24}O_{10}$:**

m/z	Intensity	nom. Intensity	formula	d(Mppm)
1	73.0265	74.53	C3H6O2	26.8
2	89.0587	334.04	C4H8O2	11.4
3	148.0227	767.04	C9H16O3	4.5
4	148.0227	767.04	C11H18O3	22.9
5	193.0470	776.90	C13H20O4	13.0
6	237.0751	306.26	C12H18O5	2.7
7	237.0751	306.26	C15H18O2	14.0
8	238.0837	113.00	C12H16O5	-6.3
9	238.0837	113.00	C15H16O2	10.9
10	281.1049	492.34	C18H20O4	-3.8
11	281.1049	492.34	C14H17O6	-10.6
12	473.1906	137.72	C24H24O10	-13.4
- Chemical Structures:** The interface shows several chemical structures, including a large central structure and smaller fragments. One structure is labeled with a score of 98.5.


Analysis of Rubber Extract

Overlay TIC with Major Components (APCI-)




Name	RT	Formula	Mass	Score	Diff (MFG. ppm)	Height	Score (MFG)
Octadecanoic acid: Stearic acid: ADEKA FATTY ACID SA-200: LUNAC S-30: F-3: VL	13.07	C18 H36 O2	284.2717	99.39	-0.73	109608	99.46
Palmitic acid: Hexadecanoic acid	12.49	C16 H32 O2	256.2404	99.76	-0.68	61915	99.78
Irganox 1010: Sumilizer BP-101: ADKSTAB AO-60: Tominox TT: Anox 20	16.787	C73 H108 O12	1176.7842	99.16	-0.09	12968	98.89
9,10-EPOXYSTEARIC ACID	11.2	C18 H34 O3	298.2508	85.92	0.06	9075	85.93
3,5-Di-tert-butylbenzoic acid	13.763	C40 H74 O9	698.534	98.11	-0.95	8080	98.11
	9.464	C15 H22 O2	234.1621	84.49	-0.57	4110	84.51
	15.673	C38 H72 O8	656.523	96.4	-0.47	3394	96.4
	12.488	C17 H34 O4	302.2454	83.29	1.17	2852	83.29
	15.812	C42 H74 O P2	656.5226	68.93	-1.75	2821	68.93

MS and MS/MS of [M-H]⁻ of Irganox 1010




Mass Profiler Foil versus Rubber

Differential Analysis Impurity Profiling




Erucamide $C_{22}H_{43}NO$ MW 337.3344

1.5 fold change between Foil and Rubber




17 only in the Foil Extracts top formula hits



Name	Formula	CAS	RT	SD	Mass	SD	Abundance	RSD
2,13-Dimethyl-1,4,12,15-tetr...	C20 H32 O8		8.380	0.000	400.2099	0.0002	35080	1.10
	C23 H38 O8		9.870	0.000	442.2568	0.0001	26551	1.10
nDOP, Chissocizer nDOP: N...	C24 H38 O4	117-84-0	12.970	0.000	390.2770	0.0002	25613	0.14
Cyclic TEDDEG	C24 H24 O10		7.000	0.000	472.1367	0.0002	24075	1.10
Cyclic TEDDEG	C24 H24 O10		7.250	0.000	472.1367	0.0002	22803	1.10
	C24 H40 O8		10.203	0.006	456.2722	0.0003	22614	1.10
	C30 H24 O12		10.230	0.000	576.1265	0.0001	18445	1.10
	C29 H36 N2 O8		9.140	0.000	540.2468	0.0003	16163	1.10
	C30 H48 O12		10.130	0.000	600.3140	0.0003	13373	1.10
	C13 H22 O		10.100	0.000	194.1670	0.0002	7889	0.02
	C22 H36 O8		9.490	0.000	428.2409	0.0003	6989	1.10
Dimethyl sebacate	C12 H22 O4	106-79-6	8.450	0.000	230.1518	0.0003	6472	1.10
	C27 H44 O2		13.083	0.006	400.3344	0.0002	5482	1.10
	C27 H44		15.137	0.006	368.3443	0.0001	4062	1.10
MEHP	C16 H22 O4	4376-20-9	10.060	0.000	278.1521	0.0005	3594	0.78
Squalene	C30 H50	111-02-4	18.950	0.000	410.3904	0.0002	3205	1.10
Sansocizer DHP, Chissocize...	C22 H34 O4	3648-21-3	12.493	0.006	362.2458	0.0008	2905	1.10

Mass Profiler Unique in Rubber Extract



Compounds in Rubber Extract (+)


Name	Formula	CAS	RT	SD	Mass	SD	Abundance
1-(Decanoyloxy)-3-(octanoyloxy)-2-propan...	C39 H70 O6		13.777	0.006	634.5173	0.0002	249490
1,3-Bis(octanoyloxy)-2-propanyl (9E,12E)...	C37 H66 O6		17.023	0.006	606.4859	0.0001	237680
1-PALMITIN-3-LINOLEIN	C37 H68 O5		15.687	0.006	592.5067	0.0001	216493
1-PALMITIN-3-LINOLEIN	C37 H68 O5		15.830	0.000	592.5066	0.0002	158559
	C21 H40 O3		12.540	0.000	340.2982	0.0000	147920
1,2,4-Benzenetricarboxylic acid, 1,2,4-trid...	C39 H66 O6	4130-35-2	14.450	0.000	630.4857	0.0001	140291
Hydroxyethyl oleate	C20 H38 O3	4500-01-0	12.317	0.006	326.2825	0.0001	127793
	C37 H64 O5		17.030	0.000	588.4749	0.0001	99203
	C49 H103 N4 ...		21.300	0.000	890.7561	0.0002	90816
	C57 H100 O9		19.860	0.000	928.7352	0.0001	90258
	C39 H64 O7		15.563	0.006	644.4645	0.0003	87914
1,2,4-Benzenetricarboxylic acid, 1,2,4-trid...	C39 H66 O6	4130-35-2	14.550	0.000	630.4855	0.0002	87035
	C59 H101 Br O...		19.493	0.006	984.6606	0.0004	74418
	C39 H72 O5		16.993	0.006	620.5372	0.0002	69718
1-(Decanoyloxy)-3-(octanoyloxy)-2-propan...	C39 H70 O6		18.110	0.000	634.5164	0.0002	69460
1-(Decanoyloxy)-3-(octanoyloxy)-2-propan...	C39 H70 O6		13.850	0.000	634.5168	0.0000	69423
	C51 H105 N4 ...		18.897	0.046	932.7666	0.0002	59555
	C50 H101 O12 P		21.210	0.010	924.7037	0.0001	59392
			0.938	0.008	120.1149	0.0001	57139
	C32 H34 O9		16.800	0.000	562.2196	0.0001	56015
2-Oxiraneoctanoic acid, 3-octyl-, octyl ester	C26 H50 O3	106-84-3	13.983	0.006	410.3758	0.0001	51641
	C39 H72 O5		17.113	0.006	620.5372	0.0001	51543




Identification of Extractable Compounds from PolyPropylene Bottle




Components using Data Mining (MFE)



MS/MS m/z 537.5350 [M+H]⁺ C₃₄H₆₈N₂O₂




MS/MS m/z 565.5669 $[M+H]^+$ $C_{36}H_{72}N_2O_2$



Agilent Technologies

MS/MS m/z 593.5977 $[M+H]^+$ $C_{38}H_{76}N_2O_2$



Agilent Technologies

Molecular Structure Correlator $C_{34}H_{68}N_2O_2$

Agilent MassHunter Molecular Structure Correlator 8.05.00 - Reflux PP AutoMMS-N,N'-Ethylenebis(amide) M=15, ce=27

File Settings Help

Compound Formula

M = 536.5277, 3 formula candidates from MFG

#	Formula	len	max diff	d(M)ppm	d(M)ppm
1	C34H68N2O2	13	14	2.7	0.7
2	C34H68N2O2	0	0	-4.3	4.3
3	C34H68N2O2	5	3	8.2	8.2

Fragment formulas for C34H68N2O2

#	m/z	intensity	nom. mass	formula	d(M)ppm
1	281.2561	3126.50	25.96	C18H37N2	12
2	282.2767	1020.91	100.00	C18H39NO	-2.0
3	283.2922	3108.53	25.94	C18H37NO	16.9
4	298.3058	3006.03	31.86	C18H39NO2	0.6
5	300.3058	753.42	6.27	C18H41NO2	12.0

3 mixed compound formulas

Fragments of structure #1 - elucidated: 80.0% ions, 95.0% Weight

#	m/z	intensity	Weight%	No. of ions	Best score
1	281.2561	3126.50	13.3	4	91.9
2	282.2767	1020.91	48.8	1	97.7
3	283.2922	3108.53	12.8	1	94.2
4	298.3058	3006.03	21.9	8	98.5
5	300.3058	753.42	4.4	0	0.0

Structure Search Parameters: Compounds/Total: 4/5, Chem/Order (Sh): 1.00, Sort by: Score, Add Structure: Add

Standard IChEster: VQZQVLOSJHROV-UPFFFAQPSAH, Score: 87.37, Misc: Save, Delete, Chem/Order: 214261

Standard IChEster: UQZQVQVQDFTH-UPFFFAQPSAH, Score: 87.37, Misc: Save, Delete, Chem/Order: 21365


Standard IChEster: XSBHFNCOJGGG-UPFFFAQPSAH, Score: 88.86, Misc: Save, Delete, Chem/Order: 255404

Parity: 5.0 dM=1.2ppm, Score=93.9, Parity: 7.0 dM=1.2ppm, Score=93.1, Parity: 7.0 dM=1.2ppm, Score=93.1, Parity: 7.0 dM=1.2ppm, Score=93.1


Agilent Technologies

Possible Structures for $C_{34}H_{68}N_2O_2$


N,N'-1,2-Ethanediybis(2-pentylundecanamide)




N,N'-ethane-1,2-diylidihexadecanamide



N,N'-1,6-Hexanediylditetradecanamide



N,N'-Dihexadecyloxamide



Agilent Technologies


Agilent Technologies

Internet Search Confirms Database Hit

KAOWAX EB-P

Kao Wax EB-P is an amide type of wax which has high melting point and special property showing the low viscosity in molten form. In molten form under high temperature, it's compatible well with resin or solvents. Kao Wax EB-P is a lubricants and/or mould releasing agent for thermoplastic resin such as ABS resin, PS, PVC and so fort.

Chemicals structure:



N,N'-Ethylene-Bis- Stearamide

Specification :

Appearance	:	White powder
Color (Gardner)	:	5 max
Acid value (mg KOH/g)	:	10 max
Total Amine value (mg KOH/g)	:	2.5 max
Melting point, °C	:	141.5 - 146.5
Volatile matter, %	:	0.5 max
Particle size	:	10 % max. on 16 mesh screen

Applications :

Lubricant for ABS resin, polystyrene and copolymers, PVC, polyolefine, phenol resin, POM, EPS, PA

Octadecanamide,N,N'-1,2-ethanediybis- (CAS:110-30-5)

Update time: Jun 04 2014
 Purity: 99 Min. Order: 1 Kilogram
 Supply Ability: 100 Year/Metric Ton Formula: C38H76N2O2
 Octadecanamide,N,N'-ethylenebis- (6CI,7CI,8CI); 1,2-Bis (octadecanamido)ethane; 1,2-Ethylenebis(stearamide); Abril 10DS; Abril wax 10DS; Accurel SL 624; Acraxax C; Acraxax C-DF 1; Acraxax CT; Acraxax HMS; Advaxax; Advaxax 240; Advaxax...

Ethylene bis(stearamide) (CAS:110-30-5)


Update time: Jun 04 2014
 Formula: C38H76N2O2
 Ethylene bis(stearamide)

Ethylene Bis Stearamide (CAS:110-30-5)


Update time: Jun 03 2014
 Purity: 98% Min. Order: 25 Kilogram
 Supply Ability: 1-500 Year/Metric Ton Formula: C38H76N2O2
 Ethylene Bis Stearamide CAS:110-30-5 EINECS NO : 203-755-6
 M.F.: C38H76N2O2 M.W.: 593.0222

 Agilent Technologies


What to learn More....



ISBN: 978-0-470-28176-5
 June 2009



ISBN: 978-0-470-17365-7
 January 2012



 Agilent Technologies

BioPharma Road Show July 2014

76

Thank you
Let's Continue the Conversation

www.agilent.com

[Access Agilent](#)

Facebook  [@Agilent.Tech](#)

Twitter  [@Agilent](#)

[LinkedIn](#) 

 Agilent Technologies

Agilent Environmental and Analytical Solutions
5/18/2015
7777

More Rapid Extractable & Leachable Analyses with Advanced Mass Spectrometers

By Cynthia Challenger, PhD

Advances in technology are improving the sensitivity and accuracy of mass spectrometry increasing its use for the analysis of extractables and leachables.

Dennis Jenke, a distinguished scientist in the Technology Resources Division of Baxter Healthcare Corporation, the coupling of high-efficiency chromatographic techniques with information-rich detectors is proving very useful for E&L analysis. In particular, he notes that the **coupling of ultra high-performance liquid chromatography (UHPLC) with high-resolution MS detectors** capable of providing **accurate mass information** enables the determination of the **empirical formulas** for otherwise unidentified E&Ls is a noteworthy recent advance.

With the sensitivity and accuracy of the instruments available today, this time can usually be significantly reduced because **structure elucidation is now possible without the need to perform laborious isolation work** for the individual unknowns," Rushing explains.

"Accurate mass is **not always able to give a single definitive empirical formula**. Oftentimes several viable candidates are identified and it is not possible, based on the accurate mass data alone, to establish which candidate formula is the correct one. There can be literally hundreds of compounds with a certain empirical formula but only a **few may be toxic**. It is **necessary to know the structure** in order to reach the ultimate end game for E&L, which is safety assessment." Dennis Jenke"

 Agilent Technologies

BioPharma Road Show July 2014

78

Concerns Associated with Leachables from Packaging

Regulators in the pharmaceutical industry are becoming increasingly aware of the need to understand whether **species that can be extracted from the primary packaging** (direct contact packaging and container closures) will **appear as leachable species across a range of dosage forms.**

Johnson & Johnson again recalls a Tylenol product for moldy smell

Updated 10/22/2010 10:50 AM | Comment | Recommend

By Linda A. Johnson, Associated Press

TRENTON, N.J. — Tylenol caplets are again being recalled by Johnson & Johnson after some consumers were sickened by ingesting the pills or **smelling a musty or moldy odor** on the bottles.

Johnson & Johnson's McNeil Consumer Healthcare unit received some complaints from people who reported **nausea, stomach pain, vomiting and diarrhea** after taking the pills, a spokesman said Monday.

This makes 13 recalls in barely a year by the health care giant, mostly for more serious problems with its nonprescription drugs, contact lenses, blood glucose test strips and hip implants.

GET MORE INFORMATION: Tylenol 8-hour caplet 50 count

It's the fifth time that the New Brunswick, N.J.-based company has recalled nonprescription medicines because of consumer complaints about an unpleasant odor.

The moldy odor is thought to be due to the presence of trace amounts of a chemical called **2,4,6-tribromoisole**.

Previously, the company said the **chemical was on wooden shipping pallets** that it has stopped using after Jan. 15.

News Releases/Statements

Kellogg Company Voluntarily Recalls Select Packages of Kellogg's® Corn Pops®, Kellogg's® Honey Smacks®, Kellogg's® Froot Loops® and Kellogg's® Apple Jacks®

Jun 25, 2010

BATTLE CREEK, Mich., June 25 /PRNewswire-FirstCall/ -- Working in consultation with the United States Food and Drug Administration (FDA), Kellogg Company (NYSE: K) is implementing a voluntary recall of certain breakfast cereals due to an **uncharacteristic off-flavor and smell** coming from the **liner in the package**.

Only products with the letters "KN" following the Better If Used Before Date are included in the recall. Products with a "KM" designation are NOT included in the recall. In addition, no products in Canada are affected.

Kellogg's® Apple Jacks®

UPC 3800039136 1: 17 ounce package with Better If Used Before Dates between APR 10 2011 and JUN 22 2011
 UPC 3800039132 3: 8.7 ounce packages with Better If Used Before Dates between JUN 03 2011 and JUN 22 2011

Leachable examples from packaging summarized in C&EN article, Issue Aug 31, 2009



C&EN CHEMICAL & ENGINEERING NEWS

Serving The Chemical, Life Sciences & Laboratory Worlds

Home Magazine News Departments Collections Blogs

Volume 87 Issue 35 (pg. 11-15)
Issue Date: August 31, 2009

COVER STORY
Chemicals Leach From Packaging
Food and drugs just can't leave their wrappings behind

By Sarah Everts

Department: Science & Technology
 News Channels: Environmental SCENE
 Keywords: leachables, food, packaging, plastic, rubber, metal, ink

UNINVITED GUESTS
Some components of food packaging end up in your food.
Credit: Shutterstock

Earlier this year, when **German** were cracking open boxes of chocolate mazel—a common **breakfast cereal** in their country—it's unlikely they thought it did to the cereal: **Smelly plastic leachables** seeped into the fact that this component of **padding ink** had slipped from the outside of the cardboard box and into the cereal. That is, until the European Food Safety Authority (EFSA) was asked to look into the matter.

Open a cereal box or a carton of juice, breathe in an asthma drug from an inhaler, or pop an arthritismine pill out of a plastic pouch's metal foil. You are probably thinking about the product you're about to consume and not about the packaging—except, perhaps, for a tinge of regret about contributing to landfill waste. But there's something to keep in mind: Even when the wrapping comes off, you inevitably ingest some of the container.

U.S. FDA 21 CFR 211.94(a) states “Drug product containers and closures shall **not be reactive, additive, or absorptive** so as to later the safety, identity, strength, quality or purity of the drug beyond the official or establishment requirements”

PQRI (Product Quality Research Institute) is a working group established to developed **regulatory guidance** for Extractable/Leachable analysis, which is also recognized by the FDA

Definitions:

Extractable:

- Chemical compounds that can be removed from a primary container or component material into the drug or biologic by exertion of an artificial or exaggerated force, potential to migrate from polymeric construct into the drug product
= What **can** come out

Leachable:

- Chemical components that migrate from a container-closure system into the final packaged drug product under normal storage conditions. Direct contact with formulation under normal conditions
- Subset of Extractable
=What **DOES** come our

Classification of Extractable Substances

- **Passive Extractable**
 - “Extracted substance that was not chemically modified as a result of the drug product – container/closure interaction. This means that passive extracted substance is also an extractable substance by definition.”
- **Reactive Extractable**
 - “An Extracted substance that was chemically modified as a results of the drug product – container/closure interaction. A reactive extracted substance is not an extractable substance but may be linked to one such as Irganox degradant.” Reaction with peroxide e.g..

Container Closure System: Drug Packaging

What Drug Packaging Should Never Do

- Contribute harmful components to the drug
- Increase toxicity of the product
- Add Genotoxic or Carcinogenic Components
- Impact the stability and Efficacy of Drug
- **Drug Manufacturing** may also contribute impurities not from packaging but from production and storage equipment such as filters, gaskets, tubing, storage bags, cleaning process

Analytical Evaluation Threshold ($\mu\text{g/g}$)

Device Component:

$$\text{AET} = \frac{\text{SCT} \times D_t}{D_d \times m}$$

D_d = Doses per day
 D_t = Total Labelled Doses
 m = Mass of Component

Drug Delivery Product:

$$\text{AET} = \frac{\text{SCT} \times D_t}{D_d}$$

Final AET:

$$\text{AET}_{\text{fin}} = \text{AET}_{\text{Est}} \{1 - (R_{f\%RSD}/100)\}$$

Extract Components Quantified Against a Single Authentic Reference Standard Are Quantified using Average response Factor ($R_{f\text{av}}$)

Risk Assessment Table by Dosage Forms

Degree of concern associated with the route of administration	Likelihood of packaging component-dosage form interaction		
	HIGH	MEDIUM	LOW
HIGHEST	Inhalation Aerosols and Solutions; Injections and Injectable Suspensions	Sterile Powders and Powders for Injections, Inhalation Powders	
HIGH	Ophthalmic Solutions and Suspensions, Transdermal Ointments and Patches, Nasal Aerosols and Sprays		
LOW	Topical Solutions and Suspensions, Topical and Lingual Aerosols, Oral Solutions, and Suspensions	Topical Powders, Oral Powders	Oral Tablets and Oral (hard and soft Gelatin) Capsules

^ Figure 1: Risk Assessment Table for various dosage forms (adopted from ref. 19)

Toxicological Threshold for OINPD from PQRI Group

- 1) The Qualification Threshold (QT) is the level (5 µg/day) below which a given leachable is not considered for safety qualifications unless it presents a structure-activity (SAR) concern.
- 2) The Safety Concern Threshold (SCT) is the dose (0.15 µg/day) below which a leachable would present negligible concern for adverse carcinogenic and noncarcinogenic effects.
- 3) Known highly toxic substances such as PAHs, nitrosamines and 2-mercaptobenzothiazole are considered to be "special case compounds" for orally inhaled and nasal drug applications and should be considered on a case-by-case basis.
- 4) By using one of these toxicological thresholds as well as dosing information, the Analytical Evaluation Threshold (AET) can be calculated. The AET is the level at or above which an OINDP pharmaceutical development team should identify and quantify a particular extractable and/or leachable and report it for potential toxicological assessment.

Spectroscopy Organic Impurity Analysis Options



Agilent Cary 630 FTIR

Identifying and confirming the structure of an impurity or degradant



1260 Infinity Purification systems

Isolation of impurities



1260 Infinity Analytical SFC system

Detection of chiral impurities



7100 CE system

Analysis of polar impurities

Container Closure System: Drug Packaging

- The sum of the packaging components that together contain and protect the dosage form. Includes both primary, secondary and tertiary packing components

What Drug Packaging Should Do:

- Storage: adequately preserve the integrity of the drug
- Deliver the drug (Inhalers, prefilled syringes, ophthalmic, patches)
- Protect the Drug during shelf life



Safety Thresholds from PQRI L/E Working Group

Safety Concern Threshold: (SCT) 0.15 μg per day, which is defined as the threshold below which an individual leachable would have a dose so low as to present **negligible safety** concerns from carcinogenic and non-carcinogenic toxic effects.

Qualification Threshold: (QT) 5 μg per day: Threshold below which a given leachable is not considered for safety qualification (toxicological assessments) unless the leachable presents structure-activity relationship (SAR) concerns

Analytical Evaluation Threshold: (AET) is determined by consideration of the SCT and the specific drug product delivery configuration (number of doses in a Drug Product vs single dose)

D. Norwood, L.M. Nagao, C.L.M. Stults; *J. Pharma Sci and Tech.*, (2013) 67(5), 413-429

Typical Extractable Leachable Project Workflow

Extraction Method and Conditions Critical

Extraction conditions should mimic as close as possible drug product formulation and contact

- Simulate conditions of contact
- Simulate product conditions with extraction solution
- Simulate exposure using temperature and time
- Multiple solvents or extracting media with varying extracting power
- Multiple and complementary extraction techniques
- extraction conditions that allow equilibrium to be achieved.

Extraction Conditions Critical

Match conditions to product conditions:

- The chemical nature of the extracting media,
- The time and duration of the extraction process,
- The temperature at which the extraction is performed,
- The stoichiometry of the extraction process (extracted surface area per unit volume of extracting solution),
- The mechanism or process by which the extraction is accomplished

Experimental Methods

Plastic IV bag System

- **Headspace GC/MS**
 - Packaging materials were cut into 1.0-cm² piece and placed into 10-mL headspace vial and purged with N₂ before analysis
 - Equilibration temperatures investigated
 - Plastics: 85, 100, 150, 200, 250 and 275 °C.
 - Full evaporation: 85 and 100 °C
- **MMI GC/MS**
 - Packaging material or liquid drug product were extracted with solvents by sonication for 5-8 hours
 - Solvent utilized: dichloromethane (DCM), hexane, ethanol, acetone
 - Organic extracts were transferred to a glass insert placed inside an autosampler vial for analysis

Instrument Parameters

Two 5977 GC/MSD Systems

- **Headspace GC/MS**
 - 7697A Headspace Sampler and a 7890A GC coupled with a 5977A MSD
- **MMI GC/MS**
 - 7693A Automatic Liquid Sampler (ALS) and a 7890A GC coupled with a 5977A MSD
 - Equipped with a Multimode Inlet (MMI) operated in solvent vent mode

Headspace GC/MS

- Liner: 0.75-mm ultra-inert, straight tapered (p/n 5190-4048)
- Column: Agilent HP-5ms, UI, 30m x 0.25mm, 0.5µm (p/n 190915-133UI)

MMI GC/MS

- Liner
 - Drug suspension: 2mm ID ultra inert, dimpled (p/n 5190-4006)
 - Patch and IV bag: 4mm ID ultra inert (p/n 5190-3162)
- Column: Agilent HP-5ms UI, 30m x 0.25mm, 0.25µm (p/n 190915-433UI)

Elemental/Inorganic Impurity Analysis ICP-MS


Agilent 7700 Series ICP-MS



- Reliable trace-level analysis of all **16 elements** whose limits are defined in USP<232>.
- **Low detection limits** ensure that all regulated elements can easily be determined at or below regulated levels, and even when large sample dilutions are required.
- Can also be used in combination with a **HPLC, GC, and CE**, providing several options for separation (or speciation) of the different chemical forms of the elements.
- ICP-MS **achieves low detection limits** for almost all elements, including those found in the more extensive analyte list proposed in the ICH Q3D, such as Au and Tl

Analytical Technologies

Potential Analytical Techniques Depending on Polarity and Molecular Weight



CE: Capillary electrophoresis
GFC: Gel-filtration chromatography
SEC: Size-exclusion chromatography
IC: Ion chromatography

Most Agilent product lines can contribute

- GC/MS, GC-QQQ, GC-QTOF, HPLC, LC/MS, ICP/MS

No one method detects all extractables


- **Headspace GC/MSD** – volatile organic compound, high migration potential species (e.g. inks, adhesives, glue, processing solvent)
- **GC-MS** – semi-volatile organic compounds, residual monomers, antioxidant, plasticizers, PAHs, slip agents
- **LC-UV-MS** – non-volatile organic compounds, large oligomers, large antioxidants, thermally labile compounds
- **Ion chromatography** – Bromide, chloride, fluoride in elastomers
- **ICP-MS** – metals from aluminum canister, glass
- **Chemiluminescence Detection** – nitrosamines

Focus on the 5977 GC/MSD

Headspace sampling
MMI – large volume liquid techniques

QTOF Extractable Leachable Data Analysis Workflow


- Combination of UHPLC separation and accurate mass TOF technology
- Effective data mining algorithms to easily FIND compounds in a sample
- Optional SW to easily COMPARE samples or sample sets to identify differences
- AMRT DBs and MS/MS Lib Search to easily IDENTIFY targeted compounds
- Several algorithms to IDENTIFY unknown compounds (MFG, MSC)
- User Interface to easily NAVIGATE RESULTS
- Custom reporting to comprehensively REPORT results
- Full AUTOMATION of data acquisition, processing and reporting



TOF/Q-TOF Extractable/Leachable

Acquisition via accurate mass LC/TOF and Q-TOF

- UHPLC for fast (throughput) highly efficient separation (minimize suppression)
- Accurate mass TOF and Q-TOF with:
 - High resolution (up to 40000 on 6540/6550)
 - Sub-ppm mass accuracy
 - Fast acquisition speed (20-50 Hz) to complement UHPLC
 - High dynamic range (10^5)
 - High sensitivity (low pg on-column)
- LC/TOF for screening and identification of TARGETED compounds
- LC/Q-TOF for confirmation of targeted and identification of UNKNOWN compounds




Unsupervised Data Mining

Molecular Feature Extractor (MFE)

- Works on the 3-dimensional data set
- Remove noise
- Group covariant mass signals
- Combine signals with chemical relations
- Create extracted compound chromatograms & spectra

=> It's crucial to find all detectable compounds in a sample efficiently




Agilent Technologies

Molecular Feature Extraction: Peak Finding Algorithm

Data is processed by proprietary feature finding algorithm

- Find chromatographic peaks
 - Find all ions that are related
 - Include any adducts, such as Na^+ or K^+
 - Include isotopes ($[\text{M}+\text{H}]^+$, $[\text{M}+\text{H}+1]^+$, ...)
 - Check for dimers
 - Create a compound chromatograms (ECC) and spectra
- Sum all ion signals into one value (Feature)
- Fully automated processing




Agilent Technologies


Mass Profiler Software:

Differential analysis

- Some applications benefit from detecting the differences in compounds between samples, e.g.
 - Comparison of batches
 - Comparison of name brand and generic



Mass Profiler




 Agilent Technologies


Identification: Accurate Mass Databases

Identification of targeted compounds via AM or AMRT databases

- Personal Compound Databases (PCD)
 - PCD for extractable and leachable with several hundred compounds
- The addition of retention time into the database search (PCD) enables unambiguous identification of isobaric compounds



Compound List						
Name	Cpd	RT	Formula (DB)	Score (DB)	D# (DB, ppm)	Mass (DB)
Cyromazine	6	0.823	C6 H10 N6	87.6	-1.08	166.0967
Cyromazine	9	1.095	C6 H10 N6	97.73	-0.65	166.0967
Cyprodinol	138	12.245	C14 H15 N3	99.32	-1.26	225.1266
Cyproconazole	149	12.88	C15 H18 Cl N3 O	95.39	-2.13	291.1138
Cymiazole	78	6.118	C12 H14 N2 S	98.04	-1.6	218.0878





 Agilent Technologies

Additional Confirmation using MS/MS Libraries

LC/MS/MS library search

- NIST LC/MS/MS library
- Agilent PCDL libraries
- Support of multiple CEs





Agilent Technologies

Identification: Correlate MS/MS Spectra

MS/MS Structural Correlation (MSC)

- Algorithm to correlate “proposed structures” with accurate mass MS/MS fragment ion spectrum.
- Favor systematic bond dissociation approach over rule based fragmentation prediction approach.



Agilent Technologies





Agilent Technologies

MassHunter MS/MS Structural Correlation (MSC)

Confirm proposed structures and aid true unknown analysis

Provides highest confidence in confirmation of proposed structures in minutes instead of hours

- Assigns fragment ions to substructures of the proposed parent structure
 - Metabolite ID
 - Metabolomics – compound identification
- Aids in the determination of true unknown compounds beyond Molecular Formula Generation via parallel MSC of all structures retrieved for a unique Molecular Formula



Agilent Technologies



Agilent Technologies