

Enhanced Helium Mode Cell Performance for Improved Interference Removal in ICP-MS

Helium Collision Mode in ICP-MS

Helium (He) mode with kinetic energy discrimination (KED) is accepted as the preferred mode for effective and reliable removal of multiple polyatomic interferences in CRC (Collision/Reaction Cell) ICP-MS. While reactive cell gases may be very effective at removing individual interferences from individual analyte isotopes, they are unable to simultaneously remove all the interferences that occur in real-world applications, where the sample composition is often unknown, complex, or variable.

However, a few elements have previously still required the use of reactive cell gases for the best detection limits. A good example is the low-level analysis of Se, where a reactive gas (typically H₂) was required in order to give single ng/L (ppt) detection limits.

Enhanced Helium Mode Performance in the ORS³

The 3rd generation Octopole Reaction System (ORS³) introduced in the Agilent 7700 Series ICP-MS has significantly improved interference removal performance in He mode, as a result of the following unique cell design innovations:

- A longer, higher frequency octopole, increasing the collision rate
- Higher cell gas flow rates for higher cell pressure
- Higher bias voltage and greater collision energy

The combined effect of these developments is illustrated in Figure 1, which compares the residual ion energy overlap for the previous design of Agilent ICP-MS CRC (ORS²) and the new cell used on the 7700 (ORS³). The ORS³ provides a much smaller overlap between the residual energies of the analyte ions (⁷⁸Se, in green) and the interfering ions (Ar₂, in red), with the result that KED at the cell exit provides more effective separation of the analyte from the polyatomic interferences.




Figure 1. Comparison of analyte/polyatomic ion separation with ORS³ cell on the 7700 (bottom) compared to previous ORS² cell (top)

An additional benefit of ORS³ He mode on the 7700 Series is that the higher collision energy provides collisional dissociation of some polyatomic ions, as the energy of the polyatomic ion's collision with the cell gas is higher than its bond energy.

Again, the Ar₂ interference on Se is a good example, as shown in Table 1, where the collision energy of Ar₂ in the ORS³ (4.88eV) is higher than the bond dissociation energy of Ar₂ (1.33eV).

	E _{cm} (1st collision) [eV]
Ar ₂ Bond Dissociation Energy	1.33eV
Normal Collision (E _i =20eV)	0.98 (<1.33)
ORS ³ Collision (E _i =100eV)	4.88 (>1.33)

Table 1. Center-of-Mass Collision Energy (E_{cm}) at normal and ORS³ ion energies, leading to collision induced dissociation (CID) with ORS³

The combination of collisional dissociation and the enhanced separation of residual kinetic energy means that the background equivalent concentration (BEC) for ⁷⁸Se is significantly reduced (to around 2 ng/L (ppt)) with the ORS³, as shown in Figure 2.




Figure 2. He cell gas optimization plot, showing Ar₂ dissociation (CID)

The He mode performance available with the ORS³ is remarkable, as illustrated for Se, where the detection limit is improved from about 150 ppt (ORS²) to <5 ppt (ORS³ – Figure 3). Detection limits for S, P and Fe are also improved with ORS³, eliminating the need for reactive cell gases in many applications.




Figure 3. Calibration for ⁷⁸Se in ORS³ He mode (DL 4.5 ppt; BEC 2.7 ppt)

For more information on the 7700 Series ICP-MS visit the Agilent Technologies web site at: www.agilent.com/chem/icpms

