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1. "Introduction of the Positive List System for Agricultural Chemical Residues in Foods," Department of Food Safety, Ministry of Health, Labour and Welfare, June 2006
<http://www.mhlw.go.jp/english/topics/food-safety/positivelist060228/introduction.html>
2. "Positive List System for Agricultural Chemical Residues in Foods," Ministry of Health, Labour and Welfare
<http://www.mhlw.go.jp/english/topics/food-safety/positivelist060228/index.html>
3. "Maximum Residue Limits (MRLs) List of Agricultural Chemicals in Foods," The Japan Food Chemical Research Foundation
<http://www.m5.ws001.squarestart.ne.jp/foundation/search.html>
4. Table in item 6 (1), Section A General Compositional Standards for Food, Part I Food: The maximum residue limits of substances used as ingredients of agricultural chemicals in foods
<http://www.mhlw.go.jp/english/topics/food-safety/positivelist060228/dl/index-1a.pdf>
5. Table in item 7 (1), Section A General Compositional Standards for Food, Part I Food: The maximum residue limits of substances used as ingredients of agricultural chemicals in foods (Provisional MRLs List)
<http://www.mhlw.go.jp/english/topics/food-safety/positivelist060228/dl/index-1b.pdf>
6. "Analytical Methods for Residual Compositional Substances of Agricultural Chemicals, Feed Additives, and Veterinary Drugs in Food", Department of Food Safety, Ministry of Health, Labour and Welfare,
<http://www.mhlw.go.jp/english/topics/food-safety/positivelist060228/dl/060526-1a.pdf>
7. S. J. Lehotay, A. de Kok, M. Hiemstra, and P. van Bodegraven (2005) *J. AOAC Int.* **88**, 595-614
8. M. Anastassiades, S. Lehotay, D. Štajnbaher, and F. J. Schenck (2003) *J. AOAC Int.*, **86**, 412-431
9. V. Giarocco, B. Quimby, and M. Klee, "Retention Time Locking: Concepts and Applications," Agilent Technologies, publication 5966-2469E, www.agilent.com/chem
10. H. Prest, P. Wylie, K. Weiner, and D. Agnew, "Efficient Screening for Pesticides and Endocrine Disruptors Using the 6890/5973 GC/MSD System," Agilent Technologies, publication 5968-4884E, www.agilent.com/chem
11. Philip L. Wylie, "Screening for 926 Pesticides and Endocrine Disruptors by GC/MS with Deconvolution Reporting Software and a New Pesticide Library," Agilent Technologies, publication 5989-5076EN, www.agilent.com/chem
12. Bruce Quimby, "Screening for Hazardous Chemicals in Homeland Security and Environmental Samples Using a GC/MS/ECD/FPD with a 731 Compound DRS Database," Agilent Technologies, publication 5989-4834EN, www.agilent.com/chem
13. Mike Szelewski and Bruce Quimby, "New Tools for Rapid Pesticide Analysis in High Matrix Samples," Agilent Technologies, publication 5989-1716EN, www.agilent.com/chem
14. Chris Sandy, "A Blind Study of Pesticide Residues in Spiked and Unspiked Fruit Extracts Using Deconvolution Reporting Software," Agilent Technologies, publication 5989-1654EN, www.agilent.com/chem
15. Philip Wylie, Michael Szelewski, Chin-Kai Meng, Christopher Sandy, "Comprehensive Pesticide Screening by GC/MSD Using Deconvolution Reporting Software," Agilent Technologies, publication 5989-1157EN, www.agilent.com/chem

16. Xiaofei Ping, Chin-Kai Meng, Michael Szelewski, "Building Agilent GC/MSD Deconvolution Reporting Libraries for Any Application," Agilent Technologies, publication 5989-2249EN, www.agilent.com/chem
17. C. Lesueur and M. Gartner, "Routine Identification and Quantification of Pesticide Multiresidues in Fruit and Vegetable Samples with Full Scan, SIM, and Deconvolution Reporting Software," *Ernahrung/Nutrition* **29** (11), 466-471 (2005)
18. Y. Hirahara, et Al., "Validation of Multiresidue Screening Methods for the Determination of 186 Pesticides in 11 Agricultural Products Using Gas Chromatography (GC)," *J. Health Science* **51** (5) 617-627 (2005)
19. "Imported Foods Monitoring Plan for 2006" Notices Nos. 0331006 and 0525001, Office of Imported Food Safety, Ministry of Health, Labour and Welfare, May 26 2006
<http://www.mhlw.go.jp/english/topics/importedfoods/dl/4.pdf>

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© 安捷伦科技有限公司，2007