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


The screenshot shows the 'MicroLab OQ Operational Qualification of System Performance' interface. At the top, it displays 'Status: Ready' with a green indicator, 'User: User', and 'Method: LaserFreqCalTest_ATR'. Below this, there are four test options, each with a checkbox and a 'Number of Runs' input field: 'Performance (Signal-to-noise)' (15 runs), 'Stability Test' (30 minutes), 'Laser Frequency Calibration Check' (5 runs), and 'Spectral Resolution Check' (5 runs). A note states 'This test is performed using: Polystyrene'. At the bottom, there are 'Home', '< Back', and 'Next >' buttons. A prompt at the bottom reads 'Please press the Next button to initiate testing'.

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The screenshot shows the 'MicroLab OQ PerformanceTest' interface. It displays 'Status: Ready' with a green indicator, 'User: User', and 'Method: PerformanceTest'. The main area contains two visual instructions: 'Prepare compartment for background' with an image of the instrument's open compartment, and 'close compartment' with an image of the closed compartment. At the bottom, there are 'Home' and 'Next >' buttons.

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The screenshot shows the 'MicroLab OQ Results Summary' interface. It displays 'Status: Ready' with a green indicator, 'User: User', and 'Method: LaserFreqCalTest_ATR'. The main area lists five test results in green boxes: 'Diagnostics Verification - PASS', 'Performance Test SNR - PASS', 'Stability Test - PASS', 'Laser Frequency Calibration Check - PASS', and 'Spectral Resolution Check - PASS'. Below these, a prompt reads 'Select the individual test above to view the results'. At the bottom, it shows 'Operation Performance Qualification: PASS' and an 'Exit' button.

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© 安捷伦科技（中国）有限公司，2022