

### Agilent Turbo Pumps

Robust, industry-leading vacuum pumps and controllers for confidence across a broad range of applications





### Table of Contents

Experience the Agilent Advantage	4	TwisTorr 804 FS	38
50 Years of expertise in vacuum science	4	Turbo-V 1K-G	40
Sustainability commitment	5	Turbo-V 1001 Navigator	42
		Turbo-V 2K-G system	44
Agilent Turbo Pumps Typical Applications	6	Turbo-V 2300 TwisTorr	46
Turbo Pumps: A Legacy of Excellence	6	Agilent Turbo Pump Controllers	48
		TwisTorr 74/84 FS on-board controller	48
Ultra-high Vacuum and Physics Research (UHV)	8	TwisTorr 74/84 FS AG rack controller	49
		TwisTorr 305 FS remote controller	50
Solutions for Mass Spectrometry	9	TwisTorr Medium-TMP on-board controller	51
Vacuum solutions for GC/MS	9	TwisTorr Medium-TMP rack controller	52
Vacuum solutions for LC/MS	10	Turbo-V 1001 Navigator controller	53
Vacuum solutions for ICP/MS	10	Turbo-V 1001 rack controller	54
Vacuum solutions for QTOF	10	Turbo-V 2300 rack controller	55
Solutions for Analytical Electron Microscopes	11	Agilent Turbo Pump Accessories	56
Modal balancing	11	Inlet screens	56
IDX double damper	11	Air cooling kits	56
Variable speed electronics	11	Water cooling kits	56
		Vibration dampers	57
Solutions for Semiconductor Electron Microscopes	13	Turbo vent valve	57
Long-standing partnership	13	Purge valve	57
Full turbomolecular pump range	13	Turbomolecular pump parameters and definitions	58
Vibration damping solutions	13	Pump selection criteria	60
Solutions for PVD and Advanced Materials	14	Agilent Molecular-Drag Technologies	61
PVD thin film deposition and optical coatings	14	TwisTorr technology	62
		MacroTorr technology	62
Solutions for e-Mobility and Other Industrial Applications	15	Agilent Floating Suspension system	63
e-Mobility batteries manufacturing	15	AI enabled adaptive modal balancing	64
Vacuum load-locks	15	, , , , , , , , , , , , , , , , , , ,	
		Agilent Vacuum Link App	65
Agilent Turbo Pump Models	16		
Technical specification comparison chart	16	Agilent Turbo Pump Technical Notes	67
TwisTorr 74 FS	20	Agilent A-PLUS software	67
TwisTorr 84 FS	22	CE/CSA, EMC Specifications Compliance	67
TwisTorr 305FS and TwisTorr 305-IC	24		
TwisTorr 305 FSQ and TwisTorr 305-ICQ	26	Agilent Services and Support	68
TwisTorr 404 FS	30	Service and Support plans	68
Turbo V-551 Navigator	32	Contact Agilent	69
TwisTorr 704 FS	36	-	



# Experience the Agilent Advantage

Imagine 1350 rotations per second, that's the incredible speed of the Agilent 80 l/s turbomolecular pump. To put it in perspective, when you're driving on a highway, your car's wheels turn about 15 times per second.

The turbo rotor spins nearly 100 times faster!

Agilent turbomolecular pumps are engineered to meet the most demanding vacuum applications across various industries, including analytical instrumentation, the semiconductor manufacturing, industrial sectors, and high-energy physics research facilities. These applications require top-tier vacuum performance, reliability, and cleanliness, supported by a service organization that offers rapid response 24/7. Achieving these standards is our daily mission.

### With over 50 years of expertise in vacuum science, Agilent is a market leader in innovative design solutions

Turbomolecular pumps are essential in high-energy physics, fusion technology, and general UHV research. Applications such as synchrotron light sources, particle accelerator rings, UHV laboratory research, and fusion reactors demand extremely clean, reliable, and cost-effective high and ultrahigh vacuum (HV and UHV) solutions. Given the inaccessibility of most pumps, maintenance-free operation is crucial.

Agilent turbomolecular pumps are engineered to deliver unparalleled reliability, performance, and cleanliness for these demanding applications. Ceramic bearing pumps, with reduced rolling friction, low stress, and high thermal stability, ensure a



### **Agilent Advantage**

Agilent turbo pumps are designed for reliability and optimum performance in real world applications to meet highest quality standard.

Find integrated pumping systems, multiflow pump solutions for scientific instruments, and dedicated solutions for vibrationsensitive applications in nanotechnologies.

https://www.agilent.com/en/product/vacuum-technologies/turbo-pumps-controllers





longer operating life compared to conventional bearings. The use of ultra-low vapor pressure solid lubricants eliminates the need for maintenance, ensuring clean operation under all conditions.

Unlike many other pumps, Agilent turbo pumps feature both upper and lower bearings located on the rough vacuum side, preventing exposure to UHV and significantly reducing contamination risks, even in cases of misuse. The patented TwisTorr stages provide the highest speed and compression ratio within the smallest footprint. Additionally, Agilent turbo pumps can be mounted in any orientation, from vertical to horizontal to upside down, facilitating system design in the most space-constrained environments.

Agilent turbo pumps can operate at higher foreline pressures, enabling the use of dry roughing pumps and offering a completely clean, oil-free, compact, and cost-effective pumping solution. For applications requiring high gas throughput, the combination of TwisTorr pumps and TriScroll dry pumps represents the state-of-the-art solution.

Agilent turbo pumps are available with integrated or onboard controllers for easy plug-and-pump operation, or with rack-mounted controllers for environments where remote electronics placement is necessary, such as radioactive settings.

For instrumentation, including electron microscopes (SEM, TEM), focused ion-beam systems (FIB), and surface analysis, Agilent offers a comprehensive range of high and ultrahigh vacuum pumps tailored to meet the stringent requirements of these applications. Our turbomolecular pumps provide fast, oil-free air evacuation of large sample chambers, a critical requirement in applications requiring modern analytical tools, like semiconductor manufacturing.

The Agilent full range of turbo pumps allows focused-beam system designers to select the optimal pump size for the best chamber evacuation time, cost of ownership, and compact size for limited space situations. We also offer customized, low-vibration turbo pumps for the most sensitive microscopy applications, along with integrated pump controllers that provide high-control flexibility with minimal electromagnetic noise generation.

Agilent application-specific SEM turbo pumps are available in various speeds, including 80, 300, 550, 700, and 1000 l/s, ensuring that all designs meet the highest standards of performance and reliability.

### **Sustainability Commitment**

At Agilent, we are committed to sustainability and environmental responsibility. Our turbomolecular pumps are designed with energy efficiency in mind, reducing power consumption and minimizing environmental impact. The use of maintenance-free, solid lubricants not only ensure clean operation but also reduces waste and the need for hazardous materials. By enabling the use of dry roughing pumps, our solutions help eliminate oil contamination, contributing to a cleaner and safer environment.

Furthermore, the long operating life and reliability of Agilent turbo pumps mean fewer replacements and less waste, supporting a more sustainable lifecycle. Our dedication to innovation and sustainability ensures that Agilent products not only meet the highest performance standards but also contribute to a greener future.



TwisTorr 305-IC and 305 FS

### Agilent Turbo Pumps: A Legacy of Excellence

With over 50 years of expertise, Agilent turbo pumps offer an unparalleled fusion of performance and features across a diverse spectrum of applications:

### Ultra-vigh Vacuum and Physics Research (UHV)

Agilent turbo pumps excel in creating and maintaining extremely low-pressure environments, essential for scientific research, precision manufacturing, and cutting-edge experiments.



### Mass Spectrometry (MS)

Turbo pumps play a pivotal role in MS instruments, ensuring accurate ionization and precise analysis. From pharmaceuticals to environmental monitoring, MS benefits from Agilent technology.



### Electron Microscopes (EM)

Agilent turbo pumps enhance imaging resolution and sample analysis in EM systems. Whether investigating materials or biological specimens or nanocircuits, these pumps elevate EM performance. Agilent turbo pumps are ideal for:

- General purpose analytical electron microscopes
- · Semiconductor electron microscopes



### PVD and Advanced Materials (PVD)

Agilent turbo pumps contribute to thin-film deposition, surface modification, and plasma processes vital for industrial coatings, and beyond.



### Automotive and Other Industrial Applications

Agilent turbo pumps offer oil-free, clean vacuum for demanding industrial applications. In summary, Agilent's Turbopumps stand as a testament to innovation, reliability, and versatility across critical scientific and industrial domains.



### Making Advanced Science and Processes Possible

For many scientific inquiries and engineered processes, it is necessary to create evacuated workspaces. Creating, measuring, and maintaining such low-pressure conditions can be challenging and require several

technologies, as well as extremely careful design, preparation, and maintenance practices.

https://www.agilent.com/en/solutions/vacuum-solutions







TwisTorr 704



FwisTorr 804



**Turbo-V 1001** 









Multi Flow Pumps









Turbo-V 551



Turbo-V 701 SEM





Double Damper ISO200





Turbo-V 1K-G





TwisTorr 305 FSQ



Turbo-V 551



Turbo-V 701



Turbo-V-1K-G





# Ultra-High Vacuum and Physics Research (UHV)

### Ultra-high Vacuum and Extreme High Vacuum. Where Advanced Research and Vacuum Technology Intersect

Whether isolating subatomic particles or simulating the conditions of space, ultrahigh and extreme high vacuum are at the center of the human search to understand the physical universe. The unique condition of such very low pressure, nearly devoid of matter, is almost never found on our planet and must be created using specialized equipment and techniques.

Mastering the creation and maintenance of vacuum is critical to the pursuit of advanced research of the natural world

### Vacuum in Particle and Plasma Physics. Understanding How the Universe is Constructed

Particle and plasma physics are branches of the larger area defined as high energy physics (HEP). The goal of HEP is to explore what the world is made of and how it works on the smallest scale. Study of matter on this scale requires isolating subatomic particles in evacuated devices such as accelerators. The challenges of engineering such unique environments are considerable, as such, systems require the breadth of vacuum technologies to create, measure, and maintain the very low pressures required.

Agilent has been supporting scientific inquiry enabled by ultrahigh (UHV) and extreme high vacuum (XHV) for many decades. As the successor to Varian, which invented the ion getter pump that made UHV possible, Agilent remains committed to supporting physics exploration that expands collective knowledge of the physical world through continued innovation of vacuum technologies, products, and services.

#### Did you Know?

Discover how Agilent enables advanced research in particle and plasma physics:

https://www.agilent.com/en/solutions/vacuum-solutions/particle-plasma-physics





### Solutions for Mass Spectrometry

Mass spectrometry (MS) is an Agilent core business and has become a fundamental analytical tool in many industries. Thanks to advances in electronics, instrument designers can implement high-performance analytical power in a cost-effective, easy-to-use system. These developments require advanced vacuum systems characterized by multichamber, high-throughput designs on high-quality instruments. These requirements, in turn, demand cost-effective, high-performance vacuum pumps.

Agilent offers a full line of pumps and controllers that meet the most challenging vacuum requirements and are optimized for the specific requirements of modern mass spectrometry systems. Agilent can further customize its pumps by providing multi-inlet pumping systems that are compact and reliable.

### Vacuum Solutions for GC/MS

Gas chromatography mass spectrometers (GC/MS) typically use one vacuum chamber in relatively low gas load environments and an intermediate vacuum interface to analyze inorganic samples.

The Agilent TwisTorr 84 with On Board Controller or printed circuit board controller (PCB) is a very cost-effective solution for this common analytical technique. The TwisTorr 305 offers a compact, cost-effective solution for larger instrument designs.

### TwisTorr 74/84 FS



TwisTorr 305 SF



### Did you Know?

Agilent has a full line of products for analytical instrumentation and Laboratories. Please visit:

https://www.agilent.com/cs/library/brochures/brochurevacuum-solutions-analytical-instrumentation-5994-0681en-agilent.pdf



### Agilent Turbo Pumps Typical Applications



#### Vacuum Solutions LC/MS

Liquid chromatography mass spectrometers typically include multichamber, high-throughput vacuum systems. The Agilent multiflow turbomolecular pumps are designed for high-throughput operation with air cooling which is an important benefit for maintaining a compact system. They are also available several multiflow versions to increase their utility and performance in this application.

The integrated controller provides a high level of control function in a small package.



#### Vacuum Solutions ICP-MS

Inductively Coupled Plasma-MS systems have a wide range of vacuum requirements. Many systems can use heavy carrier gases such as argon, while collision cell designs use helium.

Agilent turbo pumps have a high efficiency motor and TwisTorr or MacroTorr drag stages to reduce heat production under gas load. These features allow the system to pump high levels of argon. Integrated or onboard controllers assure a compact package.



#### Vacuum Solutions Q-TOF

Time-of-flight (TOF) systems are becoming very important analytical tools in drug discovery and proteomics. Vacuum requirements vary, although small size is often an important consideration.

Agilent turbo pump and controller package provides the TOF designer with a high degree of flexibility with regard to high throughput, efficient heat dissipation and compact size.


### **Multi Flow Pumps**





### Did you Know?

Agilent can offer you a complete set of vacuum solutions for Analytical Labs, including Oil-Free primary pumps, suitable for the most popular analytical instruments? Check our Analytical V4L Vacuum for Lab solutions catalog.



### Solutions for Analytical Electron Microscopes

Agilent Optimized Turbomolecular Pumping Solutions for Analytical and General Purpose Scanning Electron Microscopes (SEMs).

With extensive experience collaborating with major electron microscope manufacturers, Agilent offers a comprehensive range of turbomolecular pumps. Regardless of your SEM model, we provide solutions based on our TwisTorr 84 FS and TwisTorr 305 FS platforms, meeting stringent vibration requirements.

Noteworthy features of Agilent turbomolecular pumps for SEM include:

### **Modal Balancing:**

A novel balancing technology minimizes vibration and noise during both full-speed operation and speed transients. Users benefit from a remarkably quiet experience and consistent vibration levels over time.

#### **IDX Double Damper:**

The Agilent patented damper optimally dampens vibrations at low and high frequencies.

Tuned to the pump mass and vibration characteristics, it ensures an ideal damping rate at any working point. Additionally, the unique design electrically insulates the pump from the system, preventing potential "ground loops" that could introduce artifacts in SEM images.

### Variable Speed Electronics:

Agilent SEM turbo pumps feature digital control units, allowing precise adjustment of turbo pump rotation speed. This flexibility helps avoid coupling between pump vibrations and system-level natural modes.

Our experts are available for troubleshooting vibration and resonance issues.

In summary, Agilent turbomolecular pumps enhance SEM performance by minimizing noise, optimizing damping, and ensuring reliable operation.

#### TwisTorr 84 FS



Part number X3502-64010

TwisTorr 305 FS



X3513-64008





### **Integrated Double Damper**

Part number 9699396





Agilent's patented Double Damper provides optimal vibration damping to the pump vibration signatures, to accommodate both low and high frequencies ranges. Additionally, the design grants electrical insulation to minimize possible electric systems ground loops through the pump, contributing to an artifact-free, better SEM image.

### Did you Know?

Agilent is offering a complete line of ion pumps and control units specific for electron microscopes columns and field emission sources.

Get in touch with our Contact Centers for more information





### Solutions for Semiconductor Electron Microscopes

Agilent provides a comprehensive range of turbomolecular pumping solutions specifically designed for electron microscopes (ems) used in the semiconductor production industry. These solutions cater to various EM types, including CD-SEM, DR-SEM, and eBeam-based wafer inspection systems.

Here are the key points about the Agilent offerings:

### Long-standing Partnership:

Agilent has been a trusted partner of SEM manufacturers since the early days of the semiconductor industry. This collaboration has allowed Agilent to develop tailored solutions that meet the stringent requirements of electron microscopes used in integrated circuits (IC) manufacturing workflows.

#### Full Turbomolecular Pumps Range:

Agilent's Turbomolecular pumps cover a wide range, from 70 to 1000 l/s (liters per second).

These pumps play a crucial role in maintaining vacuum conditions within the SEM chambers, ensuring optimal imaging and analysis.

### **Vibration Damping Solutions:**

Agilent unique design solutions focus on minimizing vibration and noise.

Notably, the Double Dampers, available for ISO100, ISO 160, and ISO 200 flange sizes, provide optimal damping at both low and high frequencies. These dampers are precisely tuned to the pump models, vibration parameters, and natural frequencies, resulting in effective vibration reduction.

In summary, Agilent specialized turbomolecular pumping solutions enhance the performance and reliability of SEMs in semiconductor manufacturing, contributing to precise imaging and analysis.

Turbo-V 551SEM



Turbo-V 701SEM



Turbo-V 1001SEM



Double Damper ISO200





### Solutions for PVD and Advanced Materials

#### Solutions for PVD Thin Film Deposition - Optical Coatings

In diverse industries, spanning from glass coating to tools sterilization, manufacturers universally seek robust and reliable vacuum technologies.

Agilent, with a keen focus on total cost of ownership, consistently pioneers innovative vacuum solutions. These cutting-edge systems not only maximize throughput but also thrive under challenging and variable operating conditions, while simultaneously simplifying maintenance to ensure optimal uptime and efficiency.



The Agilent K-G turbo pumps, purpose-built for heavy industrial use, exhibit remarkable resilience.


They effortlessly handle high gas loads and swiftly recover from accidental air in-rushes.

Notably, Agilent offers a distinctive line of turbo pumping systems meticulously tailored for argon pumping and plasma applications with the Turbo-V 1K-G and Turbo-V 2K-G.

These specialized pumps are optimized for argon and other heavy molecules, commonly encountered in PVD-based processes. Their standout features include superior throughput, precise thermal control, substantial power capacity, and air-inrush protection.

But these are no ordinary turbo pumps. The Agilent K-G line seamlessly integrates critical components: The electronics, vent valve, and purge valves are all encased in rugged aluminum for maximum durability and uptime.

This comprehensive package makes the K-G line the ideal pumping solution for lens coaters, roll-to-roll coaters, and various other PVD industrial applications.



Turbo-V 1K-G



Turbo-V 2K-G





### Solutions for e-Mobility and Other Industrial Applications

### e-Mobility Batteries Manufacturing

Roll-to-roll processing is a sophisticated fabrication technique employed in manufacturing. It involves continuous handling of flexible substrates, such as thin films or sheets, as they move through a series of precisely orchestrated steps. Roll-to-roll coaters handle diverse materials; however, their most recent application is in energy storage where they are a fundamental enabling technology for manufacturing battery electrodes. Agilent has a complete portfolio of turbomolecular pumps for roll-to-roll applications including the TwisTorr 305FSQ, Turbo-V line and the Turbo K-G line.



Note: Agilent offers an e-Mobility application, including leak detectors, diffusion pumps and more.



### Vacuum Load-Locks

A wafer cassette load-lock serves as a critical interface within semiconductor manufacturing environments. Its purpose is to facilitate the transfer of semiconductor wafers, loaded in cassettes, between different process chambers while maintaining a controlled environment. Load-locks play a crucial role in wafer processing steps, such as deposition, etching, and lithography. Agilent has a broad range of turbomolecular pumps designed for rapid cycling and fast pump down time.

Looking to reduce load lock pump down time? Please check the Agilent TwisTorr 305 line and Turbo-V line

Turbo-V 551



Turbo-V 701



Turbo-V 1001



TwisTorr 305 FSQ



		Twis	Torr 74 FS	TwisT	orr 84 FS	TwisTorr 305 FS
Flange Size		KF 40	ISO 63	KF 40	ISO 63	ISO 100 K / CFF 6" ISO 160 K / CFF 8"
Pumping Speed, L/s	H <sub>2</sub> He N <sub>2</sub> Ar	28 33 44 40	42 54 60 59	36 38 49 44	53 63 67 66	220 L/s 255 L/s 250 L/s 250 L/s
Compression Ratio	H <sub>2</sub> He N <sub>2</sub> Ar	$1.0 \times 10^{9}$ $2.0 \times 10^{5}$ $1.0 \times 10^{4}$ $> 1.0 \times 10^{9}$	9	$\ge 1 \times 10^{11}$ $2 \times 10^{6}$ $5 \times 10^{4}$ $> 1 \times 10^{11}$		1.5 x 10 <sup>6</sup> > 1 x 10 <sup>8</sup> > 1 x 10 <sup>11</sup> > 1 x 10 <sup>11</sup>
Base Pressure,	With recommended mechanical pump	< 5 x 10 <sup>-10</sup>		<5 x 10 <sup>-10</sup>		<1 x 10 <sup>-10</sup> mbar
mbar	With recommended dry pump	< 3.75 x 10	O <sup>-10</sup>	<5 x 10 <sup>-10</sup>		(<1 x 10 <sup>-10</sup> Torr)
Startup Time, min		< 2		< 2		< 3 min (longer when soft start is used)
Recommended	Mechanical pump	DS 40M, D	OS 102	DS 40M, D	S 102	DS102, DS302
Forepump	Dry pump	IDP-3, IDP	-7	IDP-3, IDP-	-7, IDP-10	IDP-3 (no gas flow), IDP-7, IDP-10,
Foreline Flange, Nominal Diameter	Klamp flange	KF16 NW		KF16 NW		KF16 NW (KF25 NW - optional)







TwisTorr 84 FS



TwisTorr 305 FS

		TwisTorr 305 FSQ	TwisTorr 305-IC	TwisTorr 305-ICQ	TwisTorr 404 FS
Flange Size		ISO 100 K	ISO 100 K / CFF 6" ISO 160 K / CFF 8"	ISO 100 K	ISO 100 KF CFF6"
Pumping Speed, L/s	H <sub>2</sub> He N <sub>2</sub> Ar	220 L/s 255 L/s 250 L/s	220 L/s 255 L/s 250 L/s 250 L/s	220 L/s 255 L/s 250 L/s	445 l/s 470 l/s 355 l/s 320 l/s
Compression Ratio	H <sub>2</sub> He N <sub>2</sub> Ar	2 x 10 <sup>4</sup> 1 x 10 <sup>5</sup> 2 x 10 <sup>8</sup>	1.5 x 10 <sup>6</sup> > 1 x 10 <sup>8</sup> > 1 x 10 <sup>11</sup> > 1 x 10 <sup>11</sup>	2 x 10 <sup>4</sup> 1 x 10 <sup>5</sup> 2 x 10 <sup>8</sup>	>10 mbar >10 mbar >4 mbar >8.5 mbar
Base Pressure With Recommended Forepump		<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)	<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)	<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)	<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)
Startup Time, min		< 3 min (longer when soft start is used)	< 3 min (longer when soft start is used)	< 3 min (longer when soft start is used)	< 5 min
D d. d	Mechanical pump	DS102, DS302	DS102, DS302	DS102, DS302	DS302
Recommended Forepump	Dry pump	IDP-3 (no gas flow), IDP-7, IDP-10,	IDP-3 (no gas flow), IDP-7, IDP-10,	IDP-3 (no gas flow), IDP-7, IDP-10,	IDP-10
Foreline Flange, Nominal Diameter	Klamp flange	KF16 NW (KF25 NW - optional)	KF16 NW (KF25 NW - optional)	KF16 NW (KF25 NW - optional)	KF25 NW (KF16 NW as optional accessory)







TwisTorr 305-IC



TwisTorr 305 ICQ



TwisTorr 404 FS

		Turbo V Naviga		Turbo-V 701 Navigator	TwisTorr 704 FS	TwisTorr 804 FS
Flange Size		ISO100K CFF6"	ISO160 K CFF8"	ISO200 K-F CFF10"	IS0100K-F CFF8"	ISO200K-F, ISO250K-F, CFF10"
Pumping Speed, L/s	H <sub>2</sub> He N <sub>2</sub> Ar	450 450 350	510 600 550	510 L/s 620 L/s 690 L/s	480 L/s 640 L/s 660 L/s 625 L/s	485 L/s 660 L/s 720 L/s 690 L/s
Compression Ratio	H <sub>2</sub> He N <sub>2</sub> Ar	>1 x 10 <sup>9</sup> 1 x 10 <sup>7</sup> 1 x 10 <sup>6</sup>		1 x 10 <sup>6</sup> 1 x 10 <sup>7</sup> 1x 10 <sup>9</sup>	>4 mbar 10 mbar 10 mbar 8.5 mbar	>4 mbar 10 mbar 10 mbar 8.5 mbar
Base Pressure with Recommended Fore- pump	With recommended mechanical pump  With recommended dry pump	<1 x 10 <sup>-10</sup>		<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)	<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)	<1 x 10 <sup>-10</sup> mbar (<1 x 10 <sup>-10</sup> Torr)
Startup Time, min		< 5		< 5 min	< 5 min	< 5 min
Recommended	Mechanical pump	DS 102		DS302	DS302	DS302
Forepump	Dry pump	IDP-7, IDP-10		IDP-15, TS300	IDP-10, IDP15, TS300	IDP-10, IDP15, TS300
Foreline Flange, Nominal Diameter	Klamp flange	KF25 NW		KF25 NW	NW25 (NW40 as optional accessory)	NW25 or NW40







Turbo-V 701 Navigator



TwisTorr 704 FS



TwisTorr 804 FS

		Turbo-	V 1K-G	-	Гurbo-V 100 Navigator	1	Turbo-V 2K-G System	Turbo-V 2300 TwissTorr
Flange Size		IS0160 F	IS0200 F	CFF8" ISO160 K	CFF10" ISO250 K-F	IS0250 K-F	ISO250 F	CFF12" ISO250 F
Pumping Speed, L/s	H <sub>2</sub> He N <sub>2</sub> Ar	680 950 810 750	730 1150 1080 1040	860 820 790	900 870 950	920 900 1050	- - 1600 -	1500 1800 2050
Compression Ratio	H <sub>2</sub> He N <sub>2</sub> Ar	>4 >5	x 10 <sup>4</sup> x 10 <sup>4</sup> x 10 <sup>7</sup> x 10 <sup>8</sup>		1 x 10 <sup>6</sup> 1 x 10 <sup>7</sup> 1 x 10 <sup>9</sup>		- - 3x10 <sup>5</sup> -	4 x 10 <sup>4</sup> 8 x 10 <sup>5</sup> >8 x 10 <sup>8</sup>
Base Pressure with Recommended Forepump			O <sup>-10</sup> mbar O <sup>-10</sup> Torr)		<1 x 10 <sup>-10</sup> mba (<1 x 10 <sup>-10</sup> Tor		<1 x 10 <sup>-8</sup> mbar (<1 x 10 <sup>-8</sup> Torr)	10 <sup>-10</sup>
Startup Time, min			< 5		< 4		< 7	< 6
Recommended	Mechanical pump	>20 m³/h	>36 m³/h		DS402		>40 m³/h	DS602
Fforepump	Dry pump	TS	6600		IDP-15			TS600
Foreline Flange, Nominal Diameter	Klamp flange	KF25 NW	KF40 NW		KF40 NW		KF40 NW	KF40 NW







Turbo-V 1001 Navigator





Turbo-V 2K-G System



Turbo-V 2300 TwissTorr








Dimensions: millimeters [inches]

### TwisTorr 74 FS

### **Technical Specification**

Dumming Chard	KF40	CFF 2.75"	ISO 63	CFF 4.5"		
Pumping Speed						
N <sub>2</sub>	44 L/s 33 L/s	50 L/s 40 L/s	60 L/s 54 L/s	60 L/s 54 L/s		
He	28 L/s	31 L/s	42 L/s	42 L/s		
H <sub>2</sub>	20 L/s 40 L/s	51 L/s	59 L/s	59 L/s		
Ar ———————	40 L/S	31 1/8	39 1/8	J9 L/S		
Man Oa - Thursdan	Air co	ooling	Water	cooling		
Max Gas Throughput	(35 °C ambien	-	(25 °C wat	ter temp./		
(Backing pump 5 m³/h)	(55 Callibleit	t temperature)	35° ambie	ent temp.)		
$\overline{N_2}$	130 s	sccm	130 s	sccm		
Ar	80 s	ccm	80 s	ccm		
Compression Ratio and Fo	reline Toleranc	e				
N <sub>2</sub>	1.0 >	( 10 <sup>9</sup>	>12 r	mbar		
He	2.0 >	∢10 <sup>5</sup>	>10 mbar			
H <sub>2</sub>	1.0 >	¢ 10⁴	>4 mbar			
Ar	> 1.0	x 10 <sup>9</sup>	>14 r	nbar		
Base Pressure with Recommended Forepump		< 5 x 10 <sup>-10</sup> mbar (<	3.75 x 10 <sup>-10</sup> Torr)			
Inlet Flange	KF 40, ISO 63, CFF 4, 5", CFF 2.75"					
Foreline Flange	KF16 NW					
Rotational Speed	70000 rpm (1167	7 Hz driving freque	ency)			
Start-Up Time	< 2 minutes					
	Mechanical: Agil	ent DS 40M / DS	102	·		
Recommended Forepump	mp Dry pump: Agilent IDP-3 / IDP-7					
Operating Position	Any					
Oper. Ambient Temperature	+5 to +35 °C					
Rel. Humidity of Air	0 - 90 % (not condensing)					
Bakeout temperature	80 °C for ISO (120 °C for CFF) at inlet flange					

Lubricant	Permanent lubrication				
Cooling Requirements	Forced air (5- 35 °C ambie	nt temperature)			
Air Cooling	Air flow temperature +5 to	+35 °C			
Water Cooling	Cooling water temperature: +15 to +25 °C Minimum flow: 65 L/h				
Noise Pressure level (at 1 mt at full speed)	40 dB(A)				
Storage Temperature	-40 to +70 °C				
Max Altitude	3000 m				
Weight kg (lbs)	Pump ISO 63: Pump CFF 4.5": Pump CFF 2.75": Pump KF 40:	2.05 kg (4.50) 3.50 kg (7.70) 3.34 kg (7.35) 2.37 kg (5.22)			
Conformity to Norms	CE, C-CSA-US, RoHS compliant as per 2011/65/UE				




#### Foreline Pressure (mbar)

Compression ratio vs foreline pressure

### **Ordering Information**

Pumps	Part Number
TwisTorr 74 FS ISO 63	X3502-64170
TwisTorr 74 FS KF 40	X3502-64171
TwisTorr 74 FS CFF 4.5"	X3502-64172
TwisTorr 74 FS CFF 2.75"	X3502-64173
Controllers	
TwisTorr 74 FS AG rack controller RS232/485	X3508-64301
TwisTorr 74 FS AG rack controller Profibus	X3508-64022
TwisTorr 74 FS onboard controller 110/220 V	X3509-64030
TwisTorr 74 FS onboard controller 24 VDC	X3509-64021
TwisTorr 74 FS PCB controller	X3510-64050
Cables	
Mains cable NEMA plug 3 m long	9699958
Mains cable EU plug 2.5 m long	9699957
Mains cable UK plug 3 m long	X3501-68005
Mains cable China plug 3 m long	X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	9699883
Power cable, 1 m long, 24 V	9699869
Extension cable*	9699942 (3 m) 9699942M007 (5 m) 9699942M006 (10 m) 9699942M005 (15 m) 9699942M004 (20 m)
Inlet Screen	
Inlet screen ISO 63	X3502-68001
Inlet screen CFF 4.5"	X3502-68000
Inlet screen, KF 40	9699309
Inlet screen CFF 2.75"	9699328
Cooling	
Metal water cooling kit	X3502-68002
Plastic water cooling kit	X3502-68003
Air cooling kit (0,5 m cable)	9699290
Air cooling kit extension cable (5 m)	9699940



### Inlet Pressure (mbar)


Pumping speed vs inlet pressure

Vibration Isolator	
Vibration isolator ISO 63	9699375
Vibration isolator CFF 4.5"	9699376
Venting	
Vent valve N.O. 0.5 mm orifice (0.7 m cable)	9699844
Vent valve extension cable	9699941 (5 m) 9699941M003 (10 m) 9699941M001 (15 m) 9699941M005 (20 m) 9699941M007 (30 m) 9699941M004 (50 m)
Vent screw M5	X3502-68005
Vent valve adapter M5-M8	X1699-64039
Purge	
Purge screw	X3502-68004
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 10 SCCM ¼ Swagelok - M12	9699240
Purge valve 20 SCCM NW16KF - M12	9699241
Purge valve 20 SCCM ¼ Swagelok - M12	9699242
Purge valve 10 SCCM ¼ Swagelok - 7/16-20 male threads	9699232
Purge valve 20 SCCM ¼ Swagelok - 7/16-20 male threads	9699236
Mounting	
Controller side mounting bracket	X3502-68006
CFF 4.5 Mounting kit	X3502-68007
Metric screws kit	X3502-68008
American screws kit	X3502-68009
Active Gauges	
FRG-700 full range gauge PVG-500 Pirani vacuum gauge PCG-750 Pirani capacitance gauge CDG-500 capacitance diaphragm gauge	Ask Agilent for details

<sup>\*</sup>For a complete list of available extension cables, see the Extension Cables for Turbo Pumps dedicated page on Agilent.com.



### TwisTorr 84 FS



Dimensions: millimeters [inches]


### **Technical Specification**


Pumping Speed	KF40	CFF 2.75"	ISO 63	CFF 4.5"		
N <sub>2</sub>	49 L/s	56 L/s	67 L/s	67 L/s		
He	38 L/s	46 L/s	63 L/s	63 L/s		
H <sub>2</sub>	36 L/s	40 L/s	53 L/s	53 L/s		
Ar	44 L/s	57 L/s	66 L/s	66 L/s		
Gas Throughput At Full Rotational Speed (with	Air coolin	a (35 °C)	Water o	_		
recomm. forepump)	7.111 0001111	g (00°0)	(25 °C,	05 L/h)		
$\overline{N_2}$	100 S	CCM	100 S	CCM		
Ār	70 S0	CCM	70 S	CCM		
Compression Ratio and For	eline Toleranc	9				
$\overline{N_2}$	70 a ≥ 1.0		>14 mbar			
He	2.0 x		>12 mbar			
H <sub>2</sub>	5.0 x		>4 mbar			
Ar	> 1.0 :	< 10''	>14 r	nbar		
Base Pressure with Re- commended Forepump (5 m³/h)	< 5 >	( 10-10 mbar (<	3.75 x 10-10 T	orr)		
Inlet Flange	CFF 4	l.5" OD		63		
	CFF 2.	75" OD	KF	40		
Foreline Flange	KF16 NW					
Rotational Speed	81.000 rpm (1350 Hz driving frequency)					
Start-Up Time	< 2 minutes					
Recommended Forepump	Mechanical: Agilent DS 40M / DS 102 Dry pump: Agilent IDP-3 / IDP-7					
Operating Position	Any					
Oper. Ambient Temperature	+5 to +35 °C					


Bakeout temperature	80 °C for ISO (120 °C for CFF) at inlet flange			
Lubricant	Permanent lubrication	า		
Cooling Requirements	Forced air (5- 35 °C ar	mbient temperature)		
Air Cooling	Air flow temperature -	+5 to +35 °C		
Water Cooling	Cooling water temperature: +15 to +25 °C Minimum flow: 65 L/h (0.30 GPM) Pressure: 2 to 4 bar (45 to 75 psi)			
Noise Pressure Level (at 1 mt at full speed)	40 dB(A) *			
Storage Temperature	-40 °C to +70 °C			
Max Altitude	3000 m			
Weight kg (lbs)	Pump ISO 63: Pump CFF 4.5: Pump CFF 2.75: Pump KF 40:	2.05 kg (4.50) 3.50 kg (7.70) 3.34 kg (7.35) 2.37 kg (5.22)		
Conformity to Norms	CE, C-CSA-US, RoHS compliant as per 2011/65/UE			
*average value ± 4 dB(A) std deviation				

Rel. Humidity Of Air

0 - 90 % (not condensing)







Pumping speed vs inlet pressure

Pumps	Part Number
TwisTorr 84 FS ISO 63	X3502-64000
TwisTorr 84 FS KF 40	X3502-64001
TwisTorr 84 FS CFF 4.5"	X3502-64002
TwisTorr 84 FS CFF 2.75"	X3502-64003
Controllers	
TwisTorr 84 FS AG rack controller RS232/485	X3508-64001
TwisTorr 84 FS AG rack controller Profibus	X3508-64002
TwisTorr 84 FS onboard controller 110/220 V	X3509-64000
TwisTorr 84 FS onboard controller 24 VDCt	X3509-64001
TwisTorr 84 FS PCB Controller	X3510-64000
Cables	
Mains cable NEMA plug 3 m long	9699958
Mains cable EU plug 3 m long	9699957
Mains cable UK plug 2.5 m long	X3501-68005
Mains cable China plug 3 m long	X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	9699883
Power cable, 1 m long, 24 V	9699869
Extension cable*	9699942 (3 m) 9699942M007 (5 m) 9699942M006 (10 m) 9699942M005 (15 m) 9699942M004 (20 m)
Inlet Screen	
Inlet screen ISO 63	X3502-68001
Inlet screen CFF 4.5"	X3502-68000
Inlet screen, KF 40	9699309
Inlet screen CFF 2.75"	9699328
Cooling	
Metal water cooling kit	X3502-68002
Plastic water cooling kit	X3502-68003

Air cooling kit (0,5 m cable)	9699290
Air cooling kit extension cable (5 m)	9699940
Vibration isolator	
Vibration isolator ISO 63	9699375
Vibration isolator CFF 4.5"	9699376
Venting	
Vent valve N.O. 0.5 mm orifice (0.5 m cable)	9699844
Vent valve extension cable	9699941 (5 m) 9699941M003 (10 m) 9699941M001 (15 m) 9699941M005 (20 m)
Vent screw M5	X3502-68005
Vent adapter kit M5-M8	X1699-64039
Purge	
Purge screw	X3502-68004
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 20 SCCM NW16KF - M12	9699241
Purge valve 20 SCCM ¼ Swagelok - M12	9699242
Purge valve 10 SCCM ¼ Swagelok - ¼ Swagelok	9699232
Purge valve 20 SCCM ¼ Swagelok - ¼ Swagelok	9699236
Mounting	
Controller side mounting bracket	X3502-68006
CFF4.5 Mounting kit	X3502-68007
Metric screws kit	X3502-68008
American screws kit	X3502-68009
Active Gauges	
FRG 700 full range gauge PVG 500 Pirani vacuum gauge PCG 750 Pirani capacitance gauge CDG-500 capacitance diaphragm gauge	

<sup>\*</sup>For a complete list of available extension cables, see the Extension Cables for Turbo Pumps dedicated page on Agilent.com.





TwisTorr 305 FS

TwisTorr 305-IC

### TwisTorr 305 FS and TwisTorr 305-IC

### **Technical Specifications**

<u>-</u>	
Pumping Speed	ISO 100 K / CFF 6"/ ISO 160 K / CFF 8"
H <sub>2</sub>	220 L/s
He	255 L/s
N <sub>2</sub>	250 L/s
Ar	250 L/s

#### Max Gas Flow Rate

Ν, 250 SCCM

- Note: values refer to water-cooling pump version with: water temperature between 15 and 20 °C (non condensing)
- backing pump with pumping speed equal or above 5  $\ensuremath{\text{m}^{3}/\text{h}}$

#### **Compression Ratio**

H <sub>2</sub>	1.5 x 10 <sup>6</sup>
He	> 1 x 10 <sup>8</sup>
$N_2$	> 1 x 10 <sup>11</sup>
Ar	> 1 x 10 <sup>11</sup>

### Max Foreline Pressure Tolerance

 $N_2$ 12 mbar

#### Note:

Foreline tolerance defined as the pressure at which the turbo pump still produce a compression of 100.

For continuous operation, water cooling is recommended (water temperature between 15 and 20 °C).


Base Pressure with	<1 x 10 <sup>-10</sup> mbar
Recommended Forepump	(<1 x 10 <sup>-10</sup> Torr)

According to standard DIN 28 428, the base pressure is that measured in a leak-free test dome, 48 hours after the completion of test dome bake-out, with a turbo pump fitted with a ConFlat flange and using the recommended prevacuum pump.

Inlet Flange	ISO 100 K, CFF 6", ISO 160 K, CFF 8"	
Foreline Flange	KF16 NW (KF25 - optional)	
Max Rotation Speed	60.600 rpm (1010 Hz driving frequency)	
Start-Up Time	< 3 minutes (longer when soft start is used)	
Recommended Forepump	Dry pumps: IDP-3 (no gas flow), IDP-7, IDP-10, Mechanical: DS102, DS302	
Operating Position	Any	
Oper. Ambient Temperature	+5 to +35 °C	
Dakasut Tampawatura	ISO flange: 75 °C at inlet flange max CFF flange: 100 °C at inlet flange max	
Bakeout Temperature	Note: Measure a point close to the sealing element.	
Lubricant	Permanent lubrication	
Cooling Requirements:		
Air Cooling	Natural convection (only with no gas load) Forced air (5 - 35 °C ambient temperature)	
Water Cooling	Minimum flow: 50 L/h (0.22 GPM) Temperature: +15 °C to +30 °C Max pressure: 5 bar (75 psi)	
Noise Pressure Level (at 1m at full speed)	41 dB(A)	

Note: mean values based on a significative sample (Ar and N $_2$  compression ratio estimated); standard deviation per test: pumping speed: below  $\pm$  7%; noise pressure level:  $\pm$  10% (only pump).

Installation Category	II	
Pollution Degree	2	
Storage Temp.	-40 to +70 °C	
Max Altitude	3000 m	
Weight kg (lbs) TwisTorr 305-IC	ISO 100 K CFF 6" ISO 160 K CFF 8"	5.74 (12.6) 8.06 (17.7) 6.18 (13.6) 10.33 (22.7)
Weight kg (lbs) TwisTorr 305 FS	ISO 100 K CFF 6" ISO 160 K CFF 8"	5.84 (12.8) 8.16 (17.9) 6.28 (13.8) 10.43 (22.9)



TwisTorr 305 FS and TwisTorr 305-IC Pumping Speed

TwisTorr 305 FS and TwisTorr 305-IC Compression Ratio

### **Technical Specifications**

Remote Controller	
Voltage Frequency Power Fuse	100 - 240 Vac (voltage fluctuation +/- 10%) 50 to 60 Hz 450 VA 2 x T4 A (slow blow) 250 V
Power Supply (24 Vdc) for IC Controller: Input voltage Max input power Stand-by power Max operating power	24 Vdc 200 W 10 W 150 W with water or air cooling
Protection Fuse	8 A
Max Operating Altitude	3000 m
USB Communication	as per USB 1.1
Power Cable	Required motor input voltage is 24V +/-10%; please dimension power cable to guarantee the minimum voltage level. I.e. for AWG 20 resistence is 33,31 mOhm/m; so for a cable of 6 meters and maximum current (7.5A) the voltage lost is 1.5V.
Compliance With:	EN 61010-1 EN 61326-1 EN 1012-2 EN 12100 EN 50581 Machinery Directive 2006/42/EC Electromagnetic Compatibility Directive 2014/30/EU Directive 2011/65/EU





TwisTorr 305 FSQ

TwisTorr 305-ICQ

## TwisTorr 305 FSQ and TwisTorr 305-ICQ

### **Technical Specifications**

Pumping Speed	ISO 100 K
H <sub>2</sub>	220 L/s
He	255 L/s
$N_2$	250 L/s

#### Max Gas Flow Rate

	TwisTorr 305 FSQ	TwisTorr 305-ICQ
H <sub>2</sub>	500 SCCM	500 SCCM
He	500 SCCM	500 SCCM
$N_2$	450 SCCM	380 SCCM
Ar	100 SCCM	

Note: values refer to water-cooling pump version with:

- water temperature between 15 and 20 °C (non condensing)
- backing pump with pumping speed equal or above 5 m<sup>3</sup>/h

Compression Ratio	ISO 100
H <sub>2</sub>	2 x 10 <sup>4</sup>
He	1 x 10 <sup>5</sup>
N <sub>2</sub>	2 x 10 <sup>8</sup>

### Max Foreline Pressure Tolerance

N <sub>a</sub>	16 mbar	

Note:

Foreline tolerance defined as the pressure at which the turbo pump still produces a compression of 100.

For continuous operation, water cooling is recommended (water temperature between 15 and 20 °C).

Base Pressure with	<1 x 10 <sup>-10</sup> mbar
Recommended Forepump	(<1 x 10 <sup>-10</sup> Torr)


According to standard DIN 28 428, the base pressure is that measured in a leak-free test dome, 48 hours after the completion of test dome bake-out, with a turbo pump fitted with a ConFlat flange and using the recommended prevaculum pump.

Inlet Flange	TwisTorr 305 FSQ and 305-ICQ: ISO 100, CFF 6", ISO 160, CFF 8"
Foreline Flange	KF16 NW (KF25 - optional)
Max Rotation Speed	60.600 rpm (1010 Hz driving frequency)
Start-Up Time	< 3 minutes (longer when soft start is used)
Recommended Forepump	Dry pumps: IDP-3 (no gas flow), IDP-7, IDP-10, Mechanical: DS102, DS302
Operating Position	Any
Oper. Ambient Temperature	+5 to +35 °C
Daksasit Tampavatuva	ISO flange: 75 °C at inlet flange max CFF flange: 100 °C at inlet flange max
Bakeout Temperature	Note: Measure a point close to the sealing element.
Lubricant	Permanent lubrication
Cooling Requirements:	
Air Cooling	Natural convection (only with no gas load) Forced air (5 - 35 °C ambient temperature)
Water Cooling	Minimum flow: 50 L/h (0.22 GPM) Maximum flow: 150 L/h (0.66 GPM) Temperature: +15 to +30 °C Max pressure: 5 bar (75 psi)
Noise Pressure Level (at 1m at full speed)	41 dB(A)


Note: mean values based on a significative sample (Ar and  $\rm N_2$  compression ratio estimated); standard deviation per test: pumping speed: below  $\pm\,7\%$ ; noise pressure level:  $\pm\,10\%$  (only pump).

Installation Category	II		
Pollution Degree	2	2	
Storage Temperature	-40 to +70 °C	-40 to +70 °C	
Weight kg (lbs)	305 FSQ	305-ICQ	
Pump ISO 100 K Pump CFF 6" Pump ISO 160 K Pump CFF 8"	5.84 (12.87) 8.16 (17.98) 6.28 (13.84) 10.43 (22.99)	5.74 (12.65) 8.06 (17.76) 6.18 (13.62) 10.33 (22.77)	

Remote Controller: Voltage Frequency Power Fuse	100 - 240 Vac (voltage fluctuation +/- 10%) 50 to 60 Hz 450 VA 2 x T4 A (slow blow) 250 V	
Power Supply (24 Vdc) for IC Controller.		
Max input power: Pump stand-by average power: Pump max operating power:	300 VA 10 W 150 W	
Max Operating Altitude	3000 m	
Max Magnetic Field	<ul> <li>50 Gauss (5 mT)         <ul> <li>in the transveral direction</li> </ul> </li> <li>100 Gauss (10 mT)         <ul> <li>in the axial direction</li> </ul> </li> </ul>	
Compliance With:	EN 61010-1 EN 61326-1 EN 1012-2 EN 12100 EN 50581 Machinery Directive 2006/42/EC Electromagnetic Compatibility Directive 2014/30/EU Directive 2011/65/EU	



TwisTorr 305 FSQ, TwisTorr 305-ICQ Pumping speed



Foreline Pressure (mbar)

Pumps	Cooling	Flange	Part Number
TwisTorr 305-IC, 485A	Air	ISO 100 K	X3513-64000
TwisTorr 305-IC, 485A	Air	CFF 6"	X3513-64001
TwisTorr 305 IC, 485A	Air	ISO 160 K	X3513-64002
TwisTorr 305 IC, 485A	Air	CFF 8"	X3513-64003
TwisTorr 305 IC, 485A	Water	ISO 100 K	X3513-64004
TwisTorr 305 IC, 485A	Water	CFF 6"	X3513-64005
TwisTorr 305 IC, 485A	Water	ISO 160 K	X3513-64006
TwisTorr 305 IC, 485A	Water	CFF 8"	X3513-64007
TwisTorr 305-IC, 485P	Air	ISO 100 K	X3513-64016
TwisTorr 305-IC, 485P	Air	CFF 6"	X3513-64017
TwisTorr 305-IC, 485P	Air	ISO 160 K	X3513-64018
TwisTorr 305-IC, 485P	Air	CFF 8"	X3513-64019
TwisTorr 305-IC, 485P	Water	ISO 100 K	X3513-64020
TwisTorr 305-IC, 485P	Water	CFF 6"	X3513-64021
TwisTorr 305-IC, 485P	Water	ISO 160 K	X3513-64022
TwisTorr 305-IC, 485P	Water	CFF 8"	X3513-64023
TwisTorr 305-IC, 232	Air	ISO 100 K	X3513-64024
TwisTorr 305-IC, 232	Air	CFF 6"	X3513-64025
TwisTorr 305-IC, 232	Air	ISO 160 K	X3513-64026
TwisTorr 305-IC, 232	Air	CFF 8"	X3513-64027
TwisTorr 305-IC, 232	W <b>a</b> ter	ISO 100 K	X3513-64028
TwisTorr 305-IC, 232	W <b>a</b> ter	CFF 6"	X3513-64029
TwisTorr 305-IC, 232	W <b>a</b> ter	ISO 160 K	X3513-64030
TwisTorr 305-IC, 232	W <b>a</b> ter	CFF 8"	X3513-64031
TwisTorr 305 FS	Air	ISO 100 K	X3513-64008
TwisTorr 305 FS	Air	CFF 6"	X3513-64009
TwisTorr 305 FS	Air	ISO 160 K	X3513-64010
TwisTorr 305 FS	Air	CFF 8"	X3513-64011
TwisTorr 305 FS	Water	ISO 100 K	X3513-64012
TwisTorr 305 FS	Water	CFF 6"	X3513-64013
TwisTorr 305 FS	Water	ISO 160 K	X3513-64014
TwisTorr 305 FS	Water	CFF 8"	X3513-64015
Control	lers		
TwisTorr 305 FS Remote Controller 232-485			X3506-64130
TwisTorr 305 FS Remote Controller Profibus			X3506-64131
Cables			
Mains cable NEMA plug, 3 m long *			9699958
Mains cable European plug, 3 m long *			9699957
Mains cable China plug, 3 m long *			8121-0723
5 m Turbopump Extension Cable *			969-9942M007
10 m Turbopump Extension Cable *			969-9942M006
15 m Turbopump Extension Cable *			969-9942M005
20 m Turbopump Extension Cable *			969-9942M004
50 m Turbopump Extension Cable *			969-9942M015
5 m Turbopump Fan Extension Cable **			9699949

Inlet Screens	Part Number
Inlet Screen ISO 100 K	X3500-68000
Inlet Screen CFF 6"	9699302
Inlet Screen ISO 160 K	X3500-68001
Inlet Screen CFF 8"	9699304
Cooling	
Water Cooling Kit	9699337
Metric Water Kit 4 x 6 mm	9699347
Air cooling kit for TwisTorr 305-IC ** (Kit X3514-68001 is required)	X3500-68010
Air cooling kit for TwisTorr 305 Remote controller *	X3500-68011
Fan extension cable for Remote Controller *	9699940
5 m Vent Valve Extension cable *	9699941
Vibration isolators	
Vibration isolator ISO 100 K	9699344
Vibration isolator CFF 6"	9699334
Vibration isolator ISO 160 K	9699345
Vibration isolator CFF 8"	9699335
Venting	
Vent Valve N.O. 1, 2 mm for TwisTorr 305-IC ** (Kit X3514-68001 is required)	9699834
Vent Valve N.O. 0,5 mm for TwisTorr 305-IC ** (Kit X3514-68001 is required)	9699834M006
DB15 Mating Connector not wired 7.5A **	X3514-68000
TwisTorr 305-IC Fan/Vent Adapter kit **	X3514-68001
Vent Valve N.O. 0,5 mm Orifice *	9699844
Vent Valve N.O. 1.2 mm Orifice *	9699845
Vent Valve N.C. 1.2 mm Orifice *	9699846
Vent Valve N.C. 0,5 mm Orifice *	9699847
Purge	
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 10 SCCM ¼ Swagelock - M12	9699240
Purge valve 20 SCCM NW16KF - M12	9699241
Purge valve 20 SCCM ¼ Swagelock - M12	9699242
Purge valve 10 SCCM ¼ Swagelock - ¼ Swagelock 9699232	
Purge valve 20 SCCM ¼ Swagelock - ¼ Swagelock	9699236
Other accessories	
Serial to Bluetooth adapter (necessary for App) *	X3514-68003
KF25 Foreline flange	X3513-68000

<sup>\*</sup> For TwisTorr 305 FS

<sup>\*\*</sup> For TwisTorr 305-IC

Pumps	Cooling	Flange	Part Number
TwisTorr 305 FSQ	Air/Water	ISO100K	X3513-64068
TwisTorr 305-ICQ, 485A	Water	ISO100K	X3513-64060
TwisTorr 305-ICQ, 485A	Air	ISO100K	X3513-64061
TwisTorr 305-ICQ, 485A	Water	CFF6"	X3513-64062
TwisTorr 305-ICQ, 485A	Air	CFF6"	X3513-64063
TwisTorr 305-ICQ, 485A	Water	ISO160K	X3513-64064
TwisTorr 305-ICQ, 485A	Water	CFF8"	X3513-64065
Controllers			
TwisTorr 305 FS Remote Con	troller 232-48	35	X3506-64130
TwisTorr 305 FS Remote Con	troller Profib	us	X3506-64131
Cables			•
Mains cable NEMA plug, 3 m	long *		9699958
Mains cable European plug, 3	3 m long *		9699957
Mains cable China plug, 3 m	long *8		121-0723
5 m Turbopump Extension Ca	able *		969-9942M007
10 m Turbopump Extension (	Cable *		969-9942M006
15 m Turbopump Extension (	Cable *		969-9942M005
20 m Turbopump Extension (	Cable *		969-9942M004
50 m Turbopump Extension (	Cable *		969-9942M015
5 m Turbopump Fan Extensio	on Cable **		9699949
Inlet Screens			Part Number
Inlet Screen ISO 100 K			X3500-68000
Inlet Screen CFF 6"			9699302
Inlet Screen ISO 160 K			X3500-68001
Inlet Screen CFF 8"			9699304
Cooling			
Water Cooling Kit			9699337
Metric Water Kit 4 x 6 mm			9699347


<sup>\*</sup> For TwisTorr 305 FSQ

Air cooling kit for TwisTorr 305-IC models ** (Kit X3514-68001 is required)	X3500-68010
Air cooling kit for TwisTorr 305 Remote controller * X3500-68011	
Fan extension cable for Remote Controller *	9699940
5 m Vent Valve Extension cable *	9699941
Vibration isolators	
Vibration isolator ISO 100 K	9699344
Vibration isolator CFF 6"	9699334
Vibration isolator ISO 160 K	9699345
Vibration isolator CFF 8"	9699335
Venting	
Vent Valve N.O. 1, 2 mm for TwisTorr 305-IC models ** (Kit X3514-68001 is required)	9699834
Vent Valve N.O. 0,5 mm for TwisTorr 305-IC models ** (Kit X3514-68001 is required)	9699834M006
DB15 Mating Connector not wired 7.5A **	X3514-68000
TwisTorr 305-IC Fan/Vent Adapter kit ** X3514-6800	
Vent Valve N.O. 0,5 mm Orifice *	9699844
Vent Valve N.O. 1.2 mm Orifice *	9699845
Vent Valve N.C. 1.2 mm Orifice *	9699846
Vent Valve N.C. 0,5 mm Orifice *	9699847
Purge	
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 10 SCCM ¼ Swagelock - M12	9699240
Purge valve 20 SCCM NW16KF - M12 9699241	
Purge valve 20 SCCM ¼ Swagelock - M12 9699242	
Purge valve 10 SCCM ¼ Swagelock - ¼ Swagelock 9699232	
Purge valve 20 SCCM ¼ Swagelock - ¼ Swagelock 9699236	
Other accessories	
Serial to Bluetooth adapter (necessary for App) *	X3514-68003
KF25 Foreline flange	X3513-68000

<sup>\*\*</sup> For TwisTorr 305-ICQ



TwisTorr 404 FS



Dimensions: millimeters [inches]

### **Technical Specification**

Pumping Speed	IS0100K-F	CFF6		
He H <sub>2</sub> Ar	355 l/s 470 l/s 445 l/s 320 l/s			
Max Gas Throughput(*)	Air cooling (25 °C Air temperature)	Water cooling (15 °C water temperature / 25 °C room temperature)		
$N_2$	4.3 mbarL/s 255 SCCM	6.2 mbarL/s 367 SCCM		
Не	7.9 mbarL/s 467 SCCM	10.4 mbarL/s 615 SCCM		
Ar	1.5 mbarL/s 89 SCCM	3.3 mbarL/s 195 SCCM		
Compression Ratio and Fore	line Tolerance			
N <sub>2</sub> He H <sub>2</sub> Ar	>1 x 10 <sup>11</sup> 2 x 10 <sup>8</sup> 3 x 10 <sup>6</sup> > 1 x 10 <sup>11</sup>	>10 mbar >10 mbar >4 mbar >8.5 mbar		
(*) Foreline Tolerance defined as the pressure at which the turbopump still produce a compression of 100 and estimated in water cooling mode.				
Inlet Flange	ISO 100K, ISO 100F, CFF 6"			
Foreline Flange	NW25 (NW16 as optional accessory)			
Rotational Speed	Auto setting from 40.800 RPM to 49.500 RPM			
Start-Up Time	< 5 minutes			


Recommended Forepump	Mechanical: Agilent DS 302, TS 300 Dry pump: Agilent IDP-10, IDP-15		
Operating Position	Any		
Oper. Ambient Temperature	+5 to +35 °C		
Rel. Humidity of Air	0 to 90% (not condensir	ng)	
Bakeout Temperature	ISO pump: 80 °C at inlet flange CFF pump: 120 °C at inlet flange		
Lubricant	Permanent grease lubri	cation	
Cooling Requirements			
Air Cooling	Air temperature from +5 to 35 °C		
Water Cooling	Water temperature from +15 to +25 °C Water flow min. 100l/h		
Noise Pressure Level (at 1 m at full speed)	43dB(A)		
Storage Temperature	-40 to +70 °C		
Max Altitude	3000 m		
Weight kg (lbs)	ISO100K ISO100F CFF 6"	20.6 Kg / 45.3 lbs. 22.1 Kg / 48.6 lbs. 22 Kg / 48.4 lbs.	
Conformity to Norms			
EMC (Control Units) Safety (CE/CSA) Machinery Directive Low Voltage Directive EMC Directive (Control Units) ROHS		61326-1 61010-1 DIR 2006/42/CE DIR 2014/35/EU DIR 2014/30/EU DIR 2011/65/EU	

- Hydrogen
- Helium
- Nitrogen
- Argon

350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 300
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 350
- 35

Inlet Pressure (mbar)

Pumping speed vs inlet pressure




Compression ratio vs foreline pressure

Pumps	Part Number
TwisTorr 404 FS ISO100K KF25 water cooling	X3511-64019
TwisTorr 404 FS ISO100F KF25 water cooling	X3511-64020
TwisTorr 404 FS CFF6" KF25 water cooling	X3511-64021
Controllers	
TwisTorr Medium-TMP rack controller	X3501-64016
TwisTorr Medium-TMP on board controller	X3512-64016
Cables	
Mains cable NEMA plug 3 m long	9699958
Mains cable EU plug 3 m long	9699957
Mains cable UK plug 2.5 m long	X3501-68005
Mains cable China plug 3 m long	X1699-64144
Replace with: RS232 serial cable and A-PLUS software, 3 m long	9699883
Extension cables*	9699948 (3 m) 9699948M001 (5 m) 9699948M003 (10 m) 9699948M004 (15 m) 9699948M002 (20 m) 9699948M005 (30 m) 9699948M009 (75 m) 9699948M010 (100 m)
Inlet Screen	
Inlet screen iso100/cff6	9699302
Cooling	
Medium-TMP air cooling kit for rack controller	X3501-68001
Medium-TMP air cooling kit for on board controller	9699297
Plastic water cooling kit	9699347
Metal water cooling kit	9699337
Air cooling kit extension cable 5m	x3501-68101
Air cooling kit extension cable 10m	x3501-68051
Air cooling kit extension cable 15m	X3501-68061
Air cooling kit extension cable 20m	X3501-68021

Air cooling kit extension cable 30m	X3501-68011
Air cooling kit extension cable 50m	X3501-68071
Air cooling kit extension cable 75m	X3501-68081
Air cooling kit extension cable 100m	X3501-68091
Venting	
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002
Vent valve kit, normally closed, with 5 m cable, for pumps driven by rack controllers	X3501-68022
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834
Vent flange, NW 10 KF / M8	9699108
Vent valve rack extension cable	X3501-68004 (5 m) X3501-68054 (10 m) X3501-68064 (15 m) X3501-68074 (20 m) X3501-68084 (30 m) X3501-68034 (50 m) X3501-68094 (75 m) X3501-68104 (100 m)
Purge	
Purge valve 10 SCCM NW16KF - M12	9699239
Purge 10SCCM M12-1/4 Swagelok	9699240
Purge 20SCCM M12-NW16KF	9699241
Purge valve 20 SCCM ¼ Swagelok - M12	9699242
Spare purge screw	X3502-68004
Mounting	
Change description of X3511-68003 to: Medium TMP onboard controller side mounting bracket	X3511-68003
Medium-TMP KF16 foreline flange	X3511-68004
Medium-TMP KF25 foreline flange	X3511-68001
Medium-TMP KF40 foreline flange	X3511-68002




### Turbo V-551 Navigator



### **Technical Specification**

Pumping Speed (with inlet screen)			
CFF 6"	N2: 350 l/s	He: 450 l/s	H <sub>2</sub> : 450 l/s
CFF 8" or ISO 160	N <sub>2</sub> : 550 l/s	He: 600 l/s	H <sub>2</sub> : 510 l/s
Compression Ratio	N <sub>2</sub> : 1x10 <sup>9</sup>	He: 1x10 <sup>7</sup>	H <sub>2</sub> : 1x10 <sup>6</sup>
Base Pressure* (with minimum recommended forepump)	<1 x 10 -10 mbar (< 1 x 10 -10 Torr)		
* According to PNEUROP 5608.			
Inlet Flange	CFF 6": ISO 160 CFF 8": ISO 160-F bolted		
Foreline Flange	KF 25		
Rotational Speed	42.000 rpm		
Start-Up Time	<5 minutes		
Recommended Forepump		gilent DS 302 lent TS300, IDI	P-15

Operating Position	Any		
Cooling Requirements	Natural air convection Water optional (use water with electrical conductivity ≤ 500 µS/cm)		
Bakeout Temperature	120 °C at inlet flange max. (CF flange) 80 °C at inlet flange max. (ISO flange)		
Vibration Level (displacement)	<0.01 µm at inlet flange		
Weight kg (lbs)	ISO flange 19.4 (43.0) CFF flange 23.4 (51.6) CFF and ISO bolted flange 23.4 (51.6)		
SEM version available on request			



Compression ratio vs foreline pressure


### Pumping speed vs inlet pressure

Description	Part Number	
Pumps		
Turbo-V 551 Navigator pump, ISO 160 flange	9698922	
Turbo-V 551 Navigator pump, ISO 160-F bolted flange	9698944	
Turbo-V 551 Navigator pump, 6" CFF flange	9698925	
Turbo-V 551 Navigator pump, 8" CFF flange	9698923	
Controllers		
Turbo-V 551 Navigator controller 120/220 V - 50/60 Hz	9698976	
Turbo-V 551 Rack controller, 100-240 V	X3501-64001	
Accessories		
Mains cable NEMA plug, 3 m long	9699958	
Mains cable European plug, 3 m long	9699957	
Mains cable UK plug 2.5 m long	X3501-68005	
Mains cable China plug 3 m long	X1699-64144	
RS232 serial cable and A-PLUS software, 3 m long	9699883	
Extension cable	9699948 (3 m) 9699948M001 (5 m) 9699948M002 (20 m) 9699948M003 (10 m) 9699948M004 (15 m) 9699948M005 (30 m) 9699948M011 (65 m) 9699948M009 (75 m) 9699948M010 (100 m)	
Inlet screen DN 100	9699302	
Inlet screen DN 160	9699304	

Water cooling kit	9699337
Plastic water cooling kit	9699347
Air cooling kit for use with Navigator controller	9699339
Air cooling kit for use with standard rack controller, 24 V	X3501-68001
Bracket for Navigator controller side mouting	9699349
Vibration damper CFF 6"	9699334
Vibration damper CFF 8"	9699335
Vent valve adapter, NW10 to M8 male threads	9699108
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 10 SCCM ¼ Swagelok - M12	9699240



Turbo-V 701 Navigator





Dimensions: millimeters [inches]

### **Technical Specification**

recimical opecinication				
Pumping Speed (with inlet screen) CFF 10" or ISO 200		N <sub>2</sub> : 690 l/s	He: 620 l/s	H <sub>2</sub> : 510 l/s
Compression Ratio		N <sub>2</sub> : 1 x 10 <sup>9</sup>	He: 1 x 10 <sup>7</sup>	H <sub>2</sub> : 1 x 10 <sup>6</sup>
Base Pressure* (with minimum recommended forepump)	<1 x 10 <sup>-10</sup> r	nbar (< 1 x 10	<sup>-10</sup> Torr)	
Inlet Flange		CFF 10"	ISO 200	ISO 200-F bolted
Foreline Flange	KF 25			
Rotational Speed	42.000 rpm	า		
Start-Up Time	<5 minutes			
Recommended Forepump	Mechanical: Agilent DS 402 Dry scroll: Agilent TS300, IDP-15			
Operating Position	Any			

Cooling Requirements	Natural air convection Water optional (use water with electrical conductivity ≤ 500 µS/cm)	
Bakeout Temperature	120 °C at inlet flange max. (CF flange) 80 °C at inlet flange max. (ISO flange)	
Vibration Level (displacement)	<0.01 μm at inlet flange	
Weight kg (lbs)	ISO flange 19.4 (43.0) CFF flange 25.5 (54.2) CFF and ISO bolted flange 25.5 (54.2	
*According to PNEUROP 5608.		
SEM version available on request		





Pumping speed vs inlet pressure


Compression ratio vs foreline pressure

Description	Part Number		
Pumps			
Turbo-V 701 Navigator pump, ISO 200-K	9698926		
Turbo-V 701 Navigator pump, 10" CFF	9698927		
Controllers			
Turbo-V 701 Navigator controller 120/220 V - 50/60 Hz	9698977		
Turbo-V 701 rack controller, 120/220 V	X3501-64002		
Accessories			
Mains cable NEMA plug, 3 m long	9699958		
Mains cable European plug, 3 m long	9699957		
Mains cable UK plug, 2.5 m long	X3501-68005		
Mains cable China plug, 3 m long	X1699-64144		
RS232 serial cable and A-PLUS software, 3 m long	9699883		
Extension cable	9699948 (3 m) 9699948M001 (5 m) 9699948M002 (20 m) 9699948M003 (10 m) 9699948M004 (15 m) 9699948M005 (30 m) 9699948M009 (75 m) 9699948M010 (100 m) 9699948M011 (65 m)		
Inlet screen DN 200	9699316		
Vibration damper, ISO 200	9699346		
Vibration damper, CFF 10"	9699336		
Water cooling kit	9699337		
Plastic water cooling kit	9699347		
Air cooling kit for use with Navigator controller	9699339		

Air cooling kit for use with standard rack controller, 24 V	X3501-68001
Bracket for Navigator controller side mounting	9699349
Vent valve adapter, NW10 to M8 male threads	9699108
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834
Purge valve 10 SCCM NW16KF - M12	9699239
Purge valve 10 SCCM ¼ Swagelok – M12	9699240



TwisTorr 704 FS




Dimensions: millimeters [inches]

### **Technical Specification**

Pumping Speed	ISO 160 / CFF 8"			
N <sub>2</sub> He H <sub>2</sub> Ar	660 L/s 640 L/s 480 L/s 625 L/s			
Gas throughput at Full Rotational Speed (with recomm. forepump)*	Ambient temp. (25 °C)	Water temperature. (25 °C, 50 L/h)		
He	7.9 mbarL/s 467 SCCM	10.4 mbarL/s 615 SCCM		
N2	4.3 mbarL/s 255 SCCM	6.2 mbarL/s 367 SCCM		
Ar	1.5 mbarL/s 89 SCCM	3.3 mbarL/s 195 SCCM		
(*) Backing pump 11.6 m3/hr				
Compression Ratio and Foreline Tolerance*				
$egin{array}{l} N_2 \\ \text{He} \\ H_2 \\ \text{Ar} \end{array}$	>1 x 10 <sup>11</sup> 2 x 10 <sup>8</sup> 3 x 10 <sup>6</sup> > 1 x 10 <sup>11</sup>	10 mbar 10 mbar >4 mbar 8.5 mbar		
(*) Foreline Tolerance defined as the pressure at which the turbopump still produce a compression of 100 and estimated in water cooling mode				
Base Pressure with Recommended Forepump	< 1 x 10 <sup>-10</sup> mbar (< 1 x 10 <sup>-10</sup> Torr)			
Inlet Flange	ISO 160K, ISO 160F, CFF 8"			
Foreline Flange	NW25 (NW40 as option	NW25 (NW40 as optional accessory)		
Rotational Speed	Auto setting from 40.800 RPM to 49.500 RPM			
Start-Up Time	< 5 minutes			

Recommended Forepump	Mechanical: Agilent DS 302, TS 300 Dry pump: Agilent IDP-10, IDP-15		
Operating Position	Any		
Oper. Ambient Temperature	+5 to +35 °C		
Rel. Humidity of Air	0 - 90 % (not condensing)		
Bakeout Temp.	ISO pump: 80 °C at inlet flange CFF pump: 120 °C at inlet flange		
Lubricant	Permanent grease lubrication		
Cooling Requirements			
Air Cooling	Air temperature from +5 to 35 °C		
Water Cooling	Water temperature from +15 to +25 °C Water flow min. 100l/h		
Noise Pressure Level (at 1m at full speed)	43dB(A)		
Storage Temperature	-40 to +70 °C		
Max Altitude	3000 m		
Weight kg (lbs)	ISO160K ISO160F CFF 8"	20.6 Kg / 45.3 lbs. 22.6 Kg / 49.7 lbs. 22 Kg / 48.4 lbs.	
Conformity to Norms			
EMC (Control Units) Safety (CE/CSA) Machinery Directive Low Voltage Directive EMC Directive (Control Units) ROHS		61326-1 61010-1 DIR 2006/42/CE DIR 2014/35/EU DIR 2014/30/EU DIR 2011/65/EU	



## optimization of the state of t


### Compression ratio vs foreline pressure


Pumps	Part Number
TwisTorr 704 FS ISO 160 K, KF 25 water cooling	X3511-64015
TwisTorr 704 FS ISO160F, KF25 water cooling	X3511-64016
TwisTorr 704 FS CFF8", KF25 water cooling	X3511-64017
TwisTorr 704 FS CFF8", Long-Neck KF25 water cooling	X3511-64018
TwisTorr 704 FS, One Speed, ISO 160 K, KF 25	X3511-64040
TwisTorr 704 FS, One Speed, ISO 160 F, KF 25	X3511-64041
TwisTorr 704 FS, One Speed, CFF 8", KF 25	X3511-64042
Controllers	
TwisTorr Medium-TMP rack controller	X3501-64016
TwisTorr Medium-TMP onboard controller	X3512-64016
TwisTorr Medium-TMP onboard controller, One Speed	X3512-64006
TwisTorr Medium-TMP onboard rugged controller	X3512-64026
Cables	
Mains cable NEMA plug 3 m long	9699958
Mains cable EU plug 3 m long	9699957
Mains cable UK plug 2.5 m long	X3501-68005
Mains cable China plug 3 m long	X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	9699883
Extension cables*	9699948 (3 m) 9699948M001 (5 m) 9699948M003 (10 m) 9699948M004 (15 m) 9699948M002 (20 m) 9699948M005 (30 m) 9699948M009 (75 m) 9699948M010 (100 m)
Inlet Screen	
ISO 160/ CFF 8	9699304
Mounting	
Medium TMP onboard controller side mounting bracket	X3511-68003
Medium-TMP KF16 foreline flange	X3511-68004
Medium-TMP KF25 foreline flange	X3511-68001


Medium-TMP KF40 foreline flange	X3511-68002
Cooling	
Medium-TMP air cooling kit for rack controller	X3501-68001
Medium-TMP air cooling kit for onboard controller	9699297
Plastic water cooling kit	9699347
Metal Water cooling Kit	9699337
	X3501-68101 (5 m)
	X3501-68051 (10 m)
	X3501-68061 (15 m)
Air agaling kit autopaian aghla	X3501-68021 (20 m)
Air cooling kit extension cable	X3501-68011 (30 m)
	X3501-68071 (50 m)
	X3501-68081 (75 m)
	X3501-68091 (100 m)
Venting	
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002
Vent valve kit, normally closed, with 5 m cable, for pumps driven by rack controllers	X3501-68022
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834
Vent flange, NW 10 KF / M8	9699108
	X3501-68004 (5 m)
	X3501-68054 (10 m)
	X3501-68064 (15 m)
	X3501-68074 (20 m)
Vent valve rack extension cable	X3501-68084 (30 m)
	X3501-68034 (50 m)
	X3501-68094 (75 m)
	X3501-68104 (100 m)
Purge	
Purge valve 10 SCCM NW16KF - M12	9699239
Purge 10SCCM M12-1/4 Swagelok	9699240
Purge 20SCCM M12-NW16KF	9699241
Purge valve 20 SCCM ¼ Swagelok - M12	9699242
Spare purge screw	X3502-68004
*For a complete list of available extension cables, see the Exte	ension Cables

<sup>\*</sup>For a complete list of available extension cables, see the Extension Cables for Turbo Pumps dedicated page on Agilent.com.










TwisTorr 804 FS

Dimensions: millimeters [inches]

### **Technical Specification**

Pumping Speed	IS0200K	ISO250K	CFF10	
	ISO200F	IS0250F		
$N_2$		720 L/s		
He		660 L/s		
H <sub>2</sub>		485 L/s		
Ar		690 L/s		
Max Gas	Air cooling	Water cooling		
Throughput(*)	(25 °C Air	(15 °C water tem	perature	
	temperature)	/25 °C room tem	perature)	
N <sub>2</sub>	4.3 mbarL/s	6.1 mbar L/s		
2	255 SCCM	367 SCCM		
He	7.9 mbar L/s	10.4 mbar L/s		
	467 SCCM	615 SCCM		
Ar	1.5 mbar L/s	3.3 mbar L/s		
	89 SCCM	195 SCCM		
(*) Backing pump 11.6 m3/hr				
Compression Ratio and Fore	line Tolerance*			
N <sub>2</sub>	>1E+11	10 mbar		
He	2E+08	10 mbar		
H <sub>2</sub>	3E+06	>4 mbar		
Ar	>1E+11	8.5 mbar		
(*) Foreline Tolerance defined as the	e pressure at which the t	urbopump still produc	e a	
compression of 100 and estimated	in water cooling mode			
Base Pressure with	< 1 x 10 <sup>-10</sup> mbar			
Recommended Forepump	(< 1 x 10 <sup>-10</sup> Torr)			
	100 0001/ 100 000	E 100 0E01/		
Inlet Flange	ISO 200K, ISO 200	r, 15U 25UK,		
	ISO 250F, CFF 10"			
Foreline Flange	NW25 or NW40	NW25 or NW40		
Rotational Speed	Auto setting from 4	10.800 RPM to 49'500	RPM	
Start-Up Time	< 5 minutes			


	<del>.</del>		
Recommended Forepump	Mechanical: Agilent DS 302, TS 300 Dry pump: Agilent IDP-10, IDP-15		
Operating Position	Any		
Oper. Ambient Temperature	+5 to +35 °C		
Rel. Humidity of Air	0 to 90% (not con	densing)	
Bakeout Temperature	ISO pump: 80 °C a CFF pump: 120 °C	2	
Lubricant	Permanent grease	e lubrication	
Cooling Requirements			
Air Cooling	Air temperature fr	rom +5 to 35 °C	
Water Cooling	Water temperature from +15 to +25 °C Water flow min. 100I/h		
Noise Pressure Level (at 1m at full speed)	43dB(A)		
Storage Temperature	-40 to +70 °C		
Max Altitude	3000 m		
Weight kg (lbs)	ISO200K ISO200F ISO250K ISO250F CFF 10"	20.7 Kg / 45.5 lbs. 23.6 Kg / 51.9 lbs. 23.3 Kg / 51.2 lbs. 27.6 Kg / 60.9 lbs. 22.1 Kg / 48.6 lbs.	
Conformity to Norms			
EMC (Control Units) Safety (CE/CSA) Machinery Directive Low Voltage Directive EMC Directive (Control Units) ROHS		61326-1 61010-1 DIR 2006/42/CE DIR 2014/35/EU DIR 2014/30/EU DIR 2011/65/EU	



Pumping speed vs inlet pressure

<u> </u>	•			
Pumps	Part Number			
TwisTorr 804 FS ISO200K KF25	X3511-64007			
TwisTorr 804 FS ISO200K KF40	air cooling	X3511-64008		
TwisTorr 804 FS ISO200F KF25	X3511-64009			
TwisTorr 804 FS ISO200F KF40	X3511-64010			
TwisTorr 804 FS CFF10" KF25 a	ir cooling	X3511-64011		
TwisTorr 804 FS CFF10" KF40 a	ir cooling	X3511-64012		
TwisTorr 804 FS ISO250K KF40	air cooling	X3511-64013		
TwisTorr 804 FS ISO250F KF40	air cooling	X3511-64014		
TwisTorr 804 FS ISO200K KF25	water cooling	X3511-64022		
TwisTorr 804 FS ISO200K KF40	water cooling	X3511-64023		
TwisTorr 804 FS ISO200F KF25	water cooling	X3511-64024		
TwisTorr 804 FS ISO200F KF40	water cooling	X3511-64025		
TwisTorr 804 FS CFF10" KF25 v	vater cooling	X3511-64026		
TwisTorr 804 FS CFF10" KF40 v	vater cooling	X3511-64027		
TwisTorr 804 FS ISO250K KF40	water cooling	X3511-64028		
TwisTorr 804 FS ISO250F KF40	water cooling	X3511-64029		
Controllers				
TwisTorr Medium-TMP rack cor	itroller	X3501-64016		
TwisTorr Medium-TMP on board controller		X3512-64016		
Cables				
Mains cable NEMA plug 3 m lon	g	9699958		
Mains cable EU plug 3 m long		9699957		
Mains cable UK plug 2.5 m long		X3501-68005		
Mains cable China plug 3 m long	9	X1699-64144		
RS232 serial cable and A-PLUS	software, 3 m long	9699883		
Extension cables*				
hlu Ooran	9699948 (3 m) 9699948M001 (5 m) 9699948M003 (10 m) 9699948M004 (15 m)	9699948M002 (20 m) 9699948M005 (30 m) 9699948M009 (75 m) 9699948M010 (100 m)		
Inlet Screen				
ISO 200/ CFF 10		9699316		
ISO 250		9699350		
Cooling				
Medium-TMP air cooling kit for	rack controller	X3501-68001		

<sup>\*</sup>For a complete list of available extension cables, see the Extension Cables for Turbo Pumps dedicated page on Agilent.com.




Compression ratio vs foreline pressure

Medium-TMP Air Cooling Kit for on boa	9699297	
Plastic Water cooling Kit	9699347	
Metal Water cooling Kit	9699337	
Air cooling kit extension cable 5m	X3501-68101	
Air cooling kit extension cable 10m		X3501-68051
Air cooling kit extension cable 15m		X3501-68061
Air cooling kit extension cable 20m		X3501-68021
Air cooling kit extension cable 30m		X3501-68011
Air cooling kit extension cable 50m		X3501-68071
Air cooling kit extension cable 75m		X3501-68081
Air cooling kit extension cable 100m		X3501-68091
Venting		
Vent valve kit, normally open, with 5 m of for pumps driven by rack controllers	cable,	X3501-68002
Vent valve kit, normally closed, with 5 m for pumps driven by rack controllers	X3501-68022	
Vent valve kit, normally open, with 0.6 n for pumps driven by onboard controller		9699834
Vent flange, NW 10 KF / M8		9699108
Vent valve rack extension cable	X3501-68004 (5 r X3501-68054 (10 X3501-68064 (15 X3501-68074 (20	m) X3501-68054 (10 m m) X3501-68064 (15 m
Purge		
Purge valve 10 SCCM NW16KF - M12		9699239
Purge 10SCCM M12-1/4 Swagelok		9699240
Purge 20SCCM M12-NW16KF		9699241
Purge valve 20 SCCM ¼ Swagelok - M1	9699242	
Spare purge screw	X3502-68004	
Mounting		
Medium TMP onboard controller side m	nounting bracket	X3511-68003
Medium-TMP KF16 foreline flange	X3511-68004	
Medium-TMP KF25 foreline flange	X3511-68001	
Medium-TMP KF40 foreline flange	X3511-68002	



### Turbo-V 1K-G



### Dimensions: millimeters [inches]

### **Technical Specification**


Pumping Speed					
ISO 160	Ar: 750 l/s	N <sub>2</sub> : 810 l/s	He: 950 l/s	H <sub>2</sub> : 680 l/s	
ISO 200	Ar: 1,040 l/s	N <sub>2</sub> : 1,080 l/s	He: 1,150 l/s	H <sub>2</sub> : 730 l/s	
Compression Ratio	Ar: > 5 x 108	$N_2$ : > 5 x 10 <sup>7</sup>	He: > 4 x 10 <sup>4</sup>	H <sub>2</sub> : 1.5 x 10 <sup>4</sup>	
Base Pressure* (with minimum recommended forepump)	< 1 x 10 <sup>-10</sup> mbar				
Inlet Flange	ISO 160 F, ISO 2	200 F			
Foreline Flange	ISO 160 F: KF:				
Nominal Rotational Speed	45.500 rpm				
Start-Up Time	<5 minutes				
Minimum Recommended Forepump	> 20 m³/h (TriSe	> 20 m³/h (TriScroll 600, DS 602)			
Operating Position	Any				
Operating Ambient Temperature	+5 to +35 °C				
Cooling Requirements	Water				
Bakeout Temperature	80 °C at inlet flange max. (ISO flange)				
Vibration Level (displacement)	< 0.01 µm at inl	et flange			
Power Supply Input voltage: Input freq.: Max input power: Stand-by power: Max operating power:	100-240 Vac 50 - 60 Hz 600 VA 30 to 35 W 400 W t				
Protection Fuse (Navigator controller)	1 x 6.3 A				


Serial Communication
(Navigator Kit)


RS232 cable with a 9-pin D type male connector and a 9-pin D type female connector, and A-PLUS software (optional)

Storage Temperature	-20 to +70 °C	
Weight	26.8 kg (59.1 lbs)	

<sup>\*</sup> According to standard DIN 28 428.








Compression ratio vs foreline pressure

Description	Part Number	
Pumping Systems		
Turbo-V 1K-G ISO160 F pump	8698961R002	
Turbo-V 1K-G ISO200 F pump	8698962R001	
Controllers		
Turbo-V 1K-G Navigator controller, 120-220 V	9698978M005	
Turbo-V 1K-G rack controller, 100-240 V	X3501-64005	
Accessories		
Mains cable NEMA plug, 3 m long	9699958	
Mains cable European plug, 3 m long	9699957	
Mains cable UK plug, 2.5 m long	X3501-68005	
Mains cable China plug, 3 m long	X1699-64144	
RS232 serial cable and A-PLUS software, 3 m long	9699883	
Inlet screen ISO 160	9699304	
Inlet screen ISO 200	9699316	
Water cooling kit (hose tail G 1/4)	9699825	
Water cooling kit (Inox G ¼)	9699826	
Vent flange, NW 10 KF / M8	9699108	
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002	
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834	

Purge valve KF16-M12 20 SCCM	9699241
Purge valve 20 SCCM ¼ Swagelok - M12	9699242
Recommended Forepump	
Rotary Vane pump DS 602, with 1 ph, worldwide motor	9499335
Rotary Vane pump DS 602, with 3 ph, worldwide motor	9499336
Dry pump TriScroll 600, with 1 ph, worldwide motor	PTS06001UNIV
Dry pump TriScroll 600, with 3 ph, worldwide motor	PTS06003UNIV
Dry pump TriScroll 600 inverter, with 1 ph, worldwide motor	PTS06001INV



Turbo-V 1001 Navigator



Dimensions: millimeters [inches]


### **Technical Specification**


Pumping Speed	ISO 160 CFF 10" or ISO 200 ISO 250	N <sub>2</sub> : 790 l/s N <sub>2</sub> : 950 l/s N <sub>2</sub> : 1,050 l/s	He: 820 l/s He: 870 l/s He: 900 l/s	H <sub>2</sub> : 860 l/s H <sub>2</sub> : 900 l/s H <sub>2</sub> : 920 l/s	
Compression ratio		N <sub>2</sub> : 1 x 10 <sup>9</sup> N <sub>2</sub> : 1 x 10 <sup>9</sup>	He: 1 x 10 <sup>7</sup> He: 1 x 10 <sup>7</sup>	H <sub>2</sub> : 1 x 10 <sup>6</sup> H <sub>2</sub> : 1 x 10 <sup>6</sup>	
Base Pressure* (with minimum recommended forepump)	<1 x 10 <sup>-10</sup> mbar (<	<1 x 10 <sup>-10</sup> mbar (< 1 x 10 <sup>-10</sup> Torr)			
Inlet Flange	ISO 160, ISO 200-k	(, ISO 200-F bolted	d, CFF 10", ISO 25	50-K, ISO 250-F	
Foreline Flange	KF 40	KF 40			
Rotational Speed	38.000 rpm	38.000 rpm			
Start-Up Time	<4 minutes	<4 minutes			
Recommended Forepump		Mechanical: Agilent DS 402 Dry scroll: Agilent TS300, IDP-15			
Operating Position	Any	Any			
Cooling Requirements	Forced air or water	Forced air or water			
Bakeout Temperature		120 °C at inlet flange max. (CFF flange) 80 °C at inlet flange max. (ISO flange)			
Vibration Level (displacement)	< 0.01 µm at inlet f	< 0.01 µm at inlet flange			
*According to standard DIN 28 4	28				

Weight kg (lbs):	ISO 160 flange ISO 200 flange CFF 10" flange ISO 250 flange	19 (41.8) 19.4 (43.0) 25.5 (54.2) 21.2 (46.6)

Available with On-board Navigator controller, ½ rack controller or PCB controller; for information on controllers see also following pages.

SEM version available on request.






Description	Part Number
Complete Systems	
Turbo-V 1001 Navigator complete system, ISO 160 flange	9698840
Turbo-V 1001 Navigator complete system, ISO 200 flange	9698838
Turbo-V 1001 Navigator complete system, ISO 200-F bolted flange	9698844
Turbo-V 1001 Navigator complete system, 10" CFF flange	9698839
Turbo-V 1001 Navigator complete system, ISO 250 flange	9698841
Turbo-V 1001 Navigator complete system, ISO 250-F bolted flange	9698845
Complete system includes: Pump, inlet screen, controller mounted on the bottom, 2 Mains cables (NEMA plug and European plug).	
Pumps	
Turbo-V 1001 Navigator Pump, ISO 160 flange	9698933
Turbo-V 1001 Navigator Pump, ISO 200 flange	9698931
Turbo-V 1001 Navigator pump, ISO 200-F bolted flange	9698946
Turbo-V 1001 Navigator Pump, 8" CFF flange	9698932M003
Turbo-V 1001 Navigator pump, 10" CFF flange	9698932
Turbo-V 1001 Navigator pump, ISO 250 flange	9698934
Turbo-V 1001 Navigator pump, ISO 250-F bolted flange	9698947
Controllers	
Turbo-V 1001 Navigator controller 120/220 V - 50/60 Hz	9698978
Turbo-V 1001 rack controller, 100-240 V	X3501-64003
Pump Accessories	
Mains cable NEMA plug, 3 m long	9699958
Mains cable European plug, 3 m long	9699957
Mains cable UK plug, 2.5 m long	X3501-68005
Mains cable China plug, 3 m long	X1699-64144

RS232 serial cable and A-PLUS software, 3 m lon	ng 9699883	
Extension cables*	9699948 (3 m) 9699948M001 (5 m) 9699948M002 (20 m) 9699948M003 (10 m) 9699948M005 (30 m) 9699948M005 (75 m) 9699948M010 (100 m) 9699948M011 (65 m)	
Inlet screen DN 160	9699304	
Inlet screen DN 200	9699316	
Inlet screen DN 250	9699350	
Water cooling kit	9699337	
Plastic water cooling kit	9699347	
Air cooling kit for use with Navigator controller	9699297	
Air cooling kit for use with standard rack controller	X3501-68001	
Vibration damper ISO 160	9699345	
Vibration damper, ISO 200	9699346	
Vibration damper, CFF 10"	9699336	
Vent flange, NW 10 KF / M8	9699108	
Vent device with adjustable delay time for standard rack controller	9699831	
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002	
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers	9699834	
Purge valve 10 SCCM NW16KF - M12	9699239	
Purge valve 10 SCCM 1/4 Swagelok - M12	9699240	
Purge valve 20 SCCM NW16KF - M12	9699241	
Purge valve 20 SCCM ¼ Swagelok – M12	9699242	
Purge valve 10 SCCM ¼ Swagelok - 7/16 Swagelok	9699232	
Purge Valve 20 SCCM ¼ Swagelok - 7/16 Swagelok	9699236	
Mounting		
Side mounting bracket for Navigator controller	9699298	




Turbo-V 2K-G System




Dimensions: millimeters [inches]

### **Technical Specification**

Pumping Speed	
$N_2$	1.600 l/s
Compression Ratio	3 x 10 <sup>5</sup>
$N_2$	3 X 10 <sup>-5</sup>
Base Pressure* (with minimum recommended forepump)	< 1 x 10 <sup>-8</sup> mbar
Inlet Flange	ISO 250 F
Foreline Flange	KF 40 NW
Rotational Speed	33.300 rpm
Start-Up Time without Gas Load and with the Recommended	< 7 minutes
Recommended Forepump	> 40 m³ /h
Operating Position	Any
Operating Ambient Temperature	+5 to +40 °C
Power Supply	100 04014
Input voltage	100 - 240 Vac 50 - 60 Hz
Input freq.	30 - 00 HZ
Communication Interface	
Analogue I/O	Standard
RS232 / RS485	Standard
Profibus	Optional
Dimensions	
Height	335 mm (13.18 in.)
Diameter	335 mm (13.18 in.)
Weight	35 kg (77 lbs)





Pumping speed vs inlet pressure

Compression ratio vs foreline pressure

### **Ordering Information**

Water hoses

Description	Part Number
Pumping Systems*	
Turbo-V 2K-G system	9698871
Turbo-V 2K-G system with Profibus	9698873
*Pumping Systems include Pump with 230 V integrated electronic inlet screen with center-ring and pre-installed IP54 mating connec both standard models can be equipped with Agilent's MoniTorr de Also 120 Vac units available on request.	tors.
Accessories	
Mains cable NEMA plug, 3 m long	9699958
Mains cable European plug, 3 m long	9699957
Mains cable UK plug, 2.5 m long	X3501-68005
Mains cable China plug, 3 m long	X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	9699883
Water cooling kit for 6x8 (IDxOD) flexible tube	9699348
Water cooling kit for % in. id flexible tube	9699338
Turbo-V remote cable	9699945
Turbo-V remote panel	9698850
Inlet screen ISO 250 with centering ring	9699138
Spare Parts	
IP54 proof mating connectors	9699958

9699825



### Turbo-V 2300 TwisTorr

# BOTTOM VIEW ISO250F AND CFF 12" FORELINE KF 40 PURGE/VENT VALVES PORT VENT HOLE ON ENVELOPE O334.90 [13.19] FLANGE ISO 250F WATER COOLING FOR RUBBER TUBE Ø 9.0 O329.00 [12.95] WATER COOLING FOR RUBBER TUBE Ø 9.0 O330.00 [12.05] FLANGE O300.00 [12.05] FLANGE O300.00 [12.05] FLANGE O300.00 [12.05] WATER COOLING FOR RUBBER TUBE Ø 9.0 O329.00 [12.95] WATER COOLING FOR RUBBER TUBE Ø 9.0 O329.00 [12.95]

Dimensions: millimeters [inches]

### **Technical Specification**

Pumping Speed	
N <sub>2</sub>	2050 l/s
He (*)	1800 l/s
H <sub>2</sub> (*)	1500 l/s
Compression Ratio	
N <sub>2</sub>	>8 x 10 <sup>8</sup>
He	8 x 10 <sup>5</sup>
H <sub>2</sub>	4 x 10 <sup>4</sup>
Base Pressure* (with minimum recommended forepump)	10 <sup>-10</sup> mbar (7.5 x 10 <sup>-11</sup> Torr) (**)
Max Foreline Pressure for N <sub>2</sub>	4 mbar
Inlet Flange	ISO 250F, CFF 12" O.D.
Foreline Flange	KF 40 NW
Rotational Speed	33.300 rpm
Start-Up Time without Gas Load and with the Recommended	< 6 minutes
Recommended Forepump	TriScroll 600
Operating Position	Vertical/Upside down
Operating Ambient Temperature	+5 °C to +35 °C
Bakeout Temperature	120 °C (CFF), 80 °C (ISO)
Max Rotor Temperature	120 °C
Vibration Level (displacement)	< 0.01 µm at inlet flange
Lubricant	Permanent lubrication


Coolant water	Recommended flow: 200 l/h Temperature: +15 to +30 °C Pressure: 3 to 5 bar (45 to 75
Noise level	<60 dB(A) at 1 meter
Storage Temperature	-20 °C to +70 °C
Environment protection	IP54
Weight	ISO 250: 54.2 kg (119.5 lbs) CFF 12": 55.3 kg (121.9 lbs)

(\*)Without inlet screen

(\*\*)According to standard DIN 28 428, the base pressure is that measured in a leak-free test dome, 48 hours after the completion of test dome bake-out, with a Turbopump fitted with a CFF flange and using the recommended pre-vacuum pump

### **Controller Specifications**

Input Voltages 1-phase	100-120 VAC, 220-240 Vac,
Input Frequency	50 - 60 Hz
Maximum Input Power	1300 VA
Output Voltage	64 Vac
Output Frequency	555 Hz
Output Power Starting	560 W maximum
Output Power Normal	450 W maximum
Weight (both models)	12.5 kg (28 lbs)
Installation Category	II
Pollution Degree	2



ordering information		
Description	Part Number	
Pump		
Agilent Turbo-V 2300 TwisTorr ISO250F	9696000	
Agilent Turbo-V 2300 TwisTorr CFF12" OD	9696001	
Controllers		
Agilent Turbo-V 2300 rack controller, 100-240 V	X3501-64004	
Agilent Turbo-V 2300 Pump-controller cable kit, 5 m, for pump and vent valve operation	X3501-68003	
(*): Order controller-to-pump cable kit separately		
Accessories*		
Mains cable NEMA plug, 3 m long	9699958	
Mains cable European plug, 3 m long	9699957	
Mains cable UK plug, 2.5 m long	X3501-68005	
Mains cable China plug, 3 m long	X1699-64144	
Centering ring, ISO 250	9699144	
Inlet screen ISO 250	9699350	
Water cooling kit for 6x8 (IDxOD) flexible tube	9699348	
Water cooling kit for 3/8 in. ID flexible tube	9699338	
Purge valve 20 SCCM ¼ Swagelok ¼ Swagelok	9699236	

<sup>\*</sup>For a complete list of available extension cables, see the Extension Cables for Turbo Pumps dedicated page on Agilent.com.



### TwisTorr 74/84 FS On-board Controller

The TwisTorr 74 ans 84 FS On-board controlle is a solid-state frequency converter with the following features:

- Drives the TwisTorr 84 FS pumps.
- Powers the pump cooling fan.
- Drives the vent valve.
- Provides and acquires the pressure of the wide range gauge.
- Remote I/Os compatible with the previous version.
- Navigator default serial compatible with the previous RS 232 and 485 version.
- Able to operate with active gauges (full range gauge FRG-700 and FRG-702).

## 21.4 (0.8) 68.5 (2.7) 4x3MA LMAX=8mm 91.0 (3.6) 3x5MA LMAX=8mm 107.0 (4.2) 107.0 (4.2)

Dimensions: millimeters [inches]

The dedicated controller is a solid-state frequency converter which is driven by a single chip microcomputer and consists of two PCBs which include power supply and 3-phase output, analog and input/output section, microprocessor and digital section.

The power supply, together with the 3-phase output, converts the single phase AC mains supply or 24Vdc supply into a 3-phase, low voltage, medium frequency output which is required to power the pump. The controller can be operated by a remote host computer via the serial connection.

Windows-based software is available (optional).

The TwisTorr 74 and 84 FS AG Navigator controller can be mounted on board, on the bottom of the TwisTorr 74 and 84 FS pumps.

### **Technical Specification**

reclinical opecinication		
Input Voltages: Navigator controller Navigator controller	24 Vdc 100/240 Vac	
Input Power. Navigator controller Navigator controller	24 Vdc 80 W 100/240 Vac 180 V	/A
Input Frequency: Navigator controller	100/240 Vac	50 to 60 Hz ±5%
Fuse	2 x T4 A 250 V	
Output Voltage	60 Vac	
Output Frequency	1167 Hz for 74 FS 1350 Hz for 84 FS	
Output Power. Navigator controller, Navigator controller, (normal and autotuning)	24 Vdc 50 W 100/240 Vac 70 W	
Operating Temperature	+5 ÷ +45 °C 0 – 90% (Not cond	lensing)
Storage Temperature	-20 to +70 °C 0 - 95%	

Description	Weight kg (lbs)	Part Number
Controllers		
TwisTorr 74 FS AG on-board AG controller *110/220 V TwisTorr 74 FS AG on-board controller 24V 84FS AG on-board controller* 110/220 V 84FS AG on-board controller 24 V * Order power cable separately	0.3 (0.66)	X 3509-64030 X3509-64021 X3509-64000 X3509-64001
Cables		
Mains cable NEMA plug, 3 m lon		9699958
Mains cable European plug, 3 m long		9699957
Serial cable and A-PLUS software		9699883
Active Gauges		
FRG-700 Pirani/IMG combination gauge, KF2	15	FRG700KF25
FRG-700 Pirani/IMG combination gauge, DN4	10 CF	FRG700CF35
Gauge connection cable to AG rack controller	- 3 m	9699960
Gauge connection cable to AG rack controller	- 5 m	9699961



### 106.4 (4.19) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.4 (3.60) 91.5 (0.06) 91.5 (0.06)

### Dimensions: millimeters [inches]

### TwisTorr 74/84 FS AG Rack Controller

The TwisTorr 74 and 84 FS-AG rack controller is a microprocessor-controlled frequency converter with new, enhanced features that allow greater control and communication capabilities.

This compact, ¼ rack unit is designed for full worldwide compatibility, for vent valve control, active gauge pressure reading and pump operation parameters control, as well as for self-diagnostic and protection features.

**Universal Voltage:** Controller able to auto set according to the input voltage, providing flexibility for easy installation worldwide

### RS-232/485 Communication Protocols and Profibus (Options)

- Allows the pump to be interfaced with the system controls.
- Enables the pump to be operated via PC with A-PLUS software.
- Helps development of own customer software.

### Stop Speed Reading (SSR):

continues pump speed reading after the stop command

- Allows monitoring the pump in slow down ramp and shut down time.

### N.O. and N.C. Vent Valve Drive:

- Valve delay and opening time adjustable.
- Vent valve driven by controller automatically or by serial line.
- Opening time control through SSR function.

### **Technical Specifications**

Input	100 ÷ 240 Vac
	50/60 Hz
Maximum Input Power	210 VA
Output Voltage	76 Vac
Output Frequency	1167 Hz for 74FS series pumps 1350 Hz for 84FS series pumps
Nominal Output Power	100 W
Operating Temperature	+5 to +45 °C
Storage Temperature	−20 to +70 °C


### Active Gauge Pressure Reading Capability:

- Able to operate with active gauges (full range gauge FRG-700 and FRG-702).
- Provides accurate pressure measurements from atm down to 10-9 mbar.
- Rack mounted control electronics no longer necessary, with simple +24 Vdc power connector.
- Proven inverted magnetron design.
- Rapid start even in high-vacuum conditions.

Description	Weight kg (lbs)	Part Number
Controller		
TwisTorr 74FS AG rack controller * RS232/485 TwisTorr 74FS AG rack controller* Profibus TwisTorr 84FS AG rack controller* RS232/485 TwisTorr 84FS AG rack controller* Profibus * Order power cable separately	1.7 (3.2)	X3508-64301 X3508-64022 X3508-64001 X3508-64002
Cables		
Mains cable NEMA plug, 3 m long		9699958
Mains cable European plug, 3 m long		9699957
Mains cable UK plug, 2.5 m long		X3501-68005
Mains cable China plug, 3 m long		X1699-64144
RS232 serial cable and A-PLUS software, 3 m long		9699883
Active Gauges		
FRG-700 Pirani/IMG combination gauge, KF25		FRG700KF25
FRG-700 Pirani/IMG combination gauge, DN40 C	F	FRG700CF35
FRG-702 Pirani/IMG combination gauge, KF25 All metal - bakeable		FRG702KF25
FRG-702 Pirani/IMG combination gauge, DN40 CF All metal - bakeable		FRG702CF35
Gauge connection cable to AG rack controller, 3 r	n	9699960
Gauge connection cable to AG rack controller, 5 r	n	9699961



### TwisTorr 305 FS Rack Controller



Dimensions: millimeters [inches]


### **Technical Specifications**

Input	100 ÷ 240 Vac 50√60 H
Maximum Input Power	300 VA
Output Voltage	75 Vac
Output Frequency	1010 Hz
Nominal Output Power	150 W
Operating Temperature	+5 to +45 °C
Storage Temperature	−20 to +70 °C

Description	Weight kg (lbs)	Part Number
Controllers		
TwisTorr 305 FS AG rack controller* with RS 232/485	1.7 (3.2)	X3506-64002
TwisTorr 305 FS AG rack controller* with Profibus	1.7 (3.2)	X3506-64003
*Order power cable separately		
Cables		
Mains cable NEMA plug, 3 m long		9699958
Mains cable European plug, 3 m long		9699957
Mains cable UK plug, 2.5 m long		X3501-68005
Mains cable China plug, 3 m long		X1699-64144
Replace with: RS232 serial cable and A-PLUS	software, 3 m long	9699883
Active Gauges		
PVG-500 Pirani, KF16		PVG500KF16
PVG-500 Pirani tungsten filament, KF16 F		PVG500KF16S
FRG-700 Pirani/IMG combination gauge,		FRG700KF25
FRG-700 Pirani/IMG combination gauge,	DN40 CFF	RG700CF35
FRG-702 Pirani/IMG combination gauge, KF25 All metal - bakeable		FRG702KF25
FRG-702 Pirani/IMG combination gauge, DN40 CF All metal - bakeable		FRG702CF35
PCG-750 Pirani/CDG combination gauge	(Tungsten), KF16	PCG750KF16
Gauge connection cable to AG rack control	oller - 3 m	9699961
Gauge connection cable to AG rack control	oller - 5 m	9699961
NOTE • For other Active Gauges models pl	lease ask Agilent fordeta	ils



### TwisTorr Medium-TMP Onboard Controller



Dimensions: millimeters [inches]

### **Technical Specifications**

Input	100-240 VAC 50/60 Hz, 1 ph
Maximum Input Power	600 VA
Output Voltage	54 Vac, 3 ph
Output Frequency	825 Hz
Operating Temperature	+5 to +40 °C
Storage Temperature	−20 to +70 °C

### **Ordering Information**


Description	Part Number
Controllers	
TwisTorr Medium-TMP onboard controller	X3512-64016
Cables	
Mains cable NEMA plug, 3 m long	9699958
Mains cable European plug, 3 m long	9699957
Mains cable UK plug 2.5 m long	X3501-68005
Mains cable China plug 2.5 m long	X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	9699883

Medium-TMP extension cable	9699948M001 (5 m)
	9699948M003 (10 m)
	9699948M004 (15 m)
	9699948M002 (20 m)
	9699948M005 (30 m)
	9699948M006 (50 m)
	9699948M009 (75 m)
	9699948M010 (100 m)

Description	Part Number
Mounting	
Medium-TMP onboard ctlr. side-mounting bracket	X3511-68003
Cooling	
Medium-TMP air cooling kit for onboard controller	9699297
Venting	
Vent valve driven by onboard controller	9699834



### TwisTorr Medium-TMP Rack Controller



Dimensions: millimeters [inches]

### **Technical Specifications**

•	
Input	100-240 Vac 50/60 Hz, 1 ph Maximum
Maximum Input Power	800 VA
Output Voltage	54 Vac, 3 ph
Output Frequency	825 Hz
Operating Temperature	+5 to +45 °C
Storage Temperature	−20 to +70 °C

### **Ordering Information**

Description	Part Number	
Controllers		
TwisTorr Medium-TMP rack controller	X3501-64016	
Cables		
Mains cable NEMA plug 3 m long	9699958	
Mains cable EU plug 3 m long	9699957	
Mains cable UK plug 2.5 m long	X3501-68005	
Mains cable China plug 2.5 m long	X1699-64144	
RS232 serial cable and A-PLUS software, 3 m long	9699883	
Medium-TMP extension cable	9699948M001 (5 m) 9699948M003 (10 m) 9699948M004 (15 m) 9699948M002 (20 m) 9699948M005 (30 m) 9699948M006 (50 m) 9699948M009 (75 m) 9699948M010 (100 m)	

Description Part Number		
Cooling		
Medium-TMP Air Cooling Kit for rack	X3501-68001	
	X3501-68101	(5 m)
	X3501-68051	(10 m)
	X3501-68061	(15 m)
Air III III III	X3501-68021	(20 m)
Air cooling kit extension cable	X3501-68011	(30 m)
	X3501-68071	(50 m)
	X3501-68081	(75 m)
	X3501-68091	(100 m)
Venting		
Vent valve kit, normally open, with 5 m cable, for pumps driven by rack controllers	X3501-68002	
Vent valve kit, normally closed, with 5 m cable, for pumps driven by rack controllers	X3501-68022	
	X3501-68004	(5 m)
	X3501-68054	(10 m)
	X3501-68064	(15 m)
V 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	X3501-68074	(20 m)
Vent valve rack extension cable	X3501-68084	(30 m)
	X3501-68034	(50 m)
	X3501-68094	(75 m)
	X3501-68104	(100 m)



### Turbo-V 1001 Navigator

Dimensions: millimeters [inches]

### Controller

The Turbo-V 1001 Navigator controllers are microprocessor-controlled frequency converters, fully controllable through PC software, with self-diagnostic and protection features that ensure the highest degree of reliability. They can be mounted on board, either on the bottom or on the side of the pump, offering outstanding flexibility and simplicity.

### **Technical Specifications**

Input voltages	100/120/220/240 VAC 50/60 Hz, 1 ph
Maximum Input Power	850 VA
Output Voltage	54 VAC, 3 ph
Output Frequency	640 Hz
Operating Temperature	0 to +40 °C
Storage Temperature	-20 to +70 °C

Description	Weight kg (lbs)	Part Number
Controllers		
Turbo-V 1001 Navigator controller 120/220 V - 50/60 Hz	3.0 (6.0)	9698978
Accessories		
Mains cable NEMA plug, 3 m long	0.5 (1.0)	9699958
Mains cable European plug, 3 m long	0.5 (1.0)	9699957
Mains cable UK plug, 2.5 m long		X3501-68005
Mains cable China plug, 3 m long		X1699-64144
RS232 serial cable and A-PLUS software, 3 m long	0.5 (1.0)	9699883
External "TMP Profibus gateway"	0.4 (0.8)	9699261
Vent valve kit, normally open, with 0.6 m cable, for pumps driven by onboard controllers		9699834



### Turbo-V 1001 Rack Controller


These controllers are microprocessor-controlled frequency converters with self diagnostic and protection features that ensure the highest degree of reliability. The compact, ½ rack unit has a multifunction alphanumeric display for pump status and error code diagnostics.

The front panel has a two-line dot matrix LCD display with back lighting. It displays rotational speed as the pump starts up and indicates when full speed is reached. At any time during the operation of the pump, the speed, current, power, and bearing temperature can be displayed.

Additionally, the microprocessor acts as a pump cycle log, and can display the number of vacuum cycles, the cycle time

### **Technical Specifications**

Input	100-240 VAC 50/60 Hz, 1 ph			
Maximum Input Power	800 VA			
Output Voltage	54 Vac, 3 ph			
Output Frequency	825 Hz			
Maximum Output Power*	450 W			
Startup Power	450 W			
Operating Temperature	+5 to +45 °C			
Storage Temperature	−20 to +70 °C			
* Data valid for nitrogen.				




Dimensions: millimeters [inches]

for the current cycle, and the total operating hours on the pump. Remote operation can be accomplished with logic level contact closures and with optional computer interfaces.

PCB controllers are available. Please contact Agilent for details.

Description	Weight kg (lbs)	Part Number
Controllers		
Turbo-V 1001 rack controller, 100-240 V	15.7 (35.0)	X3501-64003
Accessories		
Mains cable NEMA plug, 3 m long		9699957
Mains cable (American plug, 120 V, 3 m lon	g)	9699958
Mains cable UK plug, 2.5 m long		X3501-68005
Mains cable China plug, 3 m long		X1699-64144
Options		
Rack adapter for controller	2.0 (4.0)	9699191
Controller to pump extension cable (5 m extension) for Turbo-V 1001*	1.0 (2.0)	9699948M001
Controller to pump extension cable (20 m extension) for Turbo-V 1001*	4.0 (8.0)	9699948M002
*For different extension cables, see medium	n rack controller extens	ion cable list
External "TMP Profibus gateway"	0.4 (0.8)	9699261
Vent valve kit, normally open, with 5 m cabl for pumps driven by rack controllers	e,	X3501-68002





### Turbo-V 2300 Rack Controller

Dimensions: millimeters [inches]

### **Technical Specifications**

Input	100-120 VAC, 220-240 VAC, 1-phase		
Maximum Input Power	1.300 VA		
Output Voltage	64 Vac		
Output Frequency	555 Hz		
Output Power Starting	560 W maximum		
Output Power Normal	450 W maximum		
Weight (both models)	12.5 kg (28 lbs)		
Installation Category	II		
Pollution Degree	2		

Description	Part Number
Controllers	
Agilent Turbo-V 2300 rack controller, 100-240 V	X3501-64004
Agilent Turbo-V 2300 pump-controller cable kit, 5 m, for pump and vent valve operation	X3501-68003
* Controller-to-pump cable kit to be ordered separately	



### Agilent Turbo Pump Accessories

### **Inlet Screens**

To prevent any possible damage to the pump blades caused by particles falling into the pump, an inlet screen is available as an accessory. The inlet screen is a curved stainless steel mesh. It provides protection against debris larger than 0.7 to 3 mm (depending on the model), while reducing the pumping speed by only about 10%.

For ordering information, see individual pumps.



Sample of inlet Screen

### **Air Cooling Kits**

A fan is provided as an option for applications requiring forced air flow. The fan is easily installed using the included installation kit. The air cooling kit can be mounted beside or underneath the pump body. The fans for all the other pump models can be installed in a side location only. The maximum allowable ambient temperature for an effective cooling action is 30 °C. To meet the air flow rate specifications, care must be taken not to restrict or cover the space around the fan so that air can flow freely. The operation of the cooling fan is controlled by the turbo controller.

For ordering information, see individual controllers.

### **Water Cooling Kits**

A water cooling kit is provided to cool the pump when operating at high inlet or high exhaust pressures. The cooling water can be supplied by an open circuit with drainage or a closed-loop refrigerated system.

Care must be taken to secure the tubes to ensure that they do not detach during operation. Four different kits with various tube sizes are available (refer to the table).


Tube Size	Material	Screw Thread	Part Number
4 x 6 mm (ID x OD)	Plastic	⅓ BSP	9699347
6 x 8 mm (ID x OD)	Plastic	1/4 BSP	9699348
1⁄4" ID	SST	⅓ BSP	9699337
3/8" ID	SST	1/4 BSP	9699338

### **Vibration Dampers**

Highly sensitive equipment such as an electron microscope and certain analytical instruments require extremely low vibration amplitudes. The dampers reduce turbopump vibration amplitude by at least a factor of 10 at their main frequencies.

They are available up to a size of 200 mm (ISO 200 or CFF 10" O.D.). To effectively decouple the equipment from the pump, the vibration damper must be installed between the connecting flanges.

For ordering information, see individual pumps.



Vibrator damper drawing

	H (mm)	H (in)	
Vibration Damper 4 ½" CF	96.5	3.80	
Vibration Damper 6" CF	101	3.97	
Vibration Damper 8" CF	110	4.33	
Vibration Damper 10" CF	113	4.45	
Vibration Damper 63 ISO	84	3.31	
Vibration Damper 100 ISO	84	3.31	
Vibration Damper 160 ISO	88	3.46	
Vibration Damper 200 ISO	88	3.46	

### **Turbo Vent Valve**

The Turbo Vent Valve, consisting of a control unit and a valve, is a complete unit for automatic venting of the turbo pump when it is switched off or during a power failure. The valve is a normally open, electromagnetically-actuated valve with a filter on the air inlet. The control unit is powered by the turbo controller and is provided with a fixed delay time of about 5 seconds to avoid undesired venting during a temporary power failure and to allow closure of the system valves before venting.


For ordering information, see individual controllers.

### **Purge Valve**

To protect the bearings of a turbomolecular pump used with aggressive process gases, a measured supply of inert gas  $(N_2, Ar)$  must flow into the pump body around the upper bearing toward the fore vacuum line. The Agilent fixed orifice purge valve is calibrated for nitrogen to provide the correct gas flow to safely operate turbo pumps. For ordering information, see individual pumps.

For ordering information, please call your local Agilent Office.

57



### Turbomolecular Pump Parameters and Definitions

### **Throughput**

"Throughput" is defined as the flow rate of pumped gas through the turbomolecular pump (and foreline pump). Throughput (Q) is measured in mbar l/s 1/60 standard cm<sup>3</sup>/min.

The maximum throughput a pumping system can handle is, in general, dependent upon the size of its foreline pump rather than the turbomolecular pump.

### **Pumping Speed**

"Pumping speed" (S) (volumetric flow rate) of a turbomolecular pump is the ratio between throughput and inlet pressure (foreline pump size must be the recommended one as a minimum).

$$S = Q/p$$

The pumping speed of a turbomolecular pump is constant over a wide pressure range and depends upon geometric factors such as diameter and rotational speed. For most turbomolecular pumps, pumping speed is nearly independent from gas species (molecular weight).

### **Compression Ratio**

"Compression Ratio" is the ratio between foreline (partial) pressure and inlet (partial) pressure for a given process gas, measured in "zero flow" conditions (performed by injecting the process gas in the pump foreline while the high vacuum port is blanked off). Compression ratio is generally indicated with a letter (K). In technical specifications of turbomolecular pumps, it is the maximum attainable value of K (at low foreline pressure).

Compression ratio is, in fact, a function of the foreline pressure as shown in Figure 1.

Compression ratio decays at high pressure depending on turbomolecular pump configuration (the number of molecular stages) and/or power limitations that slow down the rotor (gas friction increases with pressure).

The maximum compression ratio is strongly influenced by gas species: it is an exponential function of the molecular weight of the pumped gas (compression ratio is considerably lower for light gases).

### **Pumping Speed and Pressure Ratio**

The pressure ratio between foreline and inlet pressures in each operational situation is indicated by "Rp". This is, in general, equal to pumping speed ratio

$$Rp = p_{foreline} / p_{inlet} = S_{eff} / S_{foreline}$$

where  $S_{eff}$  is the effective Pumping Speed, and Sforeline is the pumping speed of the foreline pump.

In fact

$$Q = S_{eff \, pinlet} = S_{foreline \, pforeline}$$

therefore

$$S_{eff}/S_{foreline} = p_{foreline}/p_{inlet}$$

The pumping speed of a turbomolecular pump is minimally affected by pressure ratio (and foreline pump size) in most common operational conditions (when pressure ratio is much smaller than K).

Generally, however, the effective pumping speed " $S_{\rm eff}$ " is a linear function of the pressure ratio "Rp" as shown in Figure 2 (and therefore is also dependent upon the size of the backing pump).

Seff reaches its maximun value "S" (nominal pumping speed) when " $R_p$ " equals unity, and it is zero when the pressure ratio "Rp" has reached its maximum value "K".

This linear dependence can be expressed by the following relationship:

$$S_{eff} = S/(1 - 1/K + S/S_{foreline}K)$$
 (1)

Therefore, when

K >> S / S<sub>foreline</sub>

and

K >> 1

then

$$S_{eff} = S$$

when

K =1

then

$$S_{eff} = S_{foreline}$$

The above formula (1) must be used to evaluate pumping speed when operating at high pressure, especially with light gases (low K).

### **Base Pressure**

The base pressure of a turbomolecular pump is the equilibrium pressure between outgassing of pump surfaces exposed to high vacuum, including test dome, and the pumping speed of the pump.

$$p_{base} = Q_{outgas} / S_{eff}$$

In the case of ultimate operational pressure, as specified by norms, the pressure is measured after 48 hours bakeout of pump and dome (provided with metal gasket); therefore the prevailing outgassing product is H2 and equilibrium is reached with hydrogen pumping speed.

$$p_{base} = QH2 / S_{effH2}$$

When foreline pumps with relatively high base pressures are used, base pressure is sometimes limited by the compression ratio for H2O (or N2).

$$p_{base} = p_{forelineH2O} / KH2O$$

### **Pump Selection**

### How to Select a Turbo Pump

The right choice of a turbomolecular pump depends on the application. As a general rule we can reduce the choice to two types of use:

UHV (no gas flow) operations and process gas flow operations.

· UHV (no gas flow) operations.

The former case includes most cases in which the turbomolecular pump is employed to create vacuum in systems where the gas load is mainly produced by outgassing. In this application the choice is typically based on the desired base pressure within a desired time as a function of the foreseen outgassing rate, i.e.

$$S_{eff} = Q/p$$

where:

**p** is the desired base pressure (mbar)

**Q** is the total outgassing rate at the desired time (mbar l/s)

Seff is the effective pumping speed

Process gas flow operations.

The second case relates to all operations where process gases must be used. The main parameters are therefore the desired operation pressure and the process gas flow

$$S_{eff} = Q'/p'$$

where Q' is the total gas flow and p' is the operating pressure.

### How to Select the Backing Pump of a Turbomolecular Pump

The selection of a backing pump should be based on analyzing two requirements of the vacuum system:

- 1. The roughing time
- 2. The minimum recommended backing pump of the turbo
- **1. Roughing:** once the desired roughing time is established, the size of the forepump can be determined through the following formula:

$$S_{foreline} = (V/t) ln (p0/p1)$$

where

S<sub>foreline</sub> is the pumping speed of the roughing pump (I/min)

V is the volume of the chamber to be evacuated (I)

t is the desired roughing time (min)

**p0** is the starting pressure (mbar)

p1 is the end pressure (mbar)

When using a foreline pump much larger than the recommended size, a bypass line might be necessary to achieve calculated roughing time.

**2. Backing:** the backing pump must be big enough to achieve an effective pumping speed as close as possible to the nominal speed.

 $p_{foreline} = Q / S_{foreline}$ 

where

 $S_{foreline}$  is the pumping speed of the foreline pump

Q is the gas load

p is the operating foreline pressure

It should be noted that Q is the total gas load on the pump and includes process gases and turbo purge gases when used.

The size of the backing pump can be calculated according to the following rule:

S<sub>foreline</sub> ≥ 20S / K

where

**S** is the pumping speed of the turbopump

 $\mathbf{S}_{ extbf{foreline}}$  is the pumping speed of the backing pump

**K** is the maximum compression ratio of the turbopump for a given gas (i.e.: process gas) at the operating foreline pressure.

The pumping speed of the backing pump should be the higher of the two values calculated as above (roughing and backing).

Finally, it is possible to use a dry pump (scroll or diaphragm) for hydrocarbon-free operation when pumps of the MacroTorr type are used.



### Agilent Molecular Drag Technologies

Pumping stage drag technology enables higher foreline pressure, higher efficiency and allows for a smaller backing pump in a very compact design. Agilent solutions are designed using state of the art, proprietary numerical modeling.

We offer optimized solutions for:

- · Pumps with high compression ratio for UHV applications requiring lowest base pressure
- · Pumps with high pressure differential for high throughput
- Pumps with high discharge pressure allowing downsizing of the complete vacuum system

Agilent developed two innovative molecular drag stage technological platforms: TwisTorr and MacroTorr.

**TwisTorr** for high compression ratio in UHV applications requiring lowest base pressure (no gas flow).



TwisTorr 305-IC cutaway


**MacroTorr** for process gas flow operations with high pressure differential for high throughput.



TwisTorr 305 FSQ cutaway

### TwisTorr technology

- Pumping effect is created by a spinning rotor disk, which transfers momentum to gas molecules.
- Gas molecules are forced to follow spiral groove design on the stator. The specific design of the channel ensures constant local pumping speed and avoids reverse pressure gradients, minimizing power consumption.
- A single TwisTorr stage can improve the compression ratio for N<sub>2</sub> by up to 100 times compared to conventional stages, providing exceptional foreline tolerance and pumping speed.



Gas flow in centripetal and centrifugal direction through TwisTorr stages

### Centripetal pumping action

The lower surface area of the rotating disk transfers momentum to gas molecules. The spiral groove design on the upper section of the TwisTorr stator causes a centripetal pumping action.



The upper surface area of the rotating disk transfers momentum to gas molecules. The spiral groove design on the lower section of the TwisTorr stator causes a centrifugal pumping action.






### MacroTorr technology

- In the Agilent MacroTorr design, molecular impeller disks replace some of the turbo bladed stages.
- The molecular impellers consist of a disk rotating in a channel where the inlet and outlet are divided by a wall.
- The cross section of the channels decreases from the top to the bottom of the pump (from high vacuum to low vacuum or from the low pressure to the high pressure zone).

Gas molecules gain momentum after each collision with the moving surface of the impeller. The gas is then forced to pass through a hole to the next stage due to the wall.



### Designed for high gas load

The MacroTorr pumping stage is designed for operation with very high gas loads. It delivers high sustainable throughput at low operating temperatures in a compact package and extends the foreline tolerance up to 16 mbar.



### Agilent Floating Suspension

### Low vibration and stability over time

The Agilent Floating Suspension (AFS) system is a cutting-edge innovation designed to enhance the performance and longevity of turbomolecular pumps. By addressing key aspects such as bearing alignment, rotor dynamics, and noise reduction, the AFS system ensures optimal operating conditions for a wide range of high-precision applications.

### **Key Features and Benefits**

### **High Geometrical Precision for Perfect Bearing Alignment**

The AFS system boasts high geometrical precision, ensuring perfect alignment of the bearings. This precision is crucial for maintaining the integrity and performance of the turbo pump, reducing wear and tear, and extending the operational life of the bearings.

### **Improved Radial and Axial Stiffness**

The AFS system enhances both radial and axial stiffness, optimizing the dynamic behavior of the rotor. This improvement leads to better stability and reduced vibrations, which are essential for maintaining the accuracy and reliability of sensitive instruments.

### **Optimized Rotor Dynamic Behavior and Acoustic Noise Reduction**

By optimizing the rotor's dynamic behavior, the AFS system minimizes acoustic noise, creating a quieter operating environment. This feature is particularly beneficial in laboratory settings where noise reduction is critical for maintaining a conducive working atmosphere.

### **Low Vibration and Low Acoustic Noise**

The AFS system is designed to minimize both vibration and acoustic noise. Low vibration levels are essential for protecting sensitive components and ensuring accurate measurements, while low acoustic noise contributes to a more comfortable and productive working environment.

### **Optimal Working Conditions for Extended Bearing Life**

By providing optimal working conditions for the bearings, the AFS system significantly extends their operating life. This longevity reduces maintenance requirements and downtime, leading to increased efficiency and cost savings over the pump's lifecycle.

### **Exceptional Stability for SEM Applications**

The AFS system offers exceptional stability, making it ideal for the demanding requirements of scanning electron microscopy (SEM) applications. The enhanced stability ensures precise imaging and analysis, which are critical for high-resolution microscopy.

### **Excellent Thermal Stability**

The AFS system provides excellent thermal stability, maintaining consistent performance even under varying temperature conditions. This stability is crucial for applications that require precise control over operating temperatures to ensure accurate results.

The Agilent Floating Suspension system represents a significant advancement in turbopump technology, offering unparalleled precision, stability, and longevity. By incorporating the AFS system, the Agilent turbo pump can achieve superior performance and reliability, ensuring the success of your high-precision applications.










Vibration damping



Thermal stability



Optimized rotor dynamic behavior

### AI Enabled Adaptive Modal Balancing


### Lower noise and cleaner ultrahigh vacuum


Turbo pumps, which operate at speeds up to 80,000 RPM, are crucial for achieving high vacuum in various applications. To ensure their efficiency and longevity, precise balancing of rotating parts is essential. There are two main balancing techniques: influence-coefficient and modal-based, with modal balancing being the most advanced and precise.

Agilent has implemented an adaptive modal balancing (AMB) technique, using a proprietary algorithm that balances relevant vibration modes across three planes, improving precision and reducing vibration. This new approach integrates key processes into automated software, combining influence-coefficient and modal balancing methods with machine learning.

The Agilent proprietary Al-enabled algorithm continuously monitors the turbo pump balancing parameters all along the manufacturing line. It learns from each successfully balanced pump, analyzing and adjusting the modal balancing parameters to ensure optimal pump performance and production efficiency.

Key advantages: reduced noise, reliable vibration stability in a wider temperature range, and improved quality of turbo pumps. Ultimately the process provides a better customer experience!







Modal balancing planes rotor locations

Modal Balancing vs Influence Coefficient Method

### Did you Know?

Turbo pump balancing is achieved through the ABM process, which involves adjusting the position of tiny screws on the balancing planes. This method ensures that no material is removed or ablated from any rotating part, eliminating the risk of generating particles or dust. This strategy, inspired by the stringent cleanliness requirements of ion pump technology, is a key part of our manufacturing process.

Agilent turbo pumps are very clean!



### The Vacuum Link App for All Turbo Controllers

Innovation in vacuum control and data sharing



### Tips and tools

Download the new Agilent Vacuum Link app on your smartphone.

Visit https://www.agilent.com/en/product/vacuum-technologies/vacuum-leak-detection-software/vacuum-link-app for more.


Compatible with Apple devices only.




### A feature-rich app to help speed up the daily tasks.

The Agilent Vacuum Link can monitor up to three TwisTorr 305 pumps at the same tim.; A customized "Favorites" page that includes the most important parameters can be created and edited as needed.

Keeping track of pump operation is easy and fast and eliminates the need to sit in front of screen of a controller or a computer for long periods of time.



Configure your system




Read the pump parameters


### Control, export, and share data

Creating log files is easy and sharing them is quick using the regular features available on virtually any smartphone.

Log files help users review the pump parameters in a spreadsheet. Plotting parameter variations can be done using the dedicated icon.



See the pump performance in real time



Export and share data

### Agilent A-PLUS

A-PLUS is a communication, control, and monitoring software for controllers associated to selected Agilent products featuring the Agilent Window Serial Protocol.

With A-PLUS, you can simultaneously drive and control one or more Turbo Pumps, connected to a PC through an RS232 or an RS485 serial communication.

Automatic identification of the connected pumps, description of each command always on screen, and user interface adaptable to the pump status are only some of the new features developed to make the approach to pump settings easy and to reduce the number of steps during pump configuration.

Moreover, special care is given to the GUI (graphical user interface), to reproduce the environment of well-known user interfaces (such as Microsoft Windows applications), to obtain a real user friendly tool, and to reduce the user learning time.

A-PLUS software features several options like Data Logging, Chart Representation and Network Configuration, to help you configure your vacuum devices quickly, and to check your vacuum system status at any time.

Exhaustive online help is also included, providing the user with a complete, easy to learn system, tailored to customer requirements.

The previous software release 2.0.2 is included in the A-PLUS CD.



### CE/CSA, EMC Electrical Specifications Compliance

Compliance with these norms guarantees that there are no limits on the use of the controllers and turbopumps in every type of ambient, and that their use doesn't create any kind of disturbance to electronic units connected to the same line.

Agilent's new generation of Turbo controllers comply with the limits given by the following norms:

EN 55011 Class A group 1

EN 61000/3/2

EN 61000/3/3

EN 61000/4/2

EN 61000/4/3

EN 61000/4/4

EN 61010-1

UNI EN 291-1

**UNI EN 292-2** 

EN 1012-2

### Agilent Turbo Pump Service and Support Plan

This guide describes service options for your Agilent turbomolecular pump. Service options include exchange, upgrade, diagnostics, rechanneling, repair.

Tables 1 and 2 summarize each service option and provide an overview of the services included. Detailed descriptions of each option are given on the following pages.

Table 1. Summary of turbo pump service options

	When to Choose	Includes	Location	Availability	Packaging	Shipping	Warranty
Onsite Diagnostics	All pumps	Vibration analysis	Onsite	FSE/AE availability	-	-	-
Diagnostics	Health check	Fourier analysis	availability				
	Application consulting	Application optimization					
Rechanneling	Unused pumps only	Controlled grease distribution	Regional service	<2 weeks	Customer supplied	Pickup and delivery included	-
	Longer idle periods < 24 months old	Proprietary algorithm for process control	center		supplied		
	months old  Monitoring of power, current, and temperature at predefined frequency steps						
		Rebalancing input, if required					
		Final test and Fourier analysis					
Fast Track	< 5 years of infrequent use	Bearing replacement only	Regional	<2 weeks	Customer	Pickup and	6 months
Service	Clean applications only		service center	service	supplied	delivery included	
	Maintenance for pumps that still operate		center			e.aaea	
	Pump is not crashed						
Repair	Nonfunctional	Full pump repair, including cleaning and	Regional	<2 weeks	Customer	Pickup and	12 months
	Pump is not crashed	pearing replacement, balancing, and service inal test center				delivery included	
Advance	Nonfunctional	Fully rebuilt	Regional	From stock		Pickup and delivery	12 months
Exchange	Exchange is available	As-new performance and aesthetics	warehouse	use			
	Downtime reduction					meradea	
	Pump is crashed						
	(plus additional fee)						
Upgrade	Technology refresh	Fully rebuilt	Regional	From stock	Included	Pickup and delivery included	12 or 24 months, model dependent
	Downtime reduction	As-new performance and aesthetics	warehouse				
	Out of support	NEW pump option				moraucu	

Table 2. Services included

	Rechanneling	Fast Track Service	Pump Repair	Pump Exchange
Bearing Conditioning	•	•	•	•
Bearing Replacement		•	•	•
Rotor Balancing	•	•	•	•
Cleaning			•	•
Full Inspection		•	•	•
Helium Leak Check		•	•	•
Vacuum Performance Check		•	•	•
Update to "As-new" Condition			•	•
Two Week Turnaround time (TAT)	•	•	•	
Two day turnaround time (TAT)				•

### Get in touch with the Vacuum Support Team

Prompt, thorough, and accessible support is a core element of the Vacuum Products Technical Support and Service Groups. Our dedicated teams strive to ensure all users get answers to their questions or challenges in a fast and helpful manner.

Contact an Agilent Vacuum Expert using info below or Live Chat.

U.S. and Canada

Toll free: +1 800 882 7426

vpl-customercare@agilent.com

Europe/Middle East/Africa

Tel: +39 011 9979 111
Toll free: 00 800 234 234 00
vpt-customercare@agilent.com

China

Toll free: 800 820 6778

contacts.vacuum@agilent.com

India

Toll free: 1 800 572 3100 (TCL)
Toll free: 1 800 102 7906 (Airtel)

vpd\_india@agilent.com

Japan

Tel: +81 120477 111

jp-vvt-sales.pdl-ext@agilent.com

Korea

Tel: +080 222 2452

vpd\_pdl-sapkbc@agilent.com

Fax +39 011-9979-330

Learn more:

www.agilent.com/en/product/vacuum-technologies/turbo-pumps-controllers

FAQ:

https://www.agilent.com/.../turbomolecular-pumps-frequently-asked-questions

For Turbo Pump Support:

https://www.agilent.com/cs/library/brochures/br-turbo-pump-service-5994-6644en-agilent.pdf

Buy online:

www.agilent.com/store/#quickOrder

U.S., Canada and South America

1-800-882-7426 (toll-free)

vpl-customercare@agilent.com

Europe

00 800 234 234 00 (toll-free)

vpt-customercare@agilent.com

China

800 820 6778 (toll free, landline) 400 820 6778 (toll free, mobile)

Asia Pacific

inquiry\_lsca@agilent.com



