

## Butter Triglycerides

**Column:** DB-5ht

30 m x 0.32 mm I.D., 0.1  $\mu$ m

**J&W P/N:** 123-5731

**Carrier:** Hydrogen at 55 cm/sec, measured at 250°C

**Oven:** 35-250°C at 70°/min

250-400°C at 5°/min

400°C for 20 min

**Injector:**

Cool On-column

1  $\mu$ L of 9  $\mu$ g/ $\mu$ L in toluene

(approx. 1% w/w solution)


**Detector:**

FID, 400°C

Nitrogen makeup gas at 30 mL/min

Baseline Corrected

T = Total number of Carbons



C079