

Separation of *E. Coli* Strains

Column: **μSIL-WAX**
100 μm I.D. x 66 cm (22 cm to detector)

J&W P/N: 197-7202

Buffer: 10 mM sodium acetate/acetic acid, pH 4.7


Voltage: 10 kV

Injector: electrokinetic, 5 kV for 5 sec

Temperature: ambient

Detector: UV, 254 nm

Sample: A mixture of *E. Coli* Strains B and K-12
Membrane dialysis for 30 min at 4°C before injection



CE025