

Drug Screen

Column: DB-1ms

30 m x 0.25 mm I.D., 0.25 µm

J&W P/N: 122-0132


Carrier: Helium at 40 cm/sec,
measured at 50°C

Oven: 50°C for 1.0 min
50-125°C at 25°/min
125-325°C at 10°/min
325°C for 5 min

Injector: Cold Splitless using Optic II injector
50-250°C at 10°/sec,
0.75 min purge activation time

Sample: 1 µL injection of 50-150 ppm standard
Detector: FID, 300°C

1. Nicotine
2. Caffeine
3. Glutethimide
4. Lidocaine
5. PCP
6. Phenobarbital
7. Methadone primary metabolite
8. Methaqualone
9. Methadone
10. Cocaine
11. Desipramine
12. Carbamazepine
13. Trimipramine
14. Heroin
15. Fentanyl
16. Ibogaine
17. Triazolam
18. LSD



C759