


Phenoxy Acid Herbicides

Column: DB-1701
30 m x 0.32 mm ID, 0.25 μ m
P/N: 123-0732
Carrier: Helium at 35 cm/sec, measured at 50°C
Oven: 50°C for 1 min
50-160°C at 10°/min
160°C for 3 min
160-260°C at 10°/min
260°C for 5 min
Injector: Splitless, 250°C
45 sec purge activation time
1 μ L of 0.1 μ g/ μ L standard in methanol
Detector: FID, 300°C
Nitrogen makeup gas at 30mL/min

1. Dalapon
2. Dicamba
3. MCPP
4. MCPA
5. Dichloroprop
6. 2,4-D
7. Pentachlorophenol
8. Dinoseb
9. 2,4,5-TP
10. 2,4,5-T
11. 2,4-DB
12. Picloram

TMS Esters



Methyl Esters

