Agilent

Trusted Answers

GPC/SEC Column User Guide

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Installation

Tubing and connectors

Stainless steel tubing of 1/16 in od and 0.12 mm id or 0.17 mm id is recommended for column connections of analytical columns, and 0.5 mm id for 25 mm preparative columns. Connecting tubing lengths between columns, detectors and injection volumes should be minimized to avoid excessive dead volume which will diminish system performance. Column connections should be made using compatible 1/16 in nuts and ferrules. The compatibility of column connectors is illustrated in Figure 1.


Figure 1. Compatible connectors

The distance "X" for the standard column end fitting is 0.090 in and a minimum male nut length of 0.210 in is required. Examples of compatible fittings are Parker and Swagelok. Some fittings from other manufacturers may not be compatible, for example, Waters and Rheodyne. If unsure, please contact Agilent Technologies.

Column connection

Connect the GPC column in the eluent flow direction indicated and tighten the 1/16 in nut and ferrule using wrenches on the 1/16 in nut and the actual end fitting.

It is recommended that several drops of eluent have been pumped before the column outlet is connected to another column or detector to clean out the end fitting of any particulate matter which may be present.

To avoid loosening the end fittings and causing leaks, wrenches must be used on the end fitting adjacent to the connecting nut and NOT on the column barrel or the opposite end fitting, see Figure 2.


Figure 2. Do not use wrenches on the flats

Eluent Flow Rate

For conventional GPC columns using 7.5 and 7.8 mm id columns, 1.0 mL/min is an optimum flow rate for most separations. When column id is increased or decreased, the volumetric flow rate should be adjusted accordingly in order to give an equivalent linear velocity through the column.

The recommended flow rates are given in Table 1, however, in order to avoid excessively high pressures, higher viscosity eluents should be used at reduced flow rates or elevated temperature. Flow rates should be changed progressively and pressure pulses limited. At no time should the maximum operating pressure of the column be exceeded (see Table 4).

Table 1. Recommended flow rate

Column	Typical flow rates (mL/min)	Recommended flow rate (mL/min)
2.1 mm id	0.04 – 0.1	0.06
4.6 mm id	0.2 – 0.5	0.3
7.5 & 7.8 mm id	0.5 – 1.5	1.0
10 mm id	2.0 – 5.0	2.0
25 mm id	5.0 – 20.0	10.0

Sample Preparation and Injection

If maximum resolution and expected column lifetime are to be achieved, care must be taken in sample preparation.

To avoid blockage of the column frits, sample filtration is recommended (0.5 – 2.0 μm porosity depending on MW). A guard column will further protect the columns with little detrimental effect on performance.

Optimum sample volumes and concentrations are best determined for each type of analysis and are dependent on sample MW. Broad distribution polymer can generally be injected at higher concentrations than lower polydispersity samples. Overloading will not damage the column, but distorted peaks and therefore spurious results will be obtained.

Excessive injector loop volume can contribute to band broadening and reduce system performance, particularly with high efficiency or narrow bore columns. Agilent's injection volume recommendations are shown in Table 2.

Table 2. Injection volume recommendations

Column	Recommended concentration (%)	Recommended injection (μL) per column
2.1 mm id	0.01 – 0.2	0.2 – 2
4.6 mm id	0.01 – 0.2	1 – 20*
7.5 & 7.8 mm id	0.05 – 0.5	20 – 50**
10 mm id	0.05 – 0.5	20
25 mm id	0.5 – 5.0	500 – 2,000

* PL Multisolvant, MesoPore and ResiPore are 1 – 10 μL .

** PL Multisolvant, MesoPore and ResiPore are 10 – 20 μL , PL Rapide/PL Rapide Aqua are 20 μL .

All eluents should be of high purity and should be filtered and degassed prior to use.

GPC columns can be transferred to other eluents. When transferring to another eluent, miscibility and viscosity of the new eluent are of primary consideration. See Figure 3 for 7.5 mm id columns. For other id columns, apply the flow rate shown in Table 1.

When heating or cooling columns in high viscosity eluents (for example, TCB, NMP, DMF) a low solvent flow rate must always be maintained. Typically, 0.2 mL/min for 7.5 mm id and 0.1 mL/min for 4.6 mm id should be used prior to raising the temperature, see Figure 3.

Organic GPC columns are compatible with an extensive range of organic solvents ranging in polarity from perfluoroalkanes to dimethylformamide. Mixed organic solvent systems can also be used, but water should not be used except at concentrations less than 10% by volume in a miscible organic eluent. Columns are normally supplied in ethylbenzene unless otherwise stated, and can be flushed directly from ethylbenzene to THF at 0.5 mL/min. Unstabilized THF (for example, HPLC grade) is not recommended as an eluent due to the attack of peroxide on the gels.

Aqueous/polar GPC columns are normally supplied in water containing 0.02% sodium azide. Buffered eluent systems within the pH range 2 – 10 (of high and low ionic strength) may be used, with no detrimental effect on the column. When transferring from one eluent to another the compatibility and solubility of any salts or additives must be checked to prevent on-column precipitation which would irreversibly damage the column. PolarGel columns are compatible with an extensive range of organic solvents and with aqueous based eluents. Mixed organic solvent systems can be used, assuming full miscibility of the components. The only organic solvent which is recommended for use with the PL aquagel-OH and PL Rapide Aqua columns is methanol (up to 50% by volume).

PL Multisolvant columns are supplied in stabilized THF, and may be transferred directly to organic solvents, or to pure water and then to buffer. Unstabilized THF (for example, HPLC grade) is not recommended as an eluent due to the attack of peroxide on the gels. Buffered eluents within the pH range 2 – 8.5 (of high and low ionic strength) may be used. When transferring from one eluent to another the compatibility and solubility of any salts or additives must be checked to prevent on-column precipitation which would irreversibly damage the column. Polar organic solvents, such as DMF, DMSO, or DMAc are not recommended as eluents.

Organic columns (PLgel; PlusPore; PL Rapide)

Column supplied in ethylbenzene

Transfer to low viscosity solvents e.g. THF, Chloroform, Dichloromethane	Transfer to medium viscosity solvents e.g. Toluene, DMF, DMSO, HFIP	Transfer to high viscosity solvents e.g. TCB, m-Cresol, NMP
Flush column with acetone at 0.5 mL/min for 2 column volumes	Flush column with acetone at 0.5 mL/min for 2 column volumes	Set column oven to 50 °C, flow at 0.1 mL/min
Flush with new eluent at 0.5 mL/min for 2 column volumes	Flush with new eluent at 0.2 mL/min for 2 column volumes	Flush column direct with new eluent** at 50 °C at 0.1 mL/min for 2 column volumes
Increase column temperature to 30 - 40 °C* as required for analysis at 1 °C/min	Increase column temperature to 50 - 80 °C* as required for analysis at 1 °C/min	Increase column temperature to 100 - 220 °C* as required for analysis at 1 °C/min
Operate column in new eluent at required flow rate	Operate column in new eluent at required flow rate	Operate column in new eluent at required flow rate

*Always ensure operating temperature is at least 10 °C below boiling point of solvent. **Always ensure miscibility. If unsure, use acetone at room temperature.

Aqueous & Polar columns (PolarGel, PL aquagel-OH, PL Rapide Aqua)

Column supplied in water containing 0.02% NaN₃

Transfer to aqueous e.g. Water, Buffer	Transfer to polar organic** e.g. DMSO, DMF	Transfer to mixed solvent systems e.g. Water/THF*, Water/methanol*
Flush column with pure water at 1.0 mL/min for 2 column volumes	Flush column with pure water at 1.0 mL/min for 2 column volumes	Flush column with pure water at 1.0 mL/min for 2 column volumes
Flush with new buffer at 1.0 mL/min for 2 column volumes	Flush with acetone at 0.5 mL/min for 2 column volumes	Flush with new, premixed eluent at 0.2 mL/min for 2 column volumes
Increase column temperature* as required for analysis at 1 °C/min	Flush with new eluent at 0.2 mL/min for 2 column volumes	Increase column temperature* as required for analysis at 1 °C/min
Operate column in new eluent at required flow rate	Increase column temperature to 50 - 80 °C* as required for analysis at 1 °C/min	Operate column in new eluent at required flow rate
	Operate column in new eluent at required flow rate	

*Always ensure operating temperature is at least 10 °C below boiling point of solvent. ** Polargel only

Multisolvent columns (PL Multisolvent)

Column supplied in THF containing 0.015% BHT

Transfer to aqueous e.g. Water, Buffer	Transfer to organic e.g. Chloroform* Dichloromethane**
Flush column with pure water at 1.0 mL/min for 2 column volumes	Flush with new eluent at 0.5 mL/min for 2 column volumes
Flush with new buffer at 1.0 mL/min for 2 column volumes	Increase column temperature as required for analysis at 1 °C/min
Increase column temperature as required for analysis at 1 °C/min	Operate column in new eluent at required flow rate
Operate column in new eluent at required flow rate	

*Always ensure operating temperature is at least 10 °C below boiling point of solvent

Figure 3. Eluent transfer guide

Column Testing and Specifications

Every column is supplied with a test certificate indicating the test conditions and the column performance. Measurements of column performance are described below:

$$\begin{array}{ll} \text{Efficiency (1/2 ht)} & N = 5.54 \left[\frac{t}{W_{1/2}} \right]^2 / L \\ \text{(Plates/m)} & \\ \text{Efficiency (5\sigma)} & N = 25 \left[\frac{t}{W_{5\sigma}} \right]^2 / L \\ \text{(Plates/m)} & \\ \text{Symmetry} & a/b \end{array}$$

Where “t” is the peak elution time, “W1/2” is the peak width at half peak height, “W5σ” is the peak width at 4.4% of peak height, “L” is the column length in meters and “a” and “b” are the peak widths either side of the perpendicular measured at 10% of peak height.

Column efficiency is dependent on many experimental factors (system dead volume, eluent, flow rate, test probe, temperature, and so forth) and test results may differ from those quoted on the column certificate due to variability in these parameters. Band broadening effects are more severe when using high efficiency and/or narrow bore GPC/SEC columns. It is vital to ensure that the system dispersion is minimized in order to obtain the full potential of Agilent columns.

We recommend measuring performance when a column set is first used and at regular intervals after that. Consider replacing columns when efficiency falls below 80% of the starting value.

Table 3. Column specifications

Column	Typical operating pressure psi (bar)	Maximum operating pressure psi (bar)	Maximum operating temp. °C ²
PLgel 3 μm	750 (50)	2,700 (180)	110
PLgel 5 μm	450 (30)	2,200 (150)	150
PLgel 9 & 10 μm	150 (10)	2,200 (150)	220
PLgel 20 μm	50 (3)	2,200 (150)	220
PLgel Olexis	90 (6)	2,250 (150)	220
PLgel Prep 10 μm	350 (25)	2,200 (150)	150
PolyPore	450 (30)	2,200 (150)	160
ResiPore	750 (50)	2,700 (180)	120
MesoPore	750 (50)	2,700 (180)	120
OligoPore	450 (30)	1,500 (100)	120
PL Rapide-F	350 (25)	1,500 (100)	150
PL Rapide-L	225 (15)	1,100 (75)	150
PL Rapide-M	75 (5)	1,100 (75)	220
PL Rapide-H	75 (5)	1,100 (75)	220
PolarGel	450 (30)	1,600 (112)	80
PL aquagel-OH 5 μm	900 (60)	1,600 (112)	90
PL aquagel-OH 8 μm	450 (30)	1,600 (112)	90
PL aquagel-OH 15 μm	150 (10)	1,600 (112)	90
PL Rapide Aqua	225 (15)	1,600 (112)	90
PL Multisolvant (Water)	1450 (100)	1,600 (112)	80
PL Multisolvant (Organic)	750 (50)	1,600 (112)	80

¹Based on THF (organic columns) or H₂O (aqueous columns) at 20 °C, 4.6 mm id at 0.3 mL/min, 7.5 mm id at 1.0 mL/min, 10 mm id at 1.8 mL/min (2 mL/min for PL Rapide 10 x 100 mm) and 25 x 300 mm at 10 mL/min.

²At very high temperatures, column lifetimes are likely to be reduced. OligoPore and PL Rapide-F may show deterioration to lifetimes at >80°C, while PL Multisolvant may degrade at > 50 °C. Amend flow rates and temperatures accordingly.

Storage

On removing the column from the system, the end plugs must be replaced to prevent the column from drying out by evaporation, since subsequent shrinkage of the gel and disruption of the packing will occur. The end plugs need only be applied finger-tight. All eluents mentioned previously are suitable for storage, but unstabilized THF and halogenated solvents should not be used.

For long term storage of aqueous columns, flush with water and store in water containing 0.02% sodium azide.

Warranty

The columns are covered by warranty for 90 days following delivery. For columns used at or above 170 °C, the warranty period is reduced to 30 days. Agilent Technologies cannot accept liability from improper handling and use of columns above their maximum operation temperature and pressure (see Table 3) and void the warranty. For a full warranty statement, please request Agilent's General Conditions of Sale.

Maintenance

Deterioration in column performance may occur as a result of damage to the packed bed or as a result of blockage in the column frits. In the case of frit blockage, the column can be reverse flushed for 1 minute at the recommended flow rate to remove loosely retained material. If a column blockage exists despite back flushing, then some column inlet frits can be replaced by following the procedure outlined below.

Frit replacement procedure

1. Seal the end fitting at the outlet end of the column. Clamp the column vertically with the inlet end uppermost. Undo the end fitting on the inlet end of the column and gently remove it. (For aqueous columns, flush eluent into water prior to this procedure.)
2. After removing any gel from the end fitting, remove the frit, seal and spreader with the frit removal tool, and wash all traces of gel from the end fitting. Then place a new frit, seal and spreader in the bottom of the end fitting and ensure that they are properly seated.
3. Some gel (up to 2 mm) may have expanded out of the column. Using a spatula, place a dome of repair gel onto the top of the column and add sufficient solvent to bind the gel particles together.
4. Hold the column at an angle of 45° pointing downwards and replace the end fitting complete with frit, seal and spreader onto the column. Tighten the end fitting one flat past finger-tight.
5. Back flush the column with eluent at 0.5 mL/min to remove trapped air. Check for gel particles in the eluent in case the frit has not seated correctly. Return the column to normal flow conditions and check the pressure drop and column efficiency.

Agilent recommends that preparative columns are returned to Agilent for fitment of replacement frits or repair, however, a similar procedure to the above can be applied if necessary. For further technical assistance, or to use the Column Repair Service, please contact Agilent Technologies.

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Agilent Ordering Information

For more information on our products, visit our website at www.agilent.com/chem/gpcsec

For technical support, visit www.agilent.com/chem/techsupport

For your nearest Agilent office or distributor, visit www.agilent.com/chem/contactus

To place an order, visit www.agilent.com/chem/store

Agilent offers a complete line of sample preparation products to support LC and LC/MS applications.

The Agilent Bond Elut and Captiva Sample Prep family of products offer the largest choice of formats and widest range of solutions available in the market today.

Learn more at www.agilent.com/chem/sampleprep

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© Agilent Technologies, Inc. 2017, 2018, 2023