

A simple and robust CE-Marked diagnostic assay for complete mutation analysis of the *CFTR* gene by Next-Generation Sequencing. Mutations of the *CFTR* gene affect functioning of the chloride ion channels in epithelial cell membranes, leading to cystic fibrosis and other *CFTR*-related diseases.



## Application

- For first or second line testing in your CFTR work-up
- For whole blood & dried blood spots derived DNA
- For complete mutational spectrum analysis



### Assay characteristics


Gene analyzed	<i>CFTR</i>
Genomic region analyzed	8.7 kb All exons ± 30 bp flanking region, selected introns, part of promoter
Number of amplicons	48 including 11 control amplicons
Amplicon length	300-450 bp
Number of plexes	2
DNA amount required	20 ng per multiplex reaction
Validated (CE-IVD) with	MiSeq® System, Illumina
<b>Mutations included</b>	
All exonic SNVs & small indels CNVs*	~30 bp of intronic flanking sequence
HP variants	intron 9 (8): c.1210-34TG(m); c.1210-12T(n)
Selected deep intronic regions:	intron 7 (6b): c.1002-1110_1113delTAAG c.1002-1111A>C
	intron 12 (11): c.1679+1634A>G (1811+1.6kbA->G) c.1680-877G>T (1811+1643G->T) c.1680-883A>G
	intron 22 (19): c.3718-2477C>T (3849+10kbC->T)
Promoter region	up to position c.-108

\*For Research Use Only. Not for use in diagnostic procedures

### CFTR MASTR Dx on MiSeq, germline DNA

### Performance Characteristics

	DNA extracted from blood	DNA extracted from dried blood spots
On target read pair counts	98.92 % [92.70 %-99.51 %]	99.90 % [99.53 %-99.98 %]
Uniformity of amplification (0.2x mean coverage)	97.83 %	100 %



Graph representing the read counts for all 48 CFTR MASTR Dx amplicons, showing their uniform representation. To allow comparison between samples, the read counts were normalized.

## CFTR MASTR Dx on MiSeq; germline DNA

### Performance Characteristics

Performance parameter	Extracted from blood		Extracted from dried blood spots	
	Observed [95 % CI] Using JSI SeqNext	Observed [95 % CI] Using Sophia DDM	Observed Using JSI SeqNext	Observed Using Sophia DDM
Sensitivity	100 % [≥ 99.217 %]	100 % [≥ 98.202 %]	64/64 = 100 %	64/64 = 100 %
Specificity	99.990 % [99.975 % -99.994 %]	100 % [≥ 99.999 %]	1438/1438 = 100 %	1438/1438 = 100 %
Accuracy	99.990 % [99.975 % -99.994 %]	100 % [≥ 99.999 %]		
Repeatability & reproducibility	99.997 % [99.988 % -99.999 %]	100 % [≥ 99.990 %]		

## Advised maximum number of samples per run

### Workflow

Sequencing System	Illumina MiSeq®		
	Reagent Kit		
Flow cell	Nano v2 2 x 251 cycles	v2 2 x 251 cycles	v3 2 x 276 cycles
SNV variant calling Minimal coverage per allele: 50	52	786	1441
SNV and CNV* variant calling Minimal coverage per amplicon: 200	26	393	720

\*For Research Use Only. Not for use in diagnostic procedures



### Publications

#### Supporting publications:

Applicability and Efficiency of NGS in Routine Diagnosis: In-Depth Performance Analysis of a Complete Workflow for CFTR Mutation Analysis. Andrien Pagin et al. (2016)

[PLOS ONE DOI: 10.1371/journal.pone.0149426](https://doi.org/10.1371/journal.pone.0149426), 22 Feb 2016\*

\*This study was done by a third party and was verified by Multiplicom.

### Order information

Cat. No.	Product Name	Reactions	PR Number
MR-2021.024	CFTR MASTR Dx	24	PR7000-1421
MR-2021.048		48	

MID (Molecular Identifiers) kits are necessary to complete the workflow.

5991-8404EN  
Printed in Belgium, September 2017