

Agilent RNA Kits for the Agilent 2100 Bioanalyzer System

Fast Quality Control of RNA with Minimal Sample Consumption

Agilent offers RNA assays that allow characterization of a variety of RNA sample types and estimation of their concentrations within minutes, with only a few nano- or picograms of sample.

Key Features

- High assay accuracy and precision: Pre-packaged reagents and standardized assay protocols yield highly accurate and reproducible data
- Minimal sample consumption: Only 1 μL of sample required
- Fast results: Automated analysis of up to 12 samples in 30 minutes
- Quick and easy sample comparison: Automated sample alignment, one-click overlay, scaling and zooming features
- Data display options: Gel-like image, electropherogram, and tabular formats
- RNase free: Uses RNase free reagents and chips to avoid sample degradation during analysis
- Clean: Minimal exposure to hazardous materials, such as ethidium bromide
- Conveniently archive, store and share digital data



Analytical specifications	RNA 6000 Nano Total RNA	RNA 6000 Nano mRNA	RNA 6000 Pico Total RNA	RNA 6000 Pico mRNA	Small RNA
Sizing range	-	-	-	-	6 – 150 nt
Sensitivity ¹	5 ng/ μL in water	25 ng/ μL in water	50 pg/ μL in water 200 pg/ μL in TE	250 pg/ μL in water 500 pg/ μL in TE	50 pg/ μL in water ³
Quantitative precision (within a chip)	10 % CV	10 % CV	20 % CV	20 % CV	25 % CV
Quantitative accuracy ²	± 20 %	± 20 %	± 30 %	± 30 %	-
Quantitative range	25 – 500 ng/ μL	25 – 250 ng/ μL	-	-	50 – 2,000 pg/ μL of purified miRNA in water
Qualitative range	5 – 500 ng/ μL	25 – 250 ng/ μL	50 – 5,000 pg/ μL in water	250 – 5,000 pg/ μL in water	50 – 2,000 pg/ μL of purified miRNA in water
Maximum buffer concentration in sample	100 mM Tris 0.1 mM EDTA or 125 mM NaCl or 15 mM MgCl_2	100 mM Tris 0.1 mM EDTA or 125 mM NaCl or 15 mM MgCl_2	50 mM Tris 0.1 mM EDTA or 50 mM NaCl or 15 mM MgCl_2	50 mM Tris 0.1 mM EDTA or 50 mM NaCl or 15 mM MgCl_2	10 mM Tris 0.1 mM EDTA
Physical specifications					
Analysis time	30 minutes	30 minutes	30 minutes	30 minutes	30 minutes
Samples per chip	12	12	11	11	11
Sample volume required	1 μL	1 μL	1 μL	1 μL	1 μL
Kit stability	4 months	4 months	4 months	4 months	4 months
Kit size	300 samples	300 samples	275 samples	275 samples	275 samples

¹Signal-to-noise >3 (single peak)

²Determined using ladder as sample

³Measured for the 40 nt fragment of the Small RNA ladder

Ordering information

Order No.	Description
5067-1511	RNA 6000 Nano kit
5067-1512	RNA 6000 Nano reagents
5067-1529	RNA Nano ladder

Order No.	Description
5067-1513	RNA 6000 Pico kit
5067-1514	RNA 6000 Pico reagents
5067-1535	RNA Pico ladder

Order No.	Description
5067-1548	Small RNA kit
5067-1549	Small RNA reagents
5067-1550	Small RNA ladder

Small RNA Assay

Example: Quality control in genome editing

- Determination of yield and sizing of guide RNAs after in vitro transcription
- Identification of impurities, truncations, secondary structure formation

In CRISPR/Cas workflows, the synthesis of guide RNA (gRNA) is a key component for successful genome editing. gRNA created by in vitro transcription and purified may be run on the 2100 Bioanalyzer system. The Small RNA kit is ideally suited for this application as it allows for high resolution separation of small RNAs from 6-150 nucleotides. Important features to look for: impurities or secondary structure, samples are of expected size and sufficient yield for follow-up studies.




RNA 6000 Nano and Pico Assay

Example: Quality control of RNA

- Identify degradation of total RNA
- Identify degradation of mRNA
- High sensitivity with minimal sample consumption
- Unparalleled sensitivity in the picogram range

RNAse degradation of RNA samples is a common reason for failed experiments. The 2100 Bioanalyzer system provides RNA quality control results in both gel-like image as well as electrophoretic data making it easy to detect even small degradation effects. Indications for RNA degradation are:

- Decreasing ratio of ribosomal bands
- Additional peaks below the ribosomal bands
- Decrease in overall RNA signal
- Shift towards shorter fragments



www.agilent.com/genomics/bioanalyzer

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019
Printed in the USA, February 1, 2019
5991-7891EN