

# Narrow Bore High Performance Solids at 500 MHz and 600 MHz

Data Sheet



## Introduction

Narrow bore magnets fit liquids needs, but have been marginal for solids. This is because as much space in the bore as possible has been needed to fit large capacitors for each channel near the sample. These large capacitors have been necessary to handle the power requirements for short pulse widths.

## The Solution

By replacing the tune and match capacitors with tuning tubes, bore size is not a problem. High performance double and triple resonance probes can be supplied for probe sizes as narrow as 40 mm. Agilent can supply narrow bore solids probes for systems ranging from 400/54 to as high as 900/54, including 500 and 600/51.



**Agilent Technologies**


## Features

Agilent offers either HX double resonance or HXY triple resonance narrow bore T3 probes, with rotor diameters from 2.5 to 6 mm. They have performance comparable to wide bore probes with the exception of a reduced VT range (in order to protect the magnet bore). Spectrometers with these magnets are the workhorse systems for protein structure determination. Now membrane proteins can be studied using T3 probes. All it takes is standard 100/300 W amplifiers.


These probes are all 3.2 mm spinning module size. Other sizes (2.5, 4, 5 or 6 mm) are available on request. Pulse widths quoted here are for systems equipped with a nominal 100 W  $^1\text{H}$  amplifier and the appropriate number of standard 300 W broadband amplifiers. S/N values are for Double Resonance mode.

## Upgrades


Agilent offers an upgrade kit for *Unityplus* and *Inova* spectrometers. This includes  $^1\text{H}$  amplifier upgrade (if needed), VT unit upgrade, upper stack and spinning controller.




REDOR at 500 MHz



$^{13}\text{C}$  -  $^{13}\text{C}$  Correlation at 500 MHz



MQMAS at 500 MHz



HETCOR at 500 MHz

## Components



As most narrow bore systems are used for liquids NMR, Agilent has designed the solids accessory and probes to be as compatible as possible with a mixed use environment. Power amplifiers are the same as those used for liquids, the VT controller is compatible with liquids probes, and the preamps are higher power versions of the standard ones.

The solids VT upper stack is inserted through the middle of the liquids upper barrel, obviating the need to move the latter. The upper stack is heavily interlocked to guard against any chance of thermally stressing the magnet bore. VT gas flow is controlled through an intelligent gas panel, which senses whether liquids or solids VT gas is available. The spinning controller is the same as used for wide bore probes.

The result, as can be judged from the  $\text{Pb}(\text{NO}_3)_2$  spectra on the left, is a well controlled sample temperature over the entire temperature range.

## Currently Available Configurations

	500 HX	500 HXY	600 HX	600 HXY
Tuning range	$^1\text{H}, ^{31}\text{P} - ^{15}\text{N}$	$^1\text{H}, ^{31}\text{P} - ^{15}\text{N}$	$^1\text{H}, ^{31}\text{P} - ^{15}\text{N}$	$^1\text{H}, ^{31}\text{P} - ^{15}\text{N}$
$^1\text{H}$ pw90	2.2 $\mu$	2.2 $\mu$	2.2 $\mu$	2.2 $\mu$
$^{13}\text{C}$ pw90	2.0 $\mu$	2.8 $\mu$	2.0 $\mu$	2.8 $\mu$
VT range	-75 °C to + 100 °C	-75 °C to + 100 °C	-75 °C to + 100 °C	-75 °C to + 100 °C
S/N	110:1	100:1	140:1	125:1



[www.agilent.com](http://www.agilent.com)

Product specifications and descriptions in this document are subject to change without notice.

© Agilent Technologies, Inc., 2011  
Published in USA, March 28, 2011  
Publication Number SI-A-1005

