

# **GC Column Technology - Advancements and Practical Applications**

## **Agilent Technologies**



**Agilent Technologies**

Dial 1-816-650-0621 for e-Seminar Audio


Slide 4

# Capillary Column Research

- **Column deactivation chemistry:**
  - extend performance envelope- general use phases
  - optimize column chemistry - specific use columns
- **Stationary phase chemistries:**
  - improved upper temp limits- increased thermal stability & lower bleed
  - improved low temperature performance- increased flexibility for volatiles
  - improved chemical & physical stability- longer life, improved resolution
  - improved selectivities/ new applications




# Capillary Column Cross Section




# Capillary Column Types

## Porous Layer Open Tube (PLOT)



## Wall Coated Open Tube (WCOT)



# Optimized Column Deactivation

- Compatible with desired stationary phase polymer proper surface energy (polymer "wets" surface) provides functional groups for covalent bonding
- **Permits good chromatography: minimal absorption, symmetric peaks for compounds of interest:**  
  
**If general-use phase:** broad analyte spectrum capability-organic acids to organic bases
- Deactivation stable to solvent injection or wash



# Optimized Deactivation

- **“Broad Range” is difficult in practice.**
- **Deactivation affects performance**



# Common Deactivation Chemistries and Column Performance

## Deactivation Method

## Column Characteristics

Chlorosilanes

Acidic

Alkoxysilanes

Slightly Acidic

Silazanes or cyclic silazanes

Basic

Hydrosilanes or hydrosiloxanes

Acidic

Siloxane

Slightly Acidic



Polyethylene glycol (Carbowax)

Acidic





# Acid Performance

5% phenyl, 30 m x 0.32 mm x 0.25  $\mu$ m







1. pentanoic acid
2. hexanoic acid
3. heptanoic acid
4. octanoic acid
5. decanoic acid
6. dodecanoic acid





# Base Performance


5% phenyl, 30 m x 0.32 mm x 0.25  $\mu$ m




1. pentylamine
2. hexylamine
3. heptylamine
4. octylamine
5. nonylamine
6. decylamine



# HP-5 Equivalent Columns: Broad Spectrum Chemical Test



# DuraGuard



# Stationary Phase Advancements

- **Improved upper temp performance- Reducing bleed- “MS” phases.**
- **Extending minimum temp limits- more flexibility for volatiles- Waxes.**
- **Improved phase stability- longer life, improved resolution, expanded compatibility -PLOTs.**
- **New selectivities / applications.**





# What is Bleed?

- **thermodynamic equilibrium process occurs to some degree in all columns**
- **polysiloxane backbone releases low molecular weight, cyclic fragments**
- **occurs at low level in low temperature, O<sub>2</sub>-free, clean system**
- **increased by increased temperature, oxygen exposure or chemical contamination**




# Bleed: Why Does It Happen?

## “Back Biting” Mechanism of Product Formation



+




Repeat


***Cyclic products are thermodynamically more stable!***



# DB-5MS Structure



**DB-5**  
**5% Phenyl**





**DB-5ms**

- 1. Increased stability**
- 2. Different selectivity**
- 3. Optimized to match DB-5**



# Low Bleed Stationary Phases

## DB-35ms vs. DB-35



CLP Pesticides Analysis

\* 5ng Decachlorobiphenyl



**Agilent Technologies**



Dial 1-816-650-0621 for e-Seminar Audio

Slide 19



# Low Bleed Stationary Phases

## Spectral Purity



# DB-1ms

## The World's #1 Phase Just Got Better

- **Virtually identical selectivity to DB-1**
- **Higher temperature limit than DB-5ms (340/360°C)**
- **Legendary column to column reproducibility**
- **Engineered to perform: never hand selected**




# DB-1ms: Optimizing an Already Robust Phase

- **Changes to fused silica surface prior to polymerization**
- **Changes to polymerization chemistry**
- **Changes to chemical manifolds used to process columns**




# Molecular Weight Distribution and Bleed in DMPS Columns



**DMPS chains of inconsistent lengths**



# Molecular Weight Distribution and Bleed in DMPS Columns



DMPS chains of **consistent** lengths



# DB-1ms

## Compare the Retention Indices

Compound	DB-1	DB-1ms	MFG A	MFG B	MFG C
3,5-dimethylpyridine	952.4	954.0	959.1	966.1	958.7
1-nitrohexane	1009.4	1009.7	1009.9	1015.1	1015.0
1,4-diisopropylbenzene	1152.1	1152.7	1153.5	1157.3	1157.0
1-nonanol	1156.9	1157.2	1156.7	1161.1	1161.0
2-decanone	1169.3	1169.7	1170.4	1174.1	1173.0


Retention Indices at 90°C



# DB-1ms Inertness


**Column:** DB-1ms  
30 m x 0.25 mm I.D., 0.25 µm  
**J&W P/N:** 122-0132  
**Carrier:** Hydrogen at 39.8 cm/sec,  
measure at 125°C  
**Oven:** 125°C for 8.5 min  
125-265°C at 20°/min  
**Injector:** Split 1:50, 250°C  
**Detector:** FID, 300°C

1. 2-Ethylhexanoic Acid\*
  2. 1,6-Hexanediol
  3. 4-Chlorophenol\*
  4. 1-Methylnaphthalene
  5. Tridecane
  6. 1-Undecanol
  7. Tetradecane
  8. Dicyclohexylamine\*
  9. 2,4-Dinitrophenol\*
  10. Acenaphthene
  11. N-Nitrosodiphenylamine\*
  12. Pentachlorophenol
  13. Phenanthrene
  14. Carbazole
- \*Active Analytes



# Improved Bleed Performance

## Documented Quality Improvement





# Extending Low Temp Limits- Goals

- **Retain peak shape, efficiency at low temperatures**
- **Retain (if possible) original phase polarity (RI)**
- **Retain or expand maximum temperature capability**



# PEG (WAX) Phases


- **Carbowax- uncrosslinked PEG (HP 20M)**
- **DB-WAX- crosslinked via free radicals**
- **HP Innowax- crosslinked via side-chain functional groups**

**Goal: Reduce crystallinity, reduce melting point, reduce lower temp limit.**



# Melting Point of PEG Phases


## N / m vs. Temperature




# PEG Column Comparison

## Solvent Analysis

30 m x 0.25 mm id x 0.25  $\mu$ m df




1. pentane
2. methyl formate
3. acetone
4. ethyl acetate
5. methyl ethyl ketone
6. methanol
7. 2-methyl-2-propanol
8. methylene chloride
9. benzene
10. ethanol
11. 2-butanol
12. toluene
13. n-propanol
14. ethyl benzene
15. p-xylene
16. m-xylene
17. 1-butanol
18. o-xylene



# Break Number 1

For Questions and Answers  
Press \*1 on Your Phone to  
Ask a Question



# Phase Stability: PLOT Column Enhancements

- thicker phase layers - improved resolution
- improved particle to particle cross-linking or binding:
- decreased particle generation
  - less detector "spiking" and noise
  - less detector contamination
- improved chemical stability- longer column life
- improved selectivities



# Bonded PLOT Columns

- **Virtually no spiking- use with valved, on-line, MS applications.**
- **No trap needed**



# Bonded Porous Polymer PLOT Columns

**General applications- hydrocarbons and light polars, including:**

- HC up to C10, C2 -C4 isomers (natural gas, refinery gas).
- Halogenated HC
- CO<sub>2</sub>, Air/CO, methane
- Polar compounds: Solvents, Alcohols, Ketones, Esters, Aldehydes, Amines, Acids





# Bonded Porous Polymer PLOT Columns

HP PLOT Q: Polystyrene divinylbenzene- based

HP PLOT U: Divinylbenzene/ethylene glycol methacrylate- based

- resistant to water & alcohols
- can quantitate water and alcohols
- rinseable



# Column Stability/ Reproducibility: Alcohol Injections

GC: HP 6890

Detector: FID

Column: HP-PLOT Q


Analyte

Oven: 50°C (2 min), to  
250°C at 15°C/min


Concentration: 100 ppm (v/v)

Carrier: Helium, 25 psig

1. Methanol
2. Ethanol
3. Propanol
4. 1-Butanol
5. 2-Methyl-1-propanol
6. 1-Pentanol
7. 4-Methyl-2-pentanol
8. 2-Ethyl-1-Butanol




# Analysis of Acetone on HP PLOT U



# Selectivity comparison: PLOT Q vs U

Column: 0.53mm x 30m HP-PLOT/Q  
Oven: 150°C  
Carrier: Hydrogen 9ml/min  
Injection: 250C Split mode(12:1)  
Detector: FID 250C  
Sample: 1ul 0.1-0.3% solvents  
in methanol



- 1. Ethanol
- 2. Acetonitrile
- 3. Acetone
- 4. Dichloromethane
- 5. Diethylether
- 6. Pentane
- 7. Ethyl acetate
- 8. Hexane
- 9. Benzene



# Selectivity comparison: PLOT Q vs U


Column: HP-PLOT U 0.32mm x 30m x 10 $\mu$ m

Oven: 150C


Inlet: 180C, Split ratio 30;1

Detector: FID 200C

Sample: 1ul liquid



# Volatile Primary Amines




**Column: HP-PLOT U 0.53mm x 30m x 20um. Oven:150C(1min) 10C/min to 190C (5min). Carrier: Hydrogen, 7ml/min. Inlet: 180C, Split ratio 30;1. Detector: TCD 200C. Sample: 2ul amines in methanol, 1-10%.**



# Analysis of Noble & Fixed Gases Using HP PLOT MoleSieve

Column: HP-PLOT/MoleSieve  
30 m x 0.53 mm x 50  $\mu$ m  
HP part no. 19095P-MS0  
Carrier: Helium, 4 ml/min  
Oven: 35°C (3min) to 120°C  
(5 min) at 25°C/min  
Sample: 250  $\mu$ l, split (ratio 50:1)



- 1. Neon
- 2. Argon
- 3. Oxygen
- 4. Nitrogen
- 5. Krypton
- 6. Xenon



# HP Plot Al2O3

Three types, "S" "M" and KCL deactivated- will resolve hydrocarbon isomers.


- "S" and "M" are excellent general use columns. Baseline resolution of C1-C6; optimized for resolution of C2-C4 isomers.
- Design goal for "S" & "M": RI for ethylene & propylene as close to midpoint of flanking peaks as possible.
- "S" best for resolving acetylene from butane & propylene from isobutane
- "M" best for resolving cyclopropane from propylene (otherwise, "S" & "M" interchangeable)
- "KCL" best for dienes





# Comparison of PLOT Al2O3 Columns


Ethylene+Propylene on Competitive Al2O3 Column




Carrier: He, constant flow  
 Oven: TP1=35°C (2 min), 5°C/min to 100°C, 10°C/min to 180°C (5 min)  
 Detector: FID (250°C)  
 Sample size: Overload Conditions

- Methane
- Ethane
- Ethylene
- Propane
- Cyclopropane
- Propylene
- iso-Butane
- n-Butane
- Propadiene
- Acetylene
- trans-2-Butene
- Butene-1
- 13. iso-Butylene
- 14. cis-2-Butene
- 15. Iso-Pentane
- 16. n-Pentane
- 17. 1,3-Butadiene
- 18. Propyne
- 19. trans-2-Pentene
- 20. 2-Methyl-2-butene
- 21. Pentene-1
- 22. cis-2-Pentene
- 23. n-Hexane

Ethylene on HP Plot Al2O3 "S" deactivation



Propylene on HP Plot Al2O3 "M" deactivation



# GS-GasPro

- **Silica-based PLOT column**
- **“Bonded” stationary phase**
- **Unique selectivity**
- **Highly inert to reactive compounds**
- **Retention stability unaffected by water**
- **Mass Spec friendly PLOT column**




## Fantastic for Freons

### Halocarbons

**Column:** GS-GasPro  
30 m x 0.32 mm I.D.  
**J&W P/N:** 113-4332  
**Carrier:** Helium at 30 cm/sec  
**Oven:** 130°C for 4 min  
130-225°C at 10°/min  
225°C Hold  
**Injector:** Split 1:67, 1 µL, 250°C  
**Detector:** FID, 250°C

1. CH<sub>4</sub>
2. CHClF<sub>2</sub> (Freon 22)
3. CCl<sub>2</sub>F<sub>2</sub> (Freon 12)
4. ClCF<sub>2</sub>CF<sub>2</sub>Cl (Freon 114)
5. CHCl<sub>2</sub>F (Freon 21)
6. CCl<sub>3</sub>F (Freon 11)
7. CF<sub>2</sub>Br<sub>2</sub> (Freon 12B2)
8. CH<sub>3</sub>I
9. CH<sub>2</sub>Cl<sub>2</sub>
10. *trans*-ClCH=CHCl
11. CF<sub>3</sub>CCl<sub>3</sub> (Freon 113)
12. *cis*-ClCH=CHCl
13. CHCl<sub>3</sub>
14. ? from CCl<sub>4</sub>
15. CCl<sub>4</sub>
16. ? from CCl<sub>4</sub>
17. CH<sub>3</sub>CH<sub>2</sub>I
18. CH<sub>2</sub>Br<sub>2</sub>
19. CHCl<sub>2</sub>Br
20. C<sub>4</sub>F<sub>9</sub>I
21. CHClBr<sub>2</sub>
22. CH<sub>3</sub>CH<sub>2</sub>CH<sub>2</sub>I




C532

Time (min)



## Excellent sensitivity for sulfur gases!




# GS-CarbonPLOT

- **Monolithic carbon molecular sieve**
- **Phase formed in situ-no particles**
- **Extended temperature limit of 360°C**
- **Unique selectivity**
- **Mass Spec friendly**
- **Better price**



# Natural Gas Mixture



Carbon-PLOT




Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 49

## Break Number 2

For Questions and Answers  
Press \*1 on Your Phone to  
Ask a Question



# New Selectivities and Application Specific Columns

- HP Fast Residual Solvent Column
- DB-ALC1 & DB-ALC2; Blood Alcohols
- DB-MTBE; Resolves MTBE from the Methylpentanes
- DB-VRX; *The Volatiles Column*
- DB-HT SimDis; “Boiling Point” phase with 430C MAOT
- DB-TPH; Total Petroleum Hydrocarbons
- DB-Petro; PONA, PIANO and PNA analysis
- DB-Dioxin; dioxin analysis, resolves toxic congeners
- Chiral Columns





# Chiral Columns

- **It is nearly impossible to predict success of separation of enantiomers on any given column**
- **CycloSil-B is a good general purpose column (“DB-5” of chiral columns)**
- **Trial and Error: We have money-back guarantee if column doesn’t do separation**




# Chiral Columns

- **We have two different, but similar, types of chiral stationary phases. Both are derivatized beta-cyclodextrins**
- **CycloSil-B: silane modified**
- **Cyclodex-B: permethylated**
- **HP Chiral -10 and -20  $\beta$ : permethylated**



# Structure and Properties of $\beta$ - Cyclodextrin



**7 Glucose molecules joined through 1,4-glycosidic linkages**  
**Primary hydroxyl groups partially block the base**  
**Cavity size is about 6-6.5 Å in diameter**  
**Hydroxy groups are peralkylated**



# CycloSil-B

## Rosemary Oil

**Column:** CycloSil-B

30 m x 0.25 mm I.D., 0.25  $\mu$ m

**J&W P/N:** 122-6632


**Carrier:** Hydrogen at 40 cm/sec,  
measured at 60°C

**Oven:** 55°C for 1 min  
55-180°C at 5°/min

**Injector:** Split 50:1, 250°C

**Detector:** FID, 340°C

1. (+/-)  $\alpha$ -Pinene
2. (+/-) Camphene
3. (+/-)  $\beta$ -Pinene
4. (+/-) Limonene
5. (+/-) Linalool
6. (+/-) Camphor
7. (+/-) Terpinene-4-ol
8. (+/-) Isoborneol
9. (+/-) Borneol
10. (+/-)  $\alpha$ -Terpineol



C509




Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 55

# Chiral Separation- Oils & Fragrances



**HP Chiral-20  $\beta$**

**30m X 0.25mm X 0.25cm**

**P/N 19091G-B233**

Components:

1. nonane
2. (+) 3,3-Dimethyl-2-butanol
3. (-) 3,3-Dimethyl-2-butanol
4. (1S)-(-)-alpha-pinene
5. (1R)-(+)-alpha-pinene
6.  $\alpha$ -Terpinene
7. (S)-(-)-Limonene
8. (R)-(+)-Limonene
9.  $\gamma$ -Terpinene
10. (+/-)-Linalool
11. (-)-Menthone
12. (+/-)- $\alpha$ -Terpinen-4-ol

**Carrier:** Hydrogen, 36 cm/sec, Constant pressure

**Oven:** 50 C (2 min) to 143 C at 5 C/min

**Injection:** Split (ratio 30:1), 1  $\mu$ l, Inj Temp 250 C

**Detection:** FID, Det Temp 300 C

**Sample:** 0.25 ng/ $\mu$ l each analyte in Hexane




**Agilent Technologies**

**Dial 1-816-650-0621 for e-Seminar Audio**

**Slide 56**

# Headspace Analysis of Common Solvents on HP-624




**Column:** HP-624, 30 m x 0.53 mm x 3.0 um  
(HP Part No. 19095V-423)  
**Carrier:** Helium, 35 cm/sec, Constant Flow  
**Oven:** 40 C (5 min) to 90 C at 2 C/min  
90 C to 250 C at 30 C/min  
**Injection:** Headspace, 180 C, Split 7/1  
**Detector:** FID, 260 C

**Headspace Conditions:**  
Carrier Gas Pressure: 3.5 psi  
Oven Temp: 85 C  
Loop Temp: 95 C  
Transfer Line Temp: 110 C  
Vial Pressure: 10 psi

**Vial Equilibration Time:** 10 min  
**Pressurization Time:** 0.2 min  
**Loop Fill Time:** 0.15 min  
**Loop Equilibration Time:** 0.05 min  
**Injection Time:** 1.00 min



# Headspace Analysis of Common Solvents on Fast Residual Solvent Column



**GC:** 6890  
**Column:** HP-Fast Residual Solvent Column (HP Part No. 19095V-420)  
**Carrier:** Helium, 30 cm/sec, Constant Flow  
**Oven:** 40 C (1.7 min) to 90 C at 6 C/min  
 90 C to 250 C at 30 C/min  
**Injection:** Headspace, 180 C, Split 7/1  
**Detector:** FID, 260 C


**Headspace Conditions:**  
 Carrier Gas Pressure: 3.5 psi  
 Oven Temp: 85 C  
 Loop Temp: 95 C  
 Transfer Line Temp: 110 C  
 Vial Pressure: 10 psi

**Vial Equilibration Time:** 10 min  
**Pressurization Time:** 0.2 min  
**Loop Fill Time:** 0.15 min  
**Loop Equilibration Time:** 0.05 min  
**Injection Time:** 1.00 min



# DB-ALC1 & DB-ALC2

## Analysis of Blood Alcohols in 2 minutes!




**DB-ALC1**  
30m x 0.53mm I.D., 3µm

**GC Conditions**  
Carrier: Helium at 80 cm/sec,  
measured at 40°C  
Oven: 40°C Isothermal  
Injector: Split 1:10, 250°C  
Detector: FID, 300°C

**Headspace Conditions**  
Oven: 70°C  
Loop: 80°C  
Transfer line: 90°C  
Vial Equil. Time: 10 min  
Pressurization Time: 0.20 min  
Loop Fill Time: 0.20 min  
Loop Equil. Time: 0.05 min  
Inject Time: 0.2 min  
Sample Loop Size: 1.0 mL  
Sample Composition: 0.1% ethanol,  
0.001% others

1. Methanol
2. Acetaldehyde
3. Ethanol
4. Isopropanol
5. Acetone
6. 1-propanol




**DB-ALC2**  
30m x 0.53mm I.D., 2µm






# DB-MTBE



**Resolves  
the Pentanes  
from  
MTBE!**



# DB-HT Sim Dis




# DB-VRX: The Volatiles Column

**Column:** DB-VRX  
**75 m x mm.45 mm I.D., 2.55 µm**  
**J&W P/N:** 124-1574  
**Carrier:** Helium at 9 mL/min, measured at 35°C  
**Oven:** 35°C for 12 min  
 35-60°C at 5°/min  
 60°C for 1 min  
 60-200°C at 17°/min  
 200°C for 5 min  
**Injector:** Purge and trap (O.I.A 4560)  
 J&W LVI, 150°C  
 20 ppb per component in 5 mL water  
**Trap:** VOCARB™3000  
**Detector:** ELCD (O.I.A. 4430), with NiCat™ reaction tube in the halogen mode, 950°C reactor temperature

1. Dichlorodifluoromethane
2. Chloromethane
3. Vinyl chloride
4. Bromomethane
5. Chloroethane
6. Trichlorofluoromethane
7. 2-Chloropropane (IS)

- |                                       |                                    |  |
|---------------------------------------|------------------------------------|--|
| 8. 1,1-Dichloroethene                 | 38. 1,2-Dibromoethane              | 63. Benzyl chloride                      |
| 9. Iodomethane                        | 39. 1-Chloro-3-fluorobenzene (IS)  | 64. <i>n</i> -Butylbenzene               |
| 10. Allyl chloride                    | 40. Chlorobenzene                  | 65. 1,2-Dichlorobenzene                  |
| 11. Methylene chloride                | 41. 1,1,1,2-Tetrachloroethane      | 66. <i>Bis</i> (2-Chloroisopropyl) ether |
| 12. <i>trans</i> -1,2-Dichloroethene  | 42. Ethylbenzene                   | 67. 1,2-Dibromo-3-chloropropane          |
| 13. 1,1-Dichloroethane                | 43. <i>m</i> -Xylene               | 68. 1,2,4-Trichlorobenzene               |
| 14. Chloroprene                       | 44. <i>p</i> -Xylene               | 69. Hexachlorobutadiene                  |
| 15. <i>cis</i> -1,2-Dichloropropane   | 45. Styrene                        | 70. Naphthalene                          |
| 16. 2,2-Dichloropropane               | 46. <i>o</i> -Xylene               | 71. 1,2,3-Trichlorobenzene               |
| 17. Bromochloromethane                | 47. Bromoform                      |  |
| 18. Chloroform                        | 48. Isopropylbenzene               |  |
| 19. 1,1,1-Trichloroethane             | 49. <i>cis</i> -1,4-Dichlorobutene |  |
| 20. Carbon tetrachloride              | 50. 1,1,2,2-Tetrachloroethane      |  |
| 21. 1,1-Dichloropropene               | 51. Bromobenzene                   |  |
| 22. Benzene                           | 52. 1,2,3-Trichloropropane         |  |
| 23. 1,2-Dichloroethane                | 53. <i>n</i> -Propylbenzene        |  |
| 24. Fluorobenzene (IS)                | 54. 2-Chlorotoluene                |  |
| 25. Trichloroethene                   | 55. 1,3,5-Trimethylbenzene         |  |
| 26. 1,2-Dichloropropane               | 56. 4-Chlorotoluene                |  |
| 27. Dibromomethane                    | 57. <i>tert</i> -Butylbenzene      |  |
| 28. Trifluorotoluene (IS)             | 58. 1,2,4-Trimethylbenzene         |  |
| 29. Bromodichloromethane              | 59. <i>sec</i> -Butylbenzene       |  |
| 30. 2-Chloroethyl vinyl ether         | 60. 1,3-Dichlorobenzene            |  |
| 31. <i>cis</i> -1,3-Dichloropropene   | 61. <i>p</i> -Isopropyltoluene     |  |
| 32. Toluene                           | 62. 1,4-Dichlorobenzene            |  |
| 33. <i>trans</i> -1,3-Dichloropropane |                                    |  |
| 34. 1,1,2-Trichloroethane             |                                    |  |
| 35. Tetrachloroethene                 |                                    |  |
| 36. 1,3-Dichloropropane               |                                    |  |
| 37. Dibromochloromethane              |                                    |  |




Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 62

# DB-Petro50

Exceeds all ASTM and AFNOR requirements for PONA analysis



# DB-TPH

## High Throughput TPH Screening

### Purge and Trap of BTEX Standard using FID

Column: **DB-TPH**

30 m x 0.45 mm I.D., 1.0 µm

J&W P/N: **124-1632**

Carrier: Helium at 67 cm/sec, measured at 40°C

Oven: 40°C for 2 min

40-260°C at 12°/min

Injector: Purge and trap (O.I.A. 4560)

J&W Low Volume Injector (LVI)

Purge: 15 min at 85°C


Trap: BTEX (Supelco) at 50°C during purge

Desorb: 2 min at 260°C

Detector A: FID, 250°C

1. 2-Methylpentane
2. Benzene
3. Isooctane
4. Trifluorotoluene (IS)
5. Toluene
6. Ethylbenzene
7. *m,p*-Xylene
8. *o*-Xylene
9. Isopropylbenzene (IS)
10. 1,3,5-Trimethylbenzene
11. 1,2,4-Trimethylbenzene
12. Naphthalene

FID



C382

### Purge and Trap of BTEX Standard using PID

Column: **DB-TPH**

30 m x 0.45 mm I.D., 1.0 µm

J&W P/N: **124-1632**

Carrier: Helium at 67 cm/sec, measured at 40°C

Oven: 40°C for 2 min

40-260°C at 12°/min

Injector: Purge and trap (O.I.A. 4560)

J&W Low Volume Injector (LVI)

Purge: 15 min at 85°C

Trap: BTEX (Supelco) at 50°C during purge


Desorb: 2 min at 260°C

Detector B: PID, 200°C

Nitrogen makeup gas  
at 30 mL/min

1. 2-Methylpentane
2. Benzene
3. Isooctane
4. Trifluorotoluene (IS)
5. Toluene
6. Ethylbenzene
7. *m,p*-Xylene
8. *o*-Xylene
9. Isopropylbenzene (IS)
10. 1,3,5-Trimethylbenzene
11. 1,2,4-Trimethylbenzene
12. Naphthalene

PID



547



# DB-Dioxin

## Dioxins and Furans

**Column:** DB-Dioxin

60 m x 0.25 mm I.D., 0.15 µm

**J&W P/N:** 122-2461

**Carrier:** Helium at 34.3 cm/sec, measured at 250°C



**Oven:** 180°C for 1 min

180-270°C at 2.5°/min

270°C for 40 min

1. 1,3,6,8-TCDD
2. 2,3,7,8-TCDD
3. 1,2,8,9-TCDD
4. 1,2,4,6,8-/1,2,4,7,9-PeCDD
5. 1,2,3,7,8 + unknown-PeCDD
6. 1,2,4,8,9-PeCDD
7. 1,2,4,6,7,9-/1,2,4,6,8,9-HeCDD
8. 1,2,3,4,7,8-HeCDD
9. 1,2,3,4,6,9-HeCDD
10. 1,2,3,6,7,8-HeCDD
11. 1,2,3,7,8,9-HeCDD
12. 1,2,3,4,6,7-HeCDD
13. 1,2,3,4,6,7,9-HpCDD
14. 1,2,3,4,6,7,8-HpCDD
15. OcCDD
16. 1,3,6,8-TCDF
17. 2,3,7,8-TCDF
18. 1,2,8,9-TCDF
19. 1,3,4,6,8-PeCDF
20. 1,2,3,4,8-PeCDF
21. 1,2,3,7,8-PeCDF
22. 1,2,3,4,6-PeCDF
23. 2,3,4,7,8-PeCDF
24. 1,2,3,8,9-PeCDF
25. 1,3,4,6,7,9-HeCDF
26. 1,2,3,4,7,8/1,2,4,6,8,9-HeCDF
27. 1,2,3,6,7,8-HeCDF
28. 2,3,4,6,7,8-HeCDF
29. 1,2,3,7,8,9-HeCDF
30. 1,2,3,4,8,9-HeCDF
31. 1,2,3,4,6,7,8-HpCDF
32. 1,2,3,4,7,8,9-HpCDF
33. OcCDF

## Tetra Isomers



## Penta Isomers





# High Resolution Megabore

- **Same Outer Diameter as the Megabore**
- **No Special Hardware Required**
- **Smaller Inner Diameter (0.45mm)**
- **Maintain Phase Ratio (Beta)**
- **Methods are easy to translate!**



# High Resolution Megabore

## Same Resolution - Faster Analysis!



**BTEX**

**Carrier: Helium**

**Oven : 40°C for 3 min, 5°/min to 100°C**



**Agilent Technologies**

**Dial 1-816-650-0621 for e-Seminar Audio**

**Slide 67**



# Summary

- **Column deactivation chemistries improved.**
  - **enhanced acidic and basic analyte sensitivities**
  - **column chemistries optimized for specific applications**
- **Stationary phases optimized to extend max temp limit & lower bleed, extend minimum temperature limits, and improve physical & chemical stability.**
- **New selectivities/ applications**



## Wrap-up E-Seminar Questions

**Thank you for attending Agilent e-Seminars.  
Our e-Seminar schedule is expanding every  
week. Please check our website frequently at:  
[www.agilent.com](http://www.agilent.com)**

