

Agilent CytoGenomics 1.5

Streamlining your cytogenetic workflow



Arne Ijpma, Ph.D. - Product Manager

arne_ijpma@agilent.com



Agilent CytoGenomics 1.5

Take home message

Agilent CytoGenomics software is a workflow based easy to use analysis tool for CGH and CGH+SNP analysis to detect CNC, LOH and UPD

Agilent CytoGenomics 1.5

Constitutional and cancer sample analysis

CGH analysis to detect CNCs:

- Constitutional samples
- Cancer samples: Liquid, Solid, FFPE, up to 10% mosaicism can be detected

SNP analysis to detect LOH/UPD and confirm CGH data:

- Constitutional samples only, no cancer samples and no mosaic detection

* Next version of CytoGenomics software will support SNP analysis for hematological cancer samples with ability to detect as low as 20% clonal fraction

Agilent CytoGenomics 1.5

Positioning

- CytoGenomics is designed for the cytogenetics researcher who needs a simple streamlined software solution for cytogenetic analysis of human samples
- Agilent Genomic Workbench (AGW, previously DNA analytics) is designed for many aspects of DNA research of both human and non human samples and for the researcher who needs the ability of on the fly tweaking of all aspects of his (cytogenetics) workflow and who is not interested in a simple streamlined solution

Agilent CytoGenomics 1.5


Result of the cytogenetic analysis: Cyto Report

Cytos Report						
Sample Information						
Array ID	:	252983010002_1_1				
Global Display Name	:	US23502418_252983010002_S01_CGH_1010_Aug_10_1_1				
Green Sample	:	European Male (NA12891_v1)				
Red Sample	:	NA09208				
Polarity	:	1				
DerivativeOfLogRatioSD	:	0.167981				
Amp/Del Intervals Table						
Location	Size	Cytoband	#Probes	Amp/Del	P-value	Annotations
chr17:24457-5901054	5,876,598	p13.3 - p13.2	340	-0.866466	0.00E00	DOC2B, RPH3AL, C17orf97, FAM101B, VPS53, FAM57A, GEMIN4, ELP2P, GLOD4, RNMTL1, NXN, TIMM22, ABR, BHLHA9, TUSC5, YWHAE, CRK, MYO1C, INPP5K, LOC100306951, PITPNA, SLC43A2, SCARF1, RILP, PRPF8...
<small>Amp=Amplification Del=Deletion</small>						
Total Amp/Del Intervals: 1						
LOH Intervals Table						
Location	#Probes	LOH Score	Annotations			
chr17:72083-5853691	84	8.03E00	RPH3AL, C17orf97, FAM101B, VPS53, FAM57A, GEMIN4, ELP2P, GLOD4, RNMTL1, NXN, TIMM22, ABR, BHLHA9, TUSC5, YWHAE, CRK, MYO1C, INPP5K, LOC100306951, PITPNA, SLC43A2, SCARF1, RILP, PRPF8, TLCD2...			
Total LOH Intervals: 1						
Analysis Settings						
Aberration Algorithm	:	ADM-2	Design	:	029830_20100921	
Threshold	:	6.0	GC Correction	:	ON	
Window Size	:	2Kb	Centralization	:	ON	
Bin Size	:	10	Centralization Threshold	:	6.0	
SNP Copy Number	:	ON	SNP CN Confidence	:	0.95	
LOH	:	ON	Level	:		
Fuzzy Zero	:	OFF	LOH Threshold	:	6.0	
Combine Replicates(Intra Array)	:	ON	Nesting Level	:	OFF	
			Genome	:	hg19	
Aberration Filters	:	minProbes = 3 AND minAvgAbsLogRatio = 0.25 AND maxAberrations = 100000 AND percentPenetrance = 0	Feature Level Filters	:	glsSaturated = true OR risSaturated = true OR glsFeatNonUnifOL = true OR risFeatNonUnifOL = true	
Design Level Filters	:	NONE	Array Level Filters	:	NONE	
Metric Set Filters	:	NONE		:		

Agilent CytoGenomics 1.5

Value Proposition

- Streamlining the Cytogenetic data analysis workflow for CNC (copy number changes) and/or cnLOH / UPD detection of samples processed on Agilent SurePrint G3 CGH and CGH+SNP Microarrays.



Sample Information

Sample ID	2018010101_1_2
Client Study Name	152201018_2018010101_001_004_010_Agile_1_3_20180101
Client Sample	CGH+SNP_001_001_01
File Name	

Test Summary Report for Sample

Location	Scale	2-Content	4-Content	8-Content	16-Content	32-Content	64-Content
152201018_2018010101_001_004_010_Agile_1_3_20180101	1.0	1.0	1.0	1.0	1.0	1.0	1.0

SNP Z-Score Value for Sample

SNP ID	SNP Name	SNP Type	SNP Class	SNP Status	SNP Z-Score	SNP P-Value	SNP Q-Value
152201018_2018010101_001_004_010_Agile_1_3_20180101	152201018_2018010101_001_004_010_Agile_1_3_20180101	SNP	SNP	SNP	SNP	SNP	SNP

Analysis Settings

Algorithm	Agilent	Threshold	0.0
GC Correction	On	Window Size	200
Normalization	On	Bin Size	50
Copy Number Threshold	0.0	SNP Copy Number	On
SNP CN Threshold	0.0	LOH	On
LOH Threshold	0.0	Flurry Size	OFF
Missing Level	0	Control Threshold	On
Missing	Agilent	Aberration Filter	Agilent
Probe Level Filter	Agilent	Agilent	Agilent
Agilent Level Filter	Agilent	Agilent	Agilent
Control Threshold	Agilent	Agilent	Agilent

To enable the Cytogenetic Researcher to analyze a normal CGH (or CGH+SNP) sample in less than 5 minutes, including report generation.


Agilent CytoGenomics 1.5

Streamlined data analysis for Cytogenetic research labs


- Simple and quick workflow for data analysis
- Convenient input and output support allowing for integration with LIMS

CGH probes
Copy Number Changes

SNP probes
LOH/UPD



Run data analysis



Triage & sign off

Sample info,
export from LIMS

Laboratory LIMS

Cyto report,
upload to LIMS

Cyto Report

Sample Information					
Array ID	252983010001_1_2				
Global Display Name	US23502418_252983010001_S01_CGH_1010_Aug10_1_2_02022011				
Green Sample	European Male (NA12891_v1)				
Red Sample	-				

Text Summary Report for Sample US23502418_252983010001_S01_CGH_1010_Aug10_1_2_02022011						
Location	Size	Cytoband	#Probes	minCN	maxCN	Annotations
chr5:68112036-70536824	1,624,294	q11.2	57	0.296018	0.00600	OC1N, GTF2HC, GTF2ID, LOC100272115, GUSBP3, SERP1A, SERP1B, SMN2, SMN1, LOC100179939, LOC553391, GTF2H2B, NAPI, GTF2C, LOC547859
chr6:259318-362290	102,972	p25.3	14	0.450728	0.00600	GUSP22
chr6:31973370-32009282	35,912	p21.33	35	0.389420	0.00600	CYP21A2, TNXA, TNXB

SNP Table View for Sample US23502418_252983010001_S01_CGH_1010_Aug10_1_2_02022011					
Location	#Probes	LDH Score	Annotations		
chr1:7626484-7970067	87	7.81E00	ACAC4, BARO1B, SNCR4AC, SNCR4DA, SNCR4DB, MSH4, ASE17, STGALNAC3, STGALNAC5, PI3K, AK5, ZZZ3, USP33, FAM73A, NEXN, FUBP1, DNAJB4, GIPC2, MGC27362, PTF6F5		
chr1:17251901-197915297	517	4.60E01	C10orf1, FASLG, TNFSF18, TNFSF4, PRDX8, SLC5A11, ANKRD45, KLHL20, CEP350, DARS2, GASS, SNCRD81, SNCRD7, SNCRD9, SNCRD79, SNCRD78, SNCRD44, SNCRD77, SNCRD76, SNCRD75...		
chr1:223484657-240963558	189	1.64E01	KALY10A, KORY1, SLC39F3, C1orf51, TARBSP1, RIF2BP2, TOMM20, SNORA148, RBM34, ARID4B, GGP51, TECE, BGGALNT2, GNS4, LYST, MRR1537, NED1, GPR137B, ERO1LB, EDARADD...		
chr2:211536392-221570900	245	2.30E01	CPS1, ERBB4, MRS48F2, IKZF2, SPAG15, VWC2L, BARD1, ABCA2, ATRC, PVI, BREG, PCCR, TMEM169, XRC5C, MARCH4, SMARCAL1, RPL37A, IGFBP2, IGFBP5, TNF1...		

Analysis Settings		
Aberation Algorithm	ADAM-2	Threshold: 6.0
GC Correction	ON	Window Size: 2kb
Centralization	ON	Bin Size: 10
Centralization Threshold	6.0	SNP Copy Number: ON
SNP CN Confidence Level	0.95	LOH: ON
LOH Threshold	6.0	Fuzzy Zero: OFF
Nesting Level	0	Combine Replicates/Arrays: ON
Genome	hg19	Array:
Feature Level Filters	glsSaturated = true OR rlsSaturated = true OR glsFeatNonUniqOL = true OR rlsFeatNonUniqOL = true	Aberation Filters: minProbes = 6 AND minAvgAbsLogRatio = 0.0 AND maxAberrations = 100 AND percentPenetration = 10
Array Level Filters	NONE	Design Level Filters: NONE
Genomic Boundaries	Not Applied	Metric Set Filters: NONE

Agilent CytoGenomics 1.5

Workflow setup and user roles

- Initial set up: Configuring the Workflow

1. Configure analysis method

2. Configure CytoReport

3. Configure Workflow

- Day to day usage:
 - **Technician (manual) / Auto-processing**
 1. Upload TIFF images (or FE files)
 2. Map array to sample attributes
 3. Engage **configured workflow** to run FE (TIFF only) and analysis
 - **Scientist**
 1. Triage the sample and sign off
 2. Generate the CytoReport

Agilent CytoGenomics 1.5



Powerful algorithms for calling CNC, LOH and UPD

Log₂ ratio of intensities to call CNCs

No. of uncut alleles to call LOH/UPD

CGH probes
Copy Number Changes

SNP probes
LOH/UPD



Feature Extraction


ADM-1
ADM-2

SNP CN & LOH LOH calling

Amp/Del calling

Agilent CytoGenomics 1.5

Using SNPs to detect LOH and UPD



- Genotype SNPs using **restriction digestion** (Alu/Rsa)
- We measure the copy number of one allele at each SNP site relative to a **known reference**
- Regions of LOH are located by finding genomic regions with a **statistically significant scarcity of heterozygous call**
- **~5-10 Mb LOH/UPD resolution** across the entire genome

Genomic status

Normal diploid genome

Diploid genome with copy-neutral LOH or UPD

Hemizygous LOH

Amplification: e.g. trisomy

Genotype

AA, AB, BB

AA, BB

A, B

AAA, AAB, ABB, BBB


No. uncut allele

0, 1, 2

0, 2


0, 1

0, 1, 2, 3




Agilent CytoGenomics 1.5

Example 2: Sample with one copy deletion on chr9 p-arm



Agilent CytoGenomics 1.5

Sample Triage – Access to external databases



Preloaded track:
DGV-CNV


Link out from aberrations and associated genes to external Databases: UCSC, DGV, OMIM, Entrez

View and compare aberrations in DGV and UCSC genome browser

Agilent CytoGenomics 1.5

Sample Triage – Using CytoGenomics Database


Auto updated tracks based on aberration classification in CytoGenomics Database



Right click on 'Classification' allows query for overlapping aberrations in CytoGenomics database: results will be shown as new track along the gene view

Agilent CytoGenomics 1.5

Sample Triage – Sample and Aberration annotation



Add notes to sample


Add notes to aberrations

Suppress aberrations

Classify aberrations

Agilent CytoGenomics 1.5

Multi Sample View



Comparing 2 or more samples side by side

Agilent CytoGenomics 1.5

Additional features

- Workflow Automation Mode for automatic processing of TIFF images and report generation.
- Data upload to ISCA DB via Cartagenia bench.

[Download CytoGenomics free trial @ https://earray.chem.agilent.com/earray/](https://earray.chem.agilent.com/earray/)

Agilent CytoGenomics 1.5

System specs

- Client/server system with database (multiple users can access database at the same time).
- Runs on PC (XP, Windows 7).
- Runs on Mac (OS X Leopard and OS X Snow Leopard).
- Supported on both 32 and 64 bit machines.
- Minimum of 4GB of memory needed for 32 bit (8GB for 64 bit).
- Feature extraction is integrated component of CytoGenomics software.

[Download CytoGenomics free trial @ https://earray.chem.agilent.com/earray/](https://earray.chem.agilent.com/earray/)