

SureSelect

The Leader in Target Enrichment

Genomic Solutions Division

Yong Yi

Marketing Director

Doug Roberts

Manager, DNA Applications

1


SureSelect Technology

2

Whole Exome Sequencing



3

Summary



Maximize Sequencing Efficiency

Required Throughput = Genome Size x Average Coverage



Target = 50Mb x 100 = 5Gb

Target = 5Mb x 100 = 500Mb

Target = 0.5Mb x 100 = 50Mb

Target = 50Kb x 100 = 5Mb

Develop designs for any sequencing capacity:
- High Throughput or Desktop



SureSelect – Pioneer in NGS Target Enrichment



ARTICLES

**nature
biotechnology**

Solution hybrid selection with ultra-long oligonucleotides for massively parallel targeted sequencing

Andreas Gnirke¹, Alexandre Melnikov¹, Jared Maguire¹, Peter Rogov¹, Emily M LeProust², William Brockman^{1,5}, Timothy Fennell¹, Georgia Giannoukos¹, Sheila Fisher¹, Carsten Russ¹, Stacey Gabriel¹, David B Jaffe¹, Eric S Lander^{1,3,4} & Chad Nusbaum¹

Targeting genomic loci by massively parallel sequencing requires new methods to enrich templates to be sequenced. We developed a capture method that uses biotinylated RNA 'baits' to fish targets out of a 'pond' of DNA fragments. The RNA is transcribed from PCR-amplified oligodeoxynucleotides originally synthesized on a microarray, generating sufficient bait for multiple captures at concentrations high enough to drive the hybridization. We tested this method with 170-mer baits that target > 15,000 coding exons (2.5 Mb) and four regions (1.7 Mb total) using Illumina sequencing as read-out. About 90% of uniquely aligning bases fell on or near bait sequence; up to 50% lay on exons proper. The uniformity was such that ~60% of target bases in the exonic 'catch', and ~80% in the regional catch, had at least half the mean coverage. One lane of Illumina sequence was sufficient to call high-confidence genotypes for 89% of the targeted exon space.

The development and commercialization of a new generation of increasingly powerful sequencing methodologies and instruments¹⁻⁴ have lowered the cost per nucleotide of sequencing data by several orders of magnitude. Within a short time, several individual human genomes have been tested on target sets complex enough to match the scale of current next-generation sequencing instruments. The first method, microarray capture^{9,12,13}, uses hybridization to arrays containing synthetic oligonucleotides that match the target

- Capture fragments with longest, most efficient 120-mer cRNA baits
- Probes can be designed to any regions of interest, samples can be multiplexed
- Easy to implement and compatible with validated automation solution

SureSelect – The Leader in Target Enrichment

COMPLETE Solution

Best **PERFORMANCE**

FLEXIBLE Designs

SureSelect - Most Complete Enrichment Solution

All Exon Designs

All Exon V5

All Exon V5+UTRs


Non Human Exomes



Custom Solutions

Custom DNA

Custom RNA



Automation

Targeted Panels

Human Methyl-Seq

Chrm. X

Human Kinome

+ Library Prep

Post-capture

SureSelect ^{XT}
Illumina


SureSelect ^{XT}
SOLiD

Pre-capture

SureSelect ^{XT2}
Illumina

SureSelect - Most Complete Workflow Solution

From Sample to Analysis



SureSelect – The Leader in Target Enrichment

COMPLETE Solution

Best **PERFORMANCE**

FLEXIBLE Designs

SureSelect Technology

The Best Performance

Core Technology	Benefits
Ultra-Long RNA Baits (120-mer) Binding strength RNA:DNA > DNA:DNA	Better Sensitivity Detect more SNP, InDels, CNV, fusions
	Better Workflow 16hr hybridization, other require up to 72hrs
	Better Allelic Balance Equal representation of both alleles

SureSelect: Enabling Scientific Discovery

>350 Publications...

Nucleic Acids Research Advance Access published June 29, 2010

Detection of inherited mutations for breast and ovarian cancer using genomic capture and massively parallel sequencing

Tom Walsh^a, Ming K. Lee^a, Silvia Casadei^a, Anne M. Thornton^a, Sunday M. Stray^a, Christopher Pennil^b, Alex S. Nord^a, Jessica B. Mandell^a, Elizabeth M. Swisher^b, and Mary-Claire King^{a,1}

Ruth C. Schierberl^{1,2},
Jenny Z. Zhang³, Maria E. Furtu⁴, Jose Delgado⁵, Chantal Fourni⁶, Bernard Brais⁶,
Michel Michaelides^{9,10,11}, Richard G. Weleber⁹ and Joseph J. Higgins^{1,*}

Michael R. Stratton^{1,5} & P. Andrew Futreal¹

Herrmann⁷, Kurt Zatloukal⁸, Hans Lehrach⁷, Michal R. Schweiger^{7*}

Frederick C. ...
Sian J. ...
Ruth C. ...
Victor ...
Ovarian ...
explore ...
purifica ...
four ge ...
kinase ...
other t ...
subunit ...
domain ...
mutati ...
of 42 ...
aberran ...

Julia T. ...
Ka ...
Jessica B. Mandell^a, Elizabeth M. Swisher^b, and Mary-Claire King^{a,1}

Roman K. Thomas^{2,7,8} and William Pao^{5,*}

8 OCTOBER 2010 VOL 330 SCIENCE www.sciencemag.org
Pin Yue, Ph.D., and Sekar Kathiresan, M.D.

Jacek Majewski^{1,4,*} and Loydie A. Jerome-Majewska^{1,5,*}

Agata Rumana², John M. Opitz², Epimat Levy-Lahad², Rachel E. Kivvit² and Mary-Claire King^{a,1}

Pierre Dechelotte³, Jacek Majewski^{1,3} and Nada Jabado^{3,4,*}

Amal Abu Rayyan², Suheir Loulus², Karen B. Avraham³, Mary-Claire King¹ and Moien Kanaan^{2,*}

Bert B A de Vries¹ & Joris A Veltman¹

F. Hansen², Stacie K. Loftus¹

... and many more to come

underlying rare mendelian disorders and will likely transform the genetic analysis of monogenic traits.

SureSelect – The Leader in Target Enrichment

COMPLETE Solution

Best **PERFORMANCE**

FLEXIBLE Designs


SureSelect - Most Flexible Enrichment Solution

Custom or Catalog Solutions for DNA and RNA targets

- Proven, Refined Design Algorithms
- ISO 13485 Reagents
- Validated, Scalable Automation Solution



SureDesign –
New, intuitive
software for
custom capture
designs




SureSelect – Flexible Solution

From High Throughput to Desktop Sequencers



SureSelect All Exon V5

- **Design**
- **Performance**
- **Cancer Applications**



Human All Exon V5 Design and Synthesis

- Updated Content: Refseq, Gencode, UCSC, TCGA, CCDS, miRBase
- Proprietary probe selection/placement algorithms based on GC, Tm, Entropy, Mapability
- Long 120mer probes printed in Agilent's microarray fabrication facility
 - V5 ~ 544K individual probes spanning 50.4 Mb
 - V5 +UTR ~ 882K individual probes spanning 74.6 Mb
- High GC content printed in higher replication to boost performance
- Matched for best performance with Agilent's library preparation
- DNA probes are converted into biotinylated RNA bait



Human All Exon V5 and V5+UTRs

**The Best Performing Exome
Ready for Sequencing the Next Day**

The Fastest Workflow


From sample to sequencing the next day




**Overnight Hybridization Enables Sequencer Ready
Samples the Next Day**

Human All Exon V5: Distribution of Read Depths


V5 (4Gb)



V5+UTR (6Gb)



Human All Exon V5 Performance Statistics




The Best Performance and Fast Workflow


Parameter	All Exon V5	All Exon V5+UTRs	Competitor N	Competitor I
Target Size	50Mb	75Mb	64Mb	62Mb
# Genes	21,522	21,522	>20,000	20,794
# Targeted Exons	357,999	359,555	-	201,121
CCDS	X	X	X	X
RefSeq	X	X	X	X
GENCODE	X	X	X	X
miRBase	X	X	X	X
TCGA	X	X		
UCSC	X	X		
Amount of Seq. (total)	4Gb	6Gb	6Gb	10Gb
Overall Workflow	1.5 days	1.5 days	4.5 days	3.5 days
Add Custom Content	Yes	No	Yes	No

The Best Performance

Outstanding coverage with less sequencing




Human All Exon V5: Coverage in Tumor Suppressor Genes



- Agilent displays higher depth of coverage with the same amount of sequencing (4Gb)

Human All Exon V5: Coverage in Tumor Suppressor Genes




- Agilent displays a higher depth of coverage with the same amount of sequencing (4Gb)

Example Studies: Target Enrichment in Cancer

- **Target Enrichment for Novel Discovery**
Diffuse Large B cell Lymphoma

- **Intratumor Heterogeneity**
Renal Cell Carcinoma



Diffuse Large B cell Lymphoma (DLBCL)

- Aggressive Non-Hodgkin lymphoma
- Affects 30,000 new patients in the US every year
- Current Standard of Care –R-CHOP
 - Rituximab (anti CD20), cyclophosphamide, doxorubicin, vincristine, and prednisone
- 3 year event-free survival ~60% , remaining 40% die of disease

Exome sequencing of Non-Hodgkin Lymphoma

PNAS PNAS PNAS

Discovery and prioritization of somatic mutations in diffuse large B-cell lymphoma (DLBCL) by whole-exome sequencing

Jens G. Lohr^{a,b}, Petar Stojanov^{a,b}, Michael S. Lawrence^a, Daniel Audlair^a, Bjoern Chapuy^b, Carrie Sougnez^a, Peter Cruz-Gordillo^a, Birgit Knoechel^{a,b,c}, Yan W. Asmann^d, Susan L. Slager^d, Anne J. Novak^d, Ahmet Dogan^d, Stephen M. Ansell^d, Brian K. Link^e, Lihua Zou^a, Joshua Gould^a, Gordon Saksena^a, Nicolas Stransky^a, Claudia Rangel-Escareño^f, Juan Carlos Fernandez-Lopez^f, Alfredo Hidalgo-Miranda^f, Jorge Melendez-Zajgla^f, Enrique Hernández-Lemus^f, Angela Schwarz-Cruz y Celis^f, Ivan Imaz-Rosshandler^f, Akinyemi I. Ojesina^a, Joonil Jung^a, Chandra S. Pedomallu^a, Eric S. Lander^{a,g,h,1}, Thomas M. Habermann^d, James R. Cerhan^d, Margaret A. Shipp^b, Gad Getz^a, and Todd R. Golub^{a,b,g,i}

^aEli and Edythe Broad Institute, Cambridge, MA 02412; ^bDana-Farber Cancer Institute, Boston, MA 02115; ^cMayo Clinic College of Medicine, Rochester, MN 55902; ^dChildren's Hospital Boston, Boston, MA 02115; ^eUniversity of Iowa College of Medicine, Iowa City, IA 52245; ^fInstituto Nacional de Medicina Genómica, 14610 Mexico DF, Mexico; ^gHarvard Medical School, Boston, MA 02115; ^hMassachusetts Institute of Technology, Cambridge, MA 02142; and ⁱHoward Hughes Medical Institute, Chevy Chase, MD 20815

Contributed by Eric S. Lander, December 29, 2011 (sent for review November 22, 2011)

To gain insight into the genomic basis of diffuse large B-cell lymphoma (DLBCL), we performed massively parallel whole-exome sequencing of 55 primary tumor samples from patients with DLBCL and matched normal tissue. We identified recurrent mutations in genes that are well known to be functionally relevant in DLBCL, including *MYD88*, *CARD11*, *EZH2*, and *CREBBP*. We also identified somatic mutations in genes for which a functional role in DLBCL has not been previously suspected. These genes include *MEF2B*, *MLL2*, *BTG1*, *GNA13*, *ACTB*, *P2RY8*, *PCL0*, and *TNFRSF14*. Further, we show that *BCL2* mutations commonly occur in patients with *BCL2/IgH* rearrangements as a result of somatic hypermutation normally

somatic mutations in oncogenes and tumor-suppressor genes in B lymphocytes.

Traditionally, DLBCL has been classified by the morphology and immunophenotype of the malignant B-cells but more recently, molecular classifications have been reported. Specifically, gene expression-based classification of DLBCL has been proposed (5, 6), and the prognostic relevance for this has been demonstrated (7). It has been suggested that distinct signal transduction pathways are affected in the subtypes that are defined in this way, and that certain genetic defects preferentially occur in specific subtypes defined by the presumed cell of origin of the tumors (8–12).

PNAS 109 3879-3884 (2012)

Exome sequencing of DLBCL and Follicular Lymphoma (FL)

Frequent mutation of histone-modifying genes in non-Hodgkin lymphoma

Ryan D. Morin^{1*}, Maria Mendez-Lago^{1*}, Andrew J. Mungall¹, Rodrigo Goya¹, Karen L. Mungall¹, Richard D. Corbett¹, Nathalie A. Johnson², Tesa M. Severson¹, Readman Chiu¹, Matthew Field¹, Shaun Jackman¹, Martin Krzywinski¹, David W. Scott², Diane L. Trinh¹, Jessica Tamura-Wells¹, Sa Li¹, Marlo R. Firme¹, Sanja Rogic², Malachi Griffith¹, Susanna Chan¹, Oleksandr Yakovenko¹, Irmtraud M. Meyer³, Eric Y. Zhao¹, Duane Smailus¹, Michelle Moksa¹, Suganthi Chittaranjan¹, Lisa Rimsza⁴, Angela Brooks-Wilson^{1,5}, John J. Spinelli^{6,7}, Susana Ben-Neriah², Barbara Meissner², Bruce Woolcock², Merrill Boyle², Helen McDonald¹, Angela Tam¹, Yongjun Zhao¹, Allen Delaney¹, Thomas Zeng¹, Kane Tse¹, Yaron Butterfield¹, Inanç Birol¹, Rob Holt¹, Jacqueline Schein¹, Douglas E. Horsman², Richard Moore¹, Steven J. M. Jones¹, Joseph M. Connors², Martin Hirst¹, Randy D. Gascoyne^{2,8} & Marco A. Marra^{1,9}


Follicular lymphoma (FL) and diffuse large B-cell lymphoma (DLBCL) are the two most common non-Hodgkin lymphomas (NHLs). Here we sequenced tumour and matched normal DNA from 13 DLBCL cases and one FL case to identify genes with mutations in B-cell NHL. We analysed RNA-seq data from these and another 113 NHLs to identify genes with candidate mutations, and then re-sequenced tumour and matched normal DNA from these cases to confirm 109 genes with multiple somatic mutations. Genes with roles in histone modification were frequent targets of somatic mutation. For example, 32% of DLBCL and 89% of FL cases had somatic mutations in *MLL2*, which encodes a histone methyltransferase, and 11.4% and 13.4% of DLBCL and FL cases, respectively, had mutations in *MEF2B*, a calcium-regulated gene that cooperates with CREBBP and EP300 in acetylating histones. Our analysis suggests a previously unappreciated disruption of chromatin biology in lymphomagenesis.

Nature 476 298-303 Aug 2011

DLBCL Exome Sequencing Details

- *Lohr et al.*
 - 55 patients paired tumor and germline(normal) DNA
 - Captured with Agilent's exome capture
 - 150-fold mean sequence coverage of targeted exons
- *Morin et al.*
 - 117 patients FL and DLBCL tumor DNA
 - Captured with Agilent's exome capture

DLBCL mutation frequency



PNAS 109 3879-3884 (2012)

Agreement between two studies

- Common Pathways and Genes Identified in Both Studies

- Known contributing genes: *CD79B*, *TP53*, *CARD11*, *MYD88*, *EZH2*
- Suggested Pathogenic role: *MLL2*, *TNFRSF14*, *BTG1*, *MEF2B*, *GNA13*
- Selected Common Novel findings: *HIST1H1C*, *PIM1*
- Selected Novel Findings (*Lohr et al*): *PCLO*, *TMSL3*,
- Selected Novel Findings (DLBCL, *Morin et al*): *SGK1*, *BCL2*

Frequent nonsynonymous mutations in chromatin remodeling genes, suggesting potential cooperation in lymphomagenesis

-HDAC inhibitors under investigation (phase I/II)

Nature 476 298-303 Aug 2011

PNAS 109 3879-3884 (2012)

Conclusions

- High Coverage of exonic content in two large cohorts of DLBCL
- Confirmed mutations in known and novel sites
- Excellent agreement between studies
- Spurred new efforts to understand potential previously unknown drivers
- Progress toward new therapeutic avenues and potentially increased survival rates

Discovery Using Target Enrichment: Renal Cell Carcinoma

Aim:



To characterize *intra*-tumor heterogeneity that may foster tumor evolution/adaptation

Method:

SureSelect Exome sequencing of:


- 4 patient samples, 7-15 biopsies
- Pre-treatment (Primary vs. Metastasis)
- Intratumor sections (R1-R9)
- Multiple metastases

Gerlinger et.al. N Eng J Med 2012 366 (883-892)



C Phylogenetic Relationships of Tumor Regions


- Ubiquitous
- Shared primary
- Shared metastasis
- Private



Discovery Using Target Enrichment: Renal Cell Carcinoma

Results/Conclusions:

- Identified driver mutations (*VHL*, *KDM5C*, *SETD2*, *MTOR*, *PTEN*)
- Only *VHL* was mutated ubiquitously in all analyzed regions
- Different *regions* of the tumor have different mutations in the very *same driver gene*, e.g. *MTOR*, *SETD2*
- Intra-tumor heterogeneity may contribute to failure of targeted chemotherapy
- Suggests potential need for multiple assessments of tumor



Summary

- SureSelect exome target enrichment
 - Optimized probe design/production for improved capture
 - High coverage of exons with only 4Gb of sequencing
 - Samples ready for sequencing the next day
- Enables large cancer sequencing studies
 - Allows in-depth determination of nucleotide variations
 - Driving discovery of new potential targets and therapies
 - Increasing understanding of tumor evolution



Thank You!