



# Filtration Approaches to Sample Preparation and Cleanup

---

Alexander Ucci  
Application Engineer  
Sample Preparation Products  
[spp-support@agilent.com](mailto:spp-support@agilent.com)

# Today's Agenda

How does filtration work – a brief overview

Why filter – benefits of filtration, including applications

Introduction to Agilent filtration products

Selecting the right filtration product

Questions and wrap up


# MECHANISMS OF FILTRATION: A BRIEF OVERVIEW

Filtration

# Filtration – How does it work?


Filtration: separating substances based on particle size

- Solid particles in liquids (suspension)
- Insoluble droplets in liquids (emulsion)



# Microfiltration

## Two types of membranes




1. Depth filters: for heavily particle-laden samples, ideal for pre-filtration

---

Glass Fiber

---



2. Microporous filters: high safety against breakthrough, for complete removal of particles

---

CA : Cellulose Acetate

RC : Regenerated Cellulose

PES : Polyethersulfone

PP: Polypropylene

PTFE : Polytetrafluorethylene


Nylon

---

# Dual-Depth Filter

## Captiva Dual-Depth Filtration

- 1) First layer is a larger-porosity membrane. It traps larger particles, which travel in non-linear path, thereby preventing clogging and creating more surface area in which to capture the particulates
- 2) The second, smaller-porosity membrane, traps smaller particles resulting in all unwanted particles being removed from sample




(Close-up of  
filter layers)

# Mechanical Combined with Functional Filtration

## Captiva ND Lipids Filtration

- Combines a depth-filter for mechanical filtration with a functional component designed to remove lipids from samples



**Dual-depth filter schematic showing the differences in sizes between the two layers. A proprietary product offers a combination of mechanical and functional filtration for effective lipid removal.**

# WHY FILTRATION?

Filtration





# Filtration and Other Sample Preparation Techniques: An Overview

		<div> <div>More Specific</div> <div>←</div> <div>Instrument Separation and Detection Specificity</div> <div>←</div> <div>Less Specific</div> </div>							
		<div> <div>Less Specific</div> <div>→</div> <div>Sample Preparation Specificity</div> <div>→</div> <div>More Specific</div> </div>							
Sample Prep Technique Interference Removed	Dilute & Shoot	Filtration	Liquid/Liquid Extractions	Supported Liquid Extractions (SLE)	Dried Matrix Spotting	Precipitation filtration	QuEChERS	Precipitation-Lipid Removal 'Hybrid' Filtration	Solid Phase Extraction
Lipids	No	No	No	Some	No	No	Yes	Yes	Yes
Oligomeric Surfactants	No	No	No	No	No	No	No	Yes	Yes
Particulates	No	Yes	No	Some	No	Yes	Yes	Yes	Yes
Pigments	No	No	No	Some	No	No	Yes	No	Yes
Polar Organic Acids	No	No	Yes	Yes	No	No	Yes	No	
Proteins	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Salts	No	No	Yes	Yes	No	No	No	No	Yes
Filtration is suggested with any LC or GC method of sample preparation									

# Why Filter?

## Benefits of Filtration

Optimal  
Instrument Performance



Less  
System Downtime



Extend  
Column Lifetime



Greater  
Sample Integrity



Achieve lower detection limits

***Agilent recommends filtering prior to chromatographic analysis to remove particulates from the sample***



Agilent Technologies

# Why Filter the Sample?


Extreme Performance Requires Better Sample “Hygiene”




- Prevents blocking of capillaries, frits, and the column inlet (especially important for UHPLC)
- Results in less wear and tear on the critical moving parts of injection valves
- Results in less downtime of the instrument for repairs
- Produces improved analytical results by removing potentially interfering contamination
- “Functionalized” filtration can improve results over mechanical filtration by removing chemical interferences in addition to particulates




# Syringe Filter Benefit: Improved sub-2-micron LC Column Lifetime – Standardized Sample



# Syringe Filter Benefits: Improved sub-2 micron LC Column Lifetime – With Human Plasma Extract




# The Problem: Lipid Build-up and Ion Suppression in LC-MS Applications



If lipids are present in the sample and the gradient does not flush out lipids, they will build up as the injection number increases.

**More lipids → ion-suppression**


# Functional and Mechanical Filtration Benefit: Reduced Ion Suppression by Removing Lipids



Post-column infusion of albuterol (top trace)

Plasma sample showing elution of phospholipids (bottom trace)

➤ Ion suppression events correspond to elution of phospholipids



Post-column infusion of albuterol (top trace)

Sample prepared with Captiva ND Lipids to remove proteins, lipids, particulates (bottom trace)

➤ Ion suppression events are nearly absent – better results




# Functional and Mechanical Filtration Benefit: Reducing Backpressure Increases Over Time

## Conditions:

- 6,000 continuous injections
- Human plasma
- UHPLC instrument and column

## Results:


Filtration with Captiva ND Lipids led to minimal system backpressure build-up (~7% for entire 6,000 injections)





# Backpressure Build-up Comparison


Captiva ND Lipids



Backpressure increase during 1800 injections was 60 psi.

→ **53% less backpressure increase than PPT-centrifugation**

PPT followed by centrifugation




- Starting backpressure is higher than Captiva ND Lipids
- Backpressure increase during 1800 injections was 92 psi.




# Simple Methodology - Captiva ND Lipids

## Operating Instructions and tips for Captiva ND<sup>Lipids</sup> 96 -Well Plates



# Filtration Benefit: Pair Filtration with Other Sample Preparation Approaches to Improve Results




5 g of animal samples analyzed for veterinary drugs using a multi-step approach to sample cleanup incorporated classical SPE with functional/mechanical filtration (Captiva ND Lipids) for improved results



# Filtration Plus SPE – Improved Results and Cleaner Samples

Without Captiva  
ND<sup>Lipids</sup>

Filtered by Captiva  
ND<sup>Lipids</sup>



Sample clarity comparison – with  
and without Captiva ND Lipid  
filtration

Adding the  
functionalized filtration  
step improves  
recoveries dramatically

Veterinary Drug	Retention Time	Milk		Pork		Fish (Yellowtail)		Fish (baked Eel)	
		Without Captiva ND <sup>Lipids</sup>	Filtered by Captiva ND <sup>Lipids</sup>	Without Captiva ND <sup>Lipids</sup>	Filtered by Captiva ND <sup>Lipids</sup>	Without Captiva ND <sup>Lipids</sup>	Filtered by Captiva ND <sup>Lipids</sup>	Without Captiva ND <sup>Lipids</sup>	Filtered by Captiva ND <sup>Lipids</sup>
Sulfadiazine	7.21	77.9	97.5	109.5	105.8	64.2	103.3	62.1	81.9
Sulfathiazole	8.88	63.5	92.6	90.9	127.9	30.0	127.3	14.9	33.8
Sulfapyridine	9.53	67.9	74.3	75.8	94.9	42.1	96.0	38.2	60.5
Sulfamerazine	9.60	76.5	89.0	91.9	103.3	46.0	99.3	37.1	62.3
Sulfachlorpyridazine	9.70	57.6	74.5	65.1	93.9	28.2	90.8	23.4	40.9
Sulfamethoxazole	9.74	70.0	85.4	83.6	94.8	34.6	87.0	25.6	53.7
Sulfamethoxypyridazine	10.32	66.0	85.1	68.4	86.5	23.7	80.0	22.3	36.8
Thiamphenicol	10.55	38.1	88.4	154.5	117.1	72.3	135.1	32.7	101.2
Sulfadoxine	10.66	71.7	85.3	79.9	97.6	34.3	87.1	32.0	48.4
Sulfadimidine	11.20	66.4	89.3	81.1	96.8	31.3	88.1	26.8	50.5
Sulfamonomethoxine	11.24	68.0	79.1	72.8	90.8	28.0	92.8	21.8	41.6
5-Propylsulfonfyl-1 H-benzimidazole-2-amine	12.73	48.6	65.7	50.9	66.3	18.9	70.8	13.1	18.2
Sulfaquinoxaline	12.83	69.9	82.2	24.8	61.4	8.7	49.5	9.2	16.8
Levamisole	12.87	54.1	87.5	82.7	108.1	37.8	104.4	34.0	66.5
Sulfadimethoxine	12.99	58.4	69.2	48.1	78.7	22.5	66.2	20.9	33.3
Trimethoprim	13.02	61.0	85.9	86.4	99.3	38.5	95.1	39.3	63.1
Ormetoprim	13.70	56.2	82.6	75.0	88.9	36.1	90.3	27.1	54.2
Sulfantran	16.10	56.3	72.6	38.5	81.9	20.9	66.5	28.2	65.7
Thiabendazole	16.19	37.7	56.5	39.5	55.3	17.5	53.4	12.1	21.1
Ethopabate	16.48	55.3	76.4	49.0	78.6	28.9	70.9	17.8	38.0
alpha-Trenbolone	19.91	17.8	33.2	11.3	23.8	6.0	17.8	9.0	19.2
beta-Trenbolone	20.39	30.6	52.2	19.5	34.7	10.0	28.8	9.8	21.2
Zeranol	20.60	31.5	61.0	22.8	50.1	14.3	38.7	21.2	43.9
Melengestrol Acetate	23.00	35.6	86.4	24.4	90.5	4.7	55.3	15.2	62.6
Average Improvement		Milk	21.5%	Pork	20.0%	Yellowtail	49.8%	Eel	22.6%

Recovery as % response versus neat standard in 10% Acetonitrile (10ppb level), n=6.

Key: 0 - 30% 30 - 70% 70 - 120% > 120%



Agilent Technologies


# Filtration Benefits – Summary

- Filtration produces improved analytical results by removing potentially interfering contamination
- Functional filtration can improve results over mechanical filtration by removing chemical interferences in addition to particulates
- Protecting the instrument and maximizing analytical performance and uptime is straightforward with both mechanical and functional filtration approaches

# FILTRATION OPTIONS

Filtration

# Agilent Captiva Filtration Products: Setting a New Standard in Filtration



# Captiva Premium Syringe Filters



- Certified to be free of UV detectable extractables on HPLC
  - PES & Glass Fiber also certified for LC/MS
- Color-coded boxes for easy identification
- Comprehensive portfolio to meet your application needs

Premium Syringe Filters						
Membrane	Diameter / Pore size					
	4 mm		15 mm		25 mm (* 28 mm)	
	0.2 µm	0.45 µm	0.2 µm	0.45 µm	0.2 µm	0.45 µm
PTFE	♦	♦	♦	♦	♦	♦
Nylon			♦	♦	♦	♦
PES	♦	♦	♦	♦	♦	♦
Regenerated Cellulose	♦	♦	♦	♦	♦	♦
Cellulose Acetate					♦	♦
Glass Microfiber			♦		♦*	
Depth filters: Glass / PTFE			♦	♦	♦	♦*
Depth filters: Glass / Nylon			♦	♦	♦	♦





# Captiva Premium Syringe Filters: Advantages




Feature	Advantages	Benefits
Unique Syringe Filter construction	Cleanliness; certified to be free of UV-detectable extractables	<ul style="list-style-type: none"> <li>Minimizes introduction of contamination by the filter</li> <li>Reduce need for retesting</li> <li>Accuracy</li> </ul>
	Fastest flow rates	<ul style="list-style-type: none"> <li>Productivity / throughput</li> </ul>
	High loading capacity	<ul style="list-style-type: none"> <li>Productivity</li> </ul>
	Minimal analyte adsorption	<ul style="list-style-type: none"> <li>Maximum target analyte recovery</li> <li>Accuracy</li> <li>Reproducibility</li> </ul>
Ultrasonic weld and reinforced membrane	Withstands high pressure and force	<ul style="list-style-type: none"> <li>Rugged</li> <li>Reliable</li> </ul>


# Captiva Premium Syringe Filters: Fastest Flow Rates



Agilent **RC**, **Nylon** and **PTFE** filters compared to competitors




# Captiva Premium Syringe Filters: High Loading Capacity



- ✓ Due to nature of membrane, there is less clogging by particulates
- ✓ Assures more sample will be filtered

Capacity (volume) of 15 mm syringe filters over time (with Particulate-Laden Samples)



# Captiva Premium Syringe Filters: Performance




Agilent Captiva Syringe Filters provide consistent and higher than 90% filtration efficiency on particulate removal

	Agilent premium 0.2µm syringe filters						Agilent premium 0.45µm syringe filters					
	Nylon	PTFE	RC	PES	GF/NY	GF/PTFE	Nylon	PTFE	PES	CA	GF/NY	GF/PTFE
1	96.0	92.3	89.8	92.1	99	99.4	95.2	97	93.6	92.4	96.8	98.4
2	95.9	91.4	90.6	91.4	99	98.9	93.2	96.5	93.6	95.0	97.1	98.8
3	94.5	93.3	90.3	89.5	99.2	99.0	95.5	97.5	93.5	96.3	96.4	97.7
4	96.6	92.3	91.7	99.0	99.6	98.6	95.4	96.6	88.5	97.2	99.3	98.8
5	95.4	91.2	92.4	96.3	98.8	98.8	94.9	96.0	88.2	96	99.0	99.7
6	95.6	91.1	90.8	99.9	99.3	98.5	95.3	95.7	92.3	95.6	100	96.8
7	99.9	91.1	98.2	99.0	99.4	99.4	99.5	95.2	94.9	96.7	98.2	97.6
8	99.8	91.2	99.0	97.8	95.0	99.0	98.0	97.8	89.4	93.8	98.9	98.5
9	99.7	90.9	96.4	95.2	95.9	99.9	97.7	94.9	87.3	92.5	100.2	98.0
10	99.2	91.3	95.7	96.1	94.7	99.6	99.7	94.8	87.5	92.8	100.5	101.3
Average Eff (%)	97.3	91.6	93.5	95.6	98.0	99.1	96.4	96.2	90.9	94.8	98.6	98.6
RSD (%)	2.2	0.8	3.7	3.7	2.0	0.5	2.2	1.1	3.3	1.9	1.5	1.3

# Cleanliness of Agilent Captiva PES Filters



## ESI-Positive Mode



# Additional Tools to Help You Select the Right Syringe Filter



Membrane compatibility: online, in each box, in brochure

## Chemical Compatibility

LEGEND	
Compatible	••
Limited compatibility	•
Not compatible	-
Not analyzed	N/A

Contact time: 24 hours at 20 °C  
Chemical compatibility can be influenced by various factors. Therefore, we recommend that you confirm compatibility with the liquid you want to filter by performing a trial filtration run before you start your actual filtration.


FILTER	PP	PES	CA	PTFE	RC	Nylon	GF	MBS	PP
<b>Housing</b>									
<b>Solvents</b>									
Acetone	••	-	-	••	••	••	••	-	••
Acetonitrile	•	-	-	••	••	N/A	••	-	••
Benzene	-	-	•	••	••	••	••	-	••
Benzyl alcohol	••	-	-	••	••	••	••	-	••
n-Butyl acetate	N/A	-	-	••	••	••	••	-	••
n-Butanol	••	•	•	••	••	••	••	••	••
Carbon tetrachloride	•	-	-	••	••	••	••	-	••
Chloroform	•	-	-	••	••	••	••	-	••
Cyclohexane	••	-	-	••	••	••	••	-	••
Diethylacetamide	••	-	-	••	••	••	••	-	••
Diethyl ether	•	-	•	••	••	••	••	-	••
Dimethyl formamide	••	-	-	••	••	•	••	-	••
Dimethylsulfoxide	••	-	-	••	••	••	••	-	••
Dioxane	•	-	-	••	••	••	••	-	••
Ethanol, 98%	••	••	•	••	••	••	••	-	••
Ethyl acetate	•	-	-	••	••	••	••	-	••
Ethylene glycol	••	••	•	••	••	••	••	••	••
Formamide	N/A	••	-	••	••	••	••	••	••
Gasoline	•	•	•	••	••	••	••	••	••
Glycerin	••	••	•	••	••	••	••	•	••
n-Heptane	-	••	•	••	••	••	••	••	••
n-Hexane	-	••	•	••	••	••	••	•	••
Isopropanol	••	••	•	••	••	••	••	-	••
Isopropyl acetate	N/A	-	-	••	••	••	••	-	••
Methanol, 30%	••	••	N/A	••	••	••	••	••	••
Methanol, 98%	••	•	-	••	••	••	••	••	•
Methyl acetate	•	-	-	••	••	••	••	-	••
Methylene chloride	•	-	-	••	••	••	••	-	••
Methyl ethyl ketone	•	-	-	••	••	••	••	-	•
Methyl isobutyl ketone	•	-	-	••	••	••	••	-	•

\*CA and GF membranes in MBS housing for 28 mm size

## Syringe Filter Selection Guide: Online, in brochure

### Agilent Captiva Filtration SELECTION GUIDE

STEP 1



STEP 2



STEP 3



### Applications

Type of Filtration	Recommended	Alternatives
HPLC • UHPLC • LC/MS • GC	RC	PTFE or Nylon
ICP-MS	PTFE	Glass Fiber/PTFE (High Particle Samples)
CE	RC	Nylon
Undiluted Organic Solvents	PTFE	Nylon
Protein Analysis • Samples with Biomolecules – Buffers	PES	RC or CA
Tissue Culture Media	PES	RC or CA
High Particle-Load Samples – Organic Solvents	Glass Fiber/PTFE	-
High Particle-Load Samples – Aqueous Solutions	Glass Fiber/Nylon	-



Agilent Technologies

# Captiva, Captiva ND, and Captiva ND Lipids Plate and Cartridge Formats



In addition to syringe filters, Captiva filtration products also include 96-well plates and 3 mL cartridge formats, to support a wide variety of filtration needs.

- Plate formats suitable for automation platforms, small sample volumes, and high throughput
- 3 mL cartridges support larger sample volumes, automation (by vacuum/positive pressure manifolds), and lower throughput
- Available in Captiva filtration and in the “ND” – Non-Drip – configurations. Non Drip configuration allows for protein precipitation inside the cartridge/well. No dripping until vacuum/positive pressure applied
- “Functional” filtration with Captiva ND lipids supports lipid removal to provide cleaner samples



# Captiva Average Processing Time:

## Plates vs Cartridges




### Conditions:

- 300  $\mu$ L ACN + 100  $\mu$ L plasma
- Vacuum step only
- ~15 in Hg

### Results:

Average processing time of the Captiva products show good correlation, and supports the benefit of fast flow.





# Captiva ND and ND Lipids - Recoveries




## Conditions:

- Captiva Plates: 100  $\mu$ L plasma
- Captiva Cartridges: 500  $\mu$ L plasma (contains 4X more lipid stripping sorbent than plate)
- Drugs: Atenolol, Propranolol, Loperamide

## Results:

Cartridges showed comparable recovery of three drugs compared to plates. This demonstrates correlation



# Captiva ND Lipids

## Phospholipid-stripping performance


### Conditions:

- Captiva NDL vs off-line protein precipitation
- Sample = plasma
- Tube tested with both 100  $\mu$ L (same sample size as plate) and 500  $\mu$ L (ideal for 3 mL tube)
- Instrument = Agilent 1260 LC module with a 6460 QQQ LC/MS/MS

### Results:

Compared to off-line PPT (no lipid stripping), we can see the effectiveness of the Captiva NDL products, and how completely they remove ion-suppressing phospholipids

Technique	% phospholipid remaining in sample
PPT	100%
NDL plate	1.24%
NDL cartridge 100 $\mu$ L	0.33%
NDL cartridge 500 $\mu$ L	1.04%



# Sample Preparation Time Comparison



## PPT (centrifugation) vs. Captiva ND Lipids



Captiva ND Lipids	Time (min)
Add 0.6 mL of MeOH and 0.2 mL of plasma sample to Captiva ND 96-well plate.	5
Mix each well with a pipette 5 times and apply vacuum for filtration.	
Directly transfer injection plate for analysis.	0
Total time required for sample preparation	<b>5</b>

This time comparison is based on the preparation of 96 samples.

# Captiva ND Lipids Application: Reproducibility



1. Atenolol
2. Pindolol
3. Propranolol

LC/MS chromatograms over 5000 injections show no retention time shift.

# Captiva Accessories



- Captiva Disposable Syringes (5, 10, 20 mL)



- CaptiVac vacuum collar, Captiva deep-well collection plate, and pierceable plate covers



- Vac Elut 20 and 12 port vacuum manifolds



- Vac Elut SPS 24 vacuum manifold

# Captiva plates and Cartridges:

## Advantages



Feature	Advantages	Benefits
Non-Drip Filtration Capabilities	In-cartridge or in-plate protein precipitation	<ul style="list-style-type: none"> <li>Automation</li> <li>Reproducible recoveries</li> <li>No dripping without vacuum or positive pressure</li> </ul>
Dual-depth Filter	Non-linear path with two layers	<ul style="list-style-type: none"> <li>High capacity for particulates</li> <li>More residence time for mechanical filtration and adsorption to occur for cleaner samples</li> </ul>
Non-Drip Lipids Capability	Phospholipid removal	<ul style="list-style-type: none"> <li>Increased cleanup of samples through removal of lipids and phospholipids</li> <li>Longer instrument uptime</li> <li>Minimized ion suppression in LC-MSMS due to phospholipid interference</li> </ul>



# FILTRATION RECOMMENDATIONS & SUMMARY

Filtration

# Filtration and Other Sample Preparation Techniques

<div> <div>More Specific</div> <div>←</div> <div>Instrument Separation and Detection Specificity</div> <div>←</div> <div>Less Specific</div> </div>									
<div> <div>Less Specific</div> <div>→</div> <div>Sample Preparation Specificity</div> <div>→</div> <div>More Specific</div> </div>									
Sample Prep Technique Interference Removed	Dilute & Shoot	Filtration	Liquid/Liquid Extractions	Supported Liquid Extractions (SLE)	Dried Matrix Spotting	Precipitation	QuEChERS	Lipid Removal 'Hybrid' Filtration	Solid Phase Extraction
Lipids	No	No	No	Some	No	No	Yes	Yes	Yes
Oligomeric Surfactants	No	No	No	No	No	No	No	Yes	Yes
Particulates	No	Yes	No	Some	No	Yes	Yes	Yes	Yes
Pigments	No	No	No	Some	No	No	Yes	No	Yes
Polar Organic Acids	No	No	Yes	Yes	No	No	Yes	No	
Proteins	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Salts	No	No	Yes	Yes	No	No	No	No	Yes
Suggested Agilent Product	Agilent Autosampler Vials	Captiva Syringe Filters		Chem Elut	Bond Elut DMS	Captiva ND	Bond Elut QuEChERS	Captiva ND LIPIDS	Bond Elut Silica and Polymeric SPE
Agilent Captiva Filtration Products are recommended for use with any LC or LC-MS method									



# Captiva Filtration Product Recommendation

Agilent Captiva Product Format	Filtration Needs/Situation						
	Low Throughput	Vacuum processing	In-well protein precipitation	Particulate filtration only	Automation Friendly	Lipid removal	Larger sample size
Single cartridge	X	X		X	X		X
Single ND cartridge	X	X	X	can be done	X		X
Single ND Lipids cartridge	X	X	X		X	X	X
96-well plate		X		X	X		
96-well ND plate		X	X	can be done	X		
96-well ND Lipids plate		X	X		X	X	
Syringe filter	X	can be done		X	can be done		X



# Conclusions and Summary

- Filtration is a straightforward means of preparing samples when minimal cleanup is needed
- Mechanical filtration to remove particulates maximizes analytical uptime
- Functionalized filtration can further enhance standard sample preparation methods
- Combining filtration with other sample preparation processes can improve result quality and maximize analytical uptime
- Agilent offers a wide range of filtration products to suit your needs, from syringe filters through specialized filtration products

# Additional Resources and Application Support

## Reference Materials and Guides:

Agilent Captiva Plates and Cartridges Guide:

<http://www.chem.agilent.com/Library/flyers/Public/5990-9061EN.pdf>

Agilent Captiva Syringe Filter Interactive Selection Tool:

[www.agilent.com/chem/SelectFilters](http://www.agilent.com/chem/SelectFilters)

Agilent Captiva Syringe Filter Selection guide:

[www.chem.agilent.com/en-US/products-services/Columns-Sample-Preparation/Sample-Preparation/Filtration/Syringe-Filters/Pages/Selectionguide.aspx](http://www.chem.agilent.com/en-US/products-services/Columns-Sample-Preparation/Sample-Preparation/Filtration/Syringe-Filters/Pages/Selectionguide.aspx)

## Agilent Sample Preparation Products Technical Support Contact Information\*:

Phone: 800-227-9770, Options 3, 3, 3

Email: [spp-support@agilent.com](mailto:spp-support@agilent.com)

\* *North America*



# Acknowledgements

Agilent applications chemists, product managers, and technical support personnel who contributed to this presentation

# Questions?

