

# Installation, Care and Maintenance of Capillary Gas Chromatography Columns

Mark Sinnott  
Chromatography Applications Engineer

# Presenter – Mark Sinnott



**Mark Sinnott works for Agilent Technologies as a Technical Support Engineer in the Consumables and Supplies Division (CSD) at the capillary column manufacturing facility (the “J&W Scientific” location). In his position at Agilent, Mark performs technical support and applications assistance to gas chromatographers worldwide. He has more than 22 years of experience in gas chromatography, including environmental analysis of compounds in air, soil and water matrices, including dissolved gas analysis for the electrical industry. Mark holds a Master’s Degree in Chemistry from California State University, Sacramento, and currently resides in Sacramento.**

# Installation, Care and Maintenance of Capillary Gas Chromatography Columns

**or....**

**"It's not what your column can do for you,  
but what you can do for your column"**

# Column Installation

**"Getting off to a good start"**

# Column Installation Procedure

- **Install the column**
- **Leak and installation check**
- **Column conditioning**
- **Setting linear velocity or flow rate**
- **Bleed profile**
- **Test mix**

# Column Installation


## What type of ferrule should I use?

- **Graphite**
- **Graphite/Vespel**


# Column Installation

## Measuring the right distance

### White out



### Septa



# Cutting The Column

**Gently scribe through the polyimide coating.  
Do not attempt to cut the glass.**

## **Recommended tools:**


**Diamond or carbide tipped pencil; or sapphire  
cleaving tool, ceramic wafer  
Ocular**

## **Do not use:**


**Scissors, file, etc.**



# Example of a Bad Cut




# Examples of Column Cuts




# Column Installation

## How tight is tight?



# Overtightened Ferrule



# Column Installation

## Leak Check

**DO NOT USE SNOOP**

**Electronic leak detector**

**IPA/Water**

**Inject a non-retained peak**

# Leak and Installation Check

Inject a non-retained compound vs DB-1

**Detector**

**Compound**

---

**FID**

**Methane or Butane**

**ECD**

**MeCl<sub>2</sub> (headspace or diluted)**

**NPD**

**CH<sub>3</sub>CN-acetonitrile (headspace or diluted)**

**TCD**


**Air**

**MS**


**Air or Butane**

**The peak should be sharp and symmetrical**

# Non-Retained Peak Shapes



**Good Installation**



**Improper Installation or  
Injector Leak**

- Check for:**
- Too low of a split ratio
  - Injector or septum leak
  - Liner problem:  
(broken, leaking, misplaced)
  - Column position in injector and detector

# Calculating Linear Velocity

**Inject a non-retained compound and obtain the retention time:**

$$\bar{\mu} = \frac{L}{t_0}$$

$\bar{\mu}$  = Average linear velocity (cm/sec)

L = Column length (cm)

$t_0$  = Retention time (sec)

He 20-40 cm/sec

H<sub>2</sub> 35-55 cm/sec

$\bar{\mu}$  is dependent on column temperature



# Calculating Flow Rate

**Inject a non-retained compound and obtain the retention time:**

$$\bar{F} = \frac{\pi r^2 L}{t_0}$$

$\bar{F}$  = Flow rate (mL/min)

$r$  = Column radius (cm)

$L$  = Column length (cm)

$t_0$  = Retention time (sec)

$\bar{F}$  is dependent on column temperature

**Measuring flow with a flow meter is often inaccurate**

# Column Conditioning

**System must be leak free before conditioning column**

**Heat the column to the lower of:**

**Isothermal maximum temperature OR**


**20° to 30°C above highest operation temperature**

**Temperature programming is not necessary**


**Stop conditioning when the stable baseline is obtained:**

**1 to 2 hours in most cases**

# Contamination of system by residue on fingers during column installation



# Contamination from French Fry Grease




## Procedure:

- (1) Held french fry for 5 seconds.
- (2) Fingertip was wiped with paper towel to remove as much of the offending material as possible.
- (3) Lightly touched the part of the column sticking up above the ferrule.
- (4) Installed column into injector.
- (5) Set oven temperature to 40C.
- (6) Started oven temperature program as soon as oven reached 40C.

# Generating a Bleed Profile

## Temperature program the column without an injection\*



**\*DB-1 30m x .32mm I.D., .25 $\mu$ m**

**Temperature program // 40°C, hold 1 min // 20°/min to 320°C, hold 10 min.**

# Test Mixes

Used to determine how "good" the column is



# Column Performance Summary

Catalog: 1225032

Stationary Phase: DB-5

Serial: US5345175H

Description: 30m x 0.252mm x 0.25µm

Temperature Limits: -60°C to 325°C (350°C Pgm)

## Performance Results


| Theoretical Plates/Meter: | <u>Spec.</u> | <u>Meas.</u> |
|---------------------------|--------------|--------------|
| PENTADECANE               | 3900         | 4533         |

|              |  |       |
|--------------|--|-------|
| <b>UTE%:</b> |  |       |
| PENTADECANE  |  | 97.8% |

| <b>Retention Index:</b> |                  |        |
|-------------------------|------------------|--------|
| 1-UNDECANOL             | 1371.0 to 1372.0 | 1371.6 |
| ACENAPHTHYLENE          | 1459.3 to 1460.3 | 1460.0 |

| Compound Identification | Retent. Time | Part. Ratio | 1/2-Width |
|-------------------------|--------------|-------------|-----------|
| 1. 1,6-HEXANEDIOL       | 2.507        | .93         | .017      |
| 2. 4-CHLOROPHENOL       | 2.945        | 1.27        | .020      |
| 3. METHYL NONANOATE     | 3.210        | 1.47        | .020      |
| 4. 4-PROPYLANILINE      | 3.803        | 1.93        | .025      |
| 5. TRIDECANE            | 4.187        | 2.22        | .026      |
| 6. 1-UNDECANOL          | 5.505        | 3.23        | .034      |
| 7. ACENAPHTHYLENE       | 7.992        | 5.15        | .053      |
| 8. PENTADECANE          | 9.557        | 6.35        | .061      |

# Chromatographic Performance





# Test Mixture Components

## Compounds

**Hydrocarbons**

**FAME's, PAH's**

**Alcohols**

**Acids**

**Bases**

## Purpose

**Efficiency**

**Retention**

**Retention**

**Activity**

**Acidic Character**

**Basic Character**

# Own Test Mixture

- **More specific to your application**
- **Selective detectors**
- **Concentrations specific to your application**
- **Use same instrument conditions**
- **Easiest to simply inject a calibration standard**
- **Store for future measure of column performance**

# An Ounce of Prevention.....


# Common Causes of Column Performance Degradation

- **Physical damage to the polyimide coating**
- **Thermal damage**
- **Oxidation (O<sub>2</sub> damage)**
- **Chemical damage by samples**
- **Contamination**

# Physical Damage to The Polyimide Coating


- **The smaller the tubing diameter, the more flexible it is.**
- **Avoid scratches and abrasions**
- **Immediate breakage does not always occur upon physical damage**

# NOT what you want your column to look like!



# Thermal Damage

Degradation of the stationary phase is increased at higher temperatures. Breakage along the polymer backbone.



# Thermal Damage

## What To Do If It Happens

- **Disconnect column from detector**
- **“Bake out” overnight at isothermal limit**
- **Remove 10-15 cm from column end**



# Thermal Damage

- **Rapid degradation of the stationary phase caused by excessively high temperatures**


**Isothermal limit = Indefinite time**

**Programmed limit = 5-10 minutes**

- **Temporary "column failure" below lower temperature limit**

# Oxidation (O2 Damage)

Oxygen in the carrier gas rapidly degrades the stationary phase. The damage is accelerated at higher temperatures. Damage along the polymer backbone is irreversible.



Dimethylpolysiloxane

# Oxygen Damage


## What To Do If It Happens

- **Rapid damage to the column**
- **Usually results in irreversible column damage**


# How to Prevent Column Damage by Oxygen

- **High quality carrier gas (4 nine's or greater)**
- **Leak free injector and carrier lines**
  - Change septa**
  - Maintain gas regulator fittings**
- **Appropriate impurity traps**

# Configurations for Carrier Gas Purifiers



# Configurations for Carrier Gas Purifiers



# Chemical Damage

Bonded and cross-linked columns have excellent chemical resistance except for inorganic acids and bases



etc.

Chemical damage will be evident by excessive bleed, lack of inertness or loss of resolution/retention.

# Chemical Damage

## What To Do If It Happens


- Remove 1/2 - 1 meter from the front of the columns
- Severe cases may require removal of up to 5 meters




# What is Normal Column Bleed

Normal background signal generated by the elution of normal degradation products of the column stationary phase

# Column Bleed is Influenced by:



# Mass Spectrum of Phenylmethylpolysiloxane Column Bleed (Normal Background)



**Mass spectral library search is not always accurate**

# What is a Bleed Problem?

An abnormal elevated baseline at high temperature

**IT IS NOT**

A high baseline at low temperature

Wandering or drifting baseline at any temperature

Discrete peaks

# Column Contamination

- **Fouling of GC and column by contaminants**
- **Mimics nearly every chromatographic problems**

# Symptoms of Contamination

- **Poor peak shape**
- **Loss of separation (resolution)**
- **Changes in retention**
- **Reduced peak size**
- **Baseline disturbances (semi-volatiles only)**

# Typical Samples That Contain a Large Amount of Residues

**Biological (Blood, Urine, Tissue, Plants)**

**Soils**

**Foods**

**Waste Water**

**Sludges**

***All samples contain residues!! (even standards!)***

# Other Sources of Contamination

- **Septum and ferrule particles**
- **Gas and trap impurities**
- **Unknown sources (vials, syringes, etc.)**



# Non-Volatile Residues

Any portion of the sample that does not elute from the column or remains in the injector.


# Semi-Volatile Residues

Any portion of the sample that elutes from the column after the current chromatographic run.

# Methods to Minimize Non-Volatile Residue Problems

- **Sample cleanup**
- **Packed injection port liners**
- **Guard columns**

# Guard Column or Retention Gap



**The guard column is 3 - 5 meters of deactivated fused silica tubing with the same diameter as the analytical column. It is connected with a zero dead volume union.**

# Non-Volatile Contamination

## What To Do If It Happens

- Do not “bake out” the column
- Front End Maintenance
  - clean or change the injector liner
  - clean the injector
  - cut off 1/2 -1 meter of the front of the column
- Turn the column around
- Solvent rinse the column
- Cut the column in half

# Semi-Volatile Contamination

## What To Do If It Happens

- **“bake out” the column**
  - **Limit to 1-2 hours**
  - **Longer times may polymerize some contamination and reduces column life**
- **Solvent rinse the column**


# Methods to Minimize Non-Volatile Residue Problems

- **Sample cleanup**
- **Packed injection port liners**
- **Guard columns**


# Why Perform Sample Preparation?

- To acquire desired sensitivity/selectivity
- To reduce contamination of inlet supplies and columns
- To protect your investment if using very sensitive and expensive instrumentation

Pesticides in Avocado without SP



Pesticides in Avocado with SP



# Column Storage

- **Place septa over the ends**
- **Return to column box**



## Always Remember to:

- **Start with a good installation**
- **Maintain an oxygen free system**
- **Avoid physical, thermal, and chemical damage**
- **Take steps to prevent contamination**
- **Use Sample Prep where possible**

# Agilent/J&W Technical Support

800-227-9770 (phone: US & Canada)\*

**302-993-5304 (phone)\***

*\* Select option 41*

**866-422-5571 (fax)**

email: [gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[www.agilent.com/chem](http://www.agilent.com/chem)

