

# Low Bleed Stationary Phases for Gas Chromatography



**Agilent Technologies**

Dial 1-816-650-0621 for e-Seminar Audio

Slide 4

# Advantages of Low Bleed Phase

**Lower bleed compared to standard phases**

**Better signal to noise (sensitivity)**

**Higher upper temperature - shorter run times**

**Lower detector maintenance**

**Improved spectral purity**

**Potential for longer column life**



# Stationary Phase Bleed

**Normal degradation of the stationary phase**

**All columns exhibit some bleed**



# WHAT IS A BLEED PROBLEM?

**An abnormal elevated baseline at high temperature**

**IT IS NOT**


**A high baseline at low temperature**

**Wandering or drifting baseline at any temperature**

**Discrete peaks**



# Stationary Phase Bleed Mechanism



# “ms” Phase and Column Criteria

**Higher upper temperature limit**

**Lower bleed than columns with “equivalent” phase**





**Superior inertness to acidic and basic compounds**

**Analogous phase selectivity**




# Capillary GC Phase Development

## Early Phases



# Phenylene Polymer Chemistry




- **Better thermal stability**
- **Lower bleed**






# DB-5MS Structure



**DB-5**  
5% Phenyl



**DB-5ms**

1. Increased stability
2. Different selectivity
3. Optimized to match DB-5



# Types of Low Bleed Columns

**Phases tailored to “mimic” currently existing polymers**

**New phases unrelated to any previously existing polymers**

**Optimized polymerization chemistry and manufacturing processes**



# Low Bleed Success Stories

## “5% Phenyl”

|         |              |
|---------|--------------|
| DB-5ms  | BPX5         |
| HP-5TA  | CP-Sil8CB/MS |
| HP-5ms  | 007-5ms      |
| Rtx-5ms | MDN-5ms      |

## “35% Phenyl”

## “50% Phenyl”

## “Other” Low Bleed

DB-XLB

HT5,HT8

## “50% Cyanopropyl-phenyl”



# Low Bleed Family

**1ms - siloxane 360 MAOT**

**5ms - arylene**

**Next generation arylene:**

**XLB, 17ms, 35ms, 225ms**

**ion trap tested:**

**DB-5ms ITD**

**DB-XLB ITD**



## Phase Selectivity

**Many “backbone stiffening” substitutes give analogous “phenyl like” characteristics**


**But**

**not so for 100% Dimethylpolysiloxane**



# Break Number 1

For Questions and Answers  
Press \*1 on Your Phone to  
Ask a Question



## Early Attempts at “100% DMPS” Low Bleed Column

| Column     | Delta Bleed (pA) |
|------------|------------------|
| DB-1       | 5.0              |
| Cyclohexyl | 5.0              |
| Butyl      | 20               |
| Ethyl      | 40               |
| Butenyl    | 40               |

All columns: 30 m x 0.25 mm I.D., 0.25  $\mu$ m  
Conditions: 135°C to 320°C



# Retention Indices for Low Density Alkane Substitutions Undecanol and Acenaphthylene

| Column     | Undecanol | Acenaphthylene |
|------------|-----------|----------------|
| DB-1       | 1355.10   | 1424.62        |
| Cyclohexyl | 1344.55   | 1436.77        |
| Butane     | 1349.85   | 1427.19        |
| Butenyl    | 1356.34   | 1432.80        |
| Ethyl      | 1356.41   | 1433.30        |

All columns: 30 m x 0.25 mm I.D., 0.25  $\mu\text{m}$

Conditions: H<sub>2</sub> carrier, Split 100:1, Isothermal 135°C





## Retention Indices for “Low Arylene” Substitutes

| Column | 3,5-DMP | 1-Nitrohexane | 1,4-DIIPB | 1-Nonanol | 2-Decanone |
|--------|---------|---------------|-----------|-----------|------------|
| DB-1   | 952.4   | 1009.4        | 1152.1    | 1156.9    | 1169.3     |
| Mfg B  | 966.1   | 1015.1        | 1157.3    | 1161.1    | 1174.1     |
| Mfg C  | 958.7   | 1015.0        | 1157.0    | 1161.0    | 1173.0     |



# **% “Phenyl” Substitution Based on Average RI Values**

---

**Mfg B**

**1.4% phenyl equivalent**

**Mfg C**

**1.0% phenyl equivalent**



# Optimizing an Already Robust Phase


**Changes to fused silica surface prior to polymerization**

**Changes to polymerization chemistry**


**Changes to chemical manifolds used to process columns**



# Column Bleed Histogram for 100% DMPS Columns




# Molecular Weight Distribution and Bleed in DMPS Columns



**DMPS chains of inconsistent lengths**




# Molecular Weight Distribution and Bleed in DMPS Columns



DMPS chains of **consistent** lengths




# Relative Improvements in Bleed for Standard versus Low Bleed DMPS



**Columns: 30 m x 0.25 mm I.D., 0.25  $\mu$ m film**



# Mass Spectral Comparison of DB-1 and DB-1ms



Columns: 30 m x 0.25 mm I.D., 0.25  $\mu$ m film, Oven Temperature 325°C



Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 27



# Retention Indices Comparison

Retention indices at 90°C

| Compound               | DB-,1  | DB-1ms |
|------------------------|--------|--------|
| 3,5-Dimethylpyridine   | 952.4  | 954.0  |
| 1-Nitrohexane          | 1009.4 | 1009.7 |
| 1,4-Diisopropylbenzene | 1152.1 | 1152.7 |
| 1-Nonanol              | 1156.9 | 1157.2 |
| 2-Decanone             | 1169.3 | 1169.7 |

Columns: 30 m x 0.25 mm I.D., 0.25 µm film



# Comparison of Selectivity

**Columns:** 30 m x 0.25 mm I.D., 0.25  $\mu$ m

**Carrier:** Helium at 42 cm/sec,  
measured at 40°C


**Oven:** 40°C for 5 min  
40-200°C at 10°/min

**Injector:** 1  $\mu$ L injection  
Split 100:1, 250°C

**Detector:** MSD

## Compound List

1. Benzyl alcohol
2. 3-Nitrotoluene
3. 4-Chlorophenol
4. Hexachlorobutadiene
5. Dicyclohexylamine
6. Acenaphthene
7. Pentadecane
8. Phenanthrene



# Inertness Testing on DB-1ms

Column: DB-1ms  
30 m x 0.25 mm I.D., 0.25  $\mu$ m

J&W P/N: 122-0132

Carrier: Hydrogen at 39.8 cm/sec,  
measured at 125°C


Oven: 125°C for 8.5 min  
125-265°C at 20°/min

Injector: Split 1:50, 250°C

Detector: FID, 300°C

1. 2-Ethylhexanoic acid\*
2. 1,6-Hexanediol\*
3. 4-Chlorophenol\*
4. 1-Methylnaphthalene
5. Tridecane
6. 1-Undecanol
7. Tetradecane
8. Dicyclohexylamine\*
9. 2,4-Dinitrophenol\*
10. Acenaphthene
11. N-Nitrosodiphenylamine\*
12. Pentachlorophenol\*
13. Phenanthrene
14. Carbazole

\* Active Analytes



# GS-CarbonPLOT

**Proprietary carbon-layer phase**

**Very high retention**

**High upper temperature limit (360°C)**

**“Bonded” phase -- cannot generate particles**

**Good for CO<sub>2</sub>, light hydrocarbons, etc.**

**Not affected by water in samples**

**Low bleed (for a PLOT column)**



# GS-GasPro

**Proprietary “silica-based” phase**

**Very inert column for reactive compounds**

**“Bonded” phase; cannot generate particles**

**Excellent for MSD and valve systems**

**Low-bleed (for a PLOT column)**

**Not affected by water in samples**

**Currently only available in 0.32 mm I.D.**



# GS-GasPro

**Applications include:**

**Hydrocarbons (C<sub>1</sub> to C<sub>10</sub>)**

**CO<sub>2</sub>**

**Sulfur gases**

**Freons**


**Hydrides**

**Permanent gases (-80°C oven)**



## Break Number 2

For Questions and Answers  
Press \*1 on Your Phone to  
Ask a Question



# Applications





# Controlled Substances & Related Samples

**GC/MS is the preferred means of identification**

**Ideal for trace level analysis -  
greater signal-to-noise ratios  
and lower detection limits**



# Drugs

**Column:** DB-1ms  
**30 m x 0.25 mm I.D., 0.25 µm**


**J&W P/N:** 122-0132  
**Carrier:** Helium at 40 cm/sec,  
measured at 50°C

**Oven:** 50°C for 1 min  
50-125°C at 25°/min  
125-325°C at 10°/min  
325° for 5 min

**Injector:** Cold Splitless using Optic II injector  
50-250°C at 10°/sec,  
0.75 min purge activation time

**Sample:** 1 µL injection of 50-150 ppm standard  
**Detector:** FID, 300°C

- |                                 |                  |
|---------------------------------|------------------|
| 1. Nicotine                     | 10. Cocaine      |
| 2. Caffeine                     | 11. Desipramine  |
| 3. Glutethimide                 | 12. Carbazepine  |
| 4. Lidocaine                    | 13. Trimipramine |
| 5. PCP                          | 14. Heroin       |
| 6. Phenobarbital                | 15. Fentanyl     |
| 7. Methadone primary metabolite | 16. Ibogaine     |
| 8. Methaqualone                 | 17. Triazolam    |
| 9. Methadone                    | 18. LSD          |



C759

**Agilent Technologies**


**Dial 1-816-650-0621 for e-Seminar Audio**

**Slide 37**

# Anabolic Steroids

**Column:** DB-1ms  
**30 m x 0.25 mm I.D., 0.25  $\mu$ m**  
**J&W P/N:** 122-0132  
**Carrier:** Helium at 40 cm/sec,  
measured at 180°C  
**Oven:** 180-320°C at 10°/min  
320° for 4 min  
**Injector:** Split 1:40,  
**Sample:** 2  $\mu$ L of 0.125  $\mu$ g/ $\mu$ L each in Methanol  
**Detector:** FID, Nitrogen makeup gas at 30 mL/min

1. Dehydroisoandrosterone (Prasterone)
2. 5 $\alpha$ -Androstan-17 $\alpha$ -ol-3-one (Stanolone)
3. 19-Nortestosterone (Nandrolone)
4. Mesterolone
5. 17 $\alpha$ -epi-testosterone
6. Testosterone
7. 1-Dehydrotestosterone (Boldenone)
8. 17 $\alpha$ -Methyltestosterone
9. 1-Dehytro-17 $\alpha$ -methyltestosterone (Methandrostenolone)
10. Norethandrolone
11. 1-Dehydrotestosterone acetate
12. Oxymetholone
13. 19-Nortestosterone-17-propionate
14. 4-Chlortestosterone-17-acetate (Clostebol)
15. 1-Dehydrotestosterone benzoate
16. 19-Nortestosterone-17-decanoate
17. 1-Dehydrotestosterone undecyclenate



**Increased upper temperature limit extends analyses to run less volatile compounds, and reduces run times through higher final temperature**




# DB-225ms

## Lower Bleed - Higher MAOT

### FAMES

**Column:** DB-225ms  
30 m x 0.25 mm I.D., 0.25  $\mu$ m  
**J&W P/N:** 122-2932  
**Carrier:** Hydrogen at 40 cm/sec  
**Oven:** 200°C for 1 min  
200-260°C at 3°/min  
**Injector:** Split 1:50, 250°C  
**Detector:** FID  
Nitrogen make-up gas at 30 mL/min



The higher isothermal upper temperature limit of DB-225ms (260°C vs. 220°C for DB-225) allows the elution of higher molecular weight FAMES (above 24:0) while maintaining a reasonable run time.



# AZO Dyes

**Column:** DB-1ms  
30 m x 0.25 mm I.D., 0.25  $\mu$ m

**J&W P/N:** 122-0132


**Carrier:** Helium at 70 cm/sec,

**Oven:** 80°C for 2 min  
80-360°C at 20°/min  
360°C for 3 min

**Injector:** Cold Splitless  
(Optic II injector)

**Injector program:** 50-375°C at 10°/sec  
1.5 min purge  
activation time

**Detector:** FID, 350°C



# Greater inertness allows analysis of active and difficult to analyze compounds



# CLP Pesticides

**Column:** DB-1ms  
30 m x 0.25 mm I.D., 0.25  $\mu$ m

**J&W P/N:** 122-0132

**Carrier:** Helium at 45 cm/sec,  
constant flow mode


**Oven:** 110°C for 0.5 min  
110-250°C at 10°/min  
250-300°C at 25°/min  
300°C for 2 min

**Injector:** Splitless, 250°C  
0.5 min purge activation time

**Sample:** 1  $\mu$ L of a 1 ppm standard in Acetone

**Detector:** MSD, 280°C transfer line  
Full scan at m/z 50-500

- |                                      |                             |
|--------------------------------------|-----------------------------|
| 1. 2,4,5,6-Tetrachloro-m-xylene (IS) | 12. 4,4'-DDE                |
| 2. $\alpha$ -BHC                     | 13. Dieldrin                |
| 3. $\beta$ -BHC                      | 14. Endrin                  |
| 4. $\gamma$ -BHC                     | 15. Endosulfan II           |
| 5. $\delta$ -BHC                     | 16. 4,4'-DDD                |
| 6. Heptachlor                        | 17. Endrin aldehyde         |
| 7. Aldrin                            | 18. Endosulfan aldehyde     |
| 8. Heptachlor epoxide                | 19. 4,4'-DDT                |
| 9. $\gamma$ -Chlordane               | 20. Endrin ketone           |
| 10. Endosulfan I                     | 21. Methoxychlor            |
| 11. $\alpha$ -Chlordane              | 22. Decachlorobiphenyl (IS) |



C758



Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 43



## Organochlorine Pesticides (EPA Method 8081A)

**Column: DB-35ms**  
30 m x 0.25 mm I.D., 0.25 µm

**J&W P/N: 122-3832**

**Carrier:** Helium at 35 cm/sec, measured at 50°C

**Oven:** 50°C for 1 min  
50-100°C at 25°/min  
100-300°C at 5°/min

**Injector:** Splitless, 250°C  
30 sec purge activation time  
1 µL of 35 µg/mL composite 8081A standards, Accustandard Inc.

**Detector:** MSD, 300°C transfer line  
full scan at m/z 50-500

1. 1,2-Dibromo-3-chloropropane
2. 4-Chloro-3-nitrobenzotrifluoride (SS)
3. Hexachloropentadiene
4. 1-Bromo-2-nitrobenzene (IS)
5. Terrazole
6. Chloroneb
7. Trifluralin
8. 2-Bromobiphenyl (SS)

**Method 8081A Standards are found on page 144-145**


9. Tetrachloro-m-xylene (SS)
10. α, α-Dibromo-m-xylene
11. Propachlor
12. Di-allate A
13. Di-allate B
14. Hexachlorobenzene
15. α-BHC
16. Pentachloronitrobenzene (IS)

17. γ-BHC
18. β-BHC
19. Heptachlor
20. Alachlor
21. δ-BHC
22. Chlorothalonil
23. Aldrin
24. Dacthal™
25. Isodrin
26. Kelthane
27. Heptachlor epoxide
28. γ-Chlordane
29. trans-Nonachlor
30. α-Chlordane
31. Endosulfan I
32. Captan
33. p,p'-DDE
34. Dieldrin
35. Chlorobenzilate™
36. Perthane™
37. Chloropropylate

38. Endrin
39. p,p'-DDD
40. Endosulfan II
41. p,p'-DDT
42. Endrin aldehyde
43. Endosulfan sulfate
44. Dibutylchlorendate (SS)
45. Captafol
46. Methoxychlor
47. Endrin ketone
48. Mirex™
49. cis-Permethrin
50. trans-Permethrin

\* Breakdown Products  
SS - Surrogate Standard  
IS - Internal Standard

*Standards used were a composite of individual solutions supplied courtesy of Accustandard Inc., 25 Science Park, New Haven, CT 06511, 800-442-5290.*



## CLP Pesticides

### Column: DB-17ms

30 m x 0.25 mm I.D., 0.25 µm

### J&W P/N: 122-4732

Carrier: Helium at 45 cm/sec,  
measured at 110°C

Oven: 110°C for 0.5 min  
110-150°C at 25°/min  
150-260°C at 12°/min  
260-320°C at 15°/min  
320°C for 2 min


Injector: Splitless, 250°C

30 sec purge  
activation time  
50 pg composite 8081  
standard, and CLP-034  
surrogate mix,  
Accustandard Inc.  
Detector: ECD, 350°C  
Nitrogen makeup  
gas at 40 mL/min

1. Tetrachloro-m-xylene (SS)
2. α-BHC
3. γ-BHC
4. β-BHC
5. Heptachlor
6. δ-BHC
7. Aldrin
8. Heptachlor epoxide
9. γ-Chlordane
10. α-Chlordane
11. Endosulfan I
12. p, p'-DDE
13. Dieldrin
14. Endrin
15. p, p'-DDD
16. Endosulfan II
17. p, p'-DDT
18. Endrin aldehyde
19. Endosulfan sulfate
20. Dibutylchlorodate (SS)
21. Methoxychlor
22. Endrin ketone

SS - Surrogate Standard

*Standards used were a composite of individual solutions supplied courtesy of Accustandard Inc., 25 Science Park, New Haven, CT 06511, 800-442-5290.*



C510



Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 45

# PAHs

**Column:** DB-1ms  
**30 m x 0.25 mm I.D., 0.25 µm**

**J&W P/N:** 122-0132

**Carrier:** Helium at 40 cm/sec,  
constant flow mode


**Oven:** 80°C for 1 min  
80-325°C at 15°/min  
325°C at 2°/min

**Injector:** Splitless, 300°C  
0.8 min purge activation time

**Sample:** 1 µL of a 5 ppm standard in Toluene

**Detector:** MSD, 325°C transfer line  
Full scan at m/z 50-450

1. Naphthalene
2. Acenaphthene
3. Acenaphthylene
4. Fluorene
5. Phenanthrene
6. Anthracene
7. Fluoranthene
8. Pyrene
9. Benzo(a)anthracene
10. Chrysene
11. Benzo(a)anthracene
12. Benzo(k)fluoranthene
13. Benzo(a)pyrene
14. Dibenz(a,h)anthracene
15. Benzo(g,h,i)perylene



C757



**Agilent Technologies**

**Dial 1-816-650-0621 for e-Seminar Audio**

**Slide 46**

# EPA Method 8082 PCB Congeners

**Column:** DB-1ms  
30 m x 0.25 mm I.D., 0.25 µm

**J&W P/N:** 122-0132  
**Carrier:** Helium at 46 cm/sec,  
measured at 100°C,  
constant pressure mode


**Oven:** 100°C for 0.8 min  
100-325°C at 15°/min  
325°C at 2°/min

**Injector:** Splitless, 275°C  
0.8 min purge activation time

**Sample:** 1 µL of a 20 ppm mixture  
**Detector:** MSD, 280°C transfer line  
Full scan at m/z 50-550

(PCB congeners are listed by their IUPAC numbers)

- |                                 |                                |
|---------------------------------|--------------------------------|
| 1. IUPAC 1                      | 11. IUPAC 110                  |
| 2. Tetrachloro-m-xylene (IS/SS) | 12. IUPAC 151                  |
| 3. IUPAC 5                      | 13. IUPAC 153                  |
| 4. IUPAC 18                     | 14. IUPAC 141                  |
| 5. IUPAC 31                     | 15. IUPAC 137                  |
| 6. IUPAC 52                     | 16. IUPAC 187                  |
| 7. IUPAC 44                     | 17. IUPAC 183                  |
| 8. IUPAC 66                     | 18. IUPAC 180                  |
| 9. IUPAC 101                    | 19. IUPAC 170                  |
| 10. IUPAC 87                    | 20. IUPAC 206                  |
|                                 | 21. Decachlorobiphenyl (IS/SS) |



C761



Agilent Technologies



Dial 1-816-650-0621 for e-Seminar Audio

Slide 47

**“MS”-grade PLOT columns allow analysis of permanent gases and other highly volatile compounds**




Column: GS-CarbonPLOT  
30 m x 0.32 mm I.D., 1.5 µm  
J&W P/N: 113-3132  
Carrier: Helium at 30 cm/sec (measured at 125°C)  
Oven: 125°C for 2 min  
125-350°C at 20°/min  
Injector: Split 1:20, 250°C  
Detector: MS, SIM mode  
280°C transfer line



## Sulfur Gases in Light Hydrocarbon Stream

**Column:** GS-GasPro  
30 m x 0.32 mm I.D.  
**J&W P/N:** 113-4332  
**Carrier:** Helium at 53 cm/sec  
**Oven:** 25°C for 4 min  
25 - 200°C at 10°/min  
200°C Hold  
**Injector:** 200 µL  
**Detector:** TCD

1. Nitrogen
2. Methane
3. Ethane
4. CO<sub>2</sub>
5. SF<sub>6</sub>
6. Ethylene
7. COS
8. Acetylene
9. Propane
10. H<sub>2</sub>S
11. Ethylene oxide/cyclopropane
12. Propylene
13. Propadiene
14. Isobutane
15. Butane
16. SO<sub>2</sub>
17. Methyl acetylene



C579



Agilent Technologies

Dial 1-816-650-0621 for e-Seminar Audio

Slide 50

Column: **GS-GasPro**


**5 m x 0.32 mm I.D.**

Injector: Gas sample valve, 1.0 mL fixed sample loop

Detector: Sievers SCD, reaction tube 800°C

Carrier: Helium, 39 cm/sec at 0°C (determined by  $T_R$  of COS)

Oven: 0°C for 1 min, 30° / min to 160°C, hold at 160°C for 4 min,  
50° / min to 260°C, hold at 260°C for 5 min




**SO<sub>2</sub> = 100 ppbV**  
**others = 89 ppbV**





**Column:** GS-GasPro  
60 m x 0.32 mm I.D.  
**J&W P/N:** 113-4362  
**Carrier:** Helium at 40 cm/sec (3.3 mL/min)  
measured at 80°C  
**Oven:** 80°C for 0.5 min  
80-175°C at 25°/min  
175°C for 2 min  
175-250°C at 25°/min  
**Injector:** Split 1:17, 250°C  
500 µL injection of 100 ppmV  
SUMMA cannister mixture  
**Detector:** FID, 275°C  
Nitrogen makeup gas at 32 mL/min

**Compound list on page 5**



## Summary

**Low bleed columns are available in a wide range of selectivities, including DB-1ms that provides virtually identical retention to 100% DMPS, and PLOT columns with remarkably low bleed.**



## Wrap-up E-Seminar Questions

**Thank you for attending Agilent e-Seminars.  
Our e-Seminar schedule is expanding every  
week. Please check our website frequently at:  
[www.agilent.com](http://www.agilent.com)**

