



Simplified Method Development

Low, Mid and High pH Recommendations



Agilent Technologies

pH – An Important Parameter for Method Development

Retention of ionizable compounds is strongly affected by pH


Ionizable compounds (acids and bases) may be analytes or matrix compounds

Accurate pH control improves method reproducibility

The pH range from 1 – 12 provides maximum method development flexibility



Change in Retention with pH for Ionizable Compounds is Compound Dependent



Change in Retention with pH for Ionizable Compounds is Key to Method Development

Non-charged analytes have better retention (i.e. acids at low pH and bases at high pH)

Silanols on silica ionize at mid-pH, increasing retention of basic analytes (i.e possible ion-exchange interactions)

Choose mobile phase pH and column type to optimize retention and selectivity during method development



Recommended Method Development Goals

Adequate resolution of all peaks, $R_s \geq 2.0$

Retention of first peak preferred to be at least $k=1$

Analysis time below 30 minutes, 20 minutes preferred

Robust and rugged methods

Use buffered mobile phases and try low pH first



Why Develop Methods at Low pH?

Acids are protonated for maximum retention

Silica silanols are protonated thereby minimizing ion-exchange interactions with basic compounds

Good peak shape


Long term reproducibility

Excellent mobile phase choices



Method Development Scheme

Start at Low pH



Recommended Starting Conditions for RP-HPLC Method Development Approach

Separation Variable

Preferred Initial Choice

Column

Primary Stationary Phase	Eclipse XDB-C18 or Eclipse XDB-C8
Secondary Stationary Phase	SB-C18 or SB-C8
Dimensions	4.6 x 75 mm or 4.6 x 150 mm
Particle Size	3.5 μm 5 μm
Pore Size	80Å: M.W. \leq 4000, 300Å: M.W. \geq 4000

Mobile Phase

Solvents A-B	Water-acetonitrile
% B solvent	Variable
Buffer	25 mM NaH_2PO_4 , pH \leq 3 or 0.1% TFA or Formic acid
Additives i.e. amines and ion-pair reagents	TEA, Hexane sulfonate as needed
Flow Rate	1-2 mL/min
Temperature	30 - 35°C



Choose Eclipse XDB First at Low pH

Wide useable pH range – pH 2- 9 – can use for low and mid pH without changing columns


Fully-hydroxylated ultra-pure silica improves peak shape

Double endcapping for good peak shape of basic compounds at low and mid pH

Three different bonded-phases for selectivity options



Method Development – Eclipse XDB-C18 at Low pH - Separation of Flavones



Conditions:

Column: 4.6x75mm, 3.5 μ m


Temp: Ambient

Eluent: 30% ACN
70% 25mM NaH₂PO₄
pH 2.5


Flow: 0.75 mL/min

Inj. Vol.: 4 μ L

Sample: 1) Luteolin
2) Apigenin



Excellent Reproducibility of ZORBAX Eclipse XDB at Low pH



3 Lots % RSD	
k'	α
1) 3.1	-
2) 2.6	0.56
3) 2.8	0.29
4) 2.5	0.34

Conditions:

Column: Eclipse XDB-C8


4.6x150mm, 5 μ m

Temp: 35°C

Eluent: 15% ACN

85% 25mM Phosphate

Flow: 1.0 mL/min



3 Lots % RSD	
k'	α
1) 5.4	-
2) 5.1	0.41
3) 4.6	0.59
4) 4.8	0.00

Sample: 1) Cefotaxime

2) Cefoxitin

3) Cefamandole

4) Cephalothin



StableBond –Secondary Choice at Low pH

Superior column lifetime at very low pH – down to pH 1 – due to patented sterically protecting bonding technology

Fully-hydroxylated ultra-pure silica improves peak shape


Six different 80Å bonded-phases – SB-C18, SB-C8, SB-CN, SB-Phenyl, SB-C3, SB-Aq – provide optimum selectivity with exceptional lifetime

Four different 300Å bonded-phases for selectivity options with protein and peptide separations



Selectivity Options at Low pH

StableBond vs. Eclipse



Conditions:

Column: 4.6x75mm, 3.5 μ m

Temp: Ambient

Eluent: 25% ACN
75% 25mM NaH₂PO₄
pH 2.5

Flow: 0.75 mL/min

Inj. Vol.: 20 μ L

Sample: 1) Caffeic Acid

2) Luteolin


3) Naringenin

4) Apigenin



Method Development – SB-C18 at Low pH Separation of Plant Extract

Flavones, Flavanones, and Phenolic Esters




- To obtain $k=1$ for caffeic acid requires 22 minute analysis time

Method Development – Change Organic Modifier

Separation of Plant Extract

Flavones, Flavanones, and Phenolic Esters



- Methanol as the organic modifier changes selectivity and increases the analysis time.




Method Development – Change Bonded-Phase

Separation of Plant Extract on SB-CN


Flavones, Flavanones, Phenolic Esters

Column: ZORBAX Rapid Resolution SB-CN, 4.6 x 75 mm, 3.5 μ m Mobile Phase: ACN: NaH₂PO₄, pH 2.5 Flow Rate: 1.0 mL/min
Temperature: RT Detection: UV 254 nm Sample: 1. Caffeic acid 2. Impurity 3. Luteolin 4. Naringenin 5. Apigenin

22% ACN: 78% Buffer



25% ACN: 75% Buffer



SB-CN with stronger mobile phase reduces analysis time by 50% and maintains retention of $k=1$ on 1st peak.



StableBond 300SB Columns Ideal for Separations of High MW Analytes

300Å pore size necessary for good peak shape and high efficiency separation of proteins and polypeptides.


Exceptional stability with low pH “TFA containing” mobile phases.

Improve peak shape for lower molecular weight analytes with large hydrodynamic volume.

Four different bonded-phases allow bonded-phase optimization of all high MW samples.



Comparison of Bonded-Phase Options – Affect on Selectivity and Retention of Polypeptides



Columns: ZORBAX 300SB, 4.6 x 150 mm, 5 μ m

Mobile Phase: Gradient, 0 - 26% B in 30 min.

A = 0.1% TFA in Water

B = 0.1% TFA in Acetonitrile

Temperature: 40°C

Sample: 2 μ g of each peptide

Flow Rate: 1.0 mL/min

Detection: UV 210nm



Why Develop RP-HPLC Methods at Mid-pH?

Compounds of interest are unstable at low pH

Improved solubility of analytes at mid pH


Increase retention of basic compounds

May have better selectivity in the pH range 3 – 8



Method Development Scheme

Mid pH Range



Recommended Conditions for Mid-pH Method Development

Separation Variable

Preferred Initial Choice

Column

Stationary Phase

Eclipse XDB-C18

Dimensions

4.6 x 75 mm or 4.6 x 150 mm

Particle Size

3.5 μm

5 μm

Mobile Phase

Solvents A-B

Water-acetonitrile

% B solvent

Variable

Buffer

25 mM Na_2HPO_4 , pH = 7

Acetate/acetic acid

Additives i.e. amines and ion-pair reagents

TEA, tetrabutylammonium as needed

Flow Rate

1-2 mL/min

Temperature

30 - 35°C



Choose Eclipse XDB for Mid-pH

Long lifetime at mid-pH with dense bonding and double endcapping

Strong silica for long lifetime

Double endcapping provides excellent peak shape

Three different bonded-phases (C18, C8, Phenyl) for selectivity optimization




Low pH Provides Insufficient Selectivity for Some Samples


Eclipse XDB-C18, 4.6 x 75 mm, 3.5 μ m

SB-C18, 4.6 x 75 mm, 3.5 μ m

pH 3



pH 3



Mobile Phase: 20% Methanol: 80% 20 mM phosphate buffer Flow Rate: 1.0 mL/min Temperature: RT
Detection: UV 254 nm Sample: 1. Nizatidine 2. Famotidine 3. Cimetidine 4. Pirenzipine



- Poor Selectivity with Eclipse at Low pH
- Selectivity Improves Somewhat with StableBond



Agilent Technologies

Improved Selectivity and Resolution with Eclipse XDB-C18 at Mid-pH

Eclipse XDB-C18, 4.6 x 75 mm, 3.5 μ m





Mobile Phase: 20% Methanol: 80% 20 mM phosphate buffer Flow Rate: 1.0 mL/min Temperature: RT
Detection: UV 254 nm Sample: 1. Nizatidine 2. Famotidine 3. Cimetidine 4. Pirenzpine


- Better selectivity and improved retention occur at pH 7
- A little TEA also improves peak shape peak #4 (di-amine)

Optimize Separations with Eclipse XDB Selectivity Options – Analysis of Sunscreens


Columns: 4.6 x 150 mm, 3.5 μm Mobile Phase: 15% H_2O : 85% MeOH Flow Rate: 1.0 mL/min Temperature: 30°C
Detection: UV 310 nm Sample: 1. Oxybenzone 2. Padimate-O 3. Ethylhexylsalicylate



1. oxybenzone



2. padimate O



3. ethylhexyl salicylate



Bonus-RP

Provides Alternate Selectivity at Mid-pH

Polar alkyl-amide bonded-phase for unique selectivity

Improves peak shape of basic compounds


Triple-endcapped for good lifetime at mid-pH

Enhanced low-pH stability (sterically protecting bonding) for alternate selectivity at low pH

Compatible with 100% aqueous mobile phases



Bonus-RP Provides Alternate Selectivity at Mid-pH



Mobile Phase: 75% mM NaCitrate, pH 6
25% MeOH
Flow Rate: 1.0 mL/min
Temperature: Ambient
Detection: UV 254 nm
Injection Vol: 3 μ L
Sample: Cephalosporins
1. Cephalexin
2. Cephaclor
3. Cephuroxime
4. Cephoxitin

Why Develop RP-HPLC Methods at High pH?

Compounds of interest not soluble at lower pH


Compounds of interest not stable at lower pH

Increase retention of basic compounds by analyzing them in non-charged form

Improve selectivity



Method Development at High pH



Recommended Conditions for High pH Method Development

Separation Variable Column

Preferred Initial Choice

Stationary Phase	Extend-C18
Dimensions	4.6 x 75 mm or 4.6 x 150 mm
Particle Size	3.5 μm 5 μm
Pore Size	80Å: M.W. \leq 4000, 300Å: M.W. \geq 4000

Mobile Phase

Solvents A-B	Water-methanol
% B solvent	Variable
Buffer	20 mM TEA pH =11 ammonium hydroxide, pH 10
Flow Rate	1-2 mL/min
Temperature	RT - 30°C



Silica-Based HPLC Columns are Now a High pH Choice

New technologies to protect silica from dissolution provide good lifetimes at high pH

Superior efficiency of silica-based columns provides high resolution

Robust methods can be established using the same parameters as used at low pH



Choose Extend-C18 for High pH


Patented bidentate C18-C18 bonding for superior high pH stability – up to pH 11.5

Improved performance over polymeric columns

Excellent peak shape with double endcapping

LC/MS at high pH (ammonium hydroxide) with high efficiency


Extend-C18
Bidentate Structure




High pH Increases Retention of Antihistamines

1. Maleate 2. Scopolamine 3. Pseudoephedrine 4. Doxylamine 5. Chlorpheniramine 6. Triprolidine 7. Diphenhydramine

pH 7
30% 20 mM Na_2HPO_4
70% MeOH



pH 11
30% 20 mM TEA
70% MeOH



The retention of this sample of basic compounds increases at high pH.




Extend-C18 Provides High Efficiency and Good Peak Shape

Mobile Phase: 65% 20 mM TEA, pH 11: 35% MeOH Temperature: RT Detection: UV 254 nm
Sample: 1. Pyridoxine 2. Pyridine 3. n-Methylbenzylamine 4. Procainamide 5. n-Acetylprocainamide


Polymeric-Based Column

4.0 x 250 mm, 5 μ m
Flow Rate: 0.5 mL/min



Extend-C18

4.6 x 250 mm, 5 μ m
Flow Rate: 1.0 mL/min



In comparison to polymeric columns, the Extend-C18 has superior efficiency and peak shape




Match Method pH and Column Choice

Choose the Best Bonded-Phase for Each pH Range

StableBond, pH 1-6

Use at Low pH


1. Uses bulky silanes
2. Non-encapped



Eclipse XDB, pH 2-9

First choice at Low and Mid pH


1. eXtra Densely Bonded dimethylalkylsilanes
2. proprietary double-encapping



Bonus-RP, pH 2-8

Use at Low and Mid pH


1. polar alkyl phase
2. triple encapped
3. uses bulky silanes



Extend-C18, pH 2-11.5

Use at High pH

1. unique bidentate structure
2. double encapped



Agilent Technologies

001241P2.PPT

Summary

This method development scheme follows an approach of trying different pH's for ionizable compounds with an optimum bonded-phase for both small molecules and large biomolecules.

Low pH	Eclipse XDB then StableBond
Mid pH	Eclipse XDB or Bonus-RP
High pH	Extend-C18




Appendix

Recommended Buffer Choices for High pH

<u>Buffer</u>	<u>pKa</u>	<u>Effective pH range</u>
Pyrrolidine	11.3	10.3 – 12.3
Triethylamine (TEA)	10.7	9.7 – 11.7
1-methyl-piperidine	10.3	9.3 – 11.3
glycine	9.8	8.8 – 10.8
TRIS	8.1	7.1 – 9.1
Borate	9.2	8.2 – 10.2
Ammonia	9.2	8.2 – 10.2
Diethylamine	10.5	9.5 – 11.5



Good Lifetime of Extend-C18 at High pH



Columns: 4.6 x 150 mm, 5 μ m

Purge: 50% ACN / 50% 0.02 M K_2HPO_4 , pH 11

Flow Rate: 1.5 mL/min

Temperature: 25°C

Detection: Silicate concentration by silicomolybdate color reaction

001369P1.PPT




Agilent Technologies

Reduce Column “Secondary Interactions” at Low pH

Column: Eclipse XDB-C8, 4.6 x 150 mm, 5 μ m

Mobile Phase: 85% 25 mM Na₂HPO₄ : 15% ACN, Flow Rate: 1.0 mL/min, Temperature: 35°C

Sample: 1. Phenylpropanolamine 2. Ephedrine 3. Amphetamine 4. Methamphetamine 5. Phenteramine




- Reducing the mobile phase pH reduces interactions with silanols that cause peak tailing.

SB-CN Optimizes Retention and Resolution

Phytoestrogens and Isoflavones


Columns: 4.6 x 75 mm, 3.5 μm Mobile Phase: 30% ACN: 70% NaH_2PO_4 , pH 2.5 Flow Rate: 1.0 mL/min
Temperature: 35°C Sample: 1. Estriol 2. Daidzen 3. Quercetin 4. Genistein 5. Diethylstilbestrol




SB-CN reduces analysis time by 50% and increases retention of early eluting peaks.
Method development procedure followed to get to this point.



Why Choose 300Å Columns?



(A,B) Enter Pores




(C) Will be Excluded

Molecules must enter pores to interact with bonded-phase.
Molecules must freely enter and exit pores to maximize efficiency.



Why Choose 300Å Columns?

Effect of Pore Size and Molecular Size on Peak Width, Gradient Separations




- Proper pore size selection results in sharper peaks for large molecules.



Improved Peak Shape for Large Molecules in Solution

Columns: 4.6 x 150 mm, 5 μm
Temperature: RT

Mobile Phase: 60% MeOH: 40% 0.1% TFA Flow Rate: 0.75 mL/min
Detection: UV 282 nm Sample: Tylosin (MW 916)



- The size of a molecule in solution determines which pore size column is best.
- The narrower peak width indicates unrestricted access to the pores.



Method Development Scheme

Start at Low pH

