


Modern Analytical Tools to Tackle an Old Problem: Mycotoxins in Food

Elisabeth Varga, Thomas Glauner, Emma Rennie, Bernhard Wuest, Michael Sulyok, Rainer Schuhmacher, Rudolf Krska, Franz Berthiller

Christian Doppler Laboratory for Mycotoxin Metabolism
& Center for Analytical Chemistry
Department of Agrobiotechnology (IFA-Tulln)

University of Natural Resources and Life Sciences, Vienna (BOKU), Austria

Mycotoxins are....



- Aspergillus spp.*
- Penicillium spp.*
- Fusarium spp.*
- Alternaria spp.*
- Claviceps spp.*
- etc.

*... **toxic** secondary fungal metabolites*



Mycotoxins

History




- Epidemics of ergotism have been reported as early as **857 A.D.** resulting from the ingestion of cereals infected with ergot
- **1960:** Turkey X Disease
Peanuts infected with *Aspergillus flavus*
- Isolation of the highly toxic fungal metabolite **afatoxin** and **intensive research** on other fungi and toxins produced by them
- Effects: acute toxic, nephrotoxic, **carcinogenic**, immunosuppressive, estrogenic, hepatotoxic,...
- Annual losses of **several hundred million tons of food and feed** worldwide
- 100+ countries have **regulations** for the control of mycotoxins in food and feed



Determination of Mycotoxins

Analytical Scheme



- Sampling is major source of error
 - regulation EC 401/2006: Methods of sampling and analysis for the official control of the levels of mycotoxins in foods
- Careful choice of extraction solvents
 - typically mixtures of water with organic solvents (MeOH, ACN, acetone)
- SPE, IAC, MycoSep[®], QuEChERS, dilute & shoot
- Rapid methods (mostly antibody-based) or
- **Chromatographic methods**


50 years of Mycotoxin Analysis

A Transition of Methods



Transition of methods over the last 50 years used in mycotoxin analysis based on the Web of Science®

Keywords: mycotox*
& TLC
& ELISA
& GC
& HPLC
& LC-MS



LC-MS Multi-Mycotoxin Methods

Reasons



- Unified methods save time and costs


- For **accurate quantification** of regulated mycotoxins
 - high degree of accuracy
 - easy to handle, cost effective

Reasons

- European Commission Regulation (EC) No 1881/2006 maximum levels in maize for
 - aflatoxins (AFB₁, AFB₂, AFG₁, AFG₂)
 - fumonisins (FB₁, FB₂)
 - deoxynivalenol
 - ochratoxin A
 - zearalenone
 - (HT-2 toxin and T-2 toxin)

BUT:


- Electrospray ionisation (ESI)
 - matrix effects hamper accurate mass spectrometric quantification



Accurate Mycotoxin Quantification Approaches




- Dilution of the sample
 - loss of sensitivity
- Matrix matched calibration
 - needed for each individual matrix
- Standard addition to each sample
 - multiple injections needed
- Internal calibration
 - ideally: internal standard (IS) behaves exactly the same as the analyte, but still distinctive
 - structurally related or similar compounds
zearalanone (ZAN) for zearalenone (ZEN)
 - [^2H]- or [^{13}C]-labelled compounds



Stable Isotope Dilution Assay (SIDA)

Multiple Mycotoxins - Sample Preparation




+IS

Stable Isotope Dilution Assay

UHPLC-MS/MS Setup




- 1290 Infinity UHPLC (Agilent Technologies)
 - C18 column
 - methanol-water gradient
- 6490 QQQ MS/MS (Agilent Technologies)
 - single run, fast polarity switching, dynamic MRM





Analyte	External calibration R_A [%] \pm RSD [%] ¹⁾
Aflatoxin B ₁	35 \pm 4
Aflatoxin B ₂	42 \pm 4
Aflatoxin G ₁	46 \pm 5
Aflatoxin G ₂	40 \pm 6
Deoxynivalenol	41 \pm 10
Fumonisin B ₁	330 \pm 6
Fumonisin B ₂	181 \pm 7
HT-2 toxin	134 \pm 6
Ochratoxin A	167 \pm 14
T-2 toxin	118 \pm 5
Zearalenone	91 \pm 6



¹⁾ Apparent recovery \pm relative standard deviation



Analyte	External calibration R_A [%] \pm RSD [%] ¹⁾	Internal calibration R_A [%] \pm RSD [%] ¹⁾
Aflatoxin B ₁	35 \pm 4	105 \pm 6
Aflatoxin B ₂	42 \pm 4	100 \pm 4
Aflatoxin G ₁	46 \pm 5	101 \pm 5
Aflatoxin G ₂	40 \pm 6	101 \pm 8
Deoxynivalenol	41 \pm 10	99 \pm 9
Fumonisin B ₁	330 \pm 6	101 \pm 10
Fumonisin B ₂	181 \pm 7	88 \pm 7
HT-2 toxin	134 \pm 6	98 \pm 7
Ochratoxin A	167 \pm 14	93 \pm 7
T-2 toxin	118 \pm 5	99 \pm 6
Zearalenone	91 \pm 6	103 \pm 11

¹⁾ Apparent recovery \pm relative standard deviation



Analyte	External calibration R_A [%] \pm RSD [%] ¹⁾	Internal calibration R_A [%] \pm RSD [%] ¹⁾	R_E [%] ²⁾
Aflatoxin B ₁	35 \pm 4	105 \pm 6	108
Aflatoxin B ₂	42 \pm 4	100 \pm 4	107
Aflatoxin G ₁	46 \pm 5	101 \pm 5	109
Aflatoxin G ₂	40 \pm 6	101 \pm 8	111
Deoxynivalenol	41 \pm 10	99 \pm 9	107
Fumonisin B ₁	330 \pm 6	101 \pm 10	103
Fumonisin B ₂	181 \pm 7	88 \pm 7	97
HT-2 toxin	134 \pm 6	98 \pm 7	109
Ochratoxin A	167 \pm 14	93 \pm 7	102
T-2 toxin	118 \pm 5	99 \pm 6	110
Zearalenone	91 \pm 6	103 \pm 11	109


¹⁾ Apparent recovery \pm relative standard deviation

²⁾ Extraction recovery

Analyte	LOQ ($\mu\text{g}/\text{kg}$) ¹⁾	MLs ($\mu\text{g}/\text{kg}$) ²⁾
Aflatoxin B ₁	0.1	2.0 – 5.0 (0.1)
Aflatoxin B ₂	0.2	sum aflatoxins 4.0 – 10.0
Aflatoxin G ₁	0.1	
Aflatoxin G ₂	0.35	
Deoxynivalenol	12	500 – 1750 (200)
Fumonisin B ₁	4.5	sum fumonisins 800 – 4000 (200)
Fumonisin B ₂	4	
HT-2 toxin	3	indicative values ³⁾ sum 15-1000
T-2 toxin	0.3	
Ochratoxin A	0.4	3.0 – 5.0 (0.5)
Zearalenone	3	50 (20)

OTA: 0.45 $\mu\text{g}/\text{kg}$


— 404.1 → 238.9 S/N 13
 ... 404.1 → 102.1 S/N 14



1) Limit of quantification (S/N = 10) for maize

2) Maximum levels for various commodities according to European Union Commission Regulation 1881/2006 and its amendments (numbers in brackets: baby food)

3) Indicative values in various cereals according to European Union Recommendation 165/2013



Validation of a Stable Isotope Dilution Assay for the Accurate Quantitation of Mycotoxins in Maize Using UHPLC/MS/MS

Application Note 5991-2808

Anal Bioanal Chem (2012) 402:2675–2686

DOI 10.1007/s00216-012-5757-5

PAPER IN FOREFRONT

Stable isotope dilution assay for the accurate determination of mycotoxins in maize by UHPLC-MS/MS

Elisabeth Varga • Thomas Glauner • Robert Köppen •
Katharina Mayer • Michael Sulyok •
Rainer Schuhmacher • Rudolf Krska • Franz Berthiller

Open access





- Unified methods save time and costs

- For **accurate quantification** of regulated mycotoxins
 - high degree of accuracy
 - easy to handle, cost effective

- For **(semi-)quantification of 200+ mycotoxins**
 - analysis of food and feed mixtures of several commodities
 - monitoring of changes in regional fungal spread (climate change)

Multi-Mycotoxin Quantification

The Need to Determine More Toxins in Food



Berthiller et al. (2005) J. Chrom. A

- **9 mycotoxins**, APCI, polarity switching
- Mycosep clean-up
- LOD deoxynivalenol: **65 pg** on column

Sulyok et al. (2006) Rapid Commun. Mass Spectrom.

- **39 mycotoxins**, ESI, 2 runs (pos. and neg.)
- no clean-up
- LOD deoxynivalenol: **10 pg** on column

extended to:



- **87 toxins** (Sulyok et al., 2007)
- **186 metabolites** (Vishwanath et al., 2009)

Malachová et al. (2014) J. Chrom. A

- **295 fungal and bacterial metabolites**
- LOD deoxynivalenol: **0.3 pg** on column

LC-Tandem-Mass Spectrometry

Multiple Reaction Monitoring (MRM)




static (classical)
MRM



dynamic
(scheduled) MRM

Multi-Mycotoxin Quantification


Method Workflow



based on: Sulyok *et al.* (2006) *Rapid Commun. Mass Spectrom.* 20(18): 2649-2659

Multi-Mycotoxin Quantification


Method Workflow



based on: Sulyok *et al.* (2006) *Rapid Commun. Mass Spectrom.* 20(18): 2649-2659

Multi-Mycotoxin Quantification

Method Workflow



LC-MS/MS (1290 series UHPLC + QQQ 6460)
eluent (acidified methanol/water mixtures)
pos. and neg. ionisation mode


Mass Hunter (manual correction);
linear, 1/x weighted calibration curves

Multi-Mycotoxin Quantification

Chromatogram for a total of 191 compounds




- 65 analytes validated for almonds, hazelnuts, peanuts and pistachios
- Semi-quantitative method for 126 analytes



Multi-Mycotoxin Quantification

Market Survey - Nuts






Anal Bioanal Chem (2013) 405:5087–5104

DOI 10.1007/s00216-013-6831-3

RESEARCH PAPER

Development and validation of a (semi-)quantitative UHPLC-MS/MS method for the determination of 191 mycotoxins and other fungal metabolites in almonds, hazelnuts, peanuts and pistachios

Elisabeth Varga • Thomas Glauner • Franz Berthiller •
Rudolf Krska • Rainer Schuhmacher • Michael Sulyok



Screening and Quantitation of 191 Mycotoxins and Other Fungal Metabolites in Almonds, Hazelnuts, Peanuts, and Pistachios Using UHPLC/MS/MS

Application Note 5991-4991



- Unified methods save time and costs

- For **accurate quantification** of regulated mycotoxins
 - high degree of accuracy
 - easy to handle, cost effective

- For **(semi-)quantification of 200+ mycotoxins**
 - analysis of food and feed mixtures of several commodities
 - monitoring of changes in regional fungal spread (climate change)

- For **screening** of a huge variety of mycotoxins
 - identification of mycotoxins in unlikely matrices
 - screening of co-occurring analytes to evaluate synergistic effects

Contaminant Screening

Using HR-MS instruments



■ Advantages

- + in MS mode -> **fast data acquisition**
- + **post-acquisition data analysis** (at least in MS mode)
- + potential screening for **unknown contaminants**
MS-mode: high mass accuracy -> sum formulas
MS/MS-mode: sum formulas of fragments
- + **unambiguous identification possible without standards**

■ Disadvantages

- higher limit of detections compared to QQQs



6550 iFunnel QTOF

Mycotoxin Screening Database



■ Contains compound information:

- name
- formula
- exact mass
- retention time
- external IDs (e.g. CAS)
- IUPAC name
- structural information

MassHunter PCDL Manager for Forensics and Toxicology - C:\MassHunter\PCDL\Mycotoxins_AM_PCDL.cdb

File Edit View PCDL Links Help

Find Compounds

Single Search Batch Search Batch Summary Edit Compounds Spectral Search Browse Spectra Edit Spectra

Mass: [M+H]⁺ Neutral [M-H]⁻

Mass tolerance: 10.0 ppm mDa

Retention time: Require

RT tolerance: 0.1 min

Ion search mode

Include neutrals

Include anions

Include cations

Formula:

Name:

Notes:

IUPAC:

CAS:

ChemSpider:

Molecule: Structure MOL Text

Notes: Fungal metabolite
Synonyms: AOH

**mycotoxins and related metabolites:
455 entries**

Print/Copy in Summary Format Single Search Results: 455 hits

Compound Name	Formula	Mass	Anion	Cation	RT (min)	CAS	ChemSpider	IUPAC Name	Spectra
▶ AOH / Alternariol	C14H10O5	258.05282	<input type="checkbox"/>	<input type="checkbox"/>		641-38-3	4514301	3,7,9-Trihydroxy-1-methyl-6H-benzo[c]chromen-6-...	6
Aphidicolin	C20H34O4	338.24571	<input type="checkbox"/>	<input type="checkbox"/>		38966-21-1	10280269	(1S,2S,5R,6R,7R,10S,12R,13R)-6,13-Bis(hydrox...	0
Apicidin	C34H49N5O6	623.36828	<input type="checkbox"/>	<input type="checkbox"/>		183506-66-3	5293532	(3S,6S,9S,15aR)-9-[(2S)-2-Butanyl]-6-[[1-methoxy...	6
Ascomycin	C43H69NO12	791.48198	<input type="checkbox"/>	<input type="checkbox"/>		104987-12-4	4445297	15,19-Epoxy-3H-pyrido[2,1-c][1,4]oxaazacyclotric...	6
Aspercolorin	C25H28N4O5	464.20597	<input type="checkbox"/>	<input type="checkbox"/>		29123-52-2		Cyclo[β-alanyl(α-prolyl-2-amino-5-methoxybenzoyl)ph...	6
Aspergillimide	C20H29N3O3	359.22089	<input type="checkbox"/>	<input type="checkbox"/>		195966-93-9		Spiro[5H,6H-5a,9a-(1minomethano)-1H-cyclopent[...	3
Asperlactone	C9H12O4	184.07356	<input type="checkbox"/>	<input type="checkbox"/>		76375-62-7	137956	(5R)-5-[(1S)-1-Hydroxyethyl]-3-[(2S,3S)-3-methyl-2-...	3
Asperoxin A	C21H19N3O5	393.13247	<input type="checkbox"/>	<input type="checkbox"/>		223130-52-7		6H,12H-Oxepino[2',3':4,5]pyrimido[1,2-a]pyrrolo[2...	3
Asperthecin	C15H10O8	318.03757	<input type="checkbox"/>	<input type="checkbox"/>		10089-00-6	2341183	1,2,5,6,8-Pentahydroxy-3-(hydroxymethyl)-9,10-an...	0
Asoinolide A	C10H16O3	184.10994	<input type="checkbox"/>	<input type="checkbox"/>		188605-14-3	10268506	(6R,7E,10R)-6-Hydroxy-10-methyl-3,4,5,6,9,10-he...	0

Screening HR-MS/MS Library



- Acquiring target HR-MS/MS spectra from standards
 - import using mass correction based on (fragment) structures
 - only explainable fragments are imported
 - exact masses for precursors and fragments

302 compounds with spectra

The screenshot displays a software interface for HR-MS/MS screening. The top section contains search parameters: Precursor ion, Tolerance (200 ppm), Collision energy (2.0 eV), Ion polarity (Any), and Ionization mode (Any). Below this is a table titled 'Spectra for compound: T2 / T-2 Toxin' with columns for Compound Name, Ion Species, Precursor Ion, CE (V), Polarity, Ionization, and Instrument. The table lists several entries for T2 / T-2 Toxin with different ion species and precursor ions. To the right of the table is a mass spectrum plot showing Abundance vs. m/z. The plot has several peaks, with the most prominent ones at m/z 215.10666 (94.89% abundance) and 305.13834 (100.00% abundance). Other labeled peaks include 105.06988 (12.05%), 385.15948 (10.08%), and 449.21698 (2.62%).

Compound Name	Ion Species	Precursor Ion	CE (V)	Polarity	Ionization	Instrument
T2 / T-2 Toxin	(M+H) ⁺	467.22756	10	Positive	ESI	QTOF
T2 / T-2 Toxin	(M+H) ⁺	467.22756	20	Positive	ESI	QTOF
T2 / T-2 Toxin	(M+H) ⁺	467.22756	40	Positive	ESI	QTOF
T2 / T-2 Toxin	(M+NH4) ⁺	484.25411	10	Positive	ESI	QTOF
T2 / T-2 Toxin	(M+NH4) ⁺	484.25411	20	Positive	ESI	QTOF
T2 / T-2 Toxin	(M+NH4) ⁺	484.25411	40	Positive	ESI	QTOF

Single Search Results: 455 hits

Compound Name	Formula	Mass	Anion	Cation	RT (min)	CA5	ChemSpider	IUPAC Name	Spectra
Sulochrin	C17H16O7	332.08960	<input type="checkbox"/>	<input type="checkbox"/>	6.210	519-57-3	141044	Methyl 2-(2,5-dihydroxy-4-methylbenzoyl)-5-hydrox...	6
T2 / T-2 Toxin	C24H34O9	466.22028	<input type="checkbox"/>	<input type="checkbox"/>	7.852	21259-20-1	21476745	(2alpha,3beta,4alpha,8alpha)-4,15-Diacetoxy-3-h...	6
T-2 Tetraol	C15H22O6	298.14164	<input type="checkbox"/>	<input type="checkbox"/>		34114-99-3	8079985	(3R,4α,8α)-12,13-Epoxytrichothec-9-ene-3,4,8,15...	0
T-2 Triol	C20H30O7	382.19915	<input type="checkbox"/>	<input type="checkbox"/>		34114-98-2	2299199	(3R,4α,8α,12E)-3,4,15-Trihydroxy-12,13-epoxytric...	3


LC-HR-MS/MS Screening

Classical Workflow





1) HR-MS full scan

total ion chromatogram
of a naturally contaminated hazelnut sample



2) Find-by-Formula (FBF)

extracted ion chromatogram &
extracted MS spectrum of mycophenolic acid



LC-HR-MS/MS Screening

Classical Workflow




3) MS/MS scan

targeted MS/MS scan
specify precursor mass, retention time (RT),
RT window and collision energy

4) Library confirmation

extracted MS/MS spectra at different
collision energies

CE 10 eV



LC-HR-MS/MS Screening

Classical Workflow



- Library hit confirming the presence of the immunosuppressive mycotoxin mycophenolic acid in a hazelnut sample

Compound Identification Results: Cpd 70: Mycophenolic acid

Automatically Show Columns

ID Techniques Applied

LibSearch-FBF

Best	Name	Formula	Mass	Mass (Tgt)	Diff (ppm)	Diff (mD)	RT	Score (Tg)	m/z	Species	Score (Lib)	ID Source
1	Mycophenolic acid	C ₁₇ H ₂₀ O ₆	320.1253	320.126	2.26	0.72	15.863	77.27	321.1349 343.11...	(M+H) ⁺ (M+Na) ⁺	95.13	LibSearch-FBF

CE	Name	Num Peaks	m/z (prec.)	Reverse Score	Score (Lib)
10	Mycophenolic acid	8	321.1333	97.41	97.41
20	Mycophenolic acid	6	321.1333	99.65	99.65
40	Mycophenolic acid	16	321.1333	92.29	92.29

Score (iso. abund)	Score (mass)	Score (MS)	Score (iso. spacing)	Abund	Species	m/z
89.64	85.77	77.27	45.42	41720.7	(M+H) ⁺	321.1349
0.02	95.59	59.55	58.92	812708.5	(M+Na) ⁺	343.1174

mycophenolic acid (C₁₇H₂₀O₆)
produced by several *Penicillium* spp.
low acute toxicity, immunosuppressive



LC-HR-MS/MS Screening


All Ions Approach



1) HR-All Ions MS/MS acquisition

Fragmentation without precursor selection


Measurement with “low energy channel” (no collision energy, MS scan) and at least one “high energy channel” (collision energy e.g. 20 eV)





2) Database search and library confirmation

- I) Database search: applying “Find by Formula” algorithm on low channel
- II) Automatic extraction of fragment ion chromatograms in high energy channel
- III) Co-elution scoring of fragment and parent EICs for qualification and library confirmation



LC-HR-MS/MS Screening

All Ions Approach



Co-elution score

Major function to qualify compounds based on similarity of peak shapes


Compound Identification Results: Cpd 3: Mycophenolic acid

Automatically Show Columns

ID Techniques Applied: FBF

Best	Name	Formula	m/z	Mass	Mass (Tgt)	Diff (ppm)	Score (Tgt)	RT
	Mycophenolic acid	C17 H20 O6	321.1342 343.1166	320.1248	320.126	3.82	87.55	15.818

m/z	Coelution Score	Flags(FIs)	Height	SNR	RT	RT Diff	CE	FV	Compound Name
207.0652	99.8	Qualified	473539		15.807	0.011	20		Mycophenolic acid
159.0441	98.2	Qualified	282553.9		15.807	0.011	40		Mycophenolic acid
303.1227	97.6	Qualified	280945.9		15.83	0.011	10		Mycophenolic acid
275.1278	97.8	Qualified	73230.9		15.83	0.011	10		Mycophenolic acid
103.0542	86	Low coelution score	52768.4		15.807	0.011			Mycophenolic acid




All Ions approach

- worse limits of detection (more noise)
- limited with collision energies
- + only one injection for each polarity
- + post-acquisition MS/MS information



- **General unknown screening** using HR-MS instrumentation
- Identification by **HR-MS/MS library spectra**
- Standards not mandatory for unambiguous identification
- **Retrospective data-analysis** also in MS/MS mode possible



**Screening and Verifying Mycotoxins
in Food with Q-TOF LC/MS and an
Accurate Mass Library**


Application Note 5991-5667

Conclusions


Summary



- **LC-MS/MS** is increasingly used for the **simultaneous quantification and identification** of mycotoxins and other food contaminants
- **Stable Isotope Dilution Analysis** is used for quantification of mycotoxins in food and feed offering **highest possible accuracy** at affordable additional costs for internal standards
- (Semi-)quantitative **LC-MS/MS methods allow** for screening and **quantification of 200+ mycotoxins**
- **LC-HR-MS screening** with exact mass MS/MS libraries allow **unambiguous identification for a practically unlimited number of compounds** even without standards and **in retrospective data analysis**



Thanks for the attention!



Christian Doppler Laboratory for Mycotoxin Metabolism
& Center for Analytical Chemistry
Department of Agrobiotechnology (IFA-Tulln)
**University of Natural Resources and Life Sciences,
Vienna (BOKU), Austria**

franz.berthiller@boku.ac.at

