

# Understanding the GC Inlet

Which one is more  
appropriate for your method?

Simon Jones  
GC Applications Engineer

# Types of Inlets

Purged Packed

Split / Splitless

Cool On Column

Programmable Temperature Vaporization

Volatiles Interface

Multi Mode Inlet

# Where to Begin???

What are the requirements of the method?

Trace level analysis?

% level analysis?

High temperature application?

Packed column??

What do you know about the sample?

Dirty or clean?

Residual solvent?


Volatility range?

# Inlet Use Guide


Inlet	Column	Mode	Sample Concentration	Comments	Sample to Column
Split / Splitless	Capillary	Split Purged Split Splitless Purged Splitless	High High Low Low	Most commonly used inlet. Very Flexible	Very Little Very Little All All
Cool-On-Column	Capillary	N/A	Low or labile	Minimal discrimination and decomposition	All
Packed	Packed Large Capillary	N/A N/A	Any Any	OK if resolution is not critical	All All
Programmed Temperature Vaporization	Capillary	Split Pulsed Split Splitless Pulsed Splitless Solvent Vent	High High Low Low Low	Not great for HOT injections.  Can concentrate analytes and vent solvent	Very Little Very Little All All Most
Volatiles Interface	Capillary	Direct Split Splitless	Low High Low	Purge & Trap / Headspace	All Very Little All
Multi-Mode	Capillary	Split Pulsed Split Splitless Pulsed Splitless Solvent Vent	High High Low Low Low	Flexibility of standard S/SL inlet and PTV	Very Little Very Little All All Most

# COC Flow diagram

## Cool-On-Column



# COLD ON-COLUMN INJECTION PORT



# COC – Mode of Operation

## Oven Track Mode

Inlet temperature stays 3°C above the oven temperature

## Temperature Programmed Mode

Can program 3 temperature ramps

# COC Benefits

Sample Discrimination does not occur

If operated correctly, accurate and precise results are obtained

Can be used to gauge liner activity

Very Gentle sample introduction – limits decomposition of analytes.  
Good for Labile compounds!

Used for high temperature applications.

Biodiesel



# COC inlet

Key parameters to be used:

\*Starting inlet temperature must be below the boiling point of the solvent being used!!!

Guard column / Retention Gap strongly recommended to help protect the analytical column, and focus the sample

# COC Troubleshooting Tips

## Bent needles


- using the wrong size needle or insert
- insert has burrs

## Plugged needles due to septum coring

## Lost peak shape


- examine inlet end of column with a magnifier and flashlight, looking for discoloration or particles

## Injection volume too large



# PTV Flow Diagram

## Programmable Temperature Vaporization



# PTV modes of operation

Split	Major component analysis
Pulsed Split	Best used with low split flows
Splitless	Trace level analysis
Pulsed Splitless	More efficient sample transfer
Solvent Vent	Large Volume injections

# PTV Inlet

Not good for Hot injections

Minimal inlet discrimination – closest to COC

Large volume injections


Solvent vent mode

Can eliminate volatile components of the sample

Rapid Heating and Cooling


Cold trapping of Gas Injection

# Volatiles Interface



# Volatiles Interface

**Small internal volume (35 ul)**  
**Inert surface (Silcosteel coated)**  
**Helium and H<sub>2</sub> carrier only – no N<sub>2</sub>**



# Volatiles Interface Modes of Operation

Split

Splitless


Direct



# Volatiles Interface

## Splitless Injection


Pre-Run to Sampling End



Switches to Forward Flow Control to provide pneumatic stability during “Desorb” - while the sampler trap is in the carrier flow path.

# Volatiles Interface


Direct Injection -- idle




# Volatiles Interface -- Uses

Used for direct connection between Headspace / Purge & Trap  
Cannot do Manual Injections!

# Purged Packed Inlet



# Purged Packed



# PP Inlet Uses

## Packed columns

Can be used with 0.53 mm , or 0.32 mm ID columns when high flows ~10 mL/min are used

When column dimensions are not defined, the inlet functions in a 'flow' mode

Packed columns best run in flow mode, capillary columns preferred to run in pressure mode.

# PP Inlet

Very small expansion volume

More active than most inlets


- Glass liner helps minimize activity

- Glass packed columns have best reproducibility

- Small surface area of the liner minimizes the amount of active sites

Not Recommended for Capillary Columns smaller than 0.53 mm

# Split/Splitless Inlet





# S/SI Modes of Operation

Split

Pulsed Split

Splitless

Pulsed Splitless



# Split Injections - Considerations

Dirty Samples are OK - backflushing

Wide Analyte Boiling Range

Solvent Properties

- Wide Boiling Point Range
- Wide Polarity Range

Discrimination can be due to liner or inlet temperature

# Split Injections - Inertness

## More inert than splitless

- Higher velocity through the inlet
- Less exposure to inlet hardware

## Glass wool is a compromise


- Exhibits some activity
- Greatly improves fluidic performance – mixing of the vaporized sample is important for uniform splitting

# Split Injections - recommended Liners

Agilent p/n 5190-2295

Wiped needle improves

- precision
- peak shape
- discrimination



# Split Injections - Maximizing Sensitivity

## Increase Injection Volume

- liner dependent (use the Pressure-Volume Calculator)
- 2  $\mu\text{L}$  maximum




## Reduce Split Ratio

- go from 50:1 to 10:1
- 5:1 practical lower limit for liquid injections (for 250 - 320  $\mu\text{m}$  i.d. columns)
- 1:1 possible for gas injections with correct liner

## Use Pulsed Injection

# Split Injector

## Injection Volume



DB-1, 15 m x 0.25 mm I.D., 0.25 µm

60°C for 1 min, 60-180°C at 20°/min; Helium at 30 cm/sec

1. n-heptane 2. toluene 3. n-decane 4. n-butylbenzene 5. n-tridecane

# Split Injections - Pulsed

May be easiest approach for active analytes  
(example: 2,4 dinitrophenol)

Using “pulsed mode” may result in peak doublets due to system ramping down at 99 psi/min

Instead, use “ramped pressure” or “ramped flow” mode to do your pulse

- set initial pressure (or flow) to 3x-5x your normal starting setpoint
- hold this higher pressure for 0.1 - 0.3 min
- ramp at 20 psi/min (or 10 mL/min/min) down to your normal starting setpoint

# Split Injections - Fast GC Considerations

Faster than splitless because you can start at a higher initial oven temp, thereby decreasing cycle time

Easiest of the injection techniques to speed up

For 100  $\mu\text{m}$  i.d. and smaller columns

- narrower i.d. liners may be necessary to maintain input peak width

Using higher flows with normal columns

- Loose some resolution
- Better inertness
- Larger injections possible



# Split Injections - Troubleshooting

Column pressures <10 psi

- The pressure pulse from evaporating solvent can cause discrimination and poor precision

Liner residence times < 0.5 sec (> 200 ml/min)

- poor mixing will cause discrimination

No glass wool

Solvents with high expansion ratio **Backflash**

Column position - top to bottom, side to side

Large bore, short columns with a high split ratio

# Splitless Injections - Considerations

Dirty samples are OK - backflushing

Analyte Boiling Range - Wide (but narrower than split)

- early eluters need bp difference vs solvent

Solvent Properties

- Wide Boiling Point Range
  - but consider bp of earliest eluting analyte
- Wide Polarity Range (but narrower than split)
  - Water and Methanol worst choices

Greater Sample Residence Time

Lower Inlet Temperatures can be used

Better for Labile Compounds

# Splitless Injections - Inertness

Less inert than COC

- liner and inlet interaction

Less inert than Split

- longer residence time in inlet and on glass wool
- used for trace analysis, so there's a greater chance of analyte loss

# Splitless Injections - Discrimination


## Improper purge time

- short purge times cause loss of late eluters
- long purge times cause solvent tail interference with early eluters


## Improper initial oven temp

- too high of a temp prevents solvent effect and a loss of early eluters
- too low of a temp extends run time

# Splitless Injections – Splitless Time (purge time on)






Purge time too long results in large solvent tail



0.75 min purge time clips solvent tail

# Splitless Injector

## Injection Volume



DB-1, 15 m x 0.25 mm I.D., 0.25 µm


60°C for 1 min, 60-180°C at 20°/min; Helium at 30 cm/sec

1. n-decane 2. n-dodecane 3. n-tetradecane 4. n-hexadecane


# Splitless Injector

## Injector Temperature

200°C



250°C



DB-1, 15 m x 0.25 mm I.D., 0.25  $\mu$ m

50°C for 0.5 min, 50-325°C at 20°/min; Helium at 30 cm/sec

Phthalates: 1. dimethyl 2. diethyl 3. dibutyl 4. benzylbutyl 5. bis(2-ethylhexyl) 6. dioctyl

# Splitless Injector

## Sample Re-focusing

Sample re-focusing improves efficiency

Use low column temperature to refocus solvent  
- called the *solvent effect*

Use cold trapping



# Splitless Injector


## Solvent Effect

Initial column temperature at least **10°C below** sample solvent boiling point

Required to obtain good peak shapes unless cold trapping occurs

Rule of thumb, if solute BP >150°C above initial column temperature, the solute will cold trap



Cold trapping has greater efficiency than solvent effect



# Splitless Injector

Initial Column Temperature

Hexane Solvent (BP = 68-69°C)



DB-1, 15 m x 0.25 mm I.D., 0.25  $\mu$ m


50°C or 70°C for 0.5 min, to 210°C at 20°/min; Helium at 30 cm/sec

1. n-decane 2. n-dodecane 3. n-tetradecane 4. n-hexadecane



# Splitless Injector

## Reverse Solvent Effect/Polarity Miss-Match




DB-1, 15 m x 0.25 mm I.D., 0.25  $\mu$ m

50°C for 1 min, 50-210°C at 20°/min; Helium at 30 cm/sec

1. 1,3-DCP 2. 3-hexanol 3. butyl acetate 4. 1-heptanol 5. 3-octanone 6. 1,2-dichlorobenzene

# Retention Gap


Also Called A Guard Column



Usually 2-10 meters long and same diameter as the column  
(or larger if needed)

# Splitless Injector

3 m x 0.25 mm I.D. Retention Gap



DB-1, 15 m x 0.25 mm I.D., 0.25  $\mu$ m

50°C for 1 min, 50-210°C at 20°/min; Helium at 30 cm/sec

1. 1,3-DCP 2. 3-hexanol 3. butyl acetate 4. 1-heptanol 5. 3-octanone 6. 1,2-dichlorobenzene

# EPC for Splitless Pulsed Injection

Pressure Pulse contains sample expansion and transfers analytes to the column faster.

## **Pulsed Splitless**

- sample containment more critical than in split injection
- sharper peaks than in traditional splitless injection
- two new parameters to set:
  - pulse pressure and pulse time

## **Typical starting point**

- Pulse pressure = double resting pressure
- Tie pulse time to purge time

# Splitless Injections – Fast GC Considerations

Slower than split because you must start at a lower initial oven temp, thereby increasing cycle time

Difficult to use with 100  $\mu\text{m}$  i.d. columns

- smaller injection size
- smaller liner volume
- retention gap

Using higher flows with normal columns

- Lose some resolution
- Better inertness
- Larger injections possible

# Splitless Injections – Starting

Injection Volume = 1  $\mu$ L

- Check the Pressure-Volume Calculator

Initial Oven Temp = 10°C < solvent boiling point

Purge Flow = 20 to 60 mL/min

Purge Time = 0.75 min

- Sweep with 2 liner volumes of carrier gas

No pulse

Try to avoid water and methanol as solvents



# Splitless Injections – Troubleshooting Tips

## Injecting too much

- column overload = poor peak shape
- inlet overload = poor reproducibility
  - ghost peaks in subsequent blanks are possible

## No glass wool

- poor mixing
- dirt on column

## Glass wool

- reacts with trace components

# Splitless Injections – Troubleshooting Tips

If you think you have an inlet issue related to splitless injections

- then

Run a 10:1 split injection

- or

Make up a standard at 10x concentration and run a 10:1 split injection

When I changed from split to splitless I didn't see an increase in response!!!

Purge Time set to '0'

# Split Vent Trap

What is it???



Where is it???



# Improved S/SI inlet Inertness



Inert Inlet Weldment



Ultra Inert Gold Seal



Intuvo GC

Guard Chip

Gold Seal  
Guard Column

# Split vs. Splitless Injection Technique - Summary

## SPLIT:

- Best Injection Efficiency
- Less sensitive
- Prone to discrimination
- Proper liner choice more important

## SPLITLESS:

- Poor Injection efficiency
  - solvent effect
  - retention gap
- Good for Trace level detection
- Solvent/column polarity match more critical

but...what if you are already running maximum injection volume, pulsed splitless and still need more sensitivity...

# MultiMode Inlet

7890 standard  
pneumatics


7890 standard  
capillary fitting



7890 turn-top

Uses 7890 S/SL  
liners, septa and  
o-rings

# MMI Inlet



# Programmable Temperature Vaporizing (MMI) Inlet – injection modes

Hot split/splitless (also pulsed)



- similar to the S/SL inlet using the **same liners**
- all previous S/SL discussions apply here

Cold split/splitless (also pulsed)

- Significantly more inert than hot splitless
- Can inject 3-5 uL with no solvent venting
- Better sensitivity than hot splitless because large vapor cloud is not formed which travels outside the liner and portions are lost

LVI-Solvent Vent

- An extension of cold splitless
- Large volume injection for maximum sensitivity

Direct Mode

Uses a Direct Connect Liner – simulates COC \* NO purge



# MultiMode (MMI) Inlet Features

## Hardware

Temperature range of -160C to 450C

Heating @ 15C/sec (900C/min)

Septum/Liner Easily Exchangeable using Turn Top Inlet

Injection Modes: Hot S/SL, Cold S/SL, all in pulsed mode, solvent vent mode, residue removal mode

Support for single stroke injections from 0.1  $\mu$ L to 250  $\mu$ L

EPC Compatible with Packed Liners

Compatible with 7890A, 5975C, 7683, CTC Combi PAL

## Software

Ten temperature ramps

Wizard for setting up large volume injections

Fully integrated into ChemStation, MSD ChemStation, EZChrom, MassHunter

# MultiMode Inlet Solves Many Problems

## **Performing large volume injection (LVI) of relatively clean samples?**

- programmable injection slows solvent evaporation and maximizes analyte transfer into the column/detector
- decrease MDL by injecting more sample

## **Injecting dirty samples?**

- matrix vent, backflush and easy liner changing minimize dirty sample affects

## **Performing analyses of high molec. wt. and/or thermally labile compounds?**



- temperature programming of Multimode inlet elutes analytes at the lowest possible temperature, minimizing breakdown and absorption
- discrimination of high molec. wt. compounds is minimal allowing HT GC

# MultiMode GC Inlet - Cold Injections

- No syringe-needle discrimination; Minimal inlet discrimination
- No special syringes, liners or consumables
- Large volume injection (5ul to 250ul) - lower detection limits
- Solvent vent/matrix vent - decrease interference / maintenance
- Flexibility (hot/cold split/splitless, temperature programmed vaporization)
- Cold trapping in liner - improves chromatographic peak shape, resolution
- Capillary column backflush with CFT - decreases cycle time, maintenance



# MMI Column Installation



- Graphite ferrules are recommended over Vespel
- No SilTite Ferrules



Thread the column into the column adapter – Stabilize the column adapter with a 5/16" wrench



Tighten the column with a 1/4" wrench – continue to hold the column adapter with a 5/16" wrench


# MMI –Intuvo GC

Still Same function









Uses Same Liners as 7890

Uses a Guard chip

Acts like a gold seal



# Inlet Column Installation Guide










Inlet	Diagram	Procedure
Split/Splitless		<p>Place a septum over the column, then the column nut and ferrule. Trim the end of the column with a column cutter.</p> <p>Pull the column back so that 4-8 mm of column is extending past the end of the ferrule.</p> <p>Thread the column nut and column into the inlet and tighten slightly past where the column grabs – retighten after heating.</p>
Purged Packed		<p>Place a septum over the column, then the column nut and ferrule. Trim the end of the column with a column cutter.</p> <p>Pull the column back so that 1-2 mm of column is extending past the end of the ferrule.</p> <p>Thread the column nut and column into the inlet and tighten slightly past where the column grabs – retighten after heating.</p>
Multimode		<p><b>NOTE:</b> Make sure the column adapter nut on the inlet base is <b>fully threaded on and spinning freely</b> – Collar Up!</p> <div style="text-align: center;">  </div> <p><b>Tighten with two wrenches</b> - 1/4" and 5/16" To prevent damage the inlet threads.</p> 
Cool On Column		<p>Insert the column all the way into the inlet until you feel the spring tension – do not withdraw. <b>The column cut is critical.</b></p> <p><b>Tighten with two wrenches</b> - 1/4" and 5/16" to avoid damaging the inlet.</p>
PTV		<p>There should be 17mm of column above the graphpak ferrule – the graphpak ferrule should be installed with the graphite end towards the inlet base. The column nut is slotted. Use a 5 mm wrench to tighten the fitting.</p>
Volatiles Interface		<p>There is a longer column nut for the VI so that you don't have to remove the inlet block. Part Number - G3504-20504</p>

## Column Installation / Pre-swaging tool



# Inlet Liners

## Split/Splitless -- MMI Liners

Description	Volume (μL)	ID (mm)
<b>Split Inlet Liners</b>		
 Low pressure drop, Ultra Inert Liner with glass wool	870	4
 Straight, Ultra Inert Liner with glass wool	990	4
<b>Splitless Inlet Liners</b>		
 Single taper, Ultra Inert Liner	900	4
 Single taper, Ultra Inert Liner with glass wool	900	4
 Splitless, double taper Ultra Inert Liner, no wool	800	4
 Dimpled, splitless, Ultra Inert Liner	200	2
 Splitless, straight, Ultra Inert Liner	250	2
 Straight, Ultra Inert Liner	60	1
 Straight Ultra Inert Liner for SPME	35	0.75

## Purged Packed Inlet liner

Disposable glass liner, 170 μL internal volume

## PTV liners

Description	ID (mm)	Volume (μL)
<b>Liners for Septumless PTV Inlet, G3501A, G3502A, G3503A</b>		
PTV liner, single baffle, glass wool, deactivated	2	180
PTV liner, single baffle, deactivated	2	200
PTV liner, multi baffled, deactivated	1.8	150
PTV liner, sintered glass, deactivated	1.5	112
<b>Liners for High Temperature PTV Inlet, G3506A</b>		
PTV liner, high temperature, quartz	3.4	713
PTV liner, high temperature, borosilicate	3.4	668


# Inlet Tools

<http://www.agilent.com/en-us/support/gas-chromatography/gccalculators>

**Calculators**

Vapor Volume Calculator | Pressure Flow Calculator | Method Translator | Solvent Vent Calculator

Liner capacity exceeded! Choose a liner of greater volume or modify method parameters.

<p><b>Solvent Properties</b></p> <p>Water</p> <p>Boiling Point (°C) : 100</p> <p>Density (g/cm<sup>3</sup>) : 0.998</p> <p>Mol Wt. (amu) : 18.02</p>	<p><b>Injection volume (µL)</b></p> <p>1.00</p> <p><b>Inlet Temperature (°C)</b></p> <p>250</p> <p><b>Inlet Pressure (gauge)</b></p> <p>14.000</p> <p><input type="radio"/> kPa <input checked="" type="radio"/> psi <input type="radio"/> bar</p>	<p><b>Estimated Volume</b></p> <p><b>1217 µL</b></p> <p><b>% Capacity</b></p> <p><b>143%</b></p> <p></p> <p><b>Solvents</b></p> <p><input type="button" value="Add"/> <input type="button" value="Remove"/> <input type="button" value="Defaults"/></p> <p><b>Liners</b></p> <p><input type="button" value="Add"/> <input type="button" value="Remove"/> <input type="button" value="Defaults"/></p>
--	--	---



# Conclusions

Try to understand the sample as much as you can.

Residues, concentrations, solvent expansion

Packed columns are used with a PP inlet only

MMI or PTV for large volume injections (trace analysis)

MMI, PTV or COC for Labile compounds, or high bp compounds

SSL inlet is the most common

MMI is a combination of the SSL and PTV gives more flexibility

does have issues with cleaning – Intuvo addresses that

# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:

Option 1 for GC/GCMS Columns and Supplies

Option 2 for LC/LCMS Columns and Supplies

Option 3 for Sample Preparation, Filtration and QuEChERS

Option 4 for Spectroscopy Supplies

Available in the USA 8-5 all time zones



[gc-column-support@Agilent.com](mailto:gc-column-support@Agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)