



Pooled CRISPR guide RNA libraries for functional genomics screening: Do you know what's in your library?


Peter Sheffield

R&D Scientific Program Director

Agilent Technologies

INTRODUCTION

CRISPR: A Programmable Tool for Genome Engineering



- CRISPR = Clustered Regularly Interspaced Short Palindromic Repeats
- A multi-component system consisting of an RNA component and a CRISPR-associated protein (Cas)
 - RNA (CRISPR) sequence defines target location
 - CRISPR-associated protein (Cas) determines functional modification made at target location



INTRODUCTION


CRISPR: A Programmable Tool for Genome Engineering

Part 1 - CRISPR-associated protein (Cas)


- CRISPR-associated protein (Cas) determines functional modification made at target location
 - wtCas9 - induces double stranded breaks - gene function knock-out
 - dCas9 - 'Dead' becomes a DNA targeting protein – CRISPRa/i
 - nCas9 - 'Nicking' induces single stranded nicks at target location, either top or bottom strand depending on Mutation

Part 2 - RNA (CRISPR) sequence defines target location

- Two RNAs
- Single guide RNA



Two guides



Single guide




Single Guide CRISPR RNAs

Can any sequence be used as a CRISPR?


Yes - as long as it is immediately upstream of a PAM site!

What is a PAM site?

PAM = Protospacer Adjacent Motif is a 2-6 base pair DNA sequence immediately adjacent to the DNA sequence targeted by the Cas9 protein, for *S. pyogenes* PAM = NGG



Single guide



Principles of Designing a Genome Wide CRISPR library

It's a numbers game!

Human Genome = 3×10^9 bp

Cas9 PAM = NGG = once every ~40 bp = ~150 Million PAM sites (2 strands)

Too many to Screen!

How do you filter the CRISPRs?

- Target certain areas of the genome
 - Exonic regions
 - Intronic regions
 - Regulatory regions
 - Gene clusters
 - Gene Families

Good CRISPR versus Bad CRISPR

- Not all CRISPRs are equal
 - Algorithms to score
 - Doench^{1, 2}
 - Zhang³ (Broad)

1. "Rational design of highly active sgRNAs for CRISPR-Cas9-mediated gene inactivation" John G Doench et al. Nature Biotechnology 32, 1262–1267 (2014)
2. "Optimized sgRNA design to maximize activity and minimize off-target effects of CRISPR-Cas9" John G Doench et al. Nature Biotechnology 34, 184–191 (2016)
3. "Genome-Scale CRISPR-Cas9 Knockout Screening in Human Cells" Ophir Shalem, Neville E. Sanjana, Ella Hartenian, Xi Shi, David A. Scott, Tarjei Mikkelsen, Dirk Heck, Benjamin L. Ebert, David E. Root, John G. Doench, and Feng Zhang. Science. 2014 Jan 3; 343(6166): 84–87.



Example: GeCKO

Genome-scale CRISPR-Cas9 Knock-Out

Designed by The Broad Institute (Zhang)

Human Genome – 150 Million CRISPRs

Criteria	Targets	GeCKO Library Design				
		n	CRISPRs	A	n	B
Coding Regions	~19,050	6	~114,300	~62,000	3	~53,000
miRNA	~1,860	4	~7,200	~3,600	3	~3,600
Non-targeting Controls	1,000	1	1,000	1,000	1	1000
Total	~21,800		~123,000	~66,000		~57,000

Pooled CRISPR Guide Libraries

Genome Wide Modification Experiments

Using a pooled CRISPR library, containing many thousands of CRISPR RNAs, multiple CRISPR/Cas reactions can be performed simultaneously on a genome wide scale.

Gene Knock-Out:

- Introduces non-sense mutations into the reading frame of any gene in the genome via Non-homologous end joining (NHEJ), a repair pathway that joins two broken ends together.
- NHEJ does not require a homologous template for repair and typically leads to the introduction of small insertions/deletions at the break site often creating frame-shifts that knockout gene function.

CRISPRa/i:

- a = gene Activation
- i = gene Inhibition
- Targets CRISPR sites within the regulatory regions of genes and using dCas9 fusions to recruit transcriptional modification machinery that either activates or inhibits gene expression

All three methods especially useful for identifying biological relevant pathways, new drug targets for disease models

Designing a Genome Wide CRISPR library: Factors to consider

CRISPR Experiment

- Gene Knock-Out
 - Preferentially target the 5' exons over the 3' exons.
 - Increases chances of creating a Null mutant.
 - Target just the genes you are interested in knocking out
- CRISPRa/i
 - Target the regulatory regions upstream of coding regions.
 - Target only the gene families of interested

Determining n CRISPR/target

- Increasing the number of CRISPRs/target builds redundancy into the design which increase the likelihood of having an affect on your target.
- Increasing the number of CRISPRs/target can also help in the analysis of the Next Gen Sequencing data by adding mathematical confidence that the changes seen in the library distribution pre- and post- data sets are statistically relevant.

Designing a Genome Wide CRISPR library: Factors to consider

Minimize bias in Design

- No CRISPR duplicates
- Ideally all CRISPR sequences in design should be unique (Ratio 1=1)
- The presence of pseudogenes and splice variants can create instances where some CRISPR sequences are duplicated. These duplications introduce bias into a design (Ratio 1≠1)
 - Either COMPLETELY remove all CRISPRs that are Duplicated and replace with a different CRISPR that is not duplicated
or
 - Eliminate all duplicates but leave one copy in design
- Both methods maintain CRISPR ratio of 1=1.
- The first method could result in the replacement of a good CRISPR with a less desirably one.
- However, the second method means you have CRISPR sequences that target multiple loci in the genome. This could lead to analysis complications downstream of your experimental data.



How do you get all the CRISPR machinery into the target cell?


For a CRISPR/Cas event to occur you need the following to be present in the cell at the same time

- CRISPR guide RNA
- Cas9 protein

Delivery Methods


- Pooled CRISPR libraries are best delivered via a DNA vector delivery method, e.g. Lentivirus
 - Lentiviral vectors are amenable to molecular biology manipulation
 - Lentiviral particles transfect many different cell lines, non-dividing cell
 - Stably integrate into the genome
 - Perfect for expression of small sgRNA's
- Cas9 protein can be delivered via multiple methods
 - Direct transfection of active protein
 - Via a DNA expression vector (viral or non-viral)
 - Direct transfection of mRNA
- Delivery of the CRISPR sgRNAs and the Cas9 protein can be done simultaneously or sequentially
- Cas9 modifications can be made in parallel on a genome-wide scale

Synthesizing a CRISPR Oligo library:



Oligo libraries:

- Agilent can print up to 240,000 user-definable sequences/chip
- Sequences up to 230 nt
- Pooled into a single tube/library




- Libraries are synthesized as ssDNA
- Libraries usually need to be amplified by PCR
- PCR acts as an enrichment step by amplifying the full length oligos
- CRISPR libraries are an ideal size for Oligo synthesis (<140 nt)



Library Construction Considerations

Example case: GeCKo library for genome-wide Cas9 mediated gene knockout

CRISPR library




GeCKo v2 libraries

species	human	mouse
genes targeted	19,050	20,611
targeting constructs/gene	6	6
miRNA targeted	1,864	1,175
targeting constructs/miRNA	4	4
control (nontargeting) sgRNAs	1,000	1,000
total sgRNA constructs	123,411	130,209
set A	66,172	
set B	57,239	
non-redundant set A	64,580	
non-redundant set B	56,869	
total non-redundant	121,449	

Sanjana et al., *Nature Methods* **11**: 783-784 (2014)

Why break a library into subsets?


Plasmid library sizes as function of OLS Complexity




Quality of Oligo libraries: representation of library members

Representation of guides analyzed using NGS

box whisker plots of guide representation in oligo libraries



Synthesis Protocol 1



Synthesis Protocol 2

of the same CRISPR design



Agilent Technologies


For Research Use Only. Not for use in diagnostic procedures. PR7000-0707

Library Construction Considerations




- Each step in the process can introduce bias into a library
- Each method needs to be as *neutral* as possible
- Each step needs to be controlled and monitored
- The Cloning and Library Expansion steps have the greatest potential to introduce bias into a library

Cloning an CRISPR oligo library into a plasmid vector




Constructing plasmid libraries from amplified OLS libraries

Agilent's SureVector Sure Guide Cloning Kit and workflow




$$\text{Ratio } 95^{\text{th}}/5^{\text{th}} \text{ Percentile} = \frac{\text{Frequency of } 95^{\text{th}} \text{ Percentile}}{\text{Frequency of } 5^{\text{th}} \text{ Percentile}}$$







The lower the 95th/5th Percentile ratio the better the library



Methods to Generate highest quality plasmid libraries from OLS designs


- Start with highest quality OLS libraries
- combine (+) strand and (-) strand designs
- control growth conditions




		pSGL-05J	pSGL-05/06J
strand		P1	Combined
90% range		4.632	 2.800
80% range		3.000	 2.200
60% range		2.069	 1.659
reads/guide		47.6	55.4



Quality of Oligo libraries: representation of library members



Synthesis Protocol 1
of the same CRISPR design




Synthesis Protocol 2










Methods to Generate highest quality plasmid libraries from OLS designs

- Start with highest quality OLS libraries
- combine strands generated by different synthesis protocols (P1 and P2)
- control growth conditions



Uncontrolled growth of bacterial cultures containing CRISPR libraries introduces bias into the library




		pSGL-05J	pSGL-05/06J
strand		P1	Combined
90% range		4.632	 2.800
80% range		3.000	 2.200
60% range		2.069	 1.659
reads/guide		47.6	55.4



So what makes a Good CRISPR Library

Three critical aspects
that make a good pooled
CRISPR Library

- Fidelity
- Customization
- Representation




Fidelity:

Are my CRISPR Sequences correct?



- A high fidelity oligo library



CRISPR Libraries

- **More correct constructs**
- 92% vs. 77%
- **Fewer avg. errors/kb**
- 4 vs. 12

- **Fewer errors** per kbase pair improves library quality
- **Reduces screening time** and false negatives



Customization: Any CRISPR Sequence you need




Easy Customization



End-user has complete control

- Any sequence
- Any species
- Any application




Agilent provides CRISPR Libraries

- Ready-to-clone
- Ready-to-amplify libraries for total flexibility




Representation:

Are all my guides present?



Improved Representation



- SureVector CRISPR Cloning Kit
- Fewer missed guides
- Less under-represented population

Libraries	Missed Guides	90 th /10 th percentile	95 th /5 th percentile	99.5 th /0.5 th percentile
pSGL-007J	1	2.32	3.04	6.72
pSGL-006J	1	2.47	3.38	11.06
pSGL-005J	1	2.38	3.19	9.83
pSGL-128J-dc	1	1.99	2.64	8.30
GeCKO (Broad)	?	8.73	16.00	NA
GeCKO (Competitor)	39	5.29	9.83	68.40
GeCKO (Competitor) expanded	204	6.00	11.95	333.00

- Plasmid libraries with 90/10 ratios < 3
- Find all the hits in your screen, not just the over-represented ones



Now You have a high Quality CRISPR Library, What Now?


Genome Wide Modification Experiments

Using a pooled CRISPR library, containing many thousands of CRISPR RNAs, multiple CRISPR/Cas reactions can be performed simultaneously on a genome wide scale.


- Gene Knock-Out
- CRISPRa/i
- Experimental Design Decisions
 - Determine Cell line
 - Number of Replicates
 - Selection protocol
 - Next Gen Sequencing protocol
 - Data Analysis pipeline



Functional screening using CRISPR/CAS




NGS of CRISPR Oligo libraries, plasmid libraries and transduced cell lines




This section is present in all CRISPR library entities Oligo library, Plasmid and transduced cell line and can be amplified using universal primers

This region of the library can be isolated via PCR

The resultant amplicon can then be used in your preferred NGS workflow.



Analysis of Next Gen Sequence Data from a CRISPR/CAS Experiment



Statistical analysis packages

- [RIGER](#)
- [HiT Select](#)
Diaz et al., *NAR* (2014)
- [MAGeCK](#)
Li et al., *Genome Biology* **15**:554 (2014)
- [casTLE](#)
Morges et al., *Nature Biotech* ePub April 11th (2016)

No Cas9 the library should not show too much bias towards any one CRISPR

All libraries have an inherent amount of background variation (noise / bias)


With Cas9 the library should start to show bias towards certain CRISPRs due to selection

However, a bad library will have more noise/bias which makes it harder to identify the significant CRISPRs. This requires your screen to encompass a greater number of cells



Agilent SureGuide Products


Ready-to-Package



Catalog gRNA Libraries

- Plasmid Library
- GeCKOv2
- Human and Mouse
- Cloned into lentivirus vector with hU6 promoter

Ready-to-Clone



Custom Human & Mouse gRNA Libraries

- Pre-amplified oligo library
- User defined subset or designed
- Mammalian systems
- Compatible with SureVector cloning

Ready-to-Amplify




Custom User Defined Libraries

- Unamplified oligo pool
- Any species, any cloning method
- Entirely custom by user design



Agilent's Commitment to Genome Engineering

- Agilent has expertise in manufacturing market-leading high fidelity, complex, ultra-long **oligonucleotide libraries**, making market leading microarrays and SureSelect NGS solutions.
- Agilent offers a portfolio of **specialized cloning Kits (SureVector)**, **competent cells**, and **high fidelity PCR enzymes** as well as market-leading **QuikChange site directed mutagenesis** technology – developed with expertise from Stratagene (now part of Agilent).
- **Broad portfolio of analytical separation and measurement tools** (e.g., LC-/GC-MS, NMR), **liquid handling automation** and **informatics solutions**.



Leading Provider of Measurement Tools & Solutions



Learn more: www.agilent.com/genomics/SureGuide

Acknowledgements:

La Jolla R&D

Peter Sheffield
Carsten Carstens
Katie Felts
Sarah Johns
Vivian Zhang



Santa Clara Marketing

Ben Borgo
Caroline Tsou

