


Your Guide for Accurate Scoring in Head and Neck Squamous Cell Carcinoma (HNSCC) Using PD-L1 IHC 22C3 pharmDx (SK006)

Use this quick scoring guide as a reference when evaluating HNSCC specimens for PD-L1 expression using PD-L1 IHC 22C3 pharmDx.

For more information on Combined Positive Score (CPS) calculation, review the HNSCC Interpretation Manual.

Steps for scoring



Note: PD-L1 expression level CPS ≥ 20 may be of interest to treating physician but does not determine eligibility for first-line treatment with KEYTRUDA® (pembrolizumab) as a single agent.

Definition of CPS and PD-L1 staining cells

CPS is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100.

$$\text{CPS} = \frac{\text{\# PD-L1 staining cells (tumor cells, lymphocytes, macrophages)}}{\text{Total \# of viable tumor cells}} \times 100$$

Note: CPS is reported as a whole number. Although the result of the calculation can exceed 100, the maximum score is defined as CPS 100.

By definition, **PD-L1 staining cells** in HNSCC are:

- **Viable tumor cells** with perceptible and convincing partial or complete linear membrane staining (at any intensity) that is perceived distinct from cytoplasmic staining
- **Lymphocytes and macrophages** (mononuclear inflammatory cells, MICs) within the tumor nests and/or adjacent supporting stroma with membrane and/or cytoplasmic staining (at any intensity). MICs must be directly associated with the response against the tumor




Figure 1 (left): HNSCC specimen stained with PD-L1 antibody exhibiting a CPS of 0.
Figure 2 (center): HNSCC specimen stained with PD-L1 antibody exhibiting a CPS of 6, however any numerical CPS between 4–8 could be assigned to this image.
Figure 3 (right): HNSCC specimen stained with PD-L1 antibody exhibiting a CPS of 100.

Intended Use

For in vitro diagnostic use.

PD-L1 IHC 22C3 pharmDx is a qualitative immunohistochemical assay using monoclonal mouse anti-PD-L1, Clone 22C3 intended for use in the detection of PD-L1 protein in formalin-fixed, paraffin-embedded (FFPE) head and neck squamous cell carcinoma (HNSCC) tissue using EnVision FLEX visualization system on Autostainer Link 48.

Head and Neck Squamous Cell Carcinoma (HNSCC)

PD-L1 protein expression in HNSCC is determined by using Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100. The specimen should be considered to have PD-L1 expression if CPS ≥ 1 .

PD-L1 IHC 22C3 pharmDx is indicated as an aid in identifying HNSCC patients for treatment with KEYTRUDA® (pembrolizumab). See the KEYTRUDA® product label for specific clinical circumstances guiding PD-L1 testing.

For descriptions of the intended use in other indications, please refer to the current version of the Instructions for Use (IFU) for PD-L1 IHC 22C3 pharmDx, Code SK006.

Tissue samples supplied by Asterand Bioscience.

KEYTRUDA is a registered trademark of Merck Sharp & Dohme Corp., a subsidiary of Merck & Co., Inc.

HNSCC CPS numerator

Tissue Elements	Included	Excluded
Tumor Cells	Convincing partial or complete linear membrane staining (at any intensity) of viable invasive tumor cells	<ul style="list-style-type: none"> Non-staining tumor cells Tumor cells with only cytoplasmic staining
Immune Cells	Membrane and/or cytoplasmic* staining (at any intensity) of mononuclear inflammatory cells (MICs) within tumor nests and adjacent supporting stroma [†] : <ul style="list-style-type: none"> Lymphocytes (including lymphocyte aggregates) Macrophages[‡] Only MICs directly associated with the response to the tumor are scored	<ul style="list-style-type: none"> Non-staining MICs MICs (including lymphoid aggregates) associated with ulcers or other inflammatory processes MICs associated with carcinoma in situ MICs associated with benign structures Neutrophils, eosinophils, and plasma cells
Other Cells	Not included	<ul style="list-style-type: none"> Carcinoma in situ Benign cells Stromal cells (including fibroblasts) Necrotic cells and/or cellular debris

HNSCC CPS denominator

Included	Excluded
All viable invasive tumor cells	Any necrotic or non-viable tumor cells
Not included	All immune cells of any type
Not included	<ul style="list-style-type: none"> Carcinoma in situ Benign cells Stromal cells (including fibroblasts) Necrotic cells and/or cellular debris

* In MICs, membrane and cytoplasmic staining are often indistinguishable due to high nuclear to cytoplasmic ratio. Therefore, membrane and/or cytoplasmic staining of MICs are included in the score; [†] Adjacent MICs are defined as being within the same 20x field as the tumor. However, MICs that are NOT directly associated with the response to the tumor should be excluded; [‡] Macrophages and histiocytes are considered the same cells

Partial and complete linear membrane staining




Figure 4: PD-L1 staining of partial (A) and complete (B) linear membrane staining of tumor cells (20x magnification).

Tumor-associated immune cells




Figure 5: When positioning the edge of the tumor mass in the approximate center of a 20x field, PD-L1 staining MICs that are present within the same field should be included in the numerator (20x magnification).




Figure 6: Simulation of a 20x microscope field showing tumor surrounded by PD-L1 staining tumor-associated MICs that should be included in the numerator.

For countries outside of the United States, see the local KEYTRUDA® (pembrolizumab) product label for approved indications and expression cutoff values to guide therapy.

D54067/02

This information is subject to change without notice.

© Agilent Technologies, Inc. 2022
 Printed in the USA, February 11, 2022
 29315 2022JAN27