


Your Guide for Accurate Scoring in Esophageal Squamous Cell Carcinoma (ESCC) Using PD-L1 IHC 22C3 pharmDx (SK006)

Use this quick scoring guide as a reference when evaluating ESCC specimens for PD-L1 expression using PD-L1 IHC 22C3 pharmDx.

For more information on Combined Positive Score (CPS) calculation, review the ESCC Interpretation Manual.

Steps for scoring



Determine specimen adequacy

Verify that the specimen has ≥ 100 viable tumor cells

Evaluate controls

Ensure that Control Cell Line Slide and lab-supplied and patient tissue controls demonstrate acceptable staining

Evaluate PD-L1 staining

Estimate the number of PD-L1 staining cells (CPS numerator) and the total number of viable tumor cells (CPS denominator)

Calculate CPS to determine PD-L1 expression level

Report numerical CPS and PD-L1 expression level:
CPS < 10 or
CPS ≥ 10

Definition of CPS and PD-L1 staining cells

CPS is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100.

$$\text{CPS} = \frac{\text{\# PD-L1 staining cells (tumor cells, lymphocytes, macrophages)}}{\text{Total \# of viable tumor cells}} \times 100$$

Note: CPS is reported as a whole number. Although the result of the calculation can exceed 100, the maximum score is defined as CPS 100.

By definition, **PD-L1 staining cells** in ESCC are:

- **Viable tumor cells** with perceptible and convincing partial or complete linear membrane staining (at any intensity) that is perceived distinct from cytoplasmic staining
- **Lymphocytes and macrophages** (mononuclear inflammatory cells, MICs) within the tumor nests and/or adjacent supporting stroma with membrane and/or cytoplasmic staining (at any intensity). MICs must be directly associated with the response against the tumor




Figure 1: ESCC specimen stained with PD-L1 antibody exhibiting a CPS of 4, however any numerical CPS between 2–6 could be assigned to this image (20 \times magnification).




Figure 2: ESCC specimen stained with PD-L1 antibody exhibiting a CPS of 45, however any numerical CPS between 40–50 could be assigned to this image (20 \times magnification).

Intended Use

For in vitro diagnostic use.

PD-L1 IHC 22C3 pharmDx is a qualitative immunohistochemical assay using monoclonal mouse anti-PD-L1, Clone 22C3 intended for use in the detection of PD-L1 protein in formalin-fixed, paraffin-embedded (FFPE) esophageal squamous cell carcinoma (ESCC) tissue using EnVision FLEX visualization system on Autostainer Link 48.

Esophageal Squamous Cell Carcinoma (ESCC)

PD-L1 protein expression in ESCC is determined by using Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100. The specimen should be considered to have PD-L1 expression if CPS ≥ 10 .

PD-L1 IHC 22C3 pharmDx is indicated as an aid in identifying ESCC patients for treatment with KEYTRUDA[®] (pembrolizumab). See the KEYTRUDA[®] product label for specific clinical circumstances guiding PD-L1 testing. For descriptions of the intended use in other indications, please refer to the current version of the Instructions for Use (IFU) for PD-L1 IHC 22C3 pharmDx, Code SK006.

CPS numerator

Tissue Elements	Included	Excluded
Tumor Cells	Convincing partial or complete linear membrane staining (at any intensity) of viable invasive tumor cells	<ul style="list-style-type: none"> – Non-staining tumor cells – Tumor cells with only cytoplasmic staining – Non-invasive neoplasia (including carcinoma in situ)
Immune Cells	Membrane and/or cytoplasmic* staining (at any intensity) of mononuclear inflammatory cells (MICs) within tumor nests and adjacent supporting stroma [†] : <ul style="list-style-type: none"> – Lymphocytes (including lymphocyte aggregates) – Macrophages[‡] Only MICs directly associated with the response to the tumor are scored	<ul style="list-style-type: none"> – Non-staining MICs – MICs associated with non-invasive neoplasia (including carcinoma in situ) – MICs associated with benign structures – MICs (including lymphoid aggregates) not directly associated with the response to the tumor – Neutrophils, eosinophils, and plasma cells
Other Cells	Not included	<ul style="list-style-type: none"> – Benign epithelial cells – Stromal cells (including fibroblasts) – Necrotic cells and/or cellular debris

CPS denominator

Included	Excluded
All viable invasive tumor cells	<ul style="list-style-type: none"> – Non-viable tumor cells – Non-invasive neoplasia (including carcinoma in situ)
Not included	All immune cells
Not included	<ul style="list-style-type: none"> – Benign cells – Stromal cells (including fibroblasts) – Necrotic cells and/or cellular debris

* In MICs, membrane and cytoplasmic staining are often indistinguishable due to high nuclear to cytoplasmic ratio. Therefore, membrane and/or cytoplasmic staining of MICs are included in the score; [†] Adjacent MICs are defined as being within the same 20x field as the tumor. However, MICs that are NOT directly associated with the response against the tumor should be excluded; [‡] Macrophages and histiocytes are considered the same cells

Partial and complete linear membrane staining




Figure 3: PD-L1 staining of partial (A) and complete (B) linear membrane staining of tumor cells (20x magnification).

Tumor-associated immune cells




Figure 4: When positioning the edge of a tumor mass in the approximate center of a 20x field, PD-L1 staining MICs (arrows) that are present within the same field should be included in the numerator.




Figure 5: Simulation of a 20x microscope field showing tumor surrounded by PD-L1 staining tumor-associated MICs that should be included in the numerator.

For countries outside of the United States, see the local KEYTRUDA® (pembrolizumab) product label for approved indications and expression cutoff values to guide therapy.

D54360/02

This information is subject to change without notice.

© Agilent Technologies, Inc. 2022
 Printed in the USA, February 08, 2022
 29319 2022JAN31

Data and biospecimens used in this project were provided by SageBio LLC, Sharon, MA, USA and by US Biolab, Rockville, MD, with appropriate ethics approval and through Trans-Hit Biomarkers Inc.