

Agilent DNA Fish ID solution Discrimination of Sturgeon and Related Species by PCR-RFLP Using the 2100 Bioanalyzer System

FOOD AUTHENTICITY



A quick, robust and easy to use protocol to identify the species based on a well accepted PCR-RFLP method allows analysis from sample (meat, fin clip or roe) to result in one working day and yields good discrimination results

The global demand for seafood has grown considerably. Limitations of the resources and the potential for increased profits lead to the problem of substitution and mislabeling for a substantial part of the market. In order to monitor shipments for customs purposes as well as in supply chain management and to protect the consumer, efficient, and cost-effective tests to identify the species are needed. In addition, regulations to protect endangered species and to fight illegal, unregulated and unreported fishing drive the need for authenticity testing.

DNA based testing methods allow sensitive detection and identification from almost all but the most heavily processed food samples.

The Agilent DNA Fish ID ensemble provides all reagents and enzymes to complete a PCR-RFLP method for a *Cytb* PCR target sequence and analysis of restriction fragment patterns on the Agilent 2100 Bioanalyzer to identify fish species based on this mitochondrial sequence. The availability of commercial screening solutions allows for more reliable and robust test results through well-matched components and facilitates testing for screening purposes by the use of mastermix formulations and streamlined protocols. The Agilent DNA Fish ID solution was evaluated for the purpose of identifying the species from fish eggs in caviar shipments. To obtain the best possible results using roe and to accommodate the high homology between sturgeon species, the protocol was modified for enhanced identification.

Species	Common name
<i>Acipenser gueldenstaedtii</i>	Danube sturgeon
<i>Acipenser medirostris</i>	Green sturgeon
<i>Acipenser nudiiventris</i>	Fringebarbel sturgeon
<i>Acipenser ruthenus</i>	Sterlet sturgeon
<i>Acipenser schrenckii</i>	Amur sturgeon
<i>Acipenser stellatus</i>	Starry sturgeon
<i>Acipenser transmontanus</i>	White sturgeon
<i>Huso dauricus</i>	Kaluga
<i>Huso huso</i>	Beluga
<i>Polyodon spathula</i>	Mississippi paddlefish
<i>Oncorhynchus keta</i>	Chum salmon

Figure 1. Species tested in the study.

Key Benefits

- Agilent 2100 Bioanalyzer has the widest range of assays available on the market covering applications for RNA, DNA or protein
- Fast, reproducible and robust method replacing gel electrophoresis
- Agilent method using convenient mastermix format reagents will have you identifying your fish within hours, no fish standard or guesswork required.
- Flexible modular system: In-house accredited methods can replace individual building blocks of the Fish ID solution
- Fifty fish species in the Agilent library, and you can easily add your own samples to the library
- Overlay to the Expert software allows rapid identification of restriction fragments and close matches.
- Agilent method is now optimized for testing roe (caviar)
- Replacing one of the supplied enzymes leads to improved discrimination abilities for sturgeon species



Discrimination of Sturgeon and Related Species by PCR-RFLP Using the 2100 Bioanalyzer System

Table 2a: Grouping of Samples According to Restriction Patterns Observed. Samples Were Grouped According to the Number of Fragments Observed With Each Restriction Enzyme. Subgrouping Was Added if There Was a Clear Difference (At Least Three Standard Deviations) in the Size of At Least One Fragment

Group	Species	Csp6 I		Dde I			Hae III		Nla III		
A	<i>Acipenser gueldenstaedtii</i>	325	122	462			180	148	291	193	
	<i>Acipenser medirostris</i>	325	122	458			180	150	293	192	
	<i>Huso dauricus</i>	326	122	463			180	151	294	193	
B ₁	<i>Acipenser transmontanus</i>	446		455			179	149	293	192	
	<i>Acipenser schrenckii</i>	445		457			180	149	292	192	
B ₂	<i>Acipenser stellatus</i>	489		458			180	149	291	193	
C	<i>Huso huso</i>	323	120	455			179	150	77	291	191
D	<i>Polyodon spathula</i>	481		480			323	150	291	165	
E	<i>Oncorhynchus keta</i>	389	86	359	352	121	433		280	192	
F ₁	<i>Acipenser nudiiventris</i>	325	93	462			180	149	293	133	52
F ₂	<i>Acipenser ruthenus</i>	349	94	455			178	149	292	134	62
G	<i>Salmo salar</i>	397	85	362	356	121	330	107	46	460	

Table 2b. Standard Deviations of Fragment Sizes. Variability of Fragment Sizing is Between 0.5% to 2% of the Average Size

Group	Species	Csp6 I		Dde I			Hae III		Nla III		
A	<i>Acipenser gueldenstaedtii</i>	1.2	1.2	2.5			0.6	0.6	2.3	2.6	
	<i>Acipenser medirostris</i>	1.4	0.9	3.7			1.1	1.3	2.4	1.5	
	<i>Huso dauricus</i>	1.4	2.1	2.8			2.8	2.1	2.1	2.1	
B ₁	<i>Acipenser transmontanus</i>	3.8		3.3			1.3	1.3	1.2	0.9	
	<i>Acipenser schrenckii</i>	2.6		3.5			1.0	0.6	0.6	0.6	
B ₂	<i>Acipenser stellatus</i>	5.3		4.5			1.7	1.4	1.0	1.7	
C	<i>Huso huso</i>	2.1	1.2	5.8			1.0	0.6	0.5	2.9	2.6
D	<i>Polyodon spathula</i>	7.7		8.9			5.4	2.2	3.2	2.4	
E	<i>Oncorhynchus keta</i>	2.2	0.6	2.6	2.6	1.0	6.9		4.1	1.7	
F ₁	<i>Acipenser nudiiventris</i>	2.0	0.6	3.2			0.6	0.6	2.0	0.6	0.5
F ₂	<i>Acipenser ruthenus</i>	1.8	1.5	2.6			1.0	1.0	0.5	0.6	0.6
G	<i>Salmo salar</i>	2.5	1.1	2.7	2.6	1.1	2.2	0.8	1.3	4.3	




Fig. 1. Sequence analysis of Cytb target region in sturgeon species. The target region for the primers supplied with the kit was analyzed using CLC Sequence Viewer (CLC bio, Denmark). The picture shows the cutting sites for Csp6 I, Dde I, Hae III and Nla III for *Acipenser dabryanus* (derived from NC_005451), *Acipenser gueldenstaedtii* (derived from NC_012576), *Acipenser sinensis* (derived from NC_012646), *Acipenser stellatus* (derived from NC_005795), *Acipenser transmontanus* (derived from NC_004743) and *Huso huso* (derived from AY_442351).

Learn more:
www.agilent.com/chem

Email:
info_agilent@agilent.com

Find a customer center in your country:
www.agilent.com/chem/contactus

This information is subject to change without notice

© Agilent Technologies, Inc., 2011
Printed in USA, June 30, 2011
5990-8455EN



Agilent Technologies