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
**Analysis Mode**  
Select appropriate Analysis Mode. Also make the following changes to radio button selections.

**Conventional Analysis**  
 Spectrum/TRA  
 Spectrum: Conventional Mass Spectra are acquired.  
 TRA: Conventional Time Resolved Data are acquired.

**Nanoparticle**  
 Single Element / Isotope  
 Single Particle Analysis using fast TRA mode.  
 Multi-Element Screening  
 Fast Time Program Analysis for multi-element nanoparticle screening.  
 2 Elements / Isotopes  
 Average element/isotope ratios of multiple particles.  
 FFF ICP-MS  
 Particles acquired using conventional TRA mode after separation.

**Single Particle Analysis Configuration**  
 Set parameters for Single Particle Analysis.

Sample Pump Tube ID: 1.02 mm  
 Sample Inlet Flow: 0.346 ml/min  
 Response Factor Calibration Solution:  
 Ionic Standard Concentration at 197 amu: 1.000 ppb  
 Reference Material: NIST RM 8012  
 Reference Element Mass: 197 amu  
 Mean Reference Particle Diameter: 27 nm  
 Reference Material Density: 19.32 g/cm<sup>3</sup>  
 Concentration of Reference Material: 5.0 ng/l  
 Unknown Sample:  
 Target Element Mass: 197 amu  
 Analyte Mass Fraction: 1.000  
 Particle Density: 19.32 g/cm<sup>3</sup>

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Offline ICP-MS MassHunter - nanotes5.b

File Edit View Instrument Hardware Startup Batch Queue Tools Help

Hardware Plasma Tune Batch Queue Data Analysis Report

Batch - nanotes5.b

Save Batch Add to Queue Validate Method Import Sample List Auto Sampler Nebulizer Pump Speed

Acq Method	Data Analysis Method	Sample List	Skip	Sample Type	Sample Name	Comment	Vial#	Replicates	Dilution
			<input type="checkbox"/>	IonicBlk	Ionic Blank		1101		
			<input type="checkbox"/>	IonicStd	Ionic Standard	10 ppt Ag	1102		
			<input type="checkbox"/>	RM	NIST RM 8012	30nm RM	1103		
			<input type="checkbox"/>	Sample	Unknown 1		1104		
			<input type="checkbox"/>	Sample	Unknown 2		1105		
			<input type="checkbox"/>	Sample	Unknown 3		1106		
			<input type="checkbox"/>	Sample	Unknown 4		1107		
			<input type="checkbox"/>	Sample	Unknown 5		1201		
			<input type="checkbox"/>	Sample	NIST RM 8012		1202		
			<input type="checkbox"/>						

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© 安捷伦科技（中国）有限公司，2017