

# Quantitative Analysis of 1,25-Dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> by LC-MS/MS Utilizing Ion Funnel Technology



Peter Christensen<sup>1</sup>, Vicky Starkie<sup>1</sup>, Kevin McCann<sup>2</sup>, <sup>1</sup>Agilent Technologies Inc. Cheadle, UK, <sup>2</sup>Agilent Technologies Inc. Santa Clara, CA, USA

## Introduction

A highly sensitive, and selective LC-MS/MS method for determination of 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> is a powerful tool for clinical researchers. While 25 hydroxy vitamin D is found in the ng/ml concentration range, 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> are typically found in the low pg/ml range, making quantitative analysis challenging except when employing highly sensitive analytical techniques. Additionally, extraction is a critical step for this analysis as removal of interfering analytes is required to quantify at the low pg/ml concentration range. Where previously published work (Casetta et al., 2010) demonstrates good sensitivity, such approaches require a complex 2D LC set-up, with post column infusion. The work presented in this poster illustrates the quantitative analysis of 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> using the Agilent 1260 UHPLC & 6490 QQQ with Ion Funnel technology coupled with ImmunoTube LC-MS/MS Kit (an extraction kit from ImmunoDiagnostic) for the extraction of 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> from plasma samples.




Figure 1. Structures and molecular weights of 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub>

## Sample Preparation

ImmunoTubes were spun down to ensure all the suspension was forced to the bottom of the tube. 500µl of calibrator / sample / control was added followed by 10µL of IS and mixed gently. ImmunoTubes were then mixed in a spiral rotator for 1hr at RT. The closed ImmunoTubes were then placed in a micro tube and centrifuged for 1min at 550 x g. Subsequently, the cover and the outlet of the ImmunoTubes were removed. The Immunotubes were then placed back into the micro tubes for centrifugation a further 2min at 550 x g. The waste collected in the micro tubes was discarded.

## Sample Preparation

Continued...  
500uL of WASHSOL was added to the ImmunoTubes and centrifuged for 2 min at 550 x g. This step was repeated twice. Each micro tube was replaced by a glass vial and 250µL of ELUREAG were added to each ImmunoTube which was centrifuged for 2min at 550 x g. The recovered eluent was evaporated under N<sub>2</sub> at 37°C. Samples were reconstituted with 165µL of activated Solution A prior to analysis.

WASHSOL and ELUREA are solutions from the ImmunoDiagnostic extraction kit.

**LC Method**  
An Agilent 1260 HPLC series binary pump with 56 vial sample tray, sampler with thermostat, temperature-controlled column compartment, 2 position/6 ports switching valve, was used.

Column : Zorbax Eclipse Plus 2.1x100 mm 1.8µm  
Column temperature: 50 °C  
Injection volume: 100 µL  
Autosampler temperature: 4 °C  
Needle wash: 3:1 MeOH:H<sub>2</sub>O, 10 seconds

**Mobile Phase**  
A: ImmunoDiagnostic Mobile Phase A  
B: ImmunoDiagnostic Mobile Phase B

Gradient	Flow	% Solvent B
0.00	0.3	0
6.00	0.3	100
6.50	0.3	100
6.51	0.3	0
8.00	0.3	0

Table 1. LC conditions

Compound	Prec Ion	Prod Ion	CE (V)
1,25(OH) <sub>2</sub> Vitamin D <sub>2</sub>	411.1	150.7	20
1,25(OH) <sub>2</sub> Vitamin D <sub>2</sub>	411.1	132.9	18
1,25(OH) <sub>2</sub> Vitamin D <sub>3</sub>	399.1	150.9	12
1,25(OH) <sub>2</sub> Vitamin D <sub>3</sub>	399.1	134.9	12
1,25(OH) <sub>2</sub> Vitamin D <sub>3</sub> - d <sub>6</sub>	405.1	150.6	12
1,25(OH) <sub>2</sub> Vitamin D <sub>3</sub> - d <sub>6</sub>	405.1	134.6	12

Table 2. MRM Parameters

## Results and Discussion




Figure 2. Lower limits of quantification (LLOQ) for neat, unextracted 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub>.

## Results and Discussion




Figure 3. Chromatograms of extracted calibrators from ImmunoDiagnostic




Figure 4. Calibration curve for 1,25-dihydroxyvitamin D<sub>2</sub>




Figure 5. Calibration curve for 1,25-dihydroxyvitamin D<sub>3</sub>




Figure 6. 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> analyzed in 2 human samples. 1,25-dihydroxyvitamin D<sub>2</sub> was not detected in either sample but 1,25-dihydroxyvitamin D<sub>3</sub> was measured at 22 and 36 pg/ml in sample 1 and 2 respectively. Quantification of sample 1 required extrapolation of the calibration curve

## Conclusion

A highly sensitive and selective method for quantifying 1,25-dihydroxyvitamin D<sub>2</sub> and D<sub>3</sub> from human plasma has been optimized. By combining the sensitivity of the 6490 QQQ with Ifunnel Technology and the ImmunoDiagnostic extraction method, quantification at low pg/ml levels has been achieved.