

Quantification by LC-MS/MS for Pharma and Biopharma Applications

Diana Uria


LC-MS Product Specialist CH




Overview

1. Introduction – QQQ for quantitation
2. Pharma applications
3. 1290 Infinity + 6490 QQQ application examples
4. Biopharma applications
5. HPLC-Chip + QQQ application examples


Triple Quadrupole: SRM




MRM (Multiple Reaction Monitoring)




Spectrum with background ions (from ESI)




Q1 lets **only** target ion 210 pass through



Collision cell breaks ion 210 apart



Q3 monitors **only** characteristic fragments 158 and 191 from ion 210 for quant and qual.



Small molecules

PHARMA APPLICATIONS



Agilent Technologies

Main points on LC-QQQ

•Sensitivity

- Extra sensitivity facilitates the sample preparation step
- DMPK, ADME... LoD is critical

•Robustness

- Source robustness: reproducibility over time
- Required for method validation and long term studies

•Software

- Ease of use – not everybody needs to be an expert
- Data review in a quick and reliable manner
- User friendly reporting – LIMS compatibility

Mass Hunter Optimizer

Compound Entry & Method Setup

Optimizer Editor Panel

New Project Load Project Save Project SaveAs Project Save Compounds Start Optimization Stop Optimization Pause Optimization

Import From DataBase Import From Excel Export To Excel Hide Instrument Status Panel Hide Plot Control Panel

Compound Setup Precursor Ion Selection Product Ion Selection Optimator Setup

Compounds with known precursor and product ions Show results (fragmentor and collision energy)

<input checked="" type="checkbox"/>	Compound Name	Vial Number	Formula	Nominal Mass	Group	Precursor Ion	Fragmentor	Product Ion	Collision Energy	Flagged	CompoundID
<input checked="" type="checkbox"/>	Alprazolam-274, 2 n	P1-A8	C17H13N4Cl	308.0828741543	Pharmaceutical	309.2	135	274.1	35	<input checked="" type="checkbox"/>	29276b8c-e2c7-4
<input checked="" type="checkbox"/>	d5-Alprazolam-210,	P1-A7	C17H8D5N4Cl	313.1142578898	Pharmaceutical	314.2	135	210.1	45	<input checked="" type="checkbox"/>	3289d3e3-d4af-4
<input checked="" type="checkbox"/>	Alprazolam-205, 2 n	P1-A8	C17H13N4Cl	308.0828741543	Pharmaceutical	309.2	135	205.1	45	<input checked="" type="checkbox"/>	5d49ca38-1293-4
<input checked="" type="checkbox"/>	d5-Alprazolam-279,	P1-A7	C17H8D5N4Cl	313.1142578898	Pharmaceutical	314.2	135	279.2	35	<input checked="" type="checkbox"/>	b767af6d-8bb2-4
<input checked="" type="checkbox"/>	Alprazolam-281, 2 n	P1-A8	C17H13N4Cl	308.0828741543	Pharmaceutical	309.2	135	281.1	30	<input checked="" type="checkbox"/>	c5b9eb6f-2aee-4
<input checked="" type="checkbox"/>	d5-Alprazolam-286,	P1-A7	C17H8D5N4Cl	313.1142578898	Pharmaceutical	314.2	135	286.2	30	<input checked="" type="checkbox"/>	fd174aea-658f-40

Compound Setup Precursor Ion Selection Product Ion Selection **Optimizer Setup**

Sample introduction

Manual infusion using syringe
 Automatic infusion using Loop injection
 Injection (with or without column)

Polarity: POSITIVE

Method name: D:\MassHunter\methods\Alprazolam-StdESI.m

Data file: D:\MassHunter\data\Optimizer\

Fragmentor

Coarse From 125 To 150 Step 5

Fine Step (This takes more time)

Collision Energy

Coarse From 20 To 50 Step 10

Fine Step 5 (This takes more time)

MRM - Database

Database

Filter Compounds

Optimized Compounds

Date From To

Group Name

Project Name

Polarity

Search Compounds

Compound Name

Formula

Method

Show All Records


Compound Information

Compound Name	Formula	Nominal Mass	Method	Precursor Ion	Fragmentor	Product Ion	Collision Energy	Abundance	Project Name
Bromuconazole-2	C13H12N3OCl2B	374.95	<input type="checkbox"/> 300 Pesticides_1	373.30	100	<input type="checkbox"/> 70.04	20	12890	Default_Database
			<input type="checkbox"/> 75 Pesticides_12	375.96	150	<input type="checkbox"/> 158.98	24	48890	Default_Database
			<input type="checkbox"/> 75 Pesticides_12	375.96	150	<input type="checkbox"/> 70.04	16	17541	Default_Database
Bupirimate	C13H24N4O3S	316.16	<input type="checkbox"/> 300 Pesticides_1	317.16	120	<input type="checkbox"/> 272	20	116145	Default_Database
			<input type="checkbox"/> 300 Pesticides_1	317.16	120	<input type="checkbox"/> 166.1	25	435977	Default_Database
			<input type="checkbox"/> 300 Pesticides_1	317.16	120	<input type="checkbox"/> 166.1	25	571341	Default_Database
Buprofezin	C16H23N3O5	305.16	<input type="checkbox"/> 300 Pesticides_1	306.16	120	<input type="checkbox"/> 201.11	10	638148	Default_Database
			<input type="checkbox"/> 300 Pesticides_1	306.16	120	<input type="checkbox"/> 116.05	15	289254	Default_Database
			<input type="checkbox"/> 300 Pesticides_1	306.16	120	<input type="checkbox"/> 201.11	10	1227690	Default_Database
Butachlor	C17H26ClNO2	311.17	<input type="checkbox"/> Chinese_AE_Pes	312.17	80	<input type="checkbox"/> 162	20		Default_Database
			<input type="checkbox"/> Chinese_AE_Pes	312.17	80	<input type="checkbox"/> 238	10		Default_Database
			<input type="checkbox"/> Japanese_AE_Pe	312.17	100	<input type="checkbox"/> 252	10		Default_Database
Butafenacil	C20H18ClF3N2O	474.08	<input type="checkbox"/> 75 Pesticides_12	492.11	90	<input type="checkbox"/> 176	20		Default_Database
			<input type="checkbox"/> 75 Pesticides_12	492.11	90	<input type="checkbox"/> 330.8	23	122757	Default_Database
			<input type="checkbox"/> 75 Pesticides_12	492.11	90	<input type="checkbox"/> 348.8	10	71265	Default_Database
Butamifos	C13H21N2O4PS	332.1	<input type="checkbox"/> Japanese_AE_Pe	333.1	100	<input type="checkbox"/> 180	5		Default_Database
			<input type="checkbox"/> Japanese_AE_Pe	333.1	100	<input type="checkbox"/> 152	10		Default_Database

Quantitation of underivatized 17β-estradiol (E2) in human plasma

- MRM (271.345.216) 10 pg/mL E2 in 100 μL of plasma


10 pg/mL underivatized E2



3 replicates
Avg. Area = 13.1
Avg. Height = 5.2
Avg. Accuracy = 108.6 %
Area precision = 3.9%

Linearity – Underderivatized E2

Estradiol - 7 Levels, 7 Levels Used, 15 Points, 15 Points Used, 0 QCs



Related literature

Determination of Fluticasone Propionate in Human Plasma –

Sensitivity of Agilent 6490 Triple Quadrupole LC/MS Saves Time by Simplifying Sample Preparation

Determination of Estradiol in Human Plasma –

Agilent 6490 Triple Quadrupole LC/MS achieves low pg/mL Detection Limits without Sample Derivatization

Robustness of the 6490 Triple Quad LC/MS with iFunnel Technology for the Analysis of Verapamil in Human Plasma




Figure 1. Absolute Area Response (2,000 injections over a period of 4 days). Area RSD = 5.3%



Peptides, proteins, oligonucleotides...

BIOPHARMA APPLICATIONS

Following from before...



The 6490 Triple Quad LC/MS
Enables the Highest Sensitivity for
Peptide Quantitation in Plasma






Figure 2: MRM chromatograms (441.3 → 621.4) for the plasminogen SIS peptide spiked at different levels in plasma matrix.

HSA Standard Peptide: LVNEVTEFAK

Log/Log Plot of Calibration Curve




Agilent Proteomics Biomarker Workflow




Method transfer QTOF → QQQ

Multiple MRMS



Multiple MRM data have similar relative abundance ratios as seen in the QTOF MS/MS scan. This provides confirmation of the peptides identity without having to generate an MS2 scan and with **greater sensitivity**

MS/MS Spectrum



MRM Builder

Generates (D)MRM method from discovery QTOF data

Spectrum Mill - MRM Builder

Spectrum Mill | P/P Summary | **MRM Builder** | Autovalidation | Build TIC | MS/MS Search | Spectrum Summary | Tool Belt | Help

Initialized MRM Builder
Results Shown Filtered by Validation Category: valid
 Data Directory: msdataSM/NT/control
 hit table read - SpecFeatures read
 valid hits read from tagSummary file - Files: 5839 Hits: 5565
Filtered to only unique peptides: 1182
 beginning to assemble proteins ... proteins assembled 0.228381 sec
 proteins filtered by unique peptides 0.128938 sec
 calculated protein coverage maps 2.876987 sec
 beginning to roll up proteins into groups ... proteins rolled up into groups 2.850592 sec
 protein groups ready for displaySBA: Score

#	Filename	z	Spectrum Intensity	Sequence	RT (min)	Peak Width (sec)	m/z Measured (Da)
1	control-005_12594_12594_3	3	6.42e+005	(R) TFGIPGYTPVVMNVEVSPFTIEMSAFGYVFPK[A]	60.63	39.32	1189.6030
	P04114_TFGIPGYTPVVMNVEVSPFTIEMSAFGYVFPK_y28(+2)	1189.6	1537.8	10;;38			
	P04114_TFGIPGYTPVVMNVEVSPFTIEMSAFGYVFPK_y23(+2)	1189.6	1279.2	10;;38			
2	control-001_12307_12305_3	3	3.34e+006	(R) VPSYTLILPSLELPLVLEIYR[N]	49.33	16.19	746.4417
	P04114_VPSYTLILPSLELPLVLEIYR_y7(+1)	746.4	817.5	10;;22			
	P04114_VPSYTLILPSLELPLVLEIYR_y18(+2)	746.4	1024.1	10;;22			

#	Filename	z	Spectrum Intensity	Sequence	RT (min)	Peak Width (sec)	m/z Measured (Da)
1	control-005_12594_12594_3	3	6.42e+005	(R) TFGIPGYTPVVMNVEVSPFTIEMSAFGYVFPK[A]	60.63	39.32	1189.6030
2	control-001_12307_12305_3	3	3.34e+006	(R) VPSYTLILPSLELPLVLEIYR[N]	49.33	16.19	746.4417

Select Results for MRM Building

Build MFJts | Save Settings | Reset

Mode: Peptide
 Format: Agilent Triple Quad DMRM

Data directories: Select...

NT/control

Search result files:
 *.spx

Validation and Sorting

Filter results by: valid

Validation preset: none

Top n peptides for MRM: Limit to 2

Rank peptides by: Score

Filter peptides by:
 Score: > 10 % SPI: > 50
 Required AAs: any Disallowed AAs: none
 Peptide pI: from 3.0 to 10.0 All
 Accession #'s:

Sort MRM List by: Precursor m/z

MRM Parameters

Destination: Text file export Screen Preview

Top n Transitions: 2 Above precursor only y-ions only

Z options: Observed precursor/fragment charge only

Dwell time (ms): 10

Decustering Potential:
 m/z breakpoint: Potential:
 Less than 300 50
 in between 70
 More than 800 100

Collision Energy: (CE= slope*m/z + intercept)
 Use actual CE if available
 z-state: slope: intercept:
 2 0.05 5
 3 0.044 3
 4 0.05 2

Enter sequences and select b/y ions for optimization

Peptide Fragment Ion Predictor

<input type="checkbox"/>	Compound Name	Nominal Mass
<input type="checkbox"/>	AAFTECCQAMDK	
<input type="checkbox"/>	YLYEIAK	
<input type="checkbox"/>	LVNEVTEFAK	
<input type="checkbox"/>	KVPQVSTPTLVEV	
<input type="checkbox"/>	RPCFSALEK	
<input type="checkbox"/>	AVMDDFAAF	
<input type="checkbox"/>	HPYFYAFEL	

Peptide Fragment Ion Predictor


<input type="checkbox"/>	Compound Name	Nominal Mass
<input type="checkbox"/>	LVNEVTEFAK	

<input type="checkbox"/>	Precursor Ion	Precursor Ion Charge
<input type="checkbox"/>	575.3116965	2

<input checked="" type="checkbox"/>	Product Ion	Product Ion Charge	Product Ion Name	AA N-term	AA C-term
<input type="checkbox"/>	114.091889	1	b1	L	V
<input type="checkbox"/>	213.160303	1	b2	V	N
<input type="checkbox"/>	327.203231	1	b3	N	E
<input type="checkbox"/>	456.245825	1	b4	E	V
<input type="checkbox"/>	555.314235	1	b5	V	T
<input type="checkbox"/>	656.361918	1	b6	T	E
<input type="checkbox"/>	785.404512	1	b7	E	F
<input type="checkbox"/>	932.472926	1	b8	F	A
<input type="checkbox"/>	1003.51004	1	b9	A	K
<input type="checkbox"/>	147.113353	1	y1	A	K
<input type="checkbox"/>	218.150467	1	y2	F	A
<input type="checkbox"/>	365.218881	1	y3	E	F
<input type="checkbox"/>	494.261475	1	y4	T	E
<input type="checkbox"/>	595.309154	1	y5	V	T
<input checked="" type="checkbox"/>	694.377568	1	y6	E	V
<input checked="" type="checkbox"/>	823.420162	1	y7	N	E
<input type="checkbox"/>	412.2139335	2	y7	N	E
<input checked="" type="checkbox"/>	937.45309	1	y8	V	N
<input type="checkbox"/>	469.2354575	2	y8	V	N

Add Cancel

Automatic optimization of collision energies for predicted MRM transitions



Results of peptide optimization

MassHunter Optimizer (HSAtest4b)

New Project Load Project Save Project SaveAs Project Save Compounds Delete Projects Extract Peptide

Import From DataBase Import From Excel Export To Excel Start Optimization Ion Breakdown Profile Stop Optimization

Compound Setup Precursor Ion Selection Product Ion Selection Optimizer Setup

Show results summary

<input checked="" type="checkbox"/>	Compound Name	Group	Formula	Nominal Mass	Vial Number
<input checked="" type="checkbox"/>	LVNEYTEFAK			1	P1-F1

Method	Polarity	Ion Source	Instrument ID	Date Optimized	Flagged
D:\MassHunter\me	Positive	HPLC-Chip		12/18/2008	<input type="checkbox"/>

Precursor Ion	Fragmentor	Abundance
575.31	130	

Product Ion	Collision Energy	Abundance
937.5	16	80394
695.3	16	27258
694.4	16	18548
456.2	16	8572
823.4	16	8508
494.3	24	5845
555.3	12	4998
656.4	12	1309
785.4	12	503
932.5	12	479


Peptide sequence

Precursor

Optimized product ions


Abundance of each transition allows customer to choose best transitions for the final method

2,000 Dynamic MRM transitions in digest of depleted human sera – no time segment boundaries




nanoLC for extra sensitivity

100 amol on-column




1 fmol on-column







Agilent HPLC-Chip/QQQ LCMS Technology

Nanospray chip configuration brings new era in high sensitivity quantitation





Sensitivity: down to low amol
Dynamic range: up to 10^5

HPLC-Chip/MS Interface: Fluid Connections to the HPLC-Chip




Chromatographic Performance

Protein Digest Mixture



Reduced MS complexity + reduced ionization competition = improved ID

Retention Time Reproducibility




Extracted ion chromatograms for 17 peaks from a BSA tryptic digest (50 fmol on-column)

	RT	SD	%RSD
EIC 487.8	3.618	0.014	0.40
EIC 752	3.788	0.011	0.29
EIC 740.6	5.018	0.010	0.20
EIC 874.4	3.968	0.012	0.31
EIC 653.6	4.289	0.012	0.28
EIC 511.7	3.681	0.012	0.31
EIC 722.7	3.547	0.012	0.35
EIC 778	4.143	0.010	0.23
EIC 526.3	4.399	0.015	0.34
EIC 547.5	4.472	0.011	0.25
EIC 746.7	5.196	0.011	0.20
EIC 519.1	4.142	0.011	0.26
EIC 508.2	4.972	0.011	0.23
EIC 582.4	4.679	0.011	0.23
EIC 461.9	3.905	0.012	0.30
EIC 474	4.759	0.011	0.22
EIC 628	4.584	0.010	0.22

RT reproducibility evaluated using 69 repeat injections

6490 QQQ sensitivity




1290/6490 \rightarrow HPLC-Chip/6460: approx. 2 – 3 x improvement in signal intensity

HSA Standard Peptide: LVNEVTEFAK

Log/Log Plot of Calibration Curve

LVNEVTEFAK - 9 Levels, 9 Levels Used, 45 Points, 45 Points Used, 0 QCs



Summary

- 1290 + 6490 Agilent system solutions offers the most robust, sensitive LC-QQQ on the market
 - 50 fg reserpine, 900:1 s/n
 - iFunnel Technology for improved robustness
 - QuanOptimizer and Database Manager for easy compound setup
- HPLC- Chip and 6400 series Triple Quadrupole MS/MS provides high sensitivity and large dynamic range
 - Robust and stable nanoflow with HPLC-Chip
 - Good retention time and MS detection reproducibility
 - Peptide Optimizer for maximum tune sensitivity
 - Dynamic MRM for max sensitivity, minimum cycle time