

Herzlich Willkommen zum

Agilent Tipps & Tricks
Seminar:
Probenvorbereitung



Teilnehmerzertifikat & Feedback Form





Teilnehmerzertifikat

Erhältlich am Ende des Seminars

Feedback Form

Ihre Meinung ist uns wichtig! Bitte lassen Sie uns wissen, ob Ihnen das Seminar zugesagt hat, ob wir etwas verbessern können und ob Sie weitere Wünsche haben.

Bitte füllen Sie den Feedback Bogen aus. Als Dankeschön erhalten Sie einen USB Stick mit wertvollen Unterlagen

Die **Vorträge** senden wir Ihnen nach der Veranstaltung in einem separaten eMail zum Download zu.






Agenda Pharma / Biolog. Flüssigkeiten nachmittags

Title	Speaker	
Session I – Optimale Ergebnisse mit Probenvorbereitung Wie lässt sich durch Probenaufreinigung der Betriebsablauf im Labor positiv beeinflussen?		30 min
Session II, Teil A – Applikative Lösungen , Methodenentwicklung & Fehlersuche		30 min
Kaffeepause		15 min
Session II, Teil B – Applikative Lösungen , Methodenentwicklung & Fehlersuche		30 min
Session III – Die beste analytische Säule für die aufgereinigte Probe		30 min
Session IV – Zu guter Letzt: Die Analysengeräte		30 min







Optimale Ergebnisse mit Probenvorbereitung

Great Results Start with Sample Preparation

Dr. Ulrike Jegle
Agilent Technologies




Tips and Tricks

Throughout this presentation will be a series of tips and tricks to help improve your results and productivity.

Watch out for the Orange Boxes!

BUT please let us if there is a TIP or TRICK that you need in your lab







Sample Prep Terminology

Analyte(s): Molecule(s) of interest

Matrix: The sample (urine, blood, saliva, tissue)

Interferences: Entities inside the sample which will inhibit analysis

of desired analyte

IS: Internal standard

LLE: Liquid-liquid extraction

SLE: Supported-liquid extraction

SPE: Solid phase extraction

LC (MS): Liquid chromatography (mass spectroscopy)

GC (MS): Gas chromatography (mass spectroscopy)


LOQ: Limit of Quantification LLOQ: Lower limit of quantification

RSD: Relative standard deviation





Why Do Sample Preparation?



- Removal of interferences which would otherwise affect detection of analyte
- Concentration of an analyte to a detectable concentration
- Solvent Switching into an analytically more compatible solvent





Effects of Endogenous Interferences

- Poor Chromatography
- Mechanical issues (particulates, blockages)
- LC Column Lifetime Issues
- Instrument Downtime
- Carry Over
- Increase in sample run time/ cost
- Ion Suppression
- Overall loss in sensitivity

Let's take a look at some common matrices





Sample Types

Plasma



Saliva

Water

Mucus

Bacteria

Water (92%) Serum albumin (7%) Blood-clotting factors Immunoglobulin (antibodies) lipoprotein particles Other proteins Various electrolytes (Na+, Cl-, Ca++, K+ etc)

Whole Blood



Salts & Electrolytes **Enzymes & Proteins** Broken down foodstuffs When saliva is sampled it is often added to a preservative solution. This contains surfactants, PEG. Tween and other salts and buffers.

Water (90% by volume)

Serum albumins

Red Cells (erythrocytes)

White cells, Platelets

Lipids

Proteins/peptides

Hormones

Salts + mineral ions

Carbohydrates

(Prescription) drugs

Water

Fats, Lipids, Proteins

Cells

Collagen

Myosin + actin

Myoglobin

Carbohydrate

Salts + electrolytes



Urine



Water (95% by volume)

Urea (deamination product)

Creatinine

Salts & Uric Acid

Pigments

Carbohydrates

Enzymes

Mucins (glycosylated proteins)

In some urine sample analysis, an enzyme called Beta-glucuronidase needs to be used. This can cause issues with analyte detection





Dilution

Sample dilution (with an internal standard added) aka "dilute and shoot" is a commonly adopted simple approach to sample analysis

Advantages

Fast and easy

- Interferences are not removed
- Analyte concentration is reduced
- Instrument and Column contamination
- Ion Suppression






Liquid-Liquid Extraction (LLE)


Technique which uses two immiscible solvents to partition analytes of different polarities

Advantages

- Low Cost
- Effective for Inorganic/Salt removal

- No concentration of analyte achieved, sample will need to be reconstituted prior to analysis
- Uses large volumes of solvents (often chlorinated)
- No selectivity for organic soluble analytes
- Matrix specific (e.g. Urine vs Tissue)
- Difficult to automate









Supported-Liquid Extraction (SLE)

Technique which uses a highly wettable inorganic sorbent (diatomaceous earth).

Advantages

- Low Cost
- Effective for Inorganic/Salt removal

- No concentration of analyte achieved, sample will need to be reconstituted prior to analysis
- Uses large volumes of solvents (often chlorinated)
- No selectivity for organic soluble analytes
- Matrix specific (e.g. Urine vs Tissue)
- Can be automated




Chem Elut Produkte





SLE Mechanism







Matrix Solid Phase Dispersion (MSPD)

Mechanical grinding/homogenisation of a functional sorbent with a solid sample followed by elution from the sorbent

Advantages


- Excellent for tissue/solid sample
- Different functional sorbents can be used

- Interferences often co-extracted
- Can use large volumes of solvents
- Labour intensive
- Difficult to automate
- Matrix specific (e.g. Solids, Tissue, Post mortem bloods)




SPE – Solid Phase Extraction The Retention/Elution Process

Sample matrix




Retention

Wash solvent



Rinsing

Elution solvent



Elution





Typical Approaches to Plasma

Dilute and Shoot: Dilution & addition of ITSD

LLE or SLE: Agitation, partition or Tube

Precipitation: Centrifugation, Filter Plate

Solid Phase Extraction: Tube or Plate

	Dilution	LLE	Precipitation (Captiva)	SPE (Plexa)
Proteins	No	Partial	Yes	Yes
Lipids	No	No	No	Yes
Salts	No	Yes	No	Yes



Protein Precipitation

Protein Precipitation (PP) is a technique primarily aimed at the removal of proteins (albumins) from plasma samples.

This is achieved by crashing protein via the addition of organic solvent. The resulting precipitate can be removed by centrifugation or filtration.

Advantages

- Effective removal of protein
- High throughput and automation friendly
- Generic methods across plasma types
- Flexibility

- Low selectivity for analytes
- Salts and Lipids not removed
- Low levels of peptides/small proteins
- Column and Instrument contamination








Sample Preparation Workflow: Protein Precipitation

- Fractionate Plasma from whole blood (centrifuge)
- Dilute plasma 3:1 with MeOH or ACN.
- Agitate the sample
 - Aspiration
 - Vortexing
 - Shaking
 - Stirring
- Filter under vacuum
- Inject sample into LCMS








Solid Phase Extraction

- Solid Phase Extraction is the most commonly used sample preparation technique.
- The use of an immobilized sorbent, pre-packed into a device.

Advantages

- Best interference removal
- Gain concentration effect
- Extract whole suites of drugs
- Flexibility and Automation

- Requires method development
- Difficult to separate analytes with similar properties/structure!







Sample Preparation Workflow: SPE

- Fractionate Plasma from whole blood (centrifuge)
- Precondition tube
- Dilute sample with appropriate buffer (add IS)
- Apply the sample
- Wash steps
- Elution (pH modified)
- Filter under vacuum
- Inject sample into LCMS

Direct Injection or Dry-down /reconstitute

	Plexa (30 mg, 1mL)
Sample	100 μL human plasma
Pretreatment	Dilute 1:3 w/ 2% NH ₄ OH (bases) 2% Phosphoric (acids)
Conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O
Washes	5 % MeOH in water
Elution	2 x 250 µL 50:50 ACN/MeOH (or 100% MeOH)





Silica and Polymer Based SPE A Plethora of Options

Bond Elut AccuCAT

Bond Elut Alumina (AL-A)

Bond Elut Alumina (AL-B)

Bond Elut Alumina (AL-N)

Bond Elut NH₂

Empore Anion-SR

Bond Elut Atrazine

BioBond Elut

Bond Elut C1

Bond Elut C2

Bond Elut C8

Elut C18

Elut C18 EWP

Bond Elut C18 INT

Bond Elut C18 OH

Bond Elut CBA

Bond Elut Certify

Bond Elut Certify II

Bond Elut CH

Bond Elut DEA

Bond Elut Diol (20H)

Bond Elut ENV

Bond Elut FL

Bond Elut NH2

Bond Elut Plexa, PCX, PAX

PAXBond Elut PBA

Bond Elut PCB

Bond Elut PH

Bond Elut PPL

Bond Elut PRS

Bond Elut PSA

Bond Elut SAX

Bond Elut SCX

Bond Elut SI



40 + Phases and over 2000 Applications





Bond Elut SPE

Non-polar

Plexa, C18, C8, C2, C1
C18 variations in carbon load
and endcapping

EnvirElut

CH - cyclohexyl

CN-E – endcapped cyano

PH – phenyl

ENV, LMS, PPL, Focus, Nexus **Plexa**

Anion Exchange

SAX – quaternary amine

PSA – primary and secondary amine

NH2 – aminopropyl

DEA – diethylaminopropyl

Polar

Plexa PAX

PSA - primary and secondary amine

NH2 - aminopropyl

DEA - diethylaminopropyl

Diol

Si - silica

Cation Exchange

Plexa PCX

SCX - benzenesulfonic acid

PRS – propylsulfonic acid

CBA – carboxylic acid

Reversible Covalent

PBA - phenylboronic acid

Specialty Phases

AccuCAT Atrazine Etc..

Mixed mode IEX/NP

Certify – SCX/C8 Certify II – SAX/C8

Plexa PCX Plexa PAX

Alumina – aluminum oxide Florisil – magnesium-silica Carbon Carbon/NH2





Formats

The most comprehensive range of formats in the industry

Tubes

1mL to 60mL Straight Barrel (50mg -10g)


Bond Elut Jr (500mg –1g)

LRC (Large Reservoir Capacity)

(100-500mg)

Mega Bond Elut

12-150mL (2g-70g)



Multi-Array

1mL, 2mL 96 Well



Automation

Hamilton

TomTEC

Gilson ASPEC

Gerstel

Spark Holland Prospekt









SPE Manifolds



96 Well

Captivac Collar: For use with Bond Elut 96 1mL and Captiva filtration plates



VacElut 12, 20 and SPS 24

vacuum manifold for SPE syringe barrels.





Basic Accessories

Bond Elut Adaptors

- Connect cartridges in series
- Fit syringe for positive pressure

Luer stopcocks

Stop flow on all Luer tip style cartridges

Frits, Reservoirs, Fritted reservoirs





SPE Methods? Where Do I Start?

Consider your analyte?

- Do the analytes dissolve into the matrix
- Does the analyte appear to be polar or non polar
- Do the analytes contain any ionic groups
- Are the compounds unstable in acid or base
- What is the method of analysis?
- What is the concentration of the analyte in the sample

Consider your matrix?

- Is the matrix polar or non polar: (olive oil or urine?)
- is the matrix high in organic or ionic content? (will dilution of the sample be necessary)
- What is the sample size?
- Does the sample require pH adjustment?
- Is an internal standard required?
- Is filtration or centrifugation necessary upfront?






The SPE Triangle

Consider LC Theory

- The same physical processes in liquid chromatography also apply in SPE
- Think of SPE as a <u>digital</u> LC process
 - Analyte is ON (retention)
 - Analyte is OFF (elution)
- We are not looking at separation
- Normal or Reverse phase?
- Ionic modifier?







TRICK: Use Google!

- **Bond Elut SPE has more** published methods than any other SPE platform
- 1000s of articles
- Official methods
- 31 years of heritage

BOND ELUT + (Analyte)

Don't re-invent the wheel!



About 14,500 results (0.34 seconds)

Everything

▼ More

Bengaluru. Karnataka

Change location

The web

Pages from India

More search tools

Retention studies of acrylamide for the design of a robust liquid ...

by J Rosén - 2007 - Cited by 8 - Related articles

Based on these investigations, a method for analysis of acrylamide in food using Bond Elut PPL, Polymer, Styrene divinyl benzene polymer, modified ... linkinghub.elsevier.com/retrieve/pii/S0021967307016275

Determination of acrylamide in Chinese traditional carbohydrate ...

by Y Zhang - 2007 - Cited by 9 - Related articles

Results showed that acrylamide was not completely adsorbed by cartridges ... linkinghub.elsevier.com/retrieve/pii/S0003267006021684

■ Show more results from elsevier.com

Analysis of Coffee for the Presence of Acrylamide by LC-MS/MS ...

by D Andrzejewski - 2004 - Cited by 71 - Related articles

10 Mar 2004 ... After elution of sample extract through the Bond Elut Accucat SPE cartridge. some acrylamide will remain on the cartridge. ... pubs.acs.org/doi/abs/10.1021/jf0349634 - Similar

New Research Developments on **Acrylamide**: Analytical Chemistry ...

by Y Zhang - 2009 - Cited by 1 - Related articles

coffee, 13C3-acrylamide, LC-MS, SIM, defatting with hexane, extraction with ... pubs.acs.org/doi/abs/10.1021/cr800318s - Similar

■ Show more results from acs.org

A Modified Sample Preparation for Acrylamide Determination in ...

by AP Arisseto - 2008 - Cited by 2 - Related articles

method used for the determination of acrylamide in potato and cereal products to achieve good.... and Bond Elut-Accucat® cartridges are used in combina- ...

www.springerlink.com/index/kôi4528374645875.pdf - Similar

IPDFI Eriksson - Analysis of acrylamide in food, air and in ...

File Format: PDF/Adobe Acrobat - View as HTML

Different limits for acrylamide contaminations. ... AccuBON I SCX. Bond Elut AccuCAT. Bond Elut C18. Bond Elut Jr-PSA. ENVI-Carb. Isolute® MF C18 ...

www.apaceurope.com/2006/seriksson.pdf





How do Sample Prep Options Compare?

Solid Phase Extraction (Bond Elut)

- Highest level of clean up
- Potential for the most selectivity

Protein precipitation (Captiva, Captiva NDLipids)

- Works for most analytes with no or very little method development
- Using Captiva can improve throughput and ease-of-use
- Captiva ND Lipids removes lipids also

Liquid/Supported Liquid Extraction (Chem Elut)

- Reasonably clean, very inexpensive
- Can be labor intensive, especially for larger sample sizes.

Filtration

Mechanical removal of particulates only

Dilute and Shoot

No clean up, no interference removal



Complexity

Cleanliness







Fragen???





Teil A – Applikative Lösungen, Methodenent-wicklung

Ion Suppression
Polymer based SPE Sorbents



Agilent Technologies

Polymeric SPE - Plexa: Technology and Applications

Advanced Polymeric SPE

- Non polar (Bond Elut Plexa)
- Ion Exchange applications (Bond Elut Plexa PCX and PAX)

Plasma based applications

- Ion suppression studies
- Multi suite drug extraction

Alternative Matrices

- Saliva (Multi-suite DoA analysis)
- Urine (Therapeutic and DoA multi-suite testing)



Plasma

Composition

- Water (92%)
- Serum albumin (7%)
- Blood-clotting factors
- Immunoglobulin (antibodies)
- lipoprotein particles
- Other proteins
- Various electrolytes (Na, Cl, Ca, K etc)





Ion Suppression

Ion Suppression is a deleterious phenomenon that can be experienced during LCMS analysis.


Causes

Endogenous matrix interferences


- Proteins and Peptides
- Lipids (phospho- and lysophospholipids)
- Salts
- Other hydrophobic species
- Dosing vehicles/preservatives (Tween, PEG)

Effect on Data

- Loss of Ion signal will inhibit quantification of analyte
- Loss of sensitivity, precision and accuracy



Ion Suppression Regions







TIP: Monitoring Phospholipid Removal

PHOSPHATIDYLCHOLINES

LYSOPHOSPHATIDYLCHOLINES

- Monitoring the 184 → 184 m/z trace in MS/MS gives a measure of the total phospholipids and lyso-phospholipids in a sample
- Measures the ionic head group common to all classes of phospholipids
- Install a new LC column first!

"We are now monitoring the 184→184 transition in all of our bioanalytical method development studies"

DMPK Group, AstraZeneca, UK





Typical Approaches to Plasma

Dilute and Shoot: Dilution & addition of ITSD

LLE: Partitioning between immiscible solvents

Precipitation: Centrifugation, Filter Plate – Unretained Method

Solid Phase Extraction: Tube or Plate – Retained Method

	Dilution	LLE	Precipitation (Captiva)	SPE (Plexa)	
Proteins	No	Partial	Yes	Yes	
Lipids	No	No	No	Yes	
Salts	No	Yes	No	Yes	



Polymeric SPE: Advantages and Benefits

Some advantages that are common to all polymeric SPE compared to silica SPE are:

- Stability in a wide pH range useful for all method options
- High capacity/loadability smaller SPE bed masses minimize solvent usage and create greater throughput
- Loading of 10-12% of bed Capacity (Interferences count!)
- Resistance to drying after conditioning creates greater errorfree ease of use and throughput for end users



Polymeric SPE: Advantage to Pharma/biofluids

The adoption of polymeric SPE for bio-fluid analysis is the industry norm:

- Simple and robust methodologies.
- It is resource heavy to method develop, unknown or novel compounds – offers different interaction possibilities.
- Multi-suite analysis large number of internal standards need extra capacity- offers different interaction possibilities.
- Very cost effective.
- Long shelf life (non-hygroscopic, stable material)





The Bond Elut Plexa Family

Bond Elut Plexa

Non-polar polymer with a hydroxylated surface functionality is the best choice for a wide range of acidic, neutral and basic analytes.

Bond Elut Plexa PCX

based on Plexa, with a strong cation exchange functionality, specifically designed for effective results for a wide range of pKa and logP basic compounds from plasma

Available SOON


Bond Elut Plexa PAX

based on Plexa, with a strong anion exchange functionality, specifically designed for effective results for a wide range of pKa and logP acidic compounds from plasma





Bond Elut Plexa - Aufbau



Amid-functionality, which could cause interferences with the matrix are not present at all.






How Bond Elut Plexa works

Apply Sample


Water Rich Hydrophilic component allows excellent phase transfer of analytes into the polymer core



Large endogenous proteins do not bind to the surface of the polymer and cannot access pore structure.

Washing


Analytes remain tightly bound in the hydrophobic core.



Interferences (lipids and proteins) washed away without leaching into the analytes of interest

Elution

Specially engineered pore structure allows excellent mass transfer out of polymer



Clean extract with high recovery





How Bond Elut Plexa PCX works


interest

Apply Sample

Water-rich Hydrophilic component allows excellent phase transfer of analytes to the cation exchange sites within the polymer.


Washing

Analytes remain tightly bound to the polymer through strong ionic interactions



Elution

Change in pH allows the elution of analytes from the strong cation



Large endogenous proteins cannot access pore structure.

Interferences can be washed away without the leaching of basic analytes Clean extract with high recovery





Plexa Recoveries Versus Other Polymer

Analyte	Plexa	Other Polymer	рКа	LogP
Albuterol	97.9	115.4	5.9	1.3
Atenolol	97.0	94.0	4.2	4.2
Loratadine	71.0	49.0	5.7	1.5
Metoprolol	92.0	74.0	5.7	1.5
Naltrexone	85.7	13.0	4.9	5.2
Pravastatin	85.0	59.0	4.9	5.2
Propranolol	55.0	35.0	4.9	5.2
Zolpidem	93.0	96.8	9.9	3.4

Plexa gives equivalent or better absolute recoveries than other polymer for a wide range of pKa and logP compounds

Conditions: Absolute Recovery from human plasma

Basic load conditions, n=6, 200 ng/ml





Absolute vs. Relative Recovery

Absolute Recoveries (AR)

Plasma + ISTD + drug

Extract

Inject

Compared to

Mobile Phase + ISTD + drug. No extraction, no loss of recoveries, no matrix effects

Relative Recoveries (RR)

Plasma + ISTD + drug

Extract

Inject

Data is then compared to a graph of extracted drug controls.

Drugs are the same and should extract at the same efficiencies.

If the standard curves were extracted at 80% and the unknowns at 80% the relative recoveries will show to be at about 100%!!

Relative recoveries (via calibration curves) can mask the affects of ion suppression, whereas absolute recoveries do not "cancel" these affects.

Minor differences in absolute vs. relative recoveries is indicative of extract cleanliness or superior extraction characteristics.






Absolute vs. Relative Recovery

Small differences between AR and RR indicate high extraction efficiency and low signal loss due to ion suppression.


Larger differences indicate sensitivity problems arising from low extraction efficiency or ion suppression








Post Column Infusion Experiments








Improved Precision - Reduced Ion Suppression







Typical PCI Experiment



Ion Suppression is a very real and concerning phenomenon in the Pharmaceutical Industry

- The loss of quantification of an analyte can create serious issues for the analytical chemist.
- The FDA recently published guidelines on how to determine the effects of ion suppression for low abundance metabolites in DMPK studies (FDA MIST regulations)
- Over 120 publications since 2005 which look at ion suppression in detail






TIP: Evaluating Ion Suppression

Absolute Recovery

Relative Recovery



Cation exchange polymer X

Simple Route:

Compare absolute and relative recoveries (if the difference is >10% after n = 3 then you have a significant matrix effect.

Comprehensive:

Run and isocratic or gradient (preferred) Post Column Infusion experiment

Extra View:

Monitor the 184→184 transition for lipids





Simple Plexa Methodology

Bond Elut Plexa

Bond Elut Plexa PCX

																			1					

Condition: 500 µL MeOH

500 µL water

Apply sample: 100 µL plasma

dil 1:3 w/ 2%

NH₄OH

Wash: 500 µL 5% MeOH

Elute: 500 µL MeOH

Bases

Condition: 500 µL MeOH

500 µL water

Apply sample: 100 µL plasma dil 1:3 w/

2% phosphoric acid

Acidic wash: 500 µL 2% formic acid

Neutral wash: 500 µL MeOH:ACN (1:1 v/v)

Elute: $500 \mu L 5\% NH_3$ in MeOH:ACN (1:1)

Methods for 30 mg Plexa in 96 well plate





PAX Method

Bond Elut Plexa PAX

<u>Acids</u>

Condition: 500 µL MeOH

500 µL water

Apply sample: 100 µL plasma

1:3 w/ 2% NH4OH

Wash: 500 µL water

500 μL 5% MeOH

Elute: 500 µL MeOH w/ 5%

Formic Acid

Methods for 30 mg Plexa in 96 well plate





TRICK: Maximizing Elution Recovery

Elution Solvent Composition

- Most methods recommend 100% MeOH or 5% NH3 in MeOH elution (CX)
- We recommend 50:50 MeOH/ACN
- Combination of polar protic and dipolar solvent breaks more analyte interactions with SPE phase

Elution Volume

- Most methods recommend 1000uL of eluting solvent
- 2 x 500uL often gives higher recovery
- For Plexa methods 2 x 250uL (30mg bed)







Multi-Suite SPE Method – Bond Elut Plexa PCX

Fractionation of Acidic, Neutral, and Basic Drugs from Plasma with Polymeric SPE cation exchange, Bond Elut Plexa PCX 10 mg

Acids:

Atorvastatin, Diclofenac, Furosemide, Pravastatin

Neutrals:

Cortisone, Cortisol

Bases:

Procainamide, Metoprolol, Paroxetine

Sample pre-treatment:

100 μL human plasma

Dilute 1:3 with 2% H₃PO₄

Condition 1. 500 µL CH₃OH

2. 500 µL DI H₂O

Load Plasma 1:3 with 2% H₃PO₄

Wash 1 500 µL 2% Formic acid

Elute 1 500 μL ACN:MeOH (1:1, v:v)

Acids, Neutrals

Elute 2 500 μL 5% NH₃ in ACN;MeOH Bases





Recovery Data: Acids & Neutrals

Acids:

Analyte	% Rec 0.5 μg/mL	% RSD	% Rec 1μg/mL	% RSD
Diclofenac	101	4	101	5
Furosemide	99	3	96	2
Pravastatin	95	4	96	6
Atorvastatin	100	4	100	5

Neutrals:

Analyte	% Rec 0.5 μg/mL	%RSD	% Rec 1 μg/mL	% RSD
Cortisone	93	4	97	6
Cortisol	101	4	101	4



Recovery Data: Bases

Bases:

Analyte	% Rec 0.5 μg/mL	% RSD	% Red 1 μg/mL	% RSD
Procainamide	100	5	98	3
Metoprolol	94	4	92	6
Paroxetine	94	5	99	4

Multi-Suite methods have excellent application opportunities

- Hospital or clinic based multi-therapeutic screening
- Drugs of Abuse confirmation
- Drug metabolism studies (where acid or neutral metabolites may need to be extracted using a different method)
- Save time, increase throughout and reduce no of SPE operations and sample injections





Urine

Composition

- Water (95% by volume)
- Urea (deamination product)
- Creatinine
- Salts & Uric Acid
- Pigments
- Carbohydrates
- Enzymes
- Mucins (glycosylated proteins)



In some urine sample analysis, analytes need to be hydrolysed from their glucuronide form using enzyme or base. This can affect extraction efficiency.



Plexa PCX Multi-Suite Urine Method: Bases

- This generic method covers a broad range of basic drug classes
- Amphetamines
- Opiates
- Benzodiazepines
- Ketamine/PCP
- Cocaine and Metabolites
- Antidepressants
- Synthetic Analgesics
- Most multi-suite SPE devices have poor recovery with THC and THC-COOH (next slide)

	Plexa PCX 96 well plate (30 mg, 1mL)
Sample	0.2 mL urine
Pretreatment	Dilute 0.6 mL 100 mM KH ₂ PO ₄
Conditioning	 0.5mL MeOH 0.5 mL H₂O
Washes	1. 0.5mL 50% MeOH in H₂ODry 5 minutes
Elution	0.5 mL 50:50:5 EtAc:MeOH:NH ₃





Plexa PCX: THC and THC-COOH

THC

Concentration	% Recovery	RSD (n=6)
1.0 ng/mL	102%	4%
50 ng/mL	107%	5%
1000 ng/mL	107%	3%

THC-COOH

Concentration	% Recovery	RSD (n=6)
1.0 ng/mL	101%	8%
50 ng/mL	97%	4%
1000 ng/mL	102%	4%

	Plexa PCX 96 well plate (30 mg, 1mL)
Sample	0.5 mL urine
Pretreatment	Dilute 0.5 mL 2% Formic Acid,
Conditioning	 0.5mL MeOH 0.5 mL H₂O
Washes	 0.5mL 2% Formic 0.5mL 5% Aqueous MeOH Dry 1 minutes
Elution	0.5 mL MeOH





Saliva (Oral Fluid)

Composition

- Water
- Mucus
- Salts & Electrolytes
- Enzymes & Proteins
- Bacteria
- Broken down foodstuffs



When saliva is sampled it is often added to a preservative solution. This contains surfactants, PEG, Tween and other salts and buffers.





Trends in Oral Fluid Anlaysis

Trends

- Growing matrix for drugs of abuse testing
- Non intrusive test, can be performed in any location
- Low incidence /chance of sample adulteration/ spoiling
- Use of surfactants/preservatives in collection device
 - Stops drug binding to plastic and cellulose bud
 - Reduces sample degradation

Challenges

- Small sample volume
- Contains LCMS ion suppressing surfactants/ oligomers





Plexa PCX for Oral Fluid Analysis

Single Drug Suite

	Plexa PCX (30 mg 3mL)
Sample	0.5 mL oral fluid
Pre-treatment	Dilute 1:2 w/ 100mM phosphate buffer
Conditioning	1. 1 mL MeOH 2. 1 mL H ₂ O
Washes	 1 mL 100mM KH₂PO₄ 1 mL 100% MeOH
Elution	1.0 mL 2% NH3 in EtAc

Multi Suite (Synergy Health)

	Plexa PCX (1mL 30mg)
Sample	600µl of buffered Saliva
Pre-treatment	500µl of 0.1M, pH6 phosphate buffer.
Conditioning	1. 1 mL MeOH 2. 1 mL H ₂ O
Washes	1. 1ml, 0.1M HCl 2. 1ml, 60:40 Methanol/0.1M HCl, then dry.
Elution	150μl, 50:50 (ACN:MeOH) 2x 150μl, 50:50:2 ACN/MeOH/ NH4OH



Recovery Data

Drug	% Recovery	%RSD
Morphine	82	15
Codeine	96	2.0
Dihydrocodeine	92	4.3
6-MAM	74	14
Amphetamine	75	5.5
Methamphetamine	93	2.8
MDA	78	8.1
MDMA	88	5.3
MDEA	97	3.3
Benzoylecgonine	95	2.8
Methadone	98	1.5
EDDP	31	7.0
Flunitrazepam	88	5.2
Aminoflunitrazepam	77	2.9
Aminonitrazepam	67	9.9
Diazepam	90	5.9
Nor-Diazepam	91	7.7
Temazepam	80	3.3





COMING SOON: Plexa PAX

Bond Elut Plexa PAX is specifically designed for effective results for a wide range of pKa and logP acidic compounds from plasma

- Made from the same base polymer technology
- Optimized loading of SAX groups ensures good retention/elution balance
- Robust, generic method, for the widest range of acidic analytes



PAX Applications

	Plexa PAX (10 mg)
Sample	100 μL human plasma
Pretreatment	Dilute 1:3 w/ 2% NH ₄ OH
Conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O
Washes	 500 μL H₂O 500 μL MeOH
Elution	500 µL 5% Formic Acid :MeOH

Compound	LogP	рКа
Atorvastatin	5.7	4.5
Diclofenac	4.2	4.2
Furosemide	1.5	4.7
Ketoprofen	3.2	5.2
Pravastatin	2.6	4.6

Mobile Phase -

A: 5mM Ammonium Formate

B: Methanol

Gradient -

t = 0 - 60% A: 40% B

t = 3:0 - 3:59 min 20% A: 80% B

t = 4:0 - 5:00 min 60% A : 40% B

Column -

Pursuit XRs Diphenyl 2.8 µm 50 x 2.0 mm

All samples evaporated and reconstituted in 100 μ L of 80:20 5mM Ammonium Formate: MeOH






Method Comparison

	Plexa PAX (10 mg)
Sample	100 μL plasma
Pretreatment	Dilute 1:3 w/ 2% NH ₄ OH
Conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O
Washes	 500 μL H₂O 500 μL MeOH
Elution	2 x 250 μL 5% Formic Acid :MeOH

Base pre-treatment further activates AX groups, improving recovery







Recovery Data

- Absolute recoveries used against other AX product (10mg Plate)
- Other vendor's method had very poor recoveries, data shown is for Plexa method on other AX
- RSDs for AX methods typically high

Absolute Recoveries (100 ng/mL)				
Compound	Bond Elut Plexa PAX		Other AX Polymer	
	% Rec	RSD(n=6)	% Rec	RSD(n=6)
Atorvastatin	87	13	37	21
Diclofenac	67	16	42	22
Furosemide	99	11	40	26
Ketoprofen	71	18	49	31
Naproxen	93	15	50	20
Pravastatin	87	14	106	10






Summary

The Bond Elut Plexa family offers the opportunity to increase throughput and improve the quality of your data


- The hydroxylated surface minimizes binding of matrix interferences offering greater accuracy and precision
- Proven minimization of ion suppression in complex matrices
- Amenable to a broad range of biological matrices
- Ultra generic methodology for urine samples
- Multi-suite methods for saliva (inc. Tween exclusion)
- Simple methodologies allow rapid method development
- Three phases offer high recoveries for most compounds and matrices simple choices for quick results





Kaffeepause 15 min





Teil B – Applikative Lösungen, Methodenent-wicklung

Precipitation – Filtration
Troubleshooting



Agilent Technologies

Typical Approaches to Plasma

RECAP:

Dilute and Shoot: Dilution & addition of ITSD

LLE or SLE: Agitation, partition or Tube

Precipitation: Centrifugation, Filter Plate

Solid Phase Extraction: Tube or Plate

	Dilution	LLE	Precipitation (Captiva)	SPE (Plexa)
Proteins	No	Partial	Yes	Yes
Lipids	No	No	No	Yes
Salts	No	Yes	No	Yes



Typical Approaches to Plasma

What happens if plasma SPE methods are not attractive :

- Percieved cost?
- Workflow change?
- Complexity?
- Laziness?
- Sample throughput


	Dilution	Precipitation (Captiva)	SPE (Plexa)
Proteins	No	Yes	Yes
Lipids	No	No	Yes
Salts	No	No	Yes



Captiva Filtration Products

They are Mechanical Filters

- Protein precipitated plasma (0.45µm)
- Fibrinogen from plasma storage (20µm)
- Tissue homogenates (0.45µm)
- Hepatocyte filtrate (liver cells) (0.45µm)
- General particulate filtration







How Captiva Improves Precipitation

Centrifugation not necessary

processing by vacuum or positive pressure

Easier automation

Captiva filter plates interface with most automation platforms (while the centrifugation method of ppt is nearly impossible to automate)

Fewer transfer steps

no pellet means dry down & reconstitution in single collection plate

- Cleaner, particulate-free samples
 - no pellet means direct trouble-free sampling from collection plate
 - filtration effectively removes more proteins from injection



TIP: Method Development in Protein Precipitation


Method development in Protein Precipitation is not complicated

- Acetonitrile good pellet "qualities" commonly used
- Methanol fine particles , needs filtration
- Organic solvent ppt typically is a 3:1 or greater (as much as 10:1)
- AVOID:
- 10% aq. TCA highly efficient, highly aqueous supernatant, destroys LC columns quickly
- metal salts highly efficient, highly aqueous supernatant, not MS compatible






Post Column Infusion: Captiva



Propranolol Post Column Infusion



IONICS

PEPTIDES

PHOSPHATIDYLCHOLINES

LIPIDS & **HYDROPHOBICS**





Advanced Precipitation Techniques

Is there a missing technology here?

- Is there a compromise between the ease of use of PP and the cleanliness of SPE?
- Salts are not a concern in reverse phase chromatography
- Lipid Removal PP?

	Precipitation (Captiva)	Lipid Removal Precipitation	SPE (Plexa)
Proteins	Yes	Yes	Yes
Lipids	No	Yes	Yes
Salts	No	No	Yes



Lipid Build-Up With a Ballistic Gradient

Sample: Protein Precipitated Porcine plasma (3:1 ACN precipitation)

Column: 40x4.0 Polar modified reverse phase

Solvents:

A: 0.1% Formic Acid


B: Acetonitrile


Gradient:

Time (min)	<u>%B</u>	Flow (µL/min)
0:00	10	500
0:40	10	500
1:40	90	500
2:20	90	500
2:21	10	500
3:00	10	500

Detection: 184→184 (nonselective Phospholipid head group channel)

Gradient Program










Lipid Removal Precipitation: Captiva NDLipids




- Captiva particulate filter → removes protein interferences
- Proprietary Lipid Stripping Media → removes lipids
- Non-Drip Membrane → ease of use






Post Column Infusion Studies: Metroprolol







Lipid Removal



Centrifuge PPT

2:1 MeOH

3:1 MeOH

2:1 MeCN

3:1 MeCN





TIP: Captiva ND Lipids Method

- Add 3:1 pH modified MeOH: plasma to the well
- Mix thoroughly
- Filter
- pH modification to charge the analyte(s)
- For most analytes, including polar basic small molecules 0.1% formic acid in MeOH is sufficient
 - For acidic analytes 5% ammonia can be used instead
- Acid concentration can be increased if necessary, but care needs to be taken to avoid sample gelation
- MeOH is the best solvent for phospholipid removal





Blood

Spotted blood samples:

 Add excised blot spot disk directly to Captiva ND^{Lipids} well and use standard method

Whole blood:


- 1. Dilute whole blood sample 1:1 with 0.1% aqueous formic. Vortex and allow to sit for 1 hr.
- 2. Add MeOH with 0.1% formic acid: diluted blood sample to Captiva ND^{Lipids} plate (3:1)
- 3. Seal the plate with DuoSeal cover and vortex for 1 minute.
- 4. Leave cover in place. Place plate on Captivac with collection plate in place. Apply mild vacuum.
- 5. Filtrate should be clear and may require dilution prior to LC/MS






Universality of CaptivaND^{Lipids}

Analytes	LogP	Relative Response vs. Protein Precipitation Alone (%)
Zolpidem	3.32	109.6
Mianserin	3.52	167.4
Tranylcypromine	1.40	251.4
Nefazodone	4.70	144.9
Amoxapine	3.10	171.1
Maprotiline	5.10	159.0
Nomifensine	2.94	162.6
Warfarin	3.51	100.5
Sulindac	3.59	136.7
Loratidine	3.65	103.4
Loperamide	6.25	106.7
Vardenafil	6.01	106.5









Technique Summary

	Filtration (Captiva)	Filtration + (Captiva Lipids)	SPE
Proteins	Yes	Yes	Yes
Lipids	No	Yes	Yes
Salts	No	No	Yes

- Simple & easy to use
- Removes proteins
- Inexpensive

- Removes proteins & lipids - reduced ion suppression
- Simple & easy to use
- 3 step process
- Low Cost

- More involved
- Removes proteins, lipids, <u>salts</u> etc
- Better for <u>lipophilic</u> drugs
- Reasonable





TIPS & TRICKS: Trouble Shooting

From our global helpdesk logs, we have identified the most common SPE issues encountered by our end users.

All of these issues can be linked to simple practical errors



- 1. Low Recovery
- 2. Poor Flow
- 3. Loss of Analyte
- 4. Dirty Extracts



TIP: Sorbent Conditioning

Sorbent conditioning is vital for good SPE performance

- An unconditioned sorbent bed can result in Poor Flow and Poor Recovery.
- Particularly for biological samples, where precipitation can occur on dry media
- Clogging and Channelling

Recommendations

Silica SPE	Polymeric SPE
>2 bed volumes of low viscosity organic solvent (e.g. MeOH)	>2 bed volumes of low viscosity organic solvent (e.g. MeOH)
Condition slowly under minimal vacuum <2" Hg	Wash with aqueous solution prior to sample load
Allow time for equilibration (30-60s)	More tolerant to accidental drying





TIP: Poor Flow - Vacuum Considerations

Vacuum control and experience comes with practice

- Flow rates should be kept on the lower side overall (2-5 mm Hg)
- 2-3 mL/min is optimal for load, wash and elution
- Ion Exchange as a primary or secondary interaction is especially affected by high flow rates (break through)
- Even C18 can show band broadening with high flow rates
- Employ variable vacuum or variable time (pulsing) depending on sample type and drying needs









TIP: Poor Recovery - Capacity

Understanding the Capacity of your SPE Phase is Critical

Sorbent Type	Capacity
Silica (Polar or Non-Polar)	1-4% of bed mass
Silica (Ion Exchange)	< 1.0 mmol/g
Polymeric	10 -12% of bed mass

- Assumes good analyte/sorbent interaction
- Capacity does not distinguish between analyte and interference!
- Do you know amounts of interferences present?
- Stacking two cartridges on top of each other will help determine if capacity is an issue






TIP: Poor Recovery - Elution

Elution volume optimization can be the simple to achieve

- Elution volumes can be minimized and recovery maximized if elution is done in aliquots.
- Four aliquots of 50uL give better results than 200uL in a single shot. 2 x 100uL is a more practical solution.
- Using mixtures of solvents can have a dramatic effect on recovery (50:50 MeOH/ACN)



A: 4 x 50 uL Aliquots B: 1 x 200 uL Aliquot





TIP: Dirty Extracts

Dirty Extracts are common in un-refined SPE methods

- Phase may be too universal (e.g. C18)
- Extraction scheme is not specific enough
- Reduce concentration of Elution solvent
- Wash step is ineffective, increase organic concentration
- Ion Suppression try alternative phase or change wash %
- Try polymeric or mixed mode SPE
- Pre-rinse cartridge with elution solvent



TRICK: If In Doubt ... Mass Balance!

If you load ???ng of a sample onto a cartridge, tracking where it goes helps determine important method alterations

Method Development

- 1. Collect "effluent" that has passed through upon sample loading
- 2. Collection all wash steps
- 3. Collect all elution steps

Easy to Quantify if you have analyte still on the phase

If the matrix contains high levels of interference, mass balance can be difficult if analyte does not retain on the load step.





Sample Prep Summary

- In the ever demanding analytical chemistry environment, sample prep is becoming an essential component of good laboratory workflow
- The advanced sample prep techniques (SPE and Precipitation) offer clear ease of use and productivity advantages
- The quantification of endogenous materials and ion suppression is vital to getting the most from your sample
- Polymeric SPE and multi-suite methods is a key platform for biological fluids
- Agilent offers the broadest platform of chemistries in the widest range of formats in the industry
- This was just a short seminar. Tell us if you have a challenge or application need in your laboratory... we probably have a solution already!





Fragen???





Vielen Dank!



Agilent Sample Prep: Accuracy Starts Here

- Options for Your Needs
- Innovative Products Designed for Lab Efficiency
- Technical Support at Every Step
- World Class Manufacturing and Quality

Provide chromatographers with highly accurate, reproducible results and enable transfer of methods from lab-to-lab and around the world





Effect of Endogenous Interferences

- Poor Chromatography
- Mechanical issues (particulates, blockages)
- LC Column Lifetime Issues
- Carry Over
- Instrument Downtime
- Ion Suppression
- Overall loss in sensitivity
- Increase in sample run time/ cost





Effect of Endogenous Interferences

- Poor Chromatography
- Mechanical issues (particulates, blockages)
- LC Column Lifetime Issues
- Carry Over
- Instrument Downtime
- Ion Suppression
- Overall loss in sensitivity
- Increase in sample run time/ cost



Protein Precipitation with Centrifugation





- 4) Remove supernatant
- 5) Analyze supernatant, often after dry down and resuspension



Where to go for more information

www.agilent.com

Home of thousands of products and hundreds of applications notes to help you improve your laboratory productivity

Contact your local sales office or representative

INSERT LOCAL CONTACT INFO HERE!

If you require application support, internal seminars or training and samples, please let us know.



