



P.C.C. Update Instructions

Purpose

The memory update now enables the user to read rate of change of Abs with time, for kinetic reactions, direct from the printout.

The printout also provides average, minimum, maximum and Standard Deviation readings for each sample. For more information refer to the UV-Vis Sales Manual Update 4/4 dated July '83.

Action

The P.C.C. memory now uses 4 x 2176 proms. Thus, an extra 24-way I.C. socket (supplied in this kit) is required to be soldered to the main board, part number 02100896 00.

The service manual indicates the position of the fourth prom, although a note indicates that it is not used. Delete this note.

Remove the early version proms (3 off) and replace with the four (4) new proms (version 03 or later).

The part number of the fourth prom is 1410513203.



This information is subject to change without notice.



8510057000

Part Number: 8510057000

Edition 01/11

Issue 4

Printed in Australia

© Agilent Technologies, Inc. 1992, 1996,
1999, 2011

Agilent Technologies
679 Springvale Road
Mulgrave, VIC 3170