
Agilent Glass Spray Chambers for ICP-OES

User Guide

The spray chamber is a crucial component of the sample introduction system for the ICP-OES spectrometer. It plays a key role in determining the sensitivity, precision and wash-out achieved during analysis.

Agilent glass spray chambers now feature a torque-controlled nebulizer adaptor to enable users to achieve a consistent seal and reproducible nebulizer positioning within the spray chamber. The compression ferrule is also attached to the nebulizer fitting for increased ease of use. This document outlines the recommended procedures for installation and maintenance of glass spray chambers.

Spray chamber installation

Our spray chambers are supplied clean and ready to use. The spray chamber is shipped sealed and ready to use for optimum precision and wash-out performance. Before first use:

- 1 Gently turn counterclockwise the torque-controlled nebulizer adaptor to remove the sealing plug (Fig. 1)
- 2 Fully insert the nebulizer into the torque-controlled adaptor and secure it in place by gently turning clockwise the nebulizer adaptor until you hear a "click" (Fig. 2).


Figure 1. Sealing plug removal


Figure 2. Nebulizer installation

CAUTION

The torque-controlled nebulizer adaptor should NOT be tightened without a nebulizer or sealing plug inserted. Failing to follow this will cause the torque-controlled adaptor to over-compress the compression ferrule, making it difficult to insert the nebulizer, which may damage the nebulizer.

- 3 Remove the top and bottom protection caps and proceed to install the spray chamber on the instrument. Use the clamp supplied to secure the ball joint socket to the ball joint on the base of the torch.
- 4 Connect the drain line by inserting the PTFE connector over the drain tube at the spray chamber base. Connect the free end of the capillary tubing to the peristaltic pump tube for the waste.

CAUTION

Hydrofluoric acid (HF) should not be used with glass or quartz. Using any amount of HF will damage the product.

Removing the glass spray chamber

WARNING

Chemical Hazard

Acids, corrosive substances or other hazardous chemical residues might be present on the spray chamber components and can cause severe burns when they come into contact with the skin. It is essential that appropriate personal protective equipment always be worn when handling the spray chamber components and the cleaning chemicals. If acid or the cleaning chemicals contact the skin, wash off with copious amounts of water and seek medical attention immediately.

- 5 Loosen the clamp securing the spray chamber to the ball joint on the base of the torch and remove the spray chamber.
- 6 Disconnect the UniFit drain connector from the drain on the spray chamber base.
- 7 Loosen the torque-controlled nebulizer adaptor and remove the nebulizer from the spray chamber.

Cleaning the glass spray chamber

All spray chambers require regular preventative maintenance to ensure optimum performance. Use the following procedure to clean and maintain your spray chamber.

When signs of degradation in precision and/or detection limits are detected proceed to aspirate a 2.5% RBS-25 solution for 15 minutes followed by a thorough rinse with blank solution. However, if this is not effective, soak the spray chamber overnight in a 25% RBS-25 solution (4x dilution).

Droplets (beading) on the internal surfaces of the spray chamber, are a clear sign that maintenance is required as they are a visible indication of spray chamber instability, and they should be removed. An overnight soak in 25% RBS-25 solution (4x dilution) is recommended.


If you have been aspirating samples with high TDS levels, clean the spray chamber with a 5% nitric acid solution before using the detergent soak procedure outlined above.

TIP


If you routinely analyze both aqueous and organic, non-miscible solutions, the use of two spray chambers - one for aqueous solutions and one for organic solutions - will significantly reduce the instrument down-time due to spray chamber maintenance.

Re-ordering information

Description	Part number
Double-pass glass cyclonic spray chamber	G8020-68102
Single-pass glass cyclonic spray chamber	G8020-68101
Torque-controlled nebulizer adaptor, 1/pk	5005-0482
Clamp for torch with ball joint socket, 1/pk	G8000-64140
UniFit drain connector, 2.0 mm od, 0.86 mm id with 480 mm capillary tubing, 3/pk	G8010-80036

Glass double pass spray chamber p/n G8020-68102

Torque-controlled nebulizer adaptor p/n 5005-0482

This information is subject to change without notice

5390-0026
DE number: DE-002487

Part Number: 5390-0026

Edition 10/24
Issue 1

© Agilent Technologies, Inc. 2024

Agilent Technologies Australia (M) Pty Ltd
679 Springvale Road
Mulgrave, VIC 3170, Australia

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