



EQP Name: DISS.02.55_PQ_R182R0_EQP

Service Type: PQ

Company Name: _____

Customer Name/Title: _____

EQP Filename: DISS.02.55_PQ_R182R0_eqp.eqp

EQP Publish Date: January 23, 2025 11:42:58 AM

Print Date: January 23, 2025 11:49:44 AM

Test program as reflected in EQP agreement.

Table of Contents

Section	Page
Scope and Purpose	3
Outside Variance Values	11
1222818-R182R0	15
DISS PQ: Standard Test Suite	21
Report and Delivery Options	22
Selected Signature Options	22
Customer Approval	23
Legal Notice	24

Scope and Purpose

Overview

An Equipment Qualification Plan (EQP) documents the work performed during qualification services provided using the Agilent Automated Compliance Engine (ACE) software. The EQP must be approved before the work is performed (see customer responsibilities) and the results are reported in an Equipment Qualification Report (EQR), which documents the following.

- Tests performed during the qualification work
- Setpoints used (for example, wavelength, temperature, and flow settings, as applicable)
- Limits (acceptance criteria) applied to tests and setpoints
- Test descriptions and rationales for tests performed (these and test parameters are also reported EQR)
- Traceability of the Agilent qualification protocol

EQP PDF review documents support customer review and approval and are digitally produced directly from the electronic EQP contained within ACE. Two different EQP approval formats can be created.

- Concise: standard format, typically used for Agilent-recommended EQPs
- Tabular: expanded format, typically used for configured EQPs

Many customers approve the Agilent-recommended EQP in the concise format. The tests, setpoints, and limits applied are typically organized in a concise table format (see information below and the Operational Qualification (OQ) test suite). The table can also display which qualification tests can be configured, which Agilent designates as variance (see Glossary).

In a concise EQP approval document used for Agilent-recommended EQPs, the test, setpoints and parameters, and limits (as applicable) are displayed in a concise abbreviated format as shown below.

NOTE This section is used by all ACE techniques, but only LC examples are provided for illustrative purposes.

Key:	Fixed HPLC setpoints/limits	UHPLC setpoints (primary module)	Variance allowed
Test	Setpoints and Parameters	Limits	
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data	
Pump Flow Accuracy and Precision (Pump)	Flow Rate 1: 0.500 ml/minute Flow Rate 2: 5.000 ml/minute Flow Rate 2: 2.000 ml/minute (G4220B) Flow Rate 2: 1.500 ml/minute (with ULD kit)	Accuracy Ô3.00% Precision RSD Ô0.50%	

Understanding the Test Specification Section in Tabular Review Documents

Note: (1) Agilent is NOT responsible for test failures for O (out of range) setpoints and limits. (2) Limit ranges must be more tightly managed than setpoint ranges because they often reflect physical measurement limits and are directly linked to the testing method. Therefore *within range* user limits are subject to best effort repairs if they cannot be met. In particular, Agilent is not be responsible for test failures for limits tighter (more demanding or challenging) than the recommended values. (3) Changes from the Agilent-recommended EQP are also displayed in the Test Summary section of the EQR. (4) For hyphenated systems (for example, LC-LCMS), an EQP for each system needs to be approved.

Pump Flow Accuracy

Analytical Agilent and non-Agilent pump

W	Setpoint:	0.500	mL/min	L:	0.100	H:	5.000
	Setpoint:	0.200	mL/min				
	Limit:	<=	5.00	%	L:	5.00	H:

Injection Precision Additional FLD

This optional test is enabled for this EQP.

Agilent or non-Agilent injector with FLD detector						
Vial verification performed?	No					
Injection Volume on Column:	5	uL		L:	1	H: 5
Area RSD Limit:	<=	1.00	%	L:	1.00	H: 5.00
Height RSD Limit:	<=	1.00	%	L:	1.00	H: 5.00

Column Temperature Accuracy

This standard default test is disabled for this EQP.

Customer Responsibilities

To perform ACE instrument qualification services, requirements associated with instrument control and use of ACE software on site need to be understood and satisfied.

- Instrument control to perform the qualification tests.
- ACE software access by the field service engineer (FSE) to load the electronic EQP, run the qualification workflow, and publish the EQR.

ACE is designed to satisfy ALCOA+ data integrity requirements and provide qualification tests that are independent of the instrument control software. This Agilent qualification design philosophy facilitates a harmonized approach to analytical instrument qualification and ensures consistent application of validated algorithms (such as those used to perform gradient pump HPLC testing), independent of instrument and software manufacturer. Additional documentation is available to provide further information on these subjects.

Agilent recommends that the software routinely used in the lab to control the instrument is used by the FSE to perform qualification tests. Therefore, appropriate access to the instrument control software, such as the customer data system (CDS), needs to be set up for the Agilent representative before they arrive on site. For ACE software access on site, Agilent recommends installation of Network ACE to load the electronic EQP, run the qualification workflow, and publish the EQR.

Customers need to address the following responsibilities.

- EQP review
 - Review the EQP supplied (standard or configured).
 - Request changes to the EQP if required (to meet user requirements such as testing range of use) or approve the supplied EQP (that satisfies user requirements).
 - Return an electronic copy of the EOP to Agilent prior to qualification delivery.

- Agree and set up instrument control access to perform the work (provide a CDS logon for the FSE).
- Agree and approve ACE usage on site (Network ACE).
- For configured EQPs, maintain change control and revision history of the specific approved EQP.
- Safely store and archive an approved copy of this EQP.
- Review and sign the EQR when the qualification work is complete and reported.

NOTE: If Agilent representatives use customer software such as a CDS account to acquire test data during the qualification, they log off from the CDS account at the end of test acquisition. Account credentials created for FSE CDS access are a customer responsibility (that is, some labs disable the account when the work is complete).

Agilent Responsibilities

Agilent representatives follow an approved structured training process, which is represented by the training certificates included in the EQR. Additionally, any customer-specific training that is required to perform the work should be specified and complete before the work is scheduled and can be optionally added as a PDF attachment to an EQR, if required. Agilent representatives do the following.

- Deliver the approved services following the test programs described in the approved EQP;
- Provide a locked and e-signed EQR upon completion of the service;
- If requested, provide an optional ink-signed EQR CD to the customer.

Statement of Intent

The qualification is delivered according to the approved EQP. The standard EQP is followed unless a configured EQP has been requested and approved. If a configured EQP has been created and approved, Agilent defines variances as changes to the default recommended values (stated in the standard EQP) that fall within a well-defined range. These changes are considered to be within the intended use range of the system under test.

Customizations are (a) values that subject the system to limits that exceed the typical operational range or (b) additional tests that are not considered part of the core program required for completion of the selected service. Because user-defined setpoints and limits may exceed the operational envelope of the equipment, Agilent reserves the right to warrant conformance only to the closest variance value. The user is notified of this stipulation at EQP setup time and the qualification report (EQR) reflects this situation.

The service is delivered according to the terms and conditions stated in the corresponding service exhibit. After approval, it is recommended that this EQP be archived with the electronic EQP file.

General Statements on the Testing Program

The hardware OQ tests described within this EQP are designed to satisfy analytical instrument qualification (AIQ) requirements for regulated laboratories. The laboratory is responsible for documenting that analytical instruments are suitable for the analytical procedures they use and the decisions they make on the reportable values / quality attributes being measured. Diligent application of an instrument life cycle framework such as USP <1058> help define and satisfy this requirement.

Agilent-recommended EQPs are derived from Agilent's interpretation of the analytical instrument life-cycle framework defined in USP <1058> and authoritative literature requirements from sources such as ISO9001, FDA, GAMP, PIC/S, ASTM 2500, and others. User requirements and intended use can vary from laboratory to laboratory. Therefore, it is recommended that prior to EQP approval, laboratories determine if additional set points and/or optional additional tests are required to test their user requirements. Additional charges may apply to any changes. This can be accommodated in a configured or customized EQP. In line with the USP <1058> requirements, Agilent qualification services include Installation Qualification (IQ) and Operational Qualification (OQ). Typically, the range of instrument use should be tested in the OQ. Laboratories are responsible for addressing PQ requirements of the USP <1058> life cycle.

Agilent CrossLab Compliance uses a balanced selection of direct metrology and chemical-based tests to directly determine the performance of the systems without unnecessary reliance on graphical extrapolation, interpolation, inferred or derived results. For example, for HPLC qualification, direct metrology measurements are used to test pump flow rate accuracy and precision, and thermal-controlled column compartment and autosampler modules. Holistic chemical testing is used for the evaluation of the following critical instrument characteristics: linearity, precision, signal to noise, and carry over.

EQP Revision History

ACE qualification protocols follow a standardized software development life cycle within Agilent. The electronic EQP files used by ACE are managed under change control within this life cycle. Changes to the core electronic EQP are recorded in the Agilent EQP Revision History and are available on request. Change control processes vary between organizations; therefore, although Agilent maintains a revision history for electronic EQPs contained within ACE, it is recommended that, for configured EQPs, laboratories maintain EQP-independent revision history information appropriate to their own policies and procedural requirements.

CDS Software Pre-requisites for Hardware Qualifications

The software used to control instruments during hardware qualification needs to be suitable. Agilent recommends that software used during the qualification has been qualified within the qualification period specified by the customer's software qualification SOP.

Glossary

Agilent EQP

ACE:

Agilent Automated Compliance Engine or Agilent software used to provide qualification services

EQP:

Equipment Qualification Plan that documents qualification work to be performed

EQR:

Equipment Qualification Report that documents qualification results

CDS:

Customer data system

Network ACE:

ACE is installed on a network node within the laboratory LAN infrastructure

Local ACE:

ACE resides on an independent external drive connected to the LAN infrastructure

Alternative Method:

ACE is installed and run from a PC not connected to the laboratory network

Standard EQP:

EQP pre-defined by Agilent that includes default tests, setpoints, and limits

Configured EQP:

EQP that has been changed based on user requirements

Variance, Range:

The range of operating values for a test parameter within ACE (flow, temperature, and so on) that were defined and tested when the test was first developed and set up

Variance, W:

Change to a test parameter that is within the Agilent variance range

Variance, O:

Change to a test parameter that is outside the Agilent variance range

Variance, DW:

Dual reporting limits within the Agilent supported range

Variance, DO:

Dual reporting limits outside the Agilent supported range

Optional Test:

An optional additional test that is not included in the standard EQP

Additional Setpoint:

Additional measurement value for a test parameter (flow, temperature, and so on)

Qualification Life Cycle Stages

URS:

User Requirement Specification that documents intended instrument use

DQ:

Design Qualification that documents suitability of instrument for intended use

IQ:

Installation Qualification that documents the installation of the instrument

OQ:

Operational Qualification that tests the range of use

PQ:

Performance Qualification that tests the system under conditions of use

RQ:

Repair Qualification that tests the instrument after an instrument repair

PM:

Planned Maintenance for scheduled preventive maintenance

Reference Sources

ASTM:

American Society for Testing and Materials

FDA:

United States Food and Drug Administration

GAMP®:

Good Automated Manufacturing Practice

ISO:

International Organization for Standardization

ISO9001:

ISO Standard for Quality Management Systems

PIC/S:

Pharmaceutical Inspection Cooperation Scheme

USP:

United States Pharmacopeia

Outside Variance Values

The changes requested in the test definition are outside the Agilent Recommended specifications for this system. Agilent will not be under any obligation regarding the OQ testing results against user-requested limits and set-points that are more stringent than the Agilent recommended ones. Agilent responsibility will be limited to meeting the closest published allowable variance value that has been exceeded. Under no conditions will Agilent be held responsible for qualification test failures that result from these user defined limits which are therefore also excluded from the coverage of service contracts. This agreement statement will be reflected in the Equipment Qualification Report (EQR).

Technique : Diss

Prednisone Qualification - Single Stage

Prednisone Qualification - Single Stage1				
Value Name	Value used	Closest Variance Value	Recommended Value	Conditions
Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
% Coefficient of Variation Limit	4.4 %	4.6 %	4.6 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
% Coefficient of Variation Limit	4.3 %	4.5 %	4.5 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
Geometric Mean Low Limit	81 %	82 %	82 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits

% Coefficient of Variation Limit	4.2 %	4.4 %	4.4 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Geometric Mean High Limit	57 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
% Coefficient of Variation Limit	5.4 %	6.2 %	6.2 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Geometric Mean High Limit	57 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
% Coefficient of Variation Limit	5.3 %	6.1 %	6.1 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Geometric Mean High Limit	57 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
% Coefficient of Variation Limit	5.2 %	6.0 %	6.0 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits

Prednisone Qualification - Two Stage

Prednisone Qualification - Two Stage1

Value Name	Value used	Closest Variance Value	Recommended Value	Conditions
------------	------------	------------------------	-------------------	------------

Stage 1 Geometric Mean High Limit	91 %	90 %	90 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	3.3 %	3.4 %	3.4 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
Stage 2 Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	4.4 %	4.5 %	4.5 %	: Apparatus: 1 Number of vessels tested: 6 User Defined USP Limits
Stage 1 Geometric Mean High Limit	91 %	90 %	90 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	3.3 %	3.4 %	3.4 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
Stage 2 Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	4.3 %	4.4 %	4.4 %	: Apparatus: 1 Number of vessels tested: 7 User Defined USP Limits
Stage 1 Geometric Mean High Limit	91 %	90 %	90 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	3.3 %	3.4 %	3.4 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Stage 2 Geometric Mean Low Limit	81 %	82 %	82 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Stage 2 Geometric Mean High Limit	94 %	92 %	92 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	4.2 %	4.3 %	4.3 %	: Apparatus: 1 Number of vessels tested: 8 User Defined USP Limits
Stage 1 Geometric Mean Low Limit	47 %	48 %	48 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Stage 1 Geometric Mean High Limit	54 %	57 %	57 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	4.1 %	4.6 %	4.6 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Stage 2 Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits

Stage 2 Geometric Mean High Limit	57 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	5.4 %	6.1 %	6.1 %	: Apparatus: 2 Number of vessels tested: 6 User Defined USP Limits
Stage 1 Geometric Mean Low Limit	47 %	48 %	48 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 1 Geometric Mean High Limit	54 %	57 %	57 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	4.1 %	4.6 %	4.6 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 2 Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 2 Geometric Mean High Limit	57 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	5.3 %	6.0 %	6.0 %	: Apparatus: 2 Number of vessels tested: 7 User Defined USP Limits
Stage 1 Geometric Mean Low Limit	47 %	48 %	48 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Stage 1 Geometric Mean High Limit	54 %	57 %	57 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Stage 1 % Coefficient of Variation Limit	4.1 %	4.6 %	4.6 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Stage 2 Geometric Mean Low Limit	45 %	46 %	46 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Stage 2 Geometric Mean High Limit	56 %	59 %	59 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits
Stage 2 % Coefficient of Variation Limit	5.2 %	5.9 %	5.9 %	: Apparatus: 2 Number of vessels tested: 8 User Defined USP Limits



USP REFERENCE STANDARD CERTIFICATE

DISSOLUTION PERFORMANCE VERIFICATION STANDARD - PREDNISONE

(10 mg nominal prednisone content per tablet)

USP Catalog No.:	1222818
USP Lot No.:	R182R0

Valid Use Date	31-JUL-2025
Storage/Handling	As per the label.
Uses	General Chapter <711> Dissolution, Performance Verification Test (PVT), Apparatus 1 and Apparatus 2

Dissolution <711>

Medium:	499 g of degassed purified water maintained at $37^{\circ} \pm 0.5^{\circ}$
Medium degassing:	Recommended degassing procedure: Heat a suitable amount of water, while stirring gently to about $41-45^{\circ}$. Filter under vacuum through a $0.45\text{-}\mu\text{m}$ -porosity filter into a suitable filtering flask equipped with a stirring device. Seal the flask and continue to apply vacuum while stirring for an additional five minutes. Measured vacuum should be less than 100 mbar. Note: Other validated degassing methods that reduce the total dissolved gas in the media can also be used
Apparatus:	Apparatus 1 (Basket) or Apparatus 2 (Paddle) at 50 RPM Note: If equipment is dedicated for use with only one apparatus (basket or paddle), then performance verification is only required for that apparatus
Time:	30 minutes
Standard Solution:	A known concentration of USP Prednisone RS in Medium. Note: An amount of alcohol not to exceed 5% of the total volume of the standard solution may be used to bring the prednisone reference standard into solution.
Sample solution:	Laboratory can choose either Single-Stage Test or Optional Two-Stage Test scheme to obtain Sample Solutions. A filtered portion of the solution under test, suitably diluted, if necessary, with Medium to obtain a concentration similar to that of the Standard solution. Note 1: The filtering method must not cause adsorptive loss of drug Note 2: Bias introduced by automated methods is to be avoided
Analysis:	UV at maximum absorbance of about 242 nm
Procedure:	Determine the quantity of prednisone, $\text{C}_{21}\text{H}_{26}\text{O}_5$, dissolved at 30 minutes in each vessel expressed as percent of the labeled amount.

Copyright 2023 The United States Pharmacopeial Convention. All rights reserved.

USP Reference Standard Certificate Template No.: FORM-00008-02

Template Effective Date: 15-DEC-2023

Page 1 of 6



USP REFERENCE STANDARD CERTIFICATE

Single-Stage Test Instructions and Acceptance Criteria

1. For each position in the assembly, test one USP Dissolution Performance Verification Standard – Prednisone (DPVS – Prednisone) RS tablet, and record the percent dissolved at the sampling time point specified. Transform the percent dissolved results to the natural log scale, determine the mean and variance. For assemblies with 12 or 14 dissolution vessels, no further testing is required
2. For assemblies with fewer than 12 positions, repeat Step 1 with an additional set of tablets. Transform the percent dissolved results to the natural log scale, determine the mean and variance.
3. Calculate the average of the two means and of the two variances obtained in Steps 1 and 2.
Use the results from Step 1 alone for assemblies that have 12 or 14 positions.
4. Convert the results of Step 3 to a geometric mean (GM) and percent coefficient of variation (%CV). See Calculation Example for details.
5. Compare the results of Step 4 to the Single-Stage acceptance criteria in Table 1. The GM must not fall outside the limits, and the %CV must not be greater than the limit. If both meet the criteria, the assembly has passed the PVT

Table 1. Performance Verification Test Acceptance Criteria for Single-Stage Test

Apparatus	No. of vessels per run	Geometric Mean, % Prednisone Dissolved	%CV
1 (Basket)	6	81-94	4.4
	7	81-94	4.3
	8	81-94	4.2
	12	81-94	4.3
	14	81-94	4.3
2 (Paddle)	6	45-57	5.4
	7	45-57	5.3
	8	45-56	5.2
	12	45-57	5.4
	14	45-57	5.3

Copyright 2023 The United States Pharmacopeial Convention. All rights reserved.

USP Reference Standard Certificate Template No.: FORM-00008-02

Template Effective Date: 15-DEC-2023

Page 2 of 6



USP REFERENCE STANDARD CERTIFICATE

Optional Two-Stage Test Instructions and Acceptance Criteria

A laboratory may choose to implement the PVT as a Two-Stage test in case of assemblies with less than 12 positions. The Two-Stage test is a statistically valid means of allowing the possibility of stopping the test at the first stage using more stringent acceptance criteria. The following are step-by-step instructions for the two-stage test

1. For each position in the assembly, test one USP DPVS – Prednisone RS tablet, and record the percent dissolved at the sampling time point specified. Transform the percent dissolved results to the natural log scale, determine the mean and variance
2. Convert the results of Step 1 to a GM and %CV and compare to the 1st Stage of Two Stages acceptance ranges in Table 2. The GM must not fall outside the limits, and the %CV must not be greater than the limit. For calculation of the GM and %CV, see Calculation Example for details
3. If results of Step 2 satisfy both acceptance criteria, the assembly has passed the PVT. Otherwise continue to Step 4. Prior to proceeding to Step 4, see Futility Factor section.
4. Repeat Step 1 with an additional set of tablets. Transform the percent dissolved results to the natural log scale determine the mean and variance for the data obtained at this step
5. Average the two means and two variances obtained in Steps 1 and 4
6. Convert the results of Step 5 to a geometric mean (GM) and percent coefficient of variation (%CV). For calculation of the GM and %CV, see Calculation Example for details
7. Compare the results of Step 6 to the 2nd Stage of Two Stages acceptance ranges in Table 2. The GM must not fall outside the limits, and the %CV must not be greater than the limit. If both meet the acceptance criteria, the assembly has passed the PVT

Table 2. Performance Verification Test Acceptance Criteria for Two-Stage Test

Apparatus	No. of vessels per run	First Stage of Two-Stage Test		Second Stage of Two-Stage Test	
		Geometric Mean, % Prednisone Dissolved	%CV	Geometric Mean, % Prednisone Dissolved	%CV
1 (Basket)	6	83-91	3.3	81-94	4.4
	7	83-91	3.3	81-94	4.3
	8	83-91	3.3	81-94	4.2
2 (Paddle)	6	47-54	4.1	45-57	5.4
	7	47-54	4.1	45-57	5.3
	8	47-54	4.1	45-56	5.2



USP REFERENCE STANDARD CERTIFICATE

Futility Factor

If optional Stage-Two test is chosen, there are circumstances when the %CV after the First Stage of Two-Stage test equals or exceeds the value in the Futility Factor table (without rounding). In such cases it is impossible to meet the %CV criterion after the Second Stage of the Two-Stage test. The lab can stop after the First Stage run. However, after any adjustments to equipment, test procedure, and so on, the PVT must be restarted with a new first run (Step 1 of the two-stage test instructions).

Futility Factor, %CV at or above value given, second stage testing will not produce passing result

Apparatus	Number of Vessels		
	6	7	8
1	6.2	6.1	5.9
2	7.6	7.5	7.4

Refer to this website for the USP Calculation Tool: <https://apps.usp.org/app/USPNF/pvtCalculationTool/>

Calculation Example (expressed as Microsoft Excel® worksheet functions):

Run 1: $x_1, x_2, \dots, x_n$ in natural log scale: $\ln x_1, \ln x_2, \dots, \ln x_n$

Run 2: $x_{n+1}, x_{n+2}, \dots, x_{2n}$ in natural log scale: $\ln x_{n+1}, \ln x_{n+2}, \dots, \ln x_{2n}$

1st Stage of Two-Stage for n=6, 7, 8 and Single-Stage for n=12, 14:

$GM1 = \exp(\text{average}(\ln x_1 : \ln x_n))$

$\%CV1 = 100 * \sqrt{\exp(\text{var}(\ln x_1 : \ln x_n)) - 1}$

Single-Stage or 2nd Stage of Two-Stage for n= 6, 7, 8:

$GM = \exp(\text{average}((\text{average}(\ln x_1 : \ln x_n)), (\text{average}(\ln x_{n+1} : \ln x_{2n})))) = \exp(\text{average}(\ln x_1 : \ln x_{2n}))$

$\%CV = 100 * \sqrt{\exp(\text{average}((\text{var}(\ln x_1 : \ln x_n)), (\text{var}(\ln x_{n+1} : \ln x_{2n})))) - 1}$

exp: exponential (e^x)

var: variance

sqrt: square root

*: multiply

100: conversion factor to percentage

For more information and guidelines about how to complete the performance verification test refer to the following website: <https://www.usp.org/small-molecules/pvt>



Copyright 2023 The United States Pharmacopeial Convention. All rights reserved.

USP Reference Standard Certificate Template No.: FORM-00008-02

Template Effective Date: 15-DEC-2023

Page 4 of 6

 USP REFERENCE STANDARD CERTIFICATE

Label

**REFERENCE STANDARD**

**DISSOLUTION PERFORMANCE VERIFICATION STANDARD -
PREDNISONE 30 Tablets**

 The nominal weight of prednisone in each tablet is 10 mg. At the time of use, open the aluminum sachet, remove the blister card, and push the tablets through the foil backing of the blister card. Use only whole tablets. Store at controlled room temperature. Keep unused or unopened blister strips in the secondary package.

Danger! Causes eye irritation. Suspected of damaging fertility or the unborn child. Causes damage to organs (endocrine system) through prolonged or repeated exposure.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe dust/fume/gas/mist/vapors/spray. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If exposed or concerned: Get medical advice/attention. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

See certificate for any additional information.

For use with specified USP compendial tests. Not for use as a drug.
See SDS prior to use at www.usp.org/sds

USP, 12601 Twinbrook Pkwy, Rockville, MD, +1-301-881-0666
Cat. No. 1222818 Material mfd. in Spain

LOT: R182R0



Danielle A. Vattimo

Quality Assurance

Certificate Version History

Version Number	Date	Reasons for Change
00 (Current)	14-MAR-2024	First issue



USP REFERENCE STANDARD CERTIFICATE

Label

Reference Standard label typically contains the name, catalog number, lot number, package size, assigned value when applicable, storage conditions, handling instructions, and country of origin information. The label may also include hazard and precautionary statements required by the Occupational Safety and Health Administration (OSHA).

Assigned Value

For USP Reference Standards with compendial quantitative use(s), an assigned value is provided on the label and/or the Certificate.

For USP Reference Standards with compendial qualitative use(s), USP may choose to provide a value, e.g., chromatographic purity, for informational purposes in the Certificate, on a case-by-case basis.

Valid Use Date

It is the responsibility of the user to ascertain that a particular lot of a USP Reference Standard has official status either as a "Current Lot" or as a "Previous Lot" within the assigned valid use date. The online USP Reference Standards Catalog and the online USP Store at www.usp.org are updated daily. USP recommends referring to one of these sources prior to using a USP Reference Standard to make sure the lot is valid for use.

Storage

Storage conditions are lot-specific and may change from one lot to another. Storage conditions on the label and/or the Certificate are valid for unopened container as received. Once the container is opened, unless otherwise specified on the label and/or the Certificate, users are responsible for storing any remaining material according to their site procedures and ensuring continued suitability for its intended use. If no specific directions or limitations are provided on the label, conditions of storage include storage at room temperature and protection from moisture, light, freezing, and excessive heat. See General Chapter <659> in the USP-NF Online for storage and handling definitions.

Instructions for Use

Follow the instructions provided on the label and/or the Certificate and in the associated USP documentary standard(s). Please refer to General Chapter <11> for additional information.

Non-USP Compendial Use

USP Reference Standards are for use in analytical or laboratory applications generally as specified in USP compendia. They are not for use in humans or animals as drugs, food, or medical devices. It may be possible to use a USP RS outside of its associated USP compendial applications; however, it is the responsibility of the user to determine the suitability of the USP RS for a non-USP use.

LEGAL NOTICE

USP WARRANTS GOOD TITLE TO USP REFERENCE STANDARDS ON DISPATCH FROM USP. THE FOREGOING WARRANTY IS IN LIEU OF ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR ANY WARRANTY THAT THE PRODUCTS, INCLUDING THIS CERTIFICATE, ARE OF MERCHANTABILITY QUALITY. USP'S LIABILITY ARISING OUT OF OR RELATING TO THE SUPPLY OF USP REFERENCE STANDARDS AND THIS CERTIFICATE SHALL IN NO EVENT INCLUDE LOSS OF PROFITS, COST OF PROCURING SUBSTITUTE GOODS OR SERVICES, OR ANY INCIDENTAL, INDIRECT, OR CONSEQUENTIAL DAMAGES OF ANY KIND, EVEN IF USP IS AWARE OF THE POSSIBILITY OF SUCH DAMAGES. WITHOUT LIMITING THE GENERALITY OF THE FOREGOING, USP DOES NOT WARRANT THAT THE USE OR RESALE OF USP REFERENCE STANDARDS, INCLUDING THEIR USE TO PERFORM TESTS AND ASSAYS PUBLISHED BY USP, WILL NOT INFRINGE UNITED STATES OR ANY OTHER PATENTS.

USP Reference Standards are not intended for use as drugs, food, or as medical devices.

This certificate may not be reproduced without the express written permission of USP.



Standard PQ Test Suite

This document describes the test program for qualifying dissolution instruments, and the following tables list all PQ tests. PQ affirms that your dissolution apparatus performs in accordance with current USP requirements. The USP Dissolution Performance Verification Test (PVT) is performed as required in the current USP General Chapters, Dissolution <711> and Drug Release <724>, in accordance with cGMPs.

Note: The actual test limits are subject to change when a new prednisone standard is released. A custom EQP is available with limits for the current USP lot.

Test	Setpoints and Parameters	Limits
Standard Preparation	N/A	% Absorbivity $\geq 99.0\%$ and $\leq 101.0\%$ for working and control standards (calculations are performed <u>only</u> if the control standard is used)
Filter Validation	N/A	Recovery $\geq 98\%$ and $\leq 102\%$ for each filtered aliquot
Prednisone Qualification (All variations of test)	Vessel temperature: 37.0°C Elapsed time: 30 minutes (target window from tablet drop to sampling)	See current USP lot's Certificate of Analysis

Consumables, Supplies, and Parts Used for Qualification

All parts, supplies, standards and consumables specified by the Agilent qualification protocol are provided by the customer. Agilent does not provide Prednisone standards due to the impact of improper storage, which can adversely affect the potency and/or purity of the standards and put the integrity of the qualification at risk. Agilent will provide equipment necessary for the measurement of physical parameters (e.g., thermometer, level, tachometer, wobble gauge, etc.). Any additional parts for maintenance or repair needed to affect qualification will be billed to customer unless otherwise covered by Agilent service and support agreement.

Test Design and Rationale

PQ service does not include physical testing (e.g., measurements of speed, wobble, centering, level, etc.). Because the physical condition of the Dissolution tester can affect the outcome of the Performance Qualification, it is recommended that an inspection of the equipment and measurements of physical parameters are performed prior to the service to ensure that the equipment conforms to pharmacopeia requirements.

Standard Preparation

Description: This test describes how to prepare standards for the PQ tests.

Procedure: Refer to prednisone certificate for details

Filter Validation

Description: This test validates the filters used for sample and standard preparation.

Procedure: This test compares absorbance readings of three filtered aliquots of working standard.

Prednisone Qualification

Description: These chemical tests verify the performance of the dissolution tester.

Procedure: Refer to prednisone certificate for details.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2024
Published in USA

Report and Delivery Options

This list below includes options selected for inclusion in the EQR. By default, repeated run logs and transaction logs are always included; all other options can be included or excluded.

(For hyphenated system types only) If different options are chosen for the primary and supported system types, the primary system options are used for both techniques.

- Show chromatograms
- Show header and footer on cover
- Include repeated run logs
- Include Transaction logs

Selected Signature Options

Status: EQP is not locked

- Reporting variance is allowed in this EQP

Customer Approval

Name:

Title:

Date:

Signature:

Name:

Title:

Date:

Signature:

Name:

Title:

Date:

Signature:

Name:

Title:

Date:

Signature:

Legal Notice

Agilent CrossLab Compliance and its primary components (ACE software tool, procedures, test design, metrology tools, chemical reference standards, and operator training materials) have been designed, tested, validated, and released for commercial use following Agilent's Life Cycle Development Quality Assurance methodology. A review of the following life cycle stages is a key part of the ACE software and associated protocol development process.

- Product Description / Specification
- Life Cycle Phases / Transition Approval
- Quality Assurance / Testing
- Documentation and Change Management
- Source Code

Validation certificates for the ACE software and associated protocols are generated automatically as part of the life cycle process, prior to software release, and are available on request. The certificate includes a declaration of intent of interest statement, including Annex 11 and 21 CFR Part 11, and is representative of the life cycle process followed.

Agilent Associate VP, R&D | ACG - Services and Support Division: Greg Fiest, Wilmington, DE USA

Agilent Quality/Regulatory Director: Dina Riley, Wilmington, DE USA

ACE software is patented. Copyright is claimed by this statement for all original work comprising Agilent CrossLab Compliance. Any unauthorized use, reproduction, or translation will be prosecuted to the maximum extent possible by law. All customer copies of EQP approval, final qualification reports, and raw data provided to customer at delivery of the service become the property of the customer.