



9ED'BUa Y.	CE 3^} dÜ^&[{ { ^} â^â
GYfj 'JW'HmdY.	UÛ
7 ca dUbmBUa Y.
7 i glca Yf'BUa Y#~]hY.
9ED': 'YbUa Y.	ŠÔT ÛÈGE GÈ~]
9ED'Di V]g\ '8 UHY.	Ù^] c{ à^! AGÍ
Df]bh8 UHY.	Ù^] c{ à^! AGÍ AGÍ Ā KÍ KÍ ÁT

HUV`Y`cZ7cbhYblg

GYWqcb	DU`Y
Ù&[]^Áæ áÁŰ`!][•^	H
Ô! [••ŠæÄÖ[{] ææ &^	FF
œÔÒÁ^••q } Áæ áÁŰ^!çæ^Ö^ qŰ!^ Á^o@ á•Áæ áÁŰq æ^ÁŰ] ! [çæ	FÍ
œÔÒÖ^ qŰ!^ Á^o@ áÁŰ^ÁÖæ^•	Fİ
ŠÔT ÛÁÛŰKŰ æ çæqÁŰ&qŰ^•c{ •ĚT ÛĚŰ] ^&æÁ^•oŰ Ź^	FJ
Û^] [!oæ áÁÖ^ qŰ!^ ÁŰ q } •	Ğ
Û^ ^&c^áÁŰq } æ !^ÁŰ] q } •	Ğ
Ô~•q { ^!ÁŰ] ! [çæ	Ğ
Š^* qŰp[çæ^	GJ

Ù^] c^{ à^! ÁĜÊĖĠ Ä ĲÍ ĲË ÁÚ

U] 𐄂 } 𐄂Á^•dK

𐄂 Á[] 𐄂 } 𐄂Áāāā } 𐄂Á•••••Á[𐄂 &` āāāÁ Á@Á 𐄂 āāāÁŬŬŬ

𐄂āāā } 𐄂Á^d[] 𐄂 dK

𐄂āāā } 𐄂Á ^•~!^ { ^} 𐄂𐄂^ ÁÁ[Á••••• 𐄂 𐄂ÁÁÁ[, 𐄂{] ^! 𐄂^! 𐄂 āÁ[Á] D

Ei U]Z]WU]cb`@Z'7mWY'GHU] Yg

WŬŬK

W•^!ÁŬ^~ ā^ { ^} 𐄂Á ^āāā } Á•••••[& { ^} •Á 𐄂} āāā•d~ { ^} 𐄂Á•^

ÖŬK

Ö••ā } ÁŬ~ 𐄂āāā } Á•••••[& { ^} •Á~ āāā Á -Á•d~ { ^} 𐄂Á[Á 𐄂} āāÁ•^

ŬK

𐄂•𐄂|𐄂 } ÁŬ~ 𐄂āāā } Á•••••[& { ^} •Á@Á••𐄂|𐄂 } Á -Á•d~ { ^} c

UŬK

U] ^!𐄂 } 𐄂ÁŬ~ 𐄂āāā } Á•••••••Á@Á 𐄂^ Á -Á•^

ÚŬK

Ú^!+!{ 𐄂 & ÁŬ~ 𐄂āāā } Á•••••••Á@Á••{ Á} ā^!Á[} āā } •Á -Á•^

ÜŬK

Ü^] 𐄂ÁŬ~ 𐄂āāā } Á•••••••Á@Á•d~ { ^} 𐄂Á•Á[Á] •d~ { ^} 𐄂Á] 𐄂

ÚT K

Ú|𐄂 } ^āÁ 𐄂 𐄂} 𐄂 & Á[Á&@ā`|āÁ!^•^ 𐄂^Á 𐄂 𐄂} 𐄂 &

FYZfYbW'Gci fWg

𐄂ŬVT K

𐄂 ^!𐄂 Á[& c Á[Á^••• * Á 𐄂Á 𐄂^! 𐄂•

ØØK

W 𐄂ÁŬ𐄂•Á[[āÁ āŬi~ * Á{ 𐄂 ā d𐄂 }

Ö𐄂 Út K

Ö[[āÁ 𐄂 { 𐄂āÁ 𐄂~ 𐄂 𐄂 * Á! 𐄂

ŬUK

𐄂 𐄂} 𐄂 } 𐄂Á!~ 𐄂 ā 𐄂 } Á[Á] 𐄂 āāā 𐄂 }

WUJ€€FK

WUÁUœ åæåÁ{ !Á~ ææ Á æ æ ^{ ^} áU•c{ •

ÚŦŦK

Úœ{ æ^ œæŦ•] ^&Ŧ } Á[[] ^!æŦ } Á&@{ ^

WÚÚK


Wj æåÁUœ•ÁÚœ{ æ[] ^æ



Agilent CrossLab Compliance Services

Agilent CrossLab Compliance Services are designed to seamlessly integrate with traditional quality systems used by firms and recognized by regulatory agencies worldwide. Analytical instruments must be suitable for their intended use. This requirement is good science in all laboratories and a regulatory requirement in pharma and biopharma laboratories. A life-cycle process for documenting and testing the suitability of laboratory instruments should be followed and Agilent recommends the life cycle framework defined in USP General Chapter <1058> on Analytical Instrument Qualification (AIQ). USP <1058> defines the governing framework and requirements that need to be satisfied, but the laboratory is responsible for how they satisfy these requirements.

- The United States Pharmacopeia (USP) is the only major pharmacopeia with a general chapter dedicated to analytical instrument qualification, making <1058> an important global regulatory reference. The information is provided in a scientific, risk-based approach to analytical instrument qualification (AIQ). However, the life-cycle framework contained within USP <1058> is not prescriptive in its implementation, making the embedded scientific and risk-based principles flexible and universally applicable.
- The scientific process followed by CrossLab uses the Agilent's Automated Compliance Engine (ACE) to deliver paperless electronic qualification. The life-cycle stages Agilent perform are highlighted in the life-cycle diagram below. As part of this life-cycle, Agilent can configure the qualification tests performed to align with user requirements.



USP <1058> AIQ Framework

NOTE: RQ services, described later in this document, can be added to standard qualification services.

ACE Workflow and Equipment Qualification Plans (EQPs)

Overview

Within the ACE workflow, the qualification tests, setpoints, and limits are defined in an EQP that can be configured to ensure that testing satisfies user requirements. When the qualification work is complete, an Equipment Qualification Report (EQR) is issued. The electronic workflow used within ACE has significant data integrity advantages over traditional paper or Excel-based qualification protocols, as validated calculations can be performed directly using electronic data such as chromatograms and metrology test values. Several of the instrument life-cycle stages are the responsibility of the laboratory, Agilent can provide compliance consultancy services and documentation which can help customers satisfy these requirements. These additional services are not included in our typical qualification offering.



High-level ACE Qualification Workflow

Standard and User-defined Limits

(Hardware qualifications only)

EQPs are available for download and approval as standard documents with Agilent recommended tests, setpoints, and limits, or they can be electronically configured by approved personnel to align with user requirements and intended range of use requirements. The degree of configuration depends on the analytical technology, but most EQPs can be configured to some degree, and one feature that can typically be changed is test limits.

EQPs are designed to be configurable (dependent on the analytical technology and standard requirements), but including additional tests or setpoints can impact the qualification time and associated cost. If a test limit is changed, ACE includes the capability to report results against the Agilent approved limit and any customer required limits (that is, both can be reported simultaneously).

If a user-defined test limit is more stringent than an Agilent recommended limit, Agilent makes no guarantee or obligation regarding the instrument passing the tighter test specification requirements. It is important to appreciate that tests performed under conditions of use (that is, to satisfy pharmaceutical monograph and application requirements) can have different limits than those defined in the QQ. It is the continuum of the combined QQ, PQ, and any point of use testing performed each time the instrument is used that together satisfy regulatory requirements.

User Requirements Specification (URS)

The purpose of user requirements is to document the intended use of the instrument within the life-cycle process and quality management system (QMS) being followed. Therefore, the URS is a customer / laboratory responsibility. Defining user requirements is often used to guide the customer in instrument selection and is stated as the first activity that should be followed in <1058>. The URS is important for two main reasons.

- It is a regulatory requirement for FDA and EU GMP that the intended use of the instrument and any software must be specified.
- Investment protection perspective means getting the right instrument for the right job.

Qualification protocols should test the instrument against any limits or specifications listed in the URS, which should document the intended range of use. Depending on the instrument complexity and how it is classified, a separate URS document may not be needed, but the URS requirements of the <1058> framework must be satisfied. A separate URS is almost always recommended for computerized systems.

An instrument performance specification is a product of the instrument development process by the supplier. It typically documents the performance the instrument can achieve. The URS should be based on intended use of the instrument and not the instrument specification. Additionally, if the intended use of a system changes, this may trigger a need to review the URS and associated qualification testing (for example, to ensure range of use is tested if used with a new analytical procedure).

Agilent offers compliance consultation services and documentation that can help customers address URS requirements.

Design Qualification (DQ)

The main function of the DQ stage of the laboratory instrument life-cycle process is to document why the selected instrument is suitable. Typically, this includes consideration of the instrument specification, how the instrument will be qualified, and the QMS followed by the instrument manufacturer. All together, these confirm that instrument performance is capable of satisfying user requirements. Depending on laboratory instrument life-cycle policy or SOPs being followed, instrument requirements and the relationship between the URS and DQ stages may vary – but as long as the <1058> framework principles are satisfied, this is not a problem, as it is left to each laboratory to justify and document its specific approaches.

The responsibility for satisfying DQ requirements primarily lies with the laboratory, with support from the supplier.

Agilent's approach to satisfying DQ requirements of USP <1058> includes the following.

- All Agilent hardware and software laboratory products, including the ACE software used to deliver qualification services, are designed, manufactured, and tested according to Agilent internal quality life-cycle development procedures.
- Certificates of Agilent testing, validation, and conformance to standards are provided with new Agilent instruments and similar certification can be provided for ACE software.
- Agilent is capable of installation, support, preventive maintenance, on-going qualification, and re-qualification after repair and user training worldwide.

Agilent offers a compliance consultation service that can help customers with DQ documentation.

Installation Qualification (IQ)

The main functions of the IQ stage are to document that laboratory is suitable (for example, critical systems typically include a site inspection / checklist), that the instrument is installed correctly in the environment, and IQ checks such as module start up are completed. IQ is provided and automated by ACE, which collects, checks, and tests Agilent hardware and software products for the following.

1. Purchase Order Details: Allows the customer to verify that the instrument being qualified matches their design requirements (if available) and purchase order.
2. Preparation and Installation Details: Gathers and records information about preparation and installation documents.
3. Documentation: Gathers and records information about reference and user manuals for initial installations.
4. Product Quality Assurance Details: Collects and records certificates and other forms that verify that the vendor has developed and built the product according to internal standards.
5. Startup: Verifies that all modules/components start up properly.
6. Installation Verification (software only): Verifies the correctness of all installation-related files.

Operational Qualification (OQ)

The main function of the OQ stage is to evaluate and document instrument performance at the intended operational range of use. OQ protocols should include a mix of metrology, functional, and operational tests. ACE qualification protocols include information about the test description and rational, setpoints, and the limits (acceptance criteria) for each technique, category, and instrument configuration.

OQ is provided and automated by ACE. ACE checks and tests for Agilent hardware and software products include the following.

- Metrological tests such as flow, temperature, pressure, and so on that ensure that the system is performing within Agilent (or user) specifications.
- Qualification results are reported in the EQR, which can include details of all test certificates, standards, and training information for the engineer performing the work. (Note that the EQR can be configured to customer requirements.)
- System or "holistic" tests verify the combined functions of the various system components
- The qualification testing can be configured to ensure URS requirements, such as range of use are tested.

For software qualification, the OQ consists of automated diagnostics regression testing and verification of the software installation. This supports continued use of the software in regulated environments (at install and as part of supporting periodic review).

In line with regulatory requirements, the EQPs should be approved before work is performed and the EQR should be reviewed and approved when the work is complete (as illustrated in Figure 2). The EQR contains all the raw data, results, and relevant information and attachments for complete compliance and traceability.

Mechanical Qualification (MQ)

(Dissolution systems only)

The main function of the MQ stage is to document that the mechanical performance of the instrument meets specifications and is functioning properly.

Performance Qualification (PQ)

The main function of the PQ stage is to document that the instrument is fit for purpose under conditions of intended use and to create an approved framework that ensures the instrument continues to perform as required. Because instrument range of use is tested within the OQ stage, it is usually not necessary to test this during PQ. It should be noted that requirements for instrument maintenance and repair fall within the PQ life cycle stage within the USP <1058> framework, as they are components of ensuring the continued performance of the instrument.

The customer is responsible for satisfying PQ requirements. (NOTE: Agilent can provide a PQ for Dissolution systems only.)

It is important to note that PQ is a lifecycle activity and not a one-time event. PQ tests may include activities such as method validation or system suitability tests (SST), but in Agilent's opinion, SSTs contribute towards ensuring continued performance of the instrument (that is, PQ testing), but do may not fully satisfy <1058> PQ requirements.

Repair Qualification (RQ)

After an instrument is repaired, tests should be performed to evaluate the effectiveness of the repair and document that repaired instrument satisfies performance requirements. Agilent offers a service called Repair Qualification (RQ), which refers to the requalification of laboratory instrument hardware after a repair. For some laboratory systems, to document the performance after repair may require a full OQ. However, for some modular or component-based systems, such as HPLC and GC for example, partial qualification testing can be justified. This is accomplished by performing the qualification tests that are applicable to only the module or system component related to the repair, reducing the time the instrument is out of service. Re-qualifying the instrument after repair is a regulatory requirement defined in USP <1058>.

Agilent offers service contracts to repair and requalify an instrument during the period between scheduled annual OQs.

The level of retesting is prescribed in the RQ section of ACE: a form is displayed for the operator showing all types of repairs possible and the retesting required. Part of an example form for an LC system is shown below.

Re-Qualification After Repair		
Pump Strategies		
Repair/Replace Strategy	Modules	OQ Testing
Internal pump head parts, active inlet valve (or AIV cartridge), (parts of) check valves, reference valves, inlet manifold or pump drive, or taking pump head apart to clean (versus repair)	Any pump	Flow Accuracy & Precision
Pulse damper, pressure transducer	Any pump	Flow Accuracy & Precision
Multi-channel gradient valve	Quaternary	Flow Accuracy & Precision Gradient Composition

The full list of RQ repair and retest guidance is available for customer review.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2025
Published in USA



Overview

Agilent recommends use of **Network ACE** for CrossLab qualification services that are enabled using the Agilent Automated Compliance Engine (ACE) software. Network ACE and Local ACE both access data directly (default methods) and are considered equivalent from a data integrity and data traceability perspective (see below). To provide additional flexibility in qualification service delivery, an alternative method is also available that accesses data indirectly. Use of the alternative method requires customer pre-approval using this form.

Available Methods

Method	Definition
Network ACE (Agilent recommended)	ACE software is installed on a network node within the laboratory LAN infrastructure. Raw data locations are always captured in the equipment qualification report (EQR), which provides end to end traceability and a fully characterized data workflow in the delivery. This method requires collaboration with the customer to load ACE behind the customer firewall.
Local ACE	ACE software resides on an independent external drive that can be driven from the system controller, where the customer data system (CDS) resides. Because the external drive is connected to the CDS, the data integrity of this method is equivalent to that of the Network ACE delivery method. Raw data is imported directly into ACE by the Data Manager tool, with the data paths always captured in the report, which provides data traceability assurance.
Alternative (Requires pre-approval)	<p>This method requires customer pre-approval due to data integrity implications. <u>Only choose this option in scenarios like the following:</u></p> <ul style="list-style-type: none"> ACE software is not run from a PC directly connected to the customer CDS, such as the FSE's laptop. System data files are transferred indirectly from the CDS to the FSE laptop instead of directly as done with Network and Local ACE methods. Data is acquired using a CDS on the FSE's laptop and transferred directly to ACE. The CDS used in this method is qualified for data collection purposes.

EQR Storage

Select the checkbox below to authorize Agilent to store copies of the EQRs generated by ACE for Agilent internal assessments. The intention of the assessment is to evaluate the delivery of the qualification service, with a focus to improve delivery and assess the appropriateness of data integrity measures. The storage is exclusively for the internal assessment by Agilent and is not shared with other organizations. It is not to be considered a backup for the EQR provided at qualification delivery.

Customer Approval of Alternative Method and EQR Storage

Authorize Agilent to use the alternative method (check for approval): ☐

Authorize Agilent to store EQRs for their internal assessment (check for approval): ☐

Approved By/Title:

Date Approved:

Comments:

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.
© Agilent Technologies, Inc. 2025
Published in USA



Introduction

With heightened scrutiny of data integrity, the Agilent Automated Compliance Engine (ACE) software must be able to access instrument-generated raw data files one of two ways: directly, using the connection between network nodes or with the server; and indirectly, through storage in a secure transfer location. (In this document, data integrity refers to the who, what, and where of data used in generating an ACE equipment qualification report, or EQR.)

ACE includes three main service delivery methods that address data integrity requirements; the rest of this document provides details to determine which one best fits a customer's needs.

Regardless of the delivery method, ACE features and delivery procedures are compatible.


Network ACE (Agilent Recommended)

Overview



ACE software is installed on a network node within the laboratory LAN infrastructure, which requires collaboration with the customer to load ACE behind their firewall. Raw data locations are always captured in the EQR, which provides end-to-end traceability and a fully characterized data workflow in the delivery.

Details



Typical Network ACE installation diagram

Installing ACE in a separate node (a.k.a. the host PC) on the same network as the system controller offers data traceability that is equivalent to an installation on the system controller itself. The system controller (where the CDS resides) and the ACE host PC are identified and seen by the server and subject to the customer's data access controls and general IT policies. The CDS's audit trail records data movements between nodes or between the client and server, and ACE's data traceability features identify the original data directory and therefore ensures end-to-end data traceability

The ACE host PC has a separate/partitioned drive for ACE software. During ACE's installation, two services are setup on the operating system (OS): one for security and the other as a watchdog. Because the ACE host PC sits on the network as a shared drive, engineers access ACE through the networked drive: ACE is not installed on ACE Virtual Viewer PCs.

Requirements

Installation

- Install on a host PC with a separate drive (different from that of the OS)
- Attach to a network that clients can access
- 500 GB
- NTFS format
- User has local administration rights
- Customer installation instruction document is available

Operational

- User has an ACE node logon with a minimum of power user rights permissions; user also has a personal ACE account and password added through the ACE licensing tool
- Up to 5 users with 3 open sessions each can access the NDA simultaneously
- Exception to ports 11121-11141 on ACE node, clients, and switch's/Smart Hubs to be open on the network

Local ACE

Overview



ACE software resides on an independent drive that can be driven from the system controller, where the CDS resides. Because the drive is connected to the CDS, this method's data integrity is equivalent to preferred 1 method's. Raw data is imported directly into ACE by ACE's Data Manager tool, and data paths are captured in reports to provide data traceability.

ACE software resides on an independent drive that can be driven from the system controller, where the CDS resides. Because the drive is connected to the CDS, this method's data integrity is equivalent to the Network ACE method. Raw data is imported directly into ACE by ACE's Data Manager tool, and data paths are captured in reports to provide data traceability.

Details

ACE is designed to run from a dedicated drive, without leaving a footprint on the host PC. Therefore, it can be connected directly to the system controller (where the CDS resides) without altering the system's qualification status. For additional protection, the drive can be driven by another host PC on the same network; also, the drive can remain on site with the customer for use by the Agilent Field Service Engineer (FSE) during service deliveries only.

Alternative Method

This method requires customer pre-approval due to data integrity implications and only applies in scenarios like the following:

- ACE software is not run from a PC directly connected to the customer CDS, such as the FSE's laptop. System data files are transferred indirectly from the CDS to the FSE laptop instead of directly as done with Network and Local ACE methods.
- Data is acquired using a CDS on the FSE's laptop and transferred directly to ACE. The CDS used in this method is qualified for data collection purposes.

www.agilent.com/chem/qualification

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2025

Published in USA



Overview

This document describes the test program for qualifying Agilent LCMS systems. The following tables list all MS source-specific OQ tests. Source abbreviations used in this document include:

AJST:	AP-ESI with Agilent Jet Stream Technology	ES:	AP-ESI
CI:	AP-CI	MM:	Multimode AP-ESI and AP-CI modes
DSAJST:	Dual spray AJST	MMCI:	Multimode AP-CI mode only
DSES:	Dual spray AP-ESI	MMES:	Multimode AP-ESI mode only

Note: Signal to Noise test is a standard test ONLY when the IDL test – included in the standard OQ for specific SQ, TQ, and Q-TOF models with some sources – is not scheduled; when IDL is scheduled, the Signal to Noise can be run as an optional EXTRA COST TEST; contact Agilent to request a custom EQP.

Key: Fixed setpoints/limits Variances allowed

Single Quadrupole: Standard OQ Test Suite

For new G6160B and G6170A models introduced in LcMs.02.62, test methodology was changed to use a C18 column instead of a restriction capillary. This new methodology includes IDL in the standard test suite.

Initial Methodology

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data (except when using customer CDS and no logon required)	Any
Scan Verification (pos. and neg. modes)	Model-specific; see Appendix: Masses for Scan Verification section	Accuracy for each mass ≤ 0.2 m/z	AJST*, CI**, ES*, MM**
Response Linearity	Injection volume on column: 5.0 μ L (ES, CI, MM)	Coeff. of det. (r2) ≥ 0.98000	Any
Injection Precision	1.0 μ L (AJST)	Area RSD $\leq 10.00\%$ Height RSD $\leq 10.00\%$	Any
Injection Carry Over	Evaluated mass is 311 m/z	Area carry over $\leq 1.00\%$ Height carry over $\leq 1.00\%$	Any
Signal to Noise		Signal to noise ≥ 200 Signal to noise ≥ 400 Signal to noise ≥ 100 Signal to noise ≥ 20 Signal to noise ≥ 100	AJST G6130B/C, G6135B/C AJST G6150B AJST Others CI, MM(CI) ES, MM(ES)

* Only sources supported for C models.

** Calibrant import restrictions may prevent test from being performed in some regions (China, for example).

New Methodology

G6160B

|

G6170A

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data (except when using customer CDS and no logon required)	Any
Scan Verification (pos. and neg. modes)	Model-specific; see Appendix: Masses for Scan Verification section	Accuracy for each mass ≤ 0.2 m/z (masses 1-2000 m/z) Accuracy for each mass ≤ 0.3 m/z (masses > 2000 m/z)	AJST G6170A ES, CI*, MM G6160B, G6170A
Response Linearity	Injection volume on column:	Coeff. of det. (r2) ≥ 0.98000	AJST G6170A ES, CI, MM G6160B, G6170A
Injection Precision	0.5 µL (AJST) 5.0 µL (ES, CI, MM-CI)	Area RSD ≤ 10.00% Height RSD ≤ 10.00%	AJST G6170A ES, CI, MM G6160B, G6170A
Injection Carry Over	1.0 µL (MM-ES)	Area carry over ≤ 1.00% Height carry over ≤ 1.00%	AJST G6170A ES, CI, MM G6160B, G6170A
Signal to Noise	Evaluated mass is 311 m/z	Signal to noise ≥ 2000 Signal to noise ≥ 1000 Signal to noise ≥ 500 Signal to noise ≥ 100	AJST G6170A** ES G6160B**, G6170A** CI, MM(ES) G6160B, G6170A MM(CI) G6160B, G6170A
Instrument Detection Limit (IDL)	Injection volume on column: 0.5 µL (AJST) 5.0 µL (ES)	Area RSD ≤ 10.00% (Equivalent to IDL ≤ 119.92 fg, derived from RSD) Area RSD ≤ 10.00% (Equivalent to IDL ≤ 1199.2 fg, derived from RSD)	AJST G6170A ES G6160B, G6170A

* Calibrant import restrictions may prevent test from being performed in some regions (China, for example).

** Optional in addition to IDL.

Tandem Quadrupole: Standard OQ Test Suite

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data (except when using customer CDS and no logon required)	Any
Scan Verification (pos. and neg. modes, Q1 and Q3)	Model-specific; see Appendix: Masses for Scan Verification section	Accuracy for each mass ≤ 0.2 m/z (masses 1-2000 m/z) Accuracy for each mass ≤ 0.3 m/z (masses > 2000 m/z)	AJST G6460A/C, G6465A/B, G6470A/B, G6475A, G6490A, G6495A/B/C/D ES G6410A/A-2K/B, G6420A, G6430A, G6460A/C, G6465A/B, G6470A/B, G6475A, G6495D CI* G6465A/B, G6475A, G6495D MM* G6475A
Response Linearity	Injection volume on column:	Coeff. of det. (r2) ≥ 0.98000	Any
Injection Precision	5.0 µL (ES, DSES, CI, MM)	Area RSD ≤ 10.00%	Any
Injection Carry Over	1.0 µL (AJST) 0.5 µL (AJST w/G6490A, G6495A/B/C/D, G6470B, G6475A)	Area carry over ≤ 1.00% Height carry over ≤ 1.00%	Any
Signal to Noise	0.5 µL (MM w/G6475A) Evaluated mass is 156 m/z	Signal to noise ≥ 2000 Signal to noise ≥ 1000 Signal to noise ≥ 100 Signal to noise ≥ 100 Signal to noise ≥ 20 Signal to noise ≥ 1000 Signal to noise ≥ 400 Signal to noise ≥ 100 Signal to noise ≥ 100	AJST G6470A, G6470B**, G6475A**, G6495A, G6495C**, G6495D** AJST G6465A/B, G6490A, G6495B AJST Others CI G6465A/B, G6470A/B, G6475A, G6490A, G6495A/B/C/D CI Others ES G6470B, G6475A, G6495D ES G6465B, G6490A, G6495A/B/C ES Others MM(CI) G6470A, G6475A, G6490A, G6495A/B/C

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
		Signal to noise ≥ 20	MM(CI) Others
		Signal to noise ≥ 400	MM(ES) G6470A, G6475A, G6490A, G6495A/B/C
		Signal to noise ≥ 100	MM(ES) Others
Instrument Detection Limit (IDL; for 2019 models and later)	Injection volume on column: 0.5 µL	Area RSD ≤ 10.00% (Equivalent to IDL ≤ 5.996 fg, derived from RSD)	AJST G6495C/D, G6470B, G6475A

* Calibrant import restrictions may prevent test from being performed in some regions (China, for example).
** Optional in addition to IDL.

TOF: Standard OQ Test Suite

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data (except when using customer CDS and no logon required)	Any
Scan Verification (pos. and neg. modes)	Model-specific; see Appendix: Masses for Scan Verification section	Accuracy for each mass ≤ 3.0 ppm	AJST, DSAJST, DSES, ES
Response Linearity	Injection volume on column: 20.0 µL (ES, DSES, CI, MM) 5.0 µL (AJST, DSAJST)	Coeff. of det. (r2) ≥ 0.98000	AJST, DSAJST, DSES, ES
Injection Precision		Area RSD ≤ 10.00%	Any
Injection Carry Over		Area carry over ≤ 1.00% Height carry over ≤ 1.00%	AJST, DSAJST, DSES, ES
Signal to Noise	Evaluated mass is 311.08 m/z	Signal to noise ≥ 100	

Q-TOF: Standard OQ Test Suite

Test	Setpoints and Parameters	Limits	Source(s) and/or Model(s)
CDS Logon Verification	N/A	Evidence of logon used to collect qualification data (except when using customer CDS and no logon required)	Any
Scan Verification (pos. and neg. modes, two zones)	Model-specific; see Appendix: Masses for Scan Verification section	<u>Quad</u> Acc. for each mass ≤ 0.5 m/z Acc. for mass > 2000 m/z ≤ 1.0 m/z (G6575A only) <u>TOF</u> Acc. for each mass ≤ 3.0 ppm	AJST, DSAJST, DSES, ES
Response Linearity	Injection volume on column: 20.0 µL (ES, DSES, CI, MM) 10.0 µL (AJST, DSAJST)	Coeff. of det. (r2) ≥ 0.98000	AJST, DSAJST, DSES, ES
Injection Precision		Area RSD: ≤ 20.00%	AJST, DSAJST, DSES, ES, CI, MM
Injection Carry Over		Area carry over ≤ 1.00% Height carry over ≤ 1.00%	AJST, DSAJST, DSES, ES
Signal to Noise	10.0 µL (ES, DSES, CI, MMCI w/G6575A) 1.0 µL (AJST, DSAJST w/G6530C, G6545A/B, G6546A, G6549A, G6550A/B, G6560C, G6575A) 1.0 µL (MMES w/G6575A) Evaluated mass is 156.08	Signal to noise ≥ 100 Signal to noise ≥ 20 Signal to noise ≥ 10	AJST, DSAJST G6550A/B AJST, DSAJST G6575A*, Others DSES, ES
Instrument Detection Limit (IDL; for 2023 models and later)	Injection volume on column: 1.0 µL	Area RSD ≤ 20.00% (Equivalent to IDL ≤ 599.600 fg, derived from RSD)	AJST, DSAJST G6575A

* Optional in addition to IDL.

Vacuum Verification (LcMs.02.59 and Earlier)

(Removed in LcMs.02.60: successful tune and Scan Verification verify correct system vacuum.)

MS Type	High Vacuum Limits
SQ	≥ 2E-6 torr and ≤ 2E-5 torr
TQ	≥ 0 torr ≤ 6.0E-5 torr
TOF	≥ 0 torr and ≤ 5E-7 torr
Q-TOF	≥ 0 torr ≤ 8.5E-5 torr (Quad), ≥ 0 torr ≤ 5E-7 torr (TOF)

Test Design and Rationale

Overview

Many GMP/GLP enforcement agency inspectors now ask firms to provide a risk assessment of their equipment and computer systems plus a science-based rationale for subsequent validation and qualification testing.

GENERAL RISK STATEMENT: Any laboratory chemical system used for raw material testing or final drug product / medical device testing in GMP or used in formal GLP studies will likely fall into a HIGH RISK category. This risk assessment will imply the need for IQ & OQ & on-going qualification. ANY USER SPECIFIC RISK ANALYSIS SUPERCEDES THIS GENERAL RISK STATEMENT.

The rest of this section outlines the science-based rationale for each test in the Agilent hardware OQ plus a brief test design and procedure description.

The recommended set of hardware OQ tests described in this EQP derives from Agilent's interpretation of FDA, USP, and GAMP guidelines and other authoritative expert literature.

OQ test design incorporates both modular and holistic testing, which is a proven and regulatory acceptable approach. Direct metrology is used to test pump flow rates and thermal-controlled column compartment and autosampler modules. Holistic chemical testing is used for the evaluation of the following critical instrument characteristics: linearity, precision, signal to noise, and carry over.

Certified reference standards and calibrated traceable thermometers and digital flowmeters are used.

Considering the number of setpoints, parameters, and conditions of each recommended OQ test, the proven concepts of worst case, range, and representative have been applied. If a property or characteristic is known to have its worst performance at one end of a range of use, this is the setpoint that should be tested and other setpoints are not required. If a property or characteristic has no known worst case, testing at the high and low points of the range of use is required. If there are too many possible use cases and conditions to realistically test (and none is a worst case), a representative sample for test is the best approach.

CDS Logon Verification

Description: To satisfy the attributable requirement of ALCOA+, evidence of the logon used to collect data must be provided unless a customer CDS is used and no logon is required.

Procedure: The test uses a screen capture to document who is logged on to the software that controls the instrument being qualified. The capture is automatically included with this test in the EQR.

Scan Verification

Description: Calibration of mass range is critical in qualitative mass spectrometry.

Procedure: The built-in Agilent autotune is performed to determine the proper calibration of the MS and ensure that masses are correctly reported across the entire mass range of the instrument.

Response Linearity

Description: Knowledge of the response curve is critical for quantitative analysis.

Procedure: A sulfa drug mix standard of four sulfonamide drugs is injected into the system at five concentrations representing a wide range for LCMS. The ions monitored are appropriate to the system type. The calculated RSQ best-fit regression line and plot of the response curve provides the statistics required to evaluate the instrument's overall response curve. This allows users to set appropriate calibration ranges and limits in their quantitative application methods.

Injection Precision

Description: System precision is critical for accuracy of quantitation. Autosampler performance and MS ionization contribute to LCMS system precision. Autosampler precision is challenged in the standard LC module tests using a UV/UV-Vis detector. A repeat precision test in MS mode further challenges the precision of source ionization and MS detection.

Procedure: A blank injection followed by six repeat injections of the sulfa drug mix followed by a final blank injection is made. The % RSD of the six injections is calculated to provide precision statistics.

Carry Over

Description: Low carry over from a previous injection is critical for accuracy of quantitative and reliability of qualitative analysis.

Autosampler performance and MS condition contribute to LCMS carry over. Autosampler carryover is challenged in the standard LC module tests using a UV/UV-Vis detector. A repeat carry over test in MS mode further challenges the full LCMS system carry over performance.

Procedure: A blank injection followed by single injection of the highest concentration standard followed by a blank injection. The last blank injection is evaluated for carry over and the result expressed as a percentage of the value for the standard injection.

Signal to Noise

Description: Sensitivity of MS detection is an important performance feature in quantitative and qualitative analysis. A signal-to-noise value of representative compounds and appropriate ions at known concentration provides sensitivity statistics.

Procedure: For all newly installed Agilent LCMS systems, a reserpine chemical standard is injected as part of the instrument checkout test to provide a starting sensitivity reading. The reserpine signal-to-noise result is provided separately to the OQ report but can be attached to the OQ report if required. For OQ at installation and ongoing OQ/recalibration, the signal-to-noise value of the sulfa drug mix is reported at the ion of interest. System performance over time can be evaluated by repeating this OQ test at suitable intervals.

Instrument Detection Limit (IDL)

Description: Injection precision and IDL are determined using a traceable standard.

Procedure: From a sequence of 16 runs, calculations are performed on eight consecutive runs (1-8, 2-9, and so on) and the best result is used for the evaluation. Calculations include the mean, standard deviation, and %RSD and, using the area RSD and the known amount injected onto the column, IDL (IDL is not set but based on the area RSD limit). This amount returns a peak in the chromatogram, which is detectable and distinguishable from the background with a 99% probability.

Vacuum Verification (LcMs.02.59 and Earlier)

(Removed in LcMs.02.60: successful tune and Scan Verification verify correct system vacuum.)

Description: A stable, high vacuum is required for high-sensitivity mass spectrometry.

Procedure: Multiple readings of the vacuum system are taken and an automated comparison of these values to the known acceptable values is made.

www.agilent.com/chem/qualification

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2025
Published in USA

Appendix: Masses for Scan Verification

Note: In this section, **ESI** refers to ES, DSES, and/or AJST, DSAJST sources, as applicable

Single Quadrupole

Mass	Positive Mode	Negative Mode	Units
G6110A, G6120A/B/C with ESI source			
1	118.09	112.99	m/z
2	622.03	601.98	m/z
3	922.01	1033.99	m/z
G6110A, G6120A/B with CI or MM source			
1	121.05	119.04	m/z
2	622.03	556.00	m/z
3	922.01	805.99	m/z
G6125B/C, G6130A/B/C, G6135B/C with ESI source (* Does not apply to G6125B/C)			
1	118.09	112.99	m/z
2	622.03	601.98	m/z
3	922.01	1033.99	m/z
4	1521.97	1633.95	m/z
5*	2121.93	2233.91	m/z
G6125B, G6130A/B, G6135B with CI or MM source (* Does not apply to G6125B)			
1	121.05	119.04	m/z
2	622.03	556.00	m/z
3	922.01	805.99	m/z
4	1521.97	1305.95	m/z
5*	2121.93	1805.92	m/z
G6140A, G6150B with ESI source			
1	118.09	112.99	m/z
2	622.03	601.98	m/z
3	922.01	1033.99	m/z
4	1221.99	N/A	m/z
G6140A, G6150B with MM source			
1	118.09	112.99	m/z
2	622.03	601.98	m/z
3	922.01	1033.99	m/z
4	1221.99	N/A	m/z
G6140A, G6150B with CI source			
1	121.05	119.04	m/z
2	622.03	655.99	m/z
3	922.01	955.97	m/z
4	1221.99	N/A	m/z
G6160A/B, G6170A with ESI source (* Only applies to G6170A)			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1221.99	1333.97	m/z
6	1521.97*	1633.95*	m/z
7	1821.95*	1933.93*	m/z
8	2121.93*	2233.91*	m/z

Tandem Quadrupole: Q1 and Q3

Mass	Positive Mode	Negative Mode	Units
G6410A with ESI source			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1221.99	1333.97	m/z
6	1521.97	1633.95	m/z
G6410A-2K, G6410B, G6420A, G6430A with ESI source			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1221.99	1333.97	m/z
6	1821.95	1933.93	m/z
G6460A/C with ESI source			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1521.97	1633.95	m/z
6	2121.93	2233.91	m/z
G6465A/B, G6490A with ESI source			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1221.99	1333.97	m/z
G6495A/B/C/D, G6475A, G6470A/B with ESI source (* Only applies to G6475A, G6495D)			
1	118.09	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5*	1221.99	1333.97	m/z
6	1521.97	1633.95	m/z
7*	1821.95	1933.93	m/z
8	2121.93	2233.91	m/z
G6475A with MM source			
1	121.05	112.99	m/z
2	322.05	302.00	m/z
3	622.03	601.98	m/z
4	922.01	1033.99	m/z
5	1221.99	1333.97	m/z
6	1521.97	1633.95	m/z
7	1821.95	1933.93	m/z
8	2121.93	2233.91	m/z
G6465A/B, G6475A, G6495D with CI source (* Only applies to G6475A, G6495D)			
1	121.05	119.04	m/z
2	322.05	316.01	m/z
3	622.03	655.99	m/z
4	922.01	955.97	m/z
5	1221.99	1255.95	m/z
6*	1521.97	1555.93	m/z
7*	1821.95	1855.91	m/z
8*	2121.93	2155.90	m/z

TOF

Mass	Positive Mode	Negative Mode	Units
ESI source			
1	118.086255	112.985587	m/z
2	322.048121	301.998139	m/z
3	622.028960	601.978977	m/z
4	922.009798	1033.988109	m/z
5	1221.990637	1333.968947	m/z
6	1521.971475	1633.949786	m/z
7	1821.952313	1933.930624	m/z
8	2121.933152	2233.911463	m/z
9	2421.913990	2533.892301	m/z
10	2721.894829	2833.873139	m/z

Q-TOF: Quad and TOF Zones

Mass	Positive Mode	Negative Mode	Units
ESI source			
1	118.086255	112.985587	m/z
2	322.048121	301.998139	m/z
3	622.028960	601.978977	m/z
4	922.009798	1033.988109	m/z
5	1221.990637	1333.968947	m/z
6	1521.971475	1633.949786	m/z
7	1821.952313	1933.930624	m/z
8	2121.933152	2233.911463	m/z
9	2421.913990	2533.892301	m/z
10	2721.894829	2833.873139	m/z

F YdcfhUbX'8 Y]j YfmCdHcbg

V@ÁāoÁ^[[, Ág &~ á^•Á] cł }•Á^|&c^áÁ!Ág &~•ā } Ág Á@ÁÔÜÛËÓ Á^~æ |cÁ^] ^æ^áÁ~ } Á[*•Ág áÁg•æcł } Á[*•Ág^Ág, æ•
ā g &~ á^áÁg|Á c@!Á] cł }•Ág Á^Ág &~ á^áÁ!Ág &~ á^áÉ

Q[!Á@] @} æ^áÁ~•c{ Ác] ^•Á} |`DÁā~!^} cł] cł }•Ág^Ág•^} Á[!Á@Á!ā æ^ Ág áÁ~]] [!c^áÁ~•c{ Ác] ^•Ág^Á!ā æ^ Á^•c{
[] cł }•Ág^Á^•áÁÁ!Á[c@Á&@ ā~^•É

ÄÜQ, Ág@[{ æ *!æ •
ÄÜQ, Á@æ^!Ág áÁ [c!Á} Ág ç!
ÄÜg &~ á^Á^] ^æ^áÁ~ } Á[*•
ÄÜg &~ á^Á!g•æcł } Á[*•

GY'WYX'G][bUhi fY'CdHcbg

Üæ •KÖÜÚÁ Á[cł] &^á
ÄÜ^ [!cł * Ágæg &^Ág Ág[[, ^áÁg Ág ÁÔÜÜ

7 i ghca Yf'5 ddfcj U

Ɓæ ^K
Vǻ^K
Öæ^K
Ùǻ}æ'!^K

Ɓæ ^K
Vǻ^K
Öæ^K
Ùǻ}æ'!^K

Ɓæ ^K
Vǻ^K
Öæ^K
Ùǻ}æ'!^K

Ɓæ ^K
Vǻ^K
Öæ^K
Ùǻ}æ'!^K

