

Gen2 Cryogenic System Site Requirements

Varian, Inc. NMR Ion Pump Cryogenic Systems

Pub. No. 01-999296-00, Rev. B 0506

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to **www.agilent.com/chem**.



Agilent Technologies

Gen2 Cryogenic System Site Requirements

Varian, Inc. NMR Ion Pump Cryogenic Systems

Pub. No. 01-999296-00, Rev. B 0506



VARIAN

Gen2 Cryogenic System Site Requirements
Varian, Inc. NMR Ion Pump Cryogenic Systems
Pub. No. 01-999296-00, Rev. B 0506

Revision History:
A 0405– Initial Release
B 0506 – updates for ion pump

Applicability of manual:
Varian NMR systems and ^{UNITY}INOVA systems

Technical contributor: Layne Howard, Susan Klein,
Judit Losonczi, Jeff Germentis
Illustrations and Graphics: Everett Schreiber
Technical editor and writer: Everett Schreiber

Copyright ©2006 by Varian, Inc.
3120 Hansen Way, Palo Alto, California 94304
1-800-356-4437
<http://www.varianinc.com>
All rights reserved. Printed in the United States.

The information in this document has been carefully checked and is believed to be entirely reliable. However, no responsibility is assumed for inaccuracies. Statements in this document are not intended to create any warranty, expressed or implied. Specifications and performance characteristics of the software described in this manual may be changed at any time without notice. Varian reserves the right to make changes in any products herein to improve reliability, function, or design. Varian does not assume any liability arising out of the application or use of any product or circuit described herein; neither does it convey any license under its patent rights nor the rights of others. Inclusion in this document does not imply that any particular feature is standard on the instrument.

VnmrS, ^{UNITY}INOVA, *MERCURY*, *UNITYplus*, UNITY, VNMR, VnmrJ, VnmrX, VnmrI, VnmrV, VnmrSGI, MAGICAL, AutoLock, AutoShim, AutoPhase, limNET, Ultra•nmr, Indirect•nmr, Auto•nmr, Triple•nmr, MagicAngle•nmr, Proton•nmr, Bioproton•nmr, ASM, and SMS are registered trademarks or trademarks of Varian, Inc. OpenWindows, Sun, Solaris, Suninstall, SPARC, and SPARCstation are registered trademarks or trademarks of Sun Microsystems, Inc. and SPARC Int. Oxford is a registered trademark of Oxford Instruments LTD. Ethernet is a registered trademark of Xerox Corporation. VxWORKS and VxWORKS POWERED are registered trademarks of WindRiver Inc. Other product names in this document are registered trademarks or trademarks of their respective holders.

Table of Contents

Chapter 1. Introduction.....	3
1.1 Varian, Inc. Ion Pump Cold Probes	3
1.2 Closed-Cycle Cryogenic System	4
Chapter 2. Site Planning.....	7
2.1 Assistance	7
2.2 General Site Planning and Spectrometer Requirements	7
Site Planning Considerations	8
Typical Installation, Probes with Ion Pump	9
Optional Installation, Probes with Ion Pump	10
2.3 Installation Supplies and Hardware	11
Customer Supplied Equipment and Hardware	11
Varian, Inc. Supplied Equipment and Hardware	11
2.4 Site Requirements	12
General Site Requirements	12
Electrical Requirements	12
Cooling Water Requirements	13
Compressed Gas Requirements	13
Component Heat Dissipation	13
Magnetic Field Considerations	14
2.5 System Component Specifications	15
Component-to-Component Distances	15
Cold Probe	15
Helium Compressor Specifications	16
CryoBay Specifications	16
Closed Cycle Chiller Specifications	16
Mini TASK System	16
Vibration Damping Pier Specifications	16
Water Chiller Specifications	17
Chapter 3. Critical Measurements and Layout Grids.....	19
3.1 Critical Measurements	19
Side View	20
Top View	20
3.2 Room Layout Grids	21
Chapter 4. Site Survey	23
4.1 Magnet	24
4.2 NMR System	24
4.3 Room and Floor	25
4.4 Utilities	25
4.5 Pit Information	25
4.6 Multi Level Room	26
4.7 Single Level Room	26
4.8 Installation and Delivery Contact	26

Chapter 1. Introduction

A Varian Cryogenic system consists of an Ion Pump Cold Probe and closed-cycle cryogenic system:

- “Varian, Inc. Ion Pump Cold Probes,” next
- “Closed-Cycle Cryogenic System” on page 4

The cryogenic system circulates cold He gas to maintain key probe components at an operating temperature of approximately 25 Kelvin.

1.1 Varian, Inc. Ion Pump Cold Probes

Varian, Inc. Ion Pump Cold Probes (Figure 1) achieve significant sensitivity gains through the application of advanced cryogenic technologies to cool key probe components. These probes also achieve significant improvement in operational and vacuum efficiency over probes that do not incorporate the advanced ion pump technologies of these probes.




Figure 1. Varian, Inc. Ion Pump Cold Probes

1.2 Closed-Cycle Cryogenic System

The Varian closed cycle cryogenic system, shown in [Figure 2](#) and described in [Table 1](#), is designed for continuous long-term operation. During normal operations, the system does not consume or require additional He after initial setup. Routine service operations uses some He gas. The system cools the probe components to their operating temperature in about 4 hours. Optimal performance is realized when the system is allowed to stabilize overnight.

Table 1. Closed Cycle Cryogenic System Components

<i>Component</i>	<i>Description</i>
Closed cycle chiller	(CCC), or helium (He) refrigerator, is shown in Figure 2 . The CCC supplies cold He gas to the heat exchanger inside the Varian, Inc. Cold Probe.
CryoBay monitor software	Runs on the control PC — automates all routine cryogenic operations.
CryoBay (cryogenics bay)	Contains the intelligent temperature controller (ITC), control pc, and He compressor, see Figure 2 and Figure 3 . The cabinet has a small foot print and reduces noise.
He compressor	The He compressor provides high-pressure He gas to the CCC.
Intelligent Temperature Controller (ITC)	The ITC regulates the temperature of the heat exchanger inside the Varian, Inc. Cold Probe.
Ion pump and Mini TASK system	Maintains the vacuum in the Varian, Inc. Cold Probe. The ion pump is connected to the probe's vacuum pump-out port.
Damping pier	Reduces the transmission of vibrations from the CCC to the probe.




Figure 2. Varian, Inc. Ion Pump Cold Probes and Cryogenic Systems




Figure 3. CryoBay - Left Side View

Chapter 2. Site Planning

Sections in this chapter

- 2.1 “Assistance,” page 7
- 2.2 “General Site Planning and Spectrometer Requirements,” page 7
- 2.3 “Installation Supplies and Hardware,” page 11
- 2.4 “Site Requirements,” page 12
- 2.5 “System Component Specifications,” page 15

2.1 Assistance

If you need assistance, contact Varian, Inc. Customer Support Center at:

Fax: 650-855-9265

Tel: 1 (800) 356-4437

E-mail: nmr.customersupport@varianinc.com

2.2 General Site Planning and Spectrometer Requirements

- “Site Planning Considerations,” page 8
- “Typical Installation, Probes with Ion Pump,” page 9
- “Optional Installation, Probes with Ion Pump,” page 10

All sites require the development of site plans that are specific to each facility. Placement of the magnet in a pit requires special planning. Use the blank grids provided in 3.2 “Room Layout Grids,” page 21, to arrange the cryogenic system.

The location of the magnet and magnetic field determines the placement of certain system components. Complete the site survey on Chapter 4 “Site Survey,” page 23 and Fax or send the survey to Varian, Inc. Customer Support. Contact Varian, Inc. Customer Support for assistance in completing this survey. Incomplete or inaccurate information can delay the installation of the Varian, Inc. Cryogenic system.

Site Planning Considerations

Consider the following when planning a site for the Varian, Inc. Cryogenic system:

- Components listed in [Table 4](#) are sensitive to stray magnetic fields, see [“Magnetic Field Considerations,”](#) page 14.
- A flexible He transfer line (minimum bending radius of 19.7 inches (50 cm)) extends from the CCC to the probe and requires an unobstructed path. The recommended installation of flexible He transfer line is a bend of 90° from the probe to the CCC. The flexible He transfer line from the CCC ends in a bayonet connector that penetrates the sidearm of the probe.
- The Mini TASK for initially pumping down the probe.
- An ion pump maintains the probe vacuum.
- The He compressor generates a significant amount of heat and requires a continuous supply of cooling water to prevent overheating. Water may be supplied from domestic water sources or by using a closed cycle water chiller see [“Cooling Water Requirements,”](#) page 13.
- The Varian, Inc. Cryogenic system requires either VNMR 6.1 C, VnmrJ 1.1 B and all current updates, VnmrJ 1.1D, or newer.
- VnmrS or ^{UNITY}INOVA NMR Spectrometer systems.
- A 28-channel room temperature shim system is the minimum shim system for which Varian, Inc. Ion Pump Cold Probe line shape specifications are guaranteed.

Typical Installation, Probes with Ion Pump

A typical site plan is shown in **Figure 4**, refer to “Component-to-Component Distances,” **page 15**, for hose lengths and maximum component-to-component distances.




Figure 4. General Floor Plan for Typical Gen 2 Cryogenic System Installation

Optional Installation, Probes with Ion Pump

An optional site plan is shown in **Figure 5** and provides a layout with service room, refer to “**Component-to-Component Distances**,” **page 15**, for hose lengths and maximum component-to-component distances. The optional water chiller (and or the CryoBay - optional location not shown in this site plan) is located in either a utilities room or service chase. A special order flexible He hose extension (65.6 feet / 20 meters) is available.




Figure 5. General Floor Plan for Gen 2 Cryogenic System with Water Chiller in a Service Hall

2.3 Installation Supplies and Hardware

- “Customer Supplied Equipment and Hardware,” [this page](#)
- “Varian, Inc. Supplied Equipment and Hardware,” [this page](#)

Customer Supplied Equipment and Hardware

Table 2 lists the equipment and supplies that the customer must provide for the installation of the cryogenic system.

Table 2. Customer Supplied Equipment and Material

<i>Requirement</i>	<i>Specification</i>
AC power connectors	220 V single phase: 600 VAC, 10 gauge, 3-conductor wire with ground (US and North America). All connectors for single phase operation that must conforming local electrical codes. For sites outside US and North America, all wires and connectors must conform to local electrical codes. 3 phase: 4-conductor wire with ground (US and North America) and connectors for three phase operation conforming to local electrical codes. For sites outside US and North America, all wires and connectors must conform to local electrical codes.
High pressure He gas	He gas (99.999% purity) — full and unopened cylinder
Compressed air or N ₂ gas	Clean, dry, compressed air or N ₂ , min. 80 psi (5.6 bar)
5 micron water filter — for compressor cooling water	Serfilco PL-P 5 x 9 3/4 pleated polyester filter cartridge and GSO 10-3/4 filter chamber (www.serfilco.com) or equivalent
Water — for cooling He compressor	See “ Cooling Water Requirements ,” page 13 .
Cooling water flow control and pressure regulation	Flow control located on the outlet side of the He compressor. Water pressure regulation located on the inlet side of the He compressor.
Pressure regulator for He gas	High purity, He cylinder regulator, 400 PSI (28 bar)— supplied to customers in North America - customers in other locations must supply the regulator and connector to 1/4” Cu tubing.
Optional chiller	See “ Water Chiller Specifications ,” page 17 for specifications in addition to “ Cooling Water Requirements ,” page 13 .

Varian, Inc. Supplied Equipment and Hardware

Varian, Inc. provides an installation kit with the system, the contents are listed in [Table 3](#). The customer must supply the electrical plugs and sockets that are in compliance with local ordinances, refer to “[Electrical Requirements](#),” [page 12](#), for requirements. An electrician and, in some cases, a plumber should be available at the beginning of the installation to make the utility connections.

Table 3. Installation Components Supplied with the System

<i>Component</i>	<i>Components (quantity and description)</i>
Quick Disconnect, Male and Female, Water with 1/2 inch barb, 1/2 inch NPT thread	2 each
Water hose, 1/2 inch ID	150 feet (45.7 meters)

Table 3. Installation Components Supplied with the System

<i>Component</i>	<i>Components (quantity and description)</i>
Copper Tubing, ASTM-B280 high purity	50 feet (15.2 meters) 1/4 inch out side diameter
Pressure Regulator (compatible with high purity, 99.999% He gas), Bottle, 400 PSI (28 bar), 1/4 inch connector for Cu tubing, 1/4 inch Swagelok	Regulator is supplied to North American sites only. Sites outside North America – the customer must supplied the regulator and connector.
Compressed air pressure gauge	

2.4 Site Requirements

The cryogenic system electrical, water cooling, compressed gas, and air conditioning requirements must be met before the installation engineer arrives on site. These requirements are in addition to those requirements of the NMR console.

- “General Site Requirements,” [this page](#)
- “Electrical Requirements,” [this page](#)
- “Cooling Water Requirements,” [page 13](#)
- “Compressed Gas Requirements,” [page 13](#)
- “Component Heat Dissipation,” [page 13](#)
- “Magnetic Field Considerations,” [page 14](#)

General Site Requirements

The site must meet all site requirements for temperature, humidity, etc. as specified in the current *Site Planning Guide*. Requirements and specifications presented here are in addition to those specified in the current *Site Planning Guide*. Where there are potential conflicts in the specifications, the more stringent specification takes precedence.

Electrical Requirements

Determine the electrical power that is provided locally. If the local power does not conform to the following options, specify the power available in the space provided in the “[Site Survey](#),” [page 23](#).

<i>Component</i>	<i>Voltage, Phase, and Line Frequency</i>	<i>Power Consumption</i>	<i>Outlets Required</i>
CryoBay (excludes He compressor)	198 to 242 Vac 10 A 1 phase 50/60 Hz	900 W	1 outlet
Mini TASK system	198 to 242 Vac 10 A 1 phase 50/60 Hz		Plugs into CryoBay
He Compressor Model 6200	200 to 230 Vac 30 A 3 phase 60 Hz 180 to 220 Vac 30 A 3 phase 50 Hz	5.9-7.8 kW	1 outlet
He Compressor Model 6000	360 to 440 Vac 20 A 3 phase 50 Hz 414 to 504 Vac 20 A 3 phase 60 Hz	5.5-7.6 kW	1 outlet

Cooling Water Requirements

Use house water supply or a recirculating water chiller, see “[Water Chiller Specifications,](#)” [page 17](#). Note that Cooling water pressure drop across compressor is typically 3 - 9 psi (0.2 to 0.6 bar).

<i>Requirement</i>	<i>Specification</i>
Maximum inlet temperature	25° C
Minimum inlet temperature	5° C
Flow rate @ 25° C	Minimum 1.3 gal (5L)/min, max.2.6 gal (10L)/min.
Inlet pressure	Minimum 20 psi (2 bar) maximum 116 psi (11.6 bar)
Antifreeze (ethylene glycol) / Water ratio	≤ 50%
Particulate concentration	<10 mg / L (requirement is met if recommended filter is installed on the supply line)
Alkalinity (calcium carbonate)	7-8.5 pH
Carbonate hardness	7 to 10 °dH
Compressor inlet and outlet hose barbs	see “ Helium Compressor Specifications, ” page 16

Compressed Gas Requirements

The cryogenic system requires compressed gas as follows:

<i>Gas</i>	<i>Supply</i>	<i>Specification</i>
Air or Nitrogen	House (independent of the magnet leg) or Cylinder	Minimum 80 psi (max 110 psi), clean and dry (dew point -40C)
He	Full Cylinder	99.999% purity

Component Heat Dissipation

The following table provides heat dissipation from cryogenic system components to help determine the air conditioning load.

<i>Component</i>	<i>watts</i>	<i>BTUs / hour</i>
Cryobay and He compressor	150	512
Mini TASK — start up (continuous)	500	1750
He compressor only	100	350
Closed Cycle Chiller (He refrigerator)	nil	nil
Water Chiller	See “ Water Chiller Specifications, ” page 17	

Magnetic Field Considerations

Detailed magnetic field data can be found in the site planning guide supplied the NMR spectrometer (the current spectrometer site planning guide is available on line, go to the NMR User Pages at www.varianinc.com) or the site planning guide supplied with the magnet.

Cryogenic systems contain one or more components that are sensitive to magnetic fields, these components are listed in [Table 4](#).

Table 4. System Components Affected by Magnetic Fields

<i>Component</i>	<i>Maximum Magnetic Field</i>
CryoBay	5 Gauss
CCC	100 Gauss
Mini TASK	10 Gauss

2.5 System Component Specifications

This section lists the specifications of the individual system components.

- “Component-to-Component Distances,” page 15
- “Cold Probe,” page 15
- “Helium Compressor Specifications,” page 16
- “CryoBay Specifications,” page 16
- “Closed Cycle Chiller Specifications,” page 16
- “Component-to-Component Distances,” page 15
- “Vibration Damping Pier Specifications,” page 16
- “Water Chiller Specifications,” page 17

The manufacturers of the individual components used in the Varian, Inc. Cryogenic System may have additional requirements not listed here. The requirements listed serve as a guide. Where the individual component manufacturer’s requirements are more stringent than those presented here, the individual component manufacturer’s requirements will take precedence.

Component-to-Component Distances

The distances provided are approximate maximum separation between components (or cabinets and components). The distances specified assume a 5 foot (1.5 meter) service or routing loss for the shortest cable, hose, or tubing that connects the components or cabinets. All distances are straight line distances. Routing around obstacles, for appearance, or traffic paths will reduce the distance between components. An extension kit is available for purchase to locate the CryoBay further from the CCC.

<i>Components or cabinets</i>	<i>Component to Component Distance for Systems by Frequency</i>	
	<i>600 MHz systems and below</i>	<i>700 MHz systems and above</i>
CryoBay to CCC	14 ft (4.3 m)	24 ft (7.3 m)
CryoBay to AC power sources	25 ft (7.6 m)	25 ft (7.6 m)
CCC to He cylinder	45 ft (13.7 m)	45 ft (13.7 m)
CryoBay to pneumatics gas supply	25 ft (7.6 m)	25 ft (7.6 m)
CryoBay to cooling water	70 ft (21.3 m)	70 ft (21.3 m)
Mini TASK to probe vacuum valve	16 ft (4.8m)	24 ft (7.3 m)
CCC to probe	5 ft (1.5 m)	8 ft (2.4 m) – 900 MHz only

Cold Probe

The probe requires a minimum clearance between the floor and bottom of the shim tube. The orientation of the probe is limited only by the position of the cryogenic sidearm of the probe in relation to the magnet legs.

<i>Probe</i>	<i>Minimum clearance</i>
500 MHz Cryogenic Probe	28 inches (71.1 cm)
600 MHz Cryogenic Probe	28 inches (71.1 cm)
*Accommodates both new and older style magnets	

Helium Compressor Specifications

Compressor is installed in the CryoBay.

	<i>Specification</i>
Weight	207 lb. (94 kg)
Height	23.6 inches (60 cm)
Length	20.2 inches (51.1 cm)
Width	17.4 inches (44 cm)
Flexible He hose diameter (supply)	1.13 inches (2.9 cm) / fitting 1.19 inches (3.07 cm)
Flexible He hose line diameter (return)	1.13 inches (2.9 cm) / fitting 1.63 inches (3.17 cm)
Cooling water supply and return	1.25 inches (3.2 cm) with fitting, 1/2 inch (1.27 cm) ID
Operating temperature	5° C to 40° C

CryoBay Specifications

	<i>Specification</i>
Height	49 inches (124.46 cm)
Depth and width	30.5 inches x 21.75 inches (77.47 cm x 55.25 cm)
Minimum clearance, front and rear	3 feet (1 meter)
Weight with He compressor	470 lb. (213.2 kg)

Closed Cycle Chiller Specifications

	<i>Specifications</i>
Height	29 inches (73.7 cm)
Depth and width	32 inches x 18 inches (45.72 cm x 81.28 cm)
Weight (with vibration plates)	267 lb. / 121 kg (498 lb. / 226 kg)

Mini TASK System

	<i>Specifications</i>
Height	16 inches (40.6 cm)
Depth and width	12 inches (30.48 cm) x 7 inches (17.8 cm)
Weight	30 lb. / 66.1 kg

Vibration Damping Pier Specifications

Vibration damping pier height is 26 inches (66.04 cm) and footprint diameter is 10 inches (25.4 cm). The pier is located between the CCC and the probe.

Water Chiller Specifications

The water chiller must meet the following specifications:

<i>Requirements</i>	<i>Specification</i>
Cooling Capacity	8 kW (27,000 BTU)
Temperature range	+5°C to +25°C
Pressure, minimum (pumping)	20 psi (1.4 bar)
Flow rate (minimum)	5 L/min

The customer is responsible for verification that the water chiller specifications meet the He compressor requirements.

Several vendors offer compatible water chillers. Haskris Co. Cooling Systems is the Varian recommended supplier for water chillers. Haskris has provided the following information about their chillers that meet the Varian requirements — these water chillers can be purchased from Varian Inc. or Haskris Co. Cooling Systems.

<i>Model</i>	<i>Description</i>	<i>Heat Dissipation</i>	<i>Site Requirements</i>	<i>Voltage/Current (phase) @ 60 Hz</i>	<i>Fusing (Time Delay)</i>
R300-R	Refrigerated water chiller with an air cooled condenser	Heat load of 10.4 kW from the Leybold compressor and Haskris chiller is transferred to room ambient air.	Room/Lab ambient air temperature min. 40°F, max. 100°F.	208/230V 18.0 A (3 phase)	25 A
R300-CR	Refrigerated water chiller with a water cooled condenser	Heat load from the Leybold compressor and Haskris chiller is transferred to building water source	Requires a source of building water 85°F or colder	208/230V 15.4 A (3 phase)	20 A
WW1-R	Non-Refrigerated, water-to-water heat exchanger	Heat load from the Leybold compressor and Haskris chiller is transferred to building water source	Requires a source of building water 55°F or colder	115 V 6.87 A(1 phase) 208/230V 3.4 A (1 phase)	10 A 5 A

Contact Haskris Co. Cooling Systems for any additional information.

Haskris Co. Cooling Systems

100 Kelly Street

Elk Grove Village, IL 60007

www.haskris.com

Contact:

Dan Falotico, Jr.

Telephone: (847)-956-6420 ext. 222

E-mail: dan@haskris.com

Chapter 3. **Critical Measurements and Layout Grids**

Sections in this chapter

- 3.1 “Critical Measurements,” page 19
- 3.2 “Room Layout Grids,” page 21

3.1 Critical Measurements

Provide answers to the following regarding your site, refer to [Figure 6](#), and [Figure 7](#) as necessary.

- “Side View,” this page
- “Top View,” page 20

Side View




Figure 6. Pit, Multi, and Single Level Lab Required Dimensions

Top View




Figure 7. Pit, Multi, and Single Level Lab Required Dimensions - Top Down View

3.2 Room Layout Grids

Use the provided grids to sketch the floor plan, [Figure 8](#), and elevation, [Figure 9](#). Include proposed location of cryogenic system components




Figure 8. Grid for Floor Plane




Figure 9. Grid for Elevation

Chapter 4. Site Survey

Fill in the site survey form and fax it to **650-855-9265**.

Use the following if you have questions regarding this form:

E-mail: nmr.customersupport@varianinc.com

tel: 1-800-356-4437

The survey contains the following sections:

- 4.1 “Magnet,” page 24
- 4.2 “NMR System,” page 24
- 4.3 “Room and Floor,” page 25
- 4.4 “Utilities,” page 25
- 4.5 “Pit Information,” page 25
- 4.6 “Multi Level Room,” page 26
- 4.7 “Single Level Room,” page 26
- 4.8 “Installation and Delivery Contact,” page 26

4.1 Magnet

Premium shielded Varian 500 and 600 and Oxford axially shielded magnets are not compatible with the ion pump.

<i>Survey Element</i>	<i>Dimension</i>	<i>Measurements, Descriptions, other information</i>
Magnet	Field	
	Bore size	
	Actively Shielded	() Yes or () No
	Distance from floor to the bottom of the shim tube	
	Flange type	
	Magnet bottom	() Flat or () Round
Manufacturer	Purchase Date	
	Serial Number	

4.2 NMR System

<i>NMR System Element</i>	<i>Component</i>
Console (Serial number:)	VnmrS ()
	Unity INOVA ()
	Unityplus ()
	Unity ()
Shim System	Varian (channels) ()
	Ultra Shim (channels) ()
	RRI shim system ()
Pulse Field Gradient Amplifier Type	Performa I ()
	Performa II () Part Number:
	Performa XYZ ()
	Other ()
Accessories	SMS ()
	Carousel ()
	NMS ()
	LC-NMR ()
	VAST ()
Software	() VNMR version: () VnmrJ version:

4.3 Room and Floor

<i>Room Elements</i>	<i>Measurements, Descriptions, other information</i>	
Magnet in Pit	() yes () No	If yes complete Pit survey
Multi Level	() yes () No	If yes complete multi level survey
Single Level	() yes () No	If yes complete single level survey
Obstructions	() yes () No	If yes indicate on survey
Floor type	() raised () solid	
Provide digital pictures of the site		

4.4 Utilities

<i>Utility</i>	<i>Measurements, Descriptions, other information</i>			
AC power	Country			
	Single phase	Voltage=	Current=	Line frequency =
	Three phase	Voltage=	Current=	Line frequency =
Cooling Water	Flow rate			
	Pressure			
Compressed Air	Pressure, at max flow			
	Flow			
Air Conditioning	BTU/hr.			

4.5 Pit Information

<i>Pit Element</i>	<i>Measurements, Descriptions, other information</i>
Rectangular Pit	Length
	Width
	Depth
Round Pit	Diameter
	Depth
Main floor elevation relative to:	Magnet Base
Location of Magnet in pit	Center line of magnet relative to pit center line
Diameter of Magnet	
Obstructions	Inside the pit
	Outside the pit
Provide digital pictures of the site	

4.6 Multi Level Room

<i>Survey Element</i>	<i>Measurements, Descriptions, other information</i>
Room upper level dimensions	Length
	Width
Lower level dimensions	Length
	Width
	Depth
Main floor elevation relative to:	Magnet base
Main floor elevation relative to:	Magnet Floor
Obstructions	Obstructions - upper level
	Obstructions - lower level
	Additional Magnets in the same room
Provide digital pictures of the site	

4.7 Single Level Room

<i>Survey Element</i>	<i>Measurements, Descriptions, other information</i>
Room dimensions	Length
	Width
Obstructions and other magnets	
Provide digital pictures of the site	

4.8 Installation and Delivery Contact

Name
E-mail address
Street address
Street address
City
State Zip Code
Country
Contact listed is for: Installation (), Delivery (), Both ()

Name (Alternate)	
E-mail address	
Street address	
Street address	
City	
State Zip Code	
Country	
Contact for: Installation (), Delivery (), Both ()	

