

Site Preparation Specification

Purpose of Procedure

To assure that the installation of instruments and systems can be completed successfully by careful preparation and evaluation of the installation site and by ensuring the availability of appropriate utilities, consumables and supplies.

Customer Responsibilities

Customers should ensure that all necessary operating supplies, consumables and usage dependent items such as columns, vials, syringes and solvents required for the successful installation of instruments and systems are available. Installation sites should be prepared in accordance with the following specifications.

Important Information

If you have problems in providing any of the following, please contact your local Agilent Technologies office for assistance. Assistance with user specific applications may be provided but should be contracted separately. Users of the instrument should be present throughout the installation and familiarization otherwise important operational, maintenance and safety information may be missed.

Procedure Checklist

Dimensions and Weight



Weight: 208 kg	Height: 100 cm
457 lbs.	39 in
Depth: 81 cm	Width: 61 cm
32 in	24 in

At least 90 cm (3 ft) on all sides of the instrument must be kept clear for maintenance.

Tick Boxes



Power Consumption



Europe & USA:
200/220/240V AC ($\pm 5\%$)
30 amps
48-66 Hz
Stability:
Sags or surges: <3 cycles duration
< $\pm 15\%$ mean voltage
Transients: spikes between 1 nano- second
and 5 milli- seconds should be
< $\pm 50\%$ mean voltage
RFI Susceptibility: =500 mV per meter



Heat Dissipation



9894 BTU / hour max




Exhaust Venting Requirements



20 liters/min
Separate exhaust for pumps and API source.



Gas Supply



Nitrogen, UHP grade (99.999%)
80-100 psi. (Two-stage regulator with at least one 1/4 inch Swagelok outlet.)
Up to 20 liters/min. (A 200 liter dewar of liquid nitrogen is recommended as source of nitrogen.)



Site Preparation Specification

Environmental Conditions:



Temperature: 15° - 30° C (59° - 86° F)
variations within $\pm 3^{\circ}\text{C}$
Rel. Humidity: 40-80% non-condensing, non-
corrosive atmosphere
Airborne dust: maximum airborne particle
density of 55 micrograms per M^3

Tick Boxes

