

# Agilent CrossLab Start Up Services

# Agilent Drivers for Waters Empower 4.1 Software Site Preparation Checklist

Thank you for purchasing software from **Agilent Technologies**. CrossLab Start Up is focused on helping customers shorten the time it takes to start realizing the full value of their instrument system investment.

Correct site preparation is the key first step in ensuring that your instruments and software systems operate reliably over an extended lifetime. This document is an **information guide and checklist** prepared for you that outlines the software requirements for the system set up in your lab.

NOTE

Waters, Empower and LAC/E are trademarks of Waters Corporation.





### Introduction

#### **Customer Information**

- The Agilent ICF Support Layer for Waters CDS (WICF) is an instrument control application for the Agilent instrument portfolio running in Waters CDS and included in the Agilent Drivers for Waters Empower package.
- If you have questions or problems in providing anything described as part of Customer Responsibilities below, please contact your local Agilent or partner support / service organization for assistance prior to delivery. In addition, Agilent and/or its partners reserve the right to reschedule the installation dependent upon the readiness of your laboratory.
- Should your site not be ready for whatever reasons, please contact Agilent as soon as possible to re-schedule any services that have been purchased.
- Other optional services such as additional training, functional verification and consultation for user-specific applications may also be provided at the time of installation when ordered with the system but should be contracted separately.
- This document covers the Site Preparation to prepare for the WICF installation and/or to connect Agilent instrumentation into an existing Empower environment. The Waters CDS must be installed and functional before applying the Site Preparation for WICF.
- Waters is responsible for the Waters CDS deployment. Agilent is responsible for the WICF development and deployment.
- To use WICF in an Empower Enterprise environment, it must be installed on every LAC/E and client (or Citrix server) interacting with the Agilent instrument, its methods, and results. WICF must not be installed on database servers or controllers with no Agilent instrument connection or interaction. For incremental upgrade scenarios the interoperability between different WICF and specific Waters ICF Support Layer versions is possible. More details are stated in the WICF Release Notes.
- The Agilent engineer will only perform an installation directly on an Empower Personal Workstation, Client or LAC/E but not do a rollout on a Terminal Server or Citrix environment. If you have such a need, please consult with your Agilent representative upfront. Server and silent installation instructions for WICF can be found in the Installation & Configuration Guide on the installation media.

## **Customer Responsibilities**

Ensure that your site meets the following specifications before the installation date. For details, see specific sections within this checklist, including





- The **computing environment** and the necessary space is made available.
- Your site meets the software, hardware, and networking specifications described later in this
  document.
- Locate your **sales order information**, software authorization codes, and/or software licenses/certificates.
- Ensure that Waters Empower CDS is properly licensed with **ICLs and System Licenses** to configure and startup Agilent instrumentation.
- If purchased, locate, and activate the paper-based Waters Instrument Control License (ICL) via the Waters Licensing Wizard.
- If Agilent is delivering only Connection or Introduction services, the Installation of WICF including a passed Software Verification must be performed upfront, according to the product specific installation guide.
- Ensure to provide the following information at the time of installation/configuration to the service engineer.
  - o Fixed IP address for the instrument (used for instrument configuration).
  - o IP address and access to the instrument controller that will be connected to the instrument.
  - o Planned Access Control Privileges for the system (used during system set up).
- The necessary software media are available.
- That a suitable **backup solution** is identified for your software.
- If Agilent is delivering **Installation, Instrument Connection or Introduction services,** users of the instrument system should be present throughout these services. Otherwise, they will miss important operational, maintenance, and safety information.
- Availability of a system/network administrator as needed to connect to your intranet and to perform installation tasks on a server environment (for example, Citrix).
- Availability of an **Empower administrator** having the privileges to connect and configure instruments.
- Please consult the **Special Requirements** section for other product-specific information.
- [Upgrades only] Ensure that the Waters ICF Support Layer 2.X, 3.X and the included components are correctly uninstalled before the installation of WICF.
- Ensure that your Agilent modules are supported by WICF and match the **minimum Firmware requirements** as well as the capacity recommendation documented in the WICF Release





Notes. If Firmware updates are necessary, please contact your Agilent sales or service representative upfront.

#### **Important Customer Web Links**

- To access Agilent training and education, visit <a href="https://www.agilent.com/chem/training">https://www.agilent.com/chem/training</a> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- To access the **Agilent Resource Center** web page, visit https://www.agilent.com/en-us/agilentresources. The following information topics are available:
  - Sample Prep and Containment
  - Chemical Standards
  - Analysis
  - Service and Support
  - Application Workflows
- To find more information on Agilent Drivers for Waters Empower (for example, Release Notes, Brochures, Technical Notes, etc.) please visit <a href="https://www.agilent.com/chem/ade">https://www.agilent.com/chem/ade</a>
- The **Agilent Community** is an excellent place to get answers, collaborate with others about applications and Agilent products, and find in-depth documents and videos relevant to Agilent technologies. Visit <a href="https://community.agilent.com/welcome">https://community.agilent.com/welcome</a>
- Videos about specific preparation requirements for your instrument can be found by searching the Agilent YouTube channel at https://www.youtube.com/user/agilent
- Need to place a service call? Flexible Repair Options | Agilent
- Links to the **Software Status Bulletin (SSB)** and the **Software Release Bulletin (SRB)** are located under User Documentation on the WICF installation media.





## Software Specifications for Workstations, LAC/Es, and Clients

#### **Special notes**

- If you have purchased hardware instrumentation in addition to the software, refer to the instrument Site Preparation Checklist for requirements regarding laboratory bench space requirements.
- The following Software Specifications apply to WICF 4.1.

Software Specification Description	Minimum	Recommended (if applicable)	Comments
Operating system name, version	Windows 10 Windows 11 Windows Server 2016 Windows Server 2019 Windows Server 2022		
Chromatographic Data System (CDS)	Waters Empower 3 FR4	until Waters Empower 3.9.0	
O/S .NET and other add-ons	.NET 4.8 or higher		For ELSD driver 1.8, Microsoft .NET 3.5 must be activated in addition.
Language settings/compatibility	English, Japanese, and Chinese (simplified)		
Regional settings/compatibility	See Empower specifications		
Account settings/privileges	Local administrator privilege required for installation		Administrative Empower privileges required for instrument connection.
Software Verification Tool		6.3.4.1	

## Computer Hardware Specifications for Workstations, Clients and Servers

#### **Special notes**

 For Computer Hardware Specifications, please refer to the Waters specifications for the Waters CDS.





#### **Networking Specifications**

#### **Special notes**

- The below Network considerations only cover the communication between controller and instrument. For Network Specifications within the Empower environment, for example, Client or database, refer to the Waters documentation.
- Only Ethernet connection between the instrument controller and the instrument is supported. The exception is SFC instruments which use serial connections. To connect an instrument controller to an instrument, follow these guidelines:
  - o Connect via an isolated switch using standard CAT-5 network cabling.
  - o LAN communication hardware should be 100/1000 mbps (or higher) speed capable.
  - LAN communication must be on the same subnet as instruments, and preferably on the same segment.

Network Specification Description	Supported
Network type, bandwidth, speed, protocol etc.	Internet Protocol Version 4 (TCP/IPv4) only Internet Protocol Version 6 (TCP/IPv6) must be deactivated
Additional network or instrument communication card requirements	Static or DHCP Reservation

### **Special Requirements**

- Please refer to the Release Notes on the installation media for supported hardware, software, instruments, and firmware requirements.
- If a firewall is in use other than Windows Defender, ensure that all required ports can be used by the Agilent drivers. For details, see the firewall section in the WICF Installation & Configuration Guide on the installation media.
- Compliance recommendation: If the lab is in a regulated environment, we recommend
  following the customer's standard operating procedures (SOP). The known and resolved
  Issues are listed in the Release Notes. Together with the Software Verification Tool, the
  Software Status and Release Bulletins, this may assist the customer to determine if any
  qualification tasks are required.
- Agilent offers the (re-)qualification of the hardware or the functional verification of the WICF driver installation/upgrade as an additional charged service. Please contact your Agilent sales representative for more information.





# **Service Engineer Review (Optional)**

## **Service Engineer Comments**

If the Service Engineer completed a review of the Site Preparation requirements with the customer, the Service Engineer should complete the following sections below.

	If there are any specific points that should be noted as part of performing the service review or other items of interest for the customer, please write in this box.		
,	Site Preparation Verification		
,	Service Request Number: Date of Review:		



Service Engineer Name:

Service Engineer Signature:

Customer Name:

Total number of pages in this document: