

Agilent CrossLab Start Up Services

Agilent OpenLab CDS v2.6 Client Software Site Preparation Checklist

Thank you for purchasing software from **Agilent Technologies**. CrossLab Start Up is focused on helping customers shorten the time it takes to start realizing the full value of their instrument system investment.

Correct site preparation is the key first step in ensuring that your instruments and software systems operate reliably over an extended lifetime. This document is an **information guide and checklist** prepared for you that outlines the software requirements for the system set up in your lab.

Introduction

Customer Information

- If you have questions or problems in providing anything described as part of Customer Responsibilities below, please contact your local Agilent or partner support / service organization for assistance prior to delivery. In addition, Agilent and/or its partners reserve the right to reschedule the installation dependent upon the readiness of your laboratory.
- Should your site not be ready for whatever reasons, please contact Agilent as soon as possible to re-schedule any services that have been purchased.
- Other optional services such as additional training, operational qualification (OQ) and consultation for user-specific applications may also be provided at the time of installation when ordered with the system but should be contracted separately.

Customer Responsibilities

Ensure that your site meets the following specifications before the installation date. For details, see specific sections within this checklist, including

- The **computing environment** and the necessary space is made available.
- The number and location of **electrical outlets** for your computer systems and peripherals are planned.
- Your site meets the **software, hardware, and networking specifications** described later in this document.
- Locate your **sales order information**, software authorization codes, and/or software licenses/certificates.
- Check for Hotfix, Microsoft updates, and patch compatibility appropriate for the customer's system. **Critical Note: Current Microsoft updates and Windows version must be installed in the customer provided PC before installation can begin.**
- The necessary **software media**, disks etc. are available including upgrade/update disks.
- That a suitable **backup solution** is identified for your software.
- If Agilent is delivering **Installation and Introduction services**, users of the instrument system should be present throughout these services. Otherwise, they will miss important operational, maintenance, and safety information.
- Availability of a **system/network administrator** as needed to connect to your intranet.

- Please consult the **Special Requirements** section for other product-specific information.
- Complete Final Check: Run the System Preparation Tool. Insert the OpenLab CDS 2.6 installation media and run the *SystemPreparationTool.exe* under \Setup\Tools\SPT to complete the final check.

The System Preparation Tool (SPT) checks and applies Windows settings on your machine. Many of the steps are automated by the tool. Applying these settings helps avoid issues during installation and use of the product. Please read the System Preparation Tool Report and manually make the mandatory and recommended adjustments that are outlined in the report.

If Agilent is delivering installation and familiarization services, users of the instrument should be present throughout these services; otherwise, they will miss important operational, maintenance and safety information.

Important Customer Web Links

- To access **Agilent University**, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- To access the **Agilent Resource Center** web page, visit <https://www.agilent.com/en-us/agilentresources>. The following information topics are available:
 - Sample Prep and Containment
 - Chemical Standards
 - Analysis
 - Service and Support
 - Application Workflows
- The **Agilent Community** is an excellent place to get answers, collaborate with others about applications and Agilent products, and find in-depth documents and videos relevant to Agilent technologies. Visit <https://community.agilent.com/welcome>
- Videos about specific preparation requirements for your instrument can be found by searching the **Agilent YouTube** channel at <https://www.youtube.com/user/agilent>
- **Need to place a service call?** [Flexible Repair Options | Agilent](#)

Site Preparation

Software Specifications for Clients

Special notes

Software Specification Description	Comments
Operating system name, version	Windows 10 Pro, 64-bit, (version 2004 or higher) Windows 10 Enterprise, 64-bit, (version 1909 or higher) Windows Server 2016 Standard or Datacenter, 64-bit Windows Server 2019 Standard or Datacenter, 64-bit
O/S .NET and other add-ons	.NET 3.5 SP1 (<i>must be enabled on Windows 10</i>) and .NET 4.x (<i>installed by OpenLab CDS v2.6 Installer</i>)
Web browser	Internet Explorer 11 Google Chrome 40 or higher Edge
Antivirus Software	Symantec Endpoint Protection Trend Micro Microsoft Security Essentials McAfee
Language settings / compatibility	English Western European Language (<i>CDS will always appear in English</i>) Chinese (<i>Non-localized instrument drivers are supported and will always appear in English</i>) Japanese (<i>Non-localized instrument drivers are supported and will always appear in English</i>) Brazilian Portuguese (<i>Non-localized instrument drivers are supported and will always appear in English</i>) Russian (<i>Non-localized instrument drivers are supported and will always appear in English</i>)
Account settings / privileges	Domain user with local administrator privilege required for installation and configuration

Computer Hardware Specifications for Workstations, Clients and Servers

Special notes

Hardware Specification Description	Comments
Processor type and speed	Intel® i5, i7, or Xeon E3 or equivalent 3.0 GHz or greater 4 Core
Memory	8 GB, Ensure that at least 4 GB is reserved for the Windows operating system.
Internal storage/devices/media	500 GB 7200 RPM SATA drives minimum or equivalent solid-state drive
USB Port	USB 2 required for installation via provided media
LAN Card	100 MB/1 GB LAN for instrument control
Graphic resolution	1600 x 900 minimum 1920 x 1080 recommended

Networking Specifications

Special notes

Domain Requirements: Domains support the flow of information and user access rights across machines in the network. This means that all machines and instruments within the networked OpenLab CDS system must reside within the same domain or have the appropriate cross domain trusts to allow name-based communications between all components in the system. In the case of a Client installation, the Client must always be able to communicate with domain components to function as expected.

The domain name server (DNS) must be able to resolve the IPv4 address of all instrument controllers and instruments. Any unresolved instrument controller or instrument will disrupt the functionality of OpenLab resulting in errors or delays. IPv6 is not supported and must be deactivated.

Network Specification Description	Comments
Network type, bandwidth, speed, protocol etc.	Internet Protocol Version 4 (TCP/IPv4) only <i>Internet Protocol Version 6 (TCP/IPv6) is not supported</i>
IP Address	Static or DHCP Reservation

Virtualization

OpenLab CDS has been tested with the following virtualization software. Note that the resource requirements are equal to those of the physical machines.

- VMware vSphere
- Terminal Server/Remote Desktop Services
- Hyper-V for Windows Server
- Citrix XenApp.

Special Requirements

- Please refer to the **CDS_Requirements.pdf** on the installation media for supported hardware, software, instruments, and firmware requirements.

Service Engineer Review (Optional)

If the Service Engineer completed a review of the Site Preparation requirements with the customer, the Service Engineer should complete the following Comments section.

Service Engineer Comments

If there are any specific points that should be noted as part of performing the site preparation review or other items of interest for the customer, please write in this box.

Site Preparation Verification

Service Request Number:

Date of Review:

Service Engineer Name:

Customer Name:

Service Engineer Signature:

Customer Signature:

Total number of pages in this document: