

# GPC on Tour



## Conventional GPC

GPC On Tour, Barcelona,  
28<sup>th</sup> February 2012

## Polymers and Molecular Weight


# What are Polymers?

Polymers are long chain molecules produced by linking small repeat units (monomers) together

There are many ways to link different types of monomer to form polymers

Polymers exhibit very different physical properties compared to the monomers, dependent on the length of the polymer chains




The presence of small amounts of very long or very short chains can have drastic effects on properties of the material



# Variations in Polymers

They can be varied in lots of ways, for example;

- Chemical Structure of Monomer Unit
- 3D Structure
- Different Monomer Units
- Length of polymer chains
- Distribution of polymer chain lengths







# Example 1 - Nitrocellulose

First synthetic polymer made in the 1890's

Hard, strong when set, durable when in moulding

Soon to be renamed gun cotton...!



# Example 2 - Nylon


New York – London (NY-Lon)

1935 – Dupont Chemical Co.

Replaced silk in military parachutes

First product was nylon fibred toothbrush

Tights came in the 1950's









# Example 4 – DNA, Deoxyribonucleic Acid

Longest natural occurring polymer

DNA in lungfish is 36 meters long per cell


DNA in humans is about 1 meter long per cell

Double helix, spiral, symmetry



# Common Polymers

Polystyrene PS



Polyethylene PE, HDPE



Polyvinylchloride PVC, UPVC



Nylon






# Molecular Weight

The molecular weight of a polymer is a way of describing how long the polymer chains are

Each monomer has a molecular weight (often called the formula weight)

Adding the monomers together to make polymers increases the molecular weight

The longer the chains, the higher the molecular weight



# Effect of Molecular Weight

For example, let's look at hydrocarbons

Very short chain hydrocarbons are the predominant component of petrol – liquid at room temperature

Longer chain hydrocarbons are present in various waxes such as candle wax – soft, pliable and easy to melt

Polythene is a very long chain hydrocarbon – tough, strong and very resistant to heat and solvents



# Polymer Molecular Weight Distributions

Samples of synthetic polymers *always* contain polymer chains with a range of chain lengths

One way to describe the length of the polymer chains is in terms of an **average** molecular weight, i.e the average of all the chain lengths in the sample

**HOWEVER....**

Different samples of the same polymer can have the same average chain length but very different distributions of chain lengths depending on the method of production

In polymer science it is the molecular weight **distribution** that is important

# Molecular Weight Averages by GPC

Number average  $M_n$

$M_n$  can be correlative with polymer colligative properties, e.g. freezing point depression

Weight average  $M_w$

$M_w$  may be correlated with properties such as melt viscosity


Z average  $M_z$

$M_z$  may be correlated with properties such as toughness

Polydispersity,  $d = \frac{M_w}{M_n}$

Polydispersity characterises the shape of the distribution

# Defining the Molecular Weight Distribution



A molecular weight distribution can be defined by a series of average values  
Except  $M_p$ , these are various moments of the average of the molecular weights of the distribution

$M_p$  is the molecular weight of the peak maxima

For any polydisperse peak:

$$M_n < M_w < M_z < M_{z+1}$$



# Shape of Distributions


Even for the same type of polymer, each of these distributions will describe a polymer that behaves differently

The red and green plots are for low and high polydispersity materials


The blue plot shows a high polydispersity material with an additional high molecular weight component

Describing these distributions is not easy, especially if they are complex

Distribution plots



# Effect of Polydispersity on a Polymer



As the broadness of the distribution decreases the strength and toughness of the polymer increases

However as the broadness of the distribution decreases the polymer becomes more difficult to process

GPC provides key information to predict the processability and material properties of a polymer

	Strength	Toughness	Brittleness	Melt viscosity	Chemical resistance	Solubility
Increasing Mw	+	+	+	+	+	-
Decreasing distribution	+	+	-	+	+	+

# Measuring Molecular Weight

There are many ways to measure molecular weights

Examples include osmometry, centrifugation and batch light scattering

Each of these methodologies gives a single measurement, and average molecular weight

For example, light scattering measures  $M_w$ , osmometry measures  $M_n$  and centrifugation measures  $M_z$

Although these methods give you a molecular weight, they do not describe a distribution


Gel permeation chromatography (sometimes called size exclusion chromatography) is a method of measuring molecular weights

The advantage of GPC is that it is a separation technique, and as such it is the only common technique that allows the measurement of the molecular weight distribution, not just a single average value

# So where in Chromatography is GPC

*Interactive* adsorption, partition, ion exchange, etc

*Non-interactive* GPC, SEC, GFC



# Nomenclature

There are many ways of measuring the molecular weight of a polymer however, ***there is only one technique that measures the molecular weight distribution.***

- **Gel Permeation Chromatography, GPC** (Polymer Industry)

One technique but multiple acronyms, also known as

- **Size Exclusion Chromatography, SEC** (Academia)  
and
- **Gel Filtration Chromatography, GFC** (Protein/Pharma)



# What is a GPC/SEC System?

- A simple isocratic LC system fitted with a GPC/SEC column is a GPC/SEC system!
- Mode of separation only difference to other HPLC methods
- Specialist detectors can be used to determine properties of the samples investigated
- Special GPC/SEC software required to perform analysis

Liquid  
Chromatograph



GPC/SEC  
Column



Detectors, Data  
Acquisition  
and Processing

Simple isocratic system  
Provides solvent flow

Column separates sample  
Into components

Detectors measure the  
properties of eluted sample,  
data system provides  
analysis



# Additional Components Used in GPC

## Concentration detectors

- Differential refractometer (RI)
- Ultraviolet absorbance (UV)
- Evaporative light scattering or mass detector (ELS, EMD)
- Infra-red (IR)

## Molecular weight sensitive detectors


- Viscometry
- Light scattering

## Additional systems

- Online degasser
- Autosampler
- Column oven
- Additional specific detectors

# Polymer Molecules in Solution


- GPC is based on the behaviour of polymer molecules in solution
- In the solid state polymers can be considered like spaghetti – a confusing mass of intertwined chains
- In solution, polymer molecules are discrete entities
- Due to entropic effects all but the most rigid of polymer chains curls up in solution to form a ball like shape




# GPC Column Packings

- GPC columns are packing with cross-linked, insoluble beads, typically copolymers of styrene and divinyl benzene for organic GPC
- These beads have a rigid pore structure that remains intact in the presence of solvent

PLgel 10  $\mu\text{m}$   $10^6\text{A}$




PLgel 10  $\mu\text{m}$   $10^3\text{A}$





# Synthesis of Porous Beads

- High cross-link content gives a rigid, low swelling product with a well-defined pore structure






# Permeation of Polymer Molecules



- Polymer coils in solution can permeate the pores on GPC packing materials
- Exclusion, partial permeation and total permeation are possible


# GPC Separation Mechanism

- Polymer is prepared as a dilute solution in the eluent and injected into the system
- The GPC column is packed with porous beads of controlled porosity and particle size
- Large molecules are not able to permeate all of the pores and have a shorter residence time in the column
- Small molecules permeate deep into the porous matrix and have a long residence time in the column
- Polymer molecules are separated according to molecular size, eluting largest first, smallest last




# Elution Profiles

- As a result of the GPC separation mechanism, polymer molecules elute from the column in order of size in solution
- Largest elute first, smallest elute last
- The separation is purely a physical partitioning, there is no interaction or binding
- The separation is isocratic
- If polymer molecules have the same molecular dimensions, they will co-elute by GPC and may not be separated by this technique
- The calibration curve describes how different size molecules elute from the column




# Typical Calibration Curves for PLgel Individual Pore Size Columns




# Determination of Polymer Molecular Weight Distribution by GPC

- Produce a GPC calibration curve for the column set relating  $\log M$  to retention time (RT)
- Chromatograph the polymer sample
- Normalise and integrate the GPC response versus retention time plot for the polymer sample
- Convert retention time to  $\log M$  via the GPC calibration curve
- Present a  $\log M$  distribution plot and calculate molecular weight averages ( $M_n$ ,  $M_w$ ) for the distribution





# EasiCal Pre-prepared Calibrants



# Polymer Calibrants for GPC

Mn - *number average molecular weight*

Mw - *weight average molecular weight*

Mv - *viscosity average molecular weight*

**Mp** - ***peak molecular weight***

Mw/Mn - *polydispersity by GPC*

*Must be extremely well characterised*

## **Most commonly used polymer calibrants**

**Polystyrene** - THF, toluene, chloroform, TCB

**Polymethyl methacrylate** - MEK, ethyl acetate, acetone, DMF

**Polyethylene oxide/glycol** - aqueous eluents, DMF, DMSO

# Example Certificates of Analysis

## Individual standard

**CERTIFICATE OF ANALYSIS**

Product: Polyethylene oxide Mp 305,500  
 Batch Number: 0001021527  
 Part Numbers: PL2083-8001, PL2083-9005, PL2083-8010


	GPC	Light scattering	Viscometry
Mp (g/mol)	305,500		
Mn (g/mol)	288,000		
Mw (g/mol)	301,300	308,900	
Mv (g/mol)	298,400		
	Mw/Mn = 1.05		[η] = 3.0162 dl/g

Batch details

Analysis details

GPC data

Results



**Analysis Conditions**

	GPC	Light scattering	Viscometry
System	Cirrus GPC / SEC	Viscotek TDA301	Viscotek TDA301
Detector	Refractive Index	Viscotek TDA301	Viscotek TDA301
Columns	4x PL aquagel-OH MIXED-H 8µm 300x7.5mm	2 x PL aquagel-OH MIXED 1 x PL aquagel-OH Guard	2 x PL aquagel-OH MIXED 1 x PL aquagel-OH Guard
Solvent	0.02% NaCl Sol	0.2M NaNO <sub>3</sub> ; 0.01M NaH <sub>2</sub> PO <sub>4</sub> @ pH 7.0	0.2M NaNO <sub>3</sub> ; 0.01M NaH <sub>2</sub> PO <sub>4</sub> @ pH 7.0
Flow rate	1.0 ml / min	1.0 ml / min	1.0 ml / min
Injection volume	20µl	100µl	100µl
Sample concentration	0.05%	1 mg / ml	1 mg / ml
Temperature	Ambient	30°C	30°C
Calibrants	PEO/PEG	Polyethylene Oxide	Polyethylene Oxide
Angle		90°	
dn/dc		0.132 q/ml	

The above characterisation data has been measured according to our Quality Control procedures.  
 Certificate of Analysis valid until expiry date – 8th October 2015  
 Agilent Manufacturing Site: Essex Road, Church Stretton, Shropshire, SY6 8AX, UK

*N.W. Titley*  
 N.W. Titley

*G.K. Harmer*  
 G.K. Harmer

Issue 1 25th May 2011

Quality Department

COA STANDARDS 1 Rev 1.04



## EasiVial

**CERTIFICATE OF ANALYSIS**

Product: Polystyrene High EasiVials (2ml)  
 Part no: PL2010-0201 & PL2010-0202  
 Batch no: VPS-H-021

**Characterisation Data Summary**

Vial Code	Batch no.	IV	Mv	Mw (Light Scattering)	Mn *	Mw *	Mw/Mn *	Mp *	Massival mg
RED	20147-12	7.0752	5,235,000	6,375,000 <sup>†</sup>	5,405,000	5,680,000	1.05	6,035,000	0.4
	20139-21	1.2525	457,000	441,000	448,000	469,000	1.05	483,000	0.8
	20131-6	0.1422	21,330	20,000	19,340	19,760	1.03	19,720	1.2
	20124-14	0.0296	1,150	1,450	1,270	1,270	1.06	1,260	1.6
YELLOW	20145-17	4.7666	3,003,000	2,748,000	2,922,000	3,001,000	1.03	3,053,000	0.4
	20137-18	0.6551	183,500	178,800	181,700	184,200	1.02	184,900	0.8
	20128-4	0.0815	8,400	8,740	8,210	8,400	1.03	8,450	1.2
	20122-21	-	-	695	544	615	1.13	590	1.6
GREEN	20141-19	2.1237	962,000	903,000	860,000	893,000	1.04	915,000	0.4
	20134-3	0.2960	59,950	61,300	59,050	59,950	1.02	60,450	0.8
	20126-15	0.0502	3,240	3,660	3,240	3,360	1.04	3,370	1.2
	20121-9	-	-	-	-	-	1.00	162	1.6

Batch details

Analysis details

\*Results of polymer characterisation by gel permeation chromatography using tetrahydrofuran as eluent at a flow rate of 1.0ml/min. PLGel GPC columns selected appropriate to the molecular weight of the polymer. Each individual polymer has been characterised using methods referred to in DIN 55672.

Light Scattering (LS) was performed on the precision detectors Inc PD2020SA light scattering detector, in conjunction with a GPC/SEC workstation using a refractive index detector. Molecular weight was determined using the batch calculation mode. The analysis was undertaken using tetrahydrofuran as an eluent with a dn/dc of 0.185g/ml.

†Light Scattering (LS) was performed on the PL-LSP spectrophotometer at 633nm, the reference value used for the Rayleigh ratio was 14.02x10<sup>-6</sup> cm<sup>-1</sup>. The analysis was undertaken in toluene at 25°C with a dn/dc value of 0.1065g/ml.

Viscometry (IV, Mv) was performed on a Schott-Gerate viscometry system using an Ubbelohde capillary viscometer. The analysis was undertaken in toluene at 30°C. A double extrapolation of the Huggins-Kraemer plot gives the intrinsic viscosity, the viscosity average molecular weight was then determined using the Mark-Houwink constants.

The above characterisation data has been measured according to  
 Varian, Inc. Quality Assurance procedures.  
 Certificate of Analysis valid until expiry date – 8th September 2013

*N.W. Titley*

N.W. Titley

*G.K. Harmer*

G.K. Harmer

Quality Department

Issue 1 9th September 2010

COA STANDARDS-2 Rev 2.04

www.agilent.com



# Calibration Methods for Conventional GPC



Aim : to produce a mathematical model for  $\log M$  versus retention time


## **Narrow standards**

## **Broad standards** (rarely used now)


- Hamielec
- Broad on Narrow
- Integral

# Calibration of GPC Columns Using Narrow Standards

- Chromatograph a series of well characterised, narrow polydispersity polymer standards
- Plot peak retention time (RT) versus peak log molecular weight (logM)
- Fit the data using a mathematical function (e.g. polynomial order 1,2,3, etc)
- The calibration curve will be characteristic of the GPC column set used



# Errors Due to Limited Calibration Region



*The column calibration should cover the full elution time region of the sample to avoid errors due to extrapolation*

# Interpreting Chromatograms

- The data obtained in a GPC experiment will be in the form of a chromatogram showing detector response as a function of retention time
- There are fundamental parameters that are present on all chromatograms


## Plasticized PVC

Columns: 3xPLgel 5 $\mu$ m MIXED-C, 300x7.5mm  
(PL1110-6500)

Eluent: THF

Flow Rate: 1.0ml/min

Detector: RI




# Peak Separation

- Peak separation in GPC is dependent upon resolution and on molecular size
- If two samples have different molecular sizes, then they will be separated to baseline assuming there is sufficient resolution
- However, if samples are the same molecular size, then they cannot be separated by GPC as the mechanism of SEC is based upon size


## Starch

Columns: 4xPLgel 20 $\mu$ m MIXED-A, 300x7.5mm (PL1110-6200)  
Eluent: DMSO + 5mM NaNO<sub>3</sub>  
Flow Rate: 1.0ml/min  
Temp: 80°C  
Detector: RI





# Excluded Peaks



# Partial Exclusion

- The dead space of the separation will be around half of the total elution volume
- Peaks eluting close to this volume may be partially excluded
- Look for sharp peaks at the front of your chromatograms

## Polyurethane

Columns: 2xPLgel 5 $\mu$ m MIXED-C, 300x7.5mm (PL1110-6500)


Eluent: DMF + 0.1% LiBr

Flow Rate: 1.0ml/min

Temp: 80°C

Detector: RI

**Partially Excluded Peak**




# Oligomeric Resolution

- In GPC, the relationship between molecular weight and retention time is logarithmic
- As a result, peaks equidistant in molecular weight elute closer together with increasing molecular weight
- This is a classic way to tell a separation is based on SEC

## Epoxy Resin

Columns: 2xMesoPore, 300x7.5mm (PL1113-6325)  
Eluent: THF  
Flow Rate: 1.0ml/min  
Inj Vol: 20 $\mu$ l  
Detector: RI

**Gap Between  
Oligomers**




# Oligomeric Resolution

- With some columns it is possible to calibrate the column using the oligomers
- The molecular weights of the initiator fragment and the repeat unit of the polymer must be known

## Polystyrene Standards


Columns: 2xOligoPore, 300x7.5mm (PL1113-6520)  
Eluent: THF  
Flow Rate: 1.0ml/min  
Detector: RI



# Interpreting Molecular Weight Distributions

## Hyaluronic Acid

Columns: PL aquagel-OH 60 15 $\mu$ m, 300x7.5mm (PL1149-6260)  
PL aquagel-OH 40 15 $\mu$ m, 300x7.5mm (PL1149-6240)  
Eluent: 0.2M NaNO<sub>3</sub>, 0.01M NaH<sub>2</sub>PO<sub>4</sub>, pH 7  
Flow Rate: 1.0ml/min  
Detector: RI



- The molecular weight distribution shows the amount of material present as a function of the molecular weight

- The MWD looks a bit like a 'mirror image' of the chromatogram

# Effect of Baseline Position


## Polyacrylamide

Columns: PL aquagel-OH 60 15 $\mu$ m, 300x7.5mm (PL1149-6260)  
PL aquagel-OH 40 15 $\mu$ m, 300x7.5mm (PL1149-6240)

Eluent: 0.2M NaNO<sub>3</sub>, 0.01M NaH<sub>2</sub>PO<sub>4</sub>, pH 7

Flow Rate: 1.0ml/min

Detector: RI




# Effect of Baseline Position

## Sodium Polyacrylate

Columns: PL aquagel-OH 40 8 $\mu$ m, 300x7.5mm (PL1149-6840)  
PL aquagel-OH 30 8 $\mu$ m, 300x7.5mm (PL1120-6830)  
Eluent: 0.2M NaNO<sub>3</sub>, 0.01M NaH<sub>2</sub>PO<sub>4</sub>, pH 7  
Flow Rate: 1.0ml/min  
Detector: RI

- The whole peak should be analysed to get a true reflection of the sample
- The peak should go down to the baseline on either side
- Leaving out components of the peak will leave an 'incomplete' MWD



# Conventional GPC

Now lets take a look at ways we can improve the quality of our  
'Conventional GPC Analysis'