



Advanced Petrochemical Applications With ICP-OES

**Tina Harville, Ph.D.
Application Engineer**



Elemental Analysis of Petrochem-Based Solutions

Refineries and Petrochemicals

Sulfur content of gasoline

Acceptance testing of imported fuels

Biodiesel analysis (ASTM D6751)

Wear metals in oil

Oil additive analysis (process control and trend)

Metals in Catalyst & Catalyst Oil

Metals in Process Solutions

Metals in Drilling Fluids


Other: Metals in Brines

Etc...!




ICP-OES Design and Performance for Petrochem


Unmatched Background Correction Techniques



- **Automatic (Polynomial Fitted Technique)**



- **FACT (Fast Automated Curve-Fitting Technique)**




- **Offpeak BC–Left & Right Together**

Continuous Wavelength Coverage



- **33,000 wavelengths ready at all times**
- **Freedom from interferences**
- **MultiCal**

Robust RF & Plasma



- **Robust, fast, frequency-matching 40MHz RF generator**

Application Considerations

The challenge:

- High vapour pressure from solvents
 - Plasma instability
 - Extinguished plasma
- Carbon build up in injector
 - Poor precision and drift
 - Down time - injector requires regular cleaning
- Nebulizer blockage
 - Poor precision
 - Down time - nebulizer needs cleaning



Application Considerations

The solution:

- Plasma stability
 - High Efficiency RF Generator
 - Reduce solvent loading
 - Double-pass spray chamber
 - Plasma torch with 1.4mm ID injector
 - Smaller id. peristaltic pump tubing
 - Reduce pump speed
 - Higher plasma power and plasma gas flows



Spray Chambers

- Double-pass glass cyclonic spray chamber
 - Double pass design
 - Reduces solvent load
 - Increased sensitivity
 - Fast washout
 - Suitable for low vapour organic solvents
 - e.g. Kerosene, ShellSol



Nebulizers

- Sample particulates may block standard concentric nebulizers
 - Dependent on sample composition
e.g. wear metal particulates
 - Poor precision
 - Down time - nebulizer needs cleaning
- Quartz nebulizer
 - Conikal, Seaspray and Slurry
 - Faster washout
- V-groove nebulizer
 - Made of PEEK polymer
 - Minimizes particulate blockage
- OneNeb Nebulizer



Quartz nebulizer



V-groove nebulizer



OneNeb nebulizer

Compatibility of Pump Tubing

Pump Tubing			
Solvent	PVC	Viton	PVC Solva
Kerosene	✗	✓	✓
Gasoline	✗	✓	✗
Fuel	✗	✓	✗
Coolant	✓	✓	N/A

✗ = unsatisfactory, ✓ = satisfactory, N/A = no data available




Typical Instrument Conditions

Plasma power	1.3 - 1.5kW
Plasma gas flow	15 – 18 L/min
Auxiliary gas flow	0.75 – 2.25 L/min (axial – radial)
Nebulizer gas flow (axial)	Optimize (0.5 – 0.8 L/min)
Nebulizer gas flow (radial)	Set bullet to top of torch (0.5 – 0.8 L/min)
AGM-1 setting	2 - 6
Stabilization delay	Optimize via timescan
Pump speed	5 – 10 rpm
Uptake Rate	< 1 mL/min
Fast Pump	Yes or No?



Elemental Analysis of Diesel Fuel by ICP-OES




Elemental Analysis of Diesel Fuel by ICP-OES

Element/Wavelength (nm)	DL (ppb)
Ag 328.068	5
Al 396.152	20
B 249.772	9
Ba 455.403	0.8
Ca 396.847	0.5
Cd 214.439	4
Cr 267.716	7
Cu 327.395	7
Fe 238.204	5
Mg 279.553	1


Element/Wavelength (nm)	DL (ppb)
Mn 257.610	1
Mo 202.032	10
Na 589.592	6
Ni 231.604	18
P 213.618	60
Si 251.611	39
Ti 336.122	2
V 392.401	9.6
Zn 213.857	4.7




Stability of Kerosene by ICP-OES



Stability – Coolant Analysis



Plasma Stability of Naphtha



Short term stability over 30 mins. for Naphtha



Organic Applications

Other challenges:


- Wear metals and additive analysis
 - Require large linear dynamic range
- Structured backgrounds
 - Molecular band observed at wavelengths $> 330\text{nm}$
- Speed of analysis
 - Throughput = time = money



Speed of Analysis

World's Most Innovative Detector

- Simultaneous ICP using custom designed and patented CCD detector
 - Full wavelength coverage 167 – 785 nm
 - Use alternate wavelengths to **extended dynamic range** and confirm results
 - True simultaneous measurement
 - Improved resolution
- Pixels arranged to exactly match 2 dimensional echelle image
- Duplex readout circuitry **halves readout time**
 - Processing electronics off-chip for high QE
 - 40 times faster due to 1 MHz processing speed
- Typical analysis times
 - **Can measure 1 sample in < 1 min.**
(all elements)



Performance Test – Linear Range

Radial, S21 Conostan Oil Standards in Kerosene

Element & Emission Line Wavelength	Maximum Concentration (mg/L)
Ag 328.068	250 +
Al 308.215	250 +
Al 396.192	100
B 249.772	250 +
Ba 455.403	100
Ba 493.408	250 +
Ca 317.933	2500
Ca 396.847	100
Cd 226.502	250 +
Cr 284.325	250 +
Cu 327.395	250 +
Fe 259.940	250 +
Fe 274.932	1000
Mg 280.270	100
Mn 257.610	250 +
Mo 202.032	250 +
Na 589.592	250 +
Ni 230.299	250 +
P 213.618	2500
Pb 220.353	1500
Si 251.608	250 +
Sn 283.998	250 +
Ti 336.122	250 +
V 311.837	250 +
Zn 206.200	2500
Zn 213.857	250



Biodiesel



Analytical Requirement for Biodiesel

ASTM requires the analysis of Biodiesel

Must satisfy all the requirements of ASTM D 6751-07be1, including:

- Sulphur and phosphorus (D6751)
- Total calcium and magnesium (D6751- 07a)
- Total sodium and potassium (D6751- 07a)
- Includes other properties such as:
 - Flash point
 - Water and sediment
 - Free glycerin



Analytical Results

– ASTM Specified Elements in Biodiesel

Element	P	S	Na	K	Ca	Mg
Wavelength (nm)	213.618	181.972	589.592	766.491	393.366	279.553
B-100a ppm	< IDL	1.695	0.162	0.507	0.013	0.011
B-100b ppm	< IDL	1.710	0.130	0.469	0.012	0.012
B-100+1.0 ppm (%R)	105	103	104	107	105	104
B-100+2.5 ppm (%R)	102	106	102	104	103	101
IDL (ppb)	19.9	76.5	4.1	21.8	0.1	0.1
Bkgd Correction	Fitted	Fitted	FACT	FACT	Fitted	Fitted



Analytical Results

– Other Elements in Biodiesel

Element	Al	Ba	Cd	Cr	Fe	Mn	Mo	Ni
Wavelength (nm)	308.215	233.527	226.502	267.716	259.940	259.372	202.032	216.555
B-100a ppm	< IDL	< IDL	< IDL	< IDL	0.032	< IDL	< IDL	< IDL
B-100b ppm	< IDL	< IDL	< IDL	< IDL	0.037	< IDL	< IDL	< IDL
B-100+1.0 ppm (% R)	102	105	103	103	104	103	103	104
B-100+2.5 ppm (% R)	103	103	101	102	103	102	103	102
IDL (ppb)	5.4	1.0	0.6	1.6	1.0	0.3	5.7	6.9
Bkgd Correction	Fitted	Fitted	Fitted	Fitted	Fitted	Fitted	Fitted	Fitted

Element	Pb	Sc	Si	Sn	Ti	V	Zn
Wavelength (nm)	220.353	361.383	288.158	283.998	334.941	309.310	213.857
B-100a ppm	< IDL	< IDL	0.141	< IDL	< IDL	< IDL	0.024
B-100b ppm	< IDL	< IDL	0.139	< IDL	< IDL	< IDL	0.024
B-100+1.0 ppm (% R)	102	104	104	102	104	105	104
B-100+2.5 ppm (% R)	101	103	103	104	102	104	102
IDL (ppb)	25.2	0.1	4.7	11.7	0.1	0.6	1.1
Bkgd Correction	Fitted	Fitted	Fitted	Fitted	Fitted	Fitted	Fitted



Naphtha



Why Determine Trace Elements in Naphtha?

- Growing need to determine trace metals in hydrocarbon samples
- The presence of trace metals can:
 - Severely hamper the catalytic reforming process
 - Poison the catalysts used e.g. S and N compounds can deactivate catalysts




Catalysts used in catalytic cracking or reforming

Photo courtesy Phillips Petroleum Company

Improved Sample Introduction Systems

- Cooled spray chamber
 - Externally jacketed design
 - Ethylene glycol mixture circulated to cool sample to -12 °C
 - Reduces solvent load
 - Made of quartz
 - Excellent for highly volatile organic solvents e.g. naphtha



Sample and Standard Preparation

Standard Preparation

- Multi-element organometallic standard (e.g. Conostan S-21)
- Select dilution ratio to achieve required concentration
- Add extra neutral base oil to ensure consistent viscosity
- Diluted with Isopropanol (iso propyl alcohol)



Detection Limit Performance (ppm)

Element	Instrument DL
Al 167.019 nm	0.41
Al 396.152 nm	0.06
Ni 231.604 nm	0.23
Ni 216.555 nm	0.61
Mg 279.553 nm	0.05
Mg 280.270 nm	0.2
V 309.310 nm	0.04
V 292.401 nm	0.07
Ca 393.366 nm	0.05
Ca 306.847 nm	0.05
Ca 422.673 nm	0.04
Cr 267.716 nm	0.05



Detection Limit Performance (ppm)

Element	Instrument DL
Cd 214.439 nm	0.09
Cd 226.502 nm	0.10
Cd 228.802 nm	0.08
Cu 324.754 nm	0.06
Cu 327.395 nm	0.06
Fe 238.204 nm	0.06
Fe 259.940 nm	0.06
Mn 257.610 nm	0.06
Na 588.995 nm	0.05
Na 589.592 nm	0.04





Wear Metals Analysis



Rapid Wear Metal Analyses


To improve sample throughput must look at:

- Sample introduction
- Instrument delays
- Rinse time
- Instrument setup and preventative maintenance time



Rapid Flow Fast Nebulization

A triple channel peristaltic pump is required




Rapid Flow Fast Nebulization

Channel 1: Pump waste away from spray chamber

Channel 2: Large ID tube creates FAST FLOW from autosampler to channel 3

From autosampler



Transfer line between large ID channel 2 tube to small ID channel 3 tube

"T" piece diverts excess sample from channel 2 to waste

To nebulizer

Sample Throughput

Throughput speed

2 orders of magnitude washout

Using Cetac ASX 500 autosampler

Sample Introduction System	Rinse Time (s)	Sample Throughput (s)
Rapid flow fast nebulisation	11	39
Unmodified sample introduction	30	85



Elemental Analysis of 25% Brine by ICP-OES

Undiluted Analysis of 25% Brine




Figure 1 Standard torch (A), high TDS torch (B) and sheath gas high solids torch (C)




Figure 2 ICP-OES Stability achieved with the sheath gas high solids torch.

Element/Wavelength (nm)	Detection Limit (µg/L)
Ag 328.068	0.7
Al 167.019	0.5
As 188.980	5
Ba 455.403	0.2
Be 313.042	0.07
Cd 214.439	0.4
Co 238.892	1
Cr 205.560	1
Fe 238.204	1.7
K 766.491	1.5
Li 670.783	0.4
Mg 279.553	0.1
Mn 257.610	0.2
Mo 202.032	2
Ni 231.604	2.5
Pb 220.353	6
Sb 206.834	5
Se 196.026	10
Sr 407.771	0.1
Ti 336.122	0.5
Tl 190.794	9
V 292.401	1
Zn 213.857	0.5

Table 1. Measured Instrument Detection Limits for a sample containing >250 g/L TDS using the Sheath Gas Torch with the Agilent 720/730 Series ICP-OES. Replicate read time was 20 s.



Application Papers

Available on the Agilent Technologies Web site

(<http://www.agilent.com>)

- SI-A-1413 Determination of metals in oils by ICP-OES
- SI-A-1417 Determination of V, Ni and Fe in crude oils and bitumen with Sc as an internal standard
- SI-A-1420 Determination of wear metals in lubricating oil with Axial ICP
- SI-A-1422 Determination of Pb in Unleaded Gasoline with Axial ICP
- SI-A-1423 Determination of trace elements in a xylene solution of oil by ICP-AES with ultrasonic nebulization and membrane desolvation
- SI-A-1427 Multi-element analysis of fuel and lubricating oils by simultaneous ICP-OES
- SI-A-1431 Improving Throughput for Oils Analysis by ICP-OES

PLUS: SI-A-1202, SI-A-1415 and SI-A-1418





THANK YOU

QUESTIONS?

