

OPTIMIZING GPC/SEC SOLUTIONS

POLYMER PEOPLE SEMINAR TOUR



SPONSORED AND HOSTED BY AGILENT TECHNOLOGIES, INC.

GPC/SEC An essential tool for polymer analysis

Ben MacCreath, PhD
Product Manager
GPC/SEC Instrumentation

26th March 2013

Introduction to Polymers

Where are they found?

Polyolefins



Analysis of polyolefins by GPC/SEC

Application Compendium

Authors
Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.



Low molecular weight resins - Analysis of low molecular weight resins and prepolymers by GPC/SEC

Application compendium

Authors
Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Resins



Engineering Polymers



Biodegradable polymers - analysis of biodegradable polymers by GPC/SEC

Application compendium

Authors
Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Excipients



Excipient analysis by GPC/SEC and other LC techniques

Application compendium



Analysis of food additives by GPC/SEC

Application compendium

Authors
Greg Saunders and Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Authors
Greg Saunders and Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Analysis of engineering polymers by GPC/SEC

Application Compendium

Authors
Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Biodegradable Polymers

Analysis of elastomers by GPC/SEC

Application compendium

Authors
Greg Saunders and Ben MacCreath
Agilent Technologies, Inc.



The Measure of Confidence

Agilent Technologies

Elastomers

Food


Molecular Weight

The molecular weight of a polymer is a way of describing how long the polymer chains are

Each monomer has a molecular weight (often called the formula weight)

Adding the monomers together to make polymers increases the molecular weight

The longer the chains, the higher the molecular weight



Effect of Molecular Weight

For example, let's look at hydrocarbons

Very short chain hydrocarbons are the predominant component of petrol – liquid at room temperature

Longer chain hydrocarbons are present in various waxes such as candle wax – soft, pliable and easy to melt

Polyethylene is a very long chain hydrocarbon – tough, strong and very resistant to heat and solvents



Polymer Molecular Weight Distributions

Samples of synthetic polymers *always* contain polymer chains with a range of chain lengths


One way to describe the length of the polymer chains is in terms of an average molecular weight, i.e the average of all the chain lengths in the sample

HOWEVER....

Different samples of the same polymer can have the same average chain length but very different distributions of chain lengths depending on the method of production

In polymer science it is the molecular weight *distribution* that is important

Molecular Weight Distribution



What is my Molecular Weight?

Molecular Weight Averages

NMR measures M_n

Light scattering measures M_w

Osmometry measures M_z

GPC measures *distribution*



Shape of Distributions

Even for the same type of polymer, each of these distributions will describe a polymer that behaves differently.


The Mn or Mw might be identical for all plots.

The red and green plots are for low and high polydispersity materials


The blue plot shows a high polydispersity material with a additional high molecular weight component

Describing these distributions is not easily, especially if they are complex

Distribution plots



Effect of Polydispersity on a Polymer



As the broadness of the distribution decreases the strength and toughness of the polymer increases

However as the broadness of the distribution decreases the polymer becomes more difficult to process




GPC provides key information to predict the processability and material properties of a polymer

	Strength	Toughness	Brittleness	Melt viscosity	Chemical resistance	Solubility
Increasing Mw	+	+	+	+	+	-
Increasing distribution	+	-	-	-	-	-

Variations in Polymers

There are a variety of ways to alter properties

- Chemical Structure of Monomer Unit
- 3D Structure
- Different Monomer Units
- Length of polymer chains
- Distribution of polymer chain lengths



Introduction to GPC/SEC

GPC Separation Mechanism


The established technique to measure the distribution of chain lengths is

Gel Permeation Chromatography, **GPC**

also known as

Size Exclusion Chromatography, **SEC**

Separation based on size of the polymer in solution (hydrodynamic volume) only.



Introduction to GPC/SEC


GPC Separation Mechanism



Types of Chromatography

Interactive adsorption, partition, ion exchange, etc

Non-interactive GPC, SEC, GFC



GPC/SEC compared to HPLC

Isocratic elution, often single solvent systems

Typically uses organic eluents, not ACN, MeOH, IPA

Typically lower sample concentrations (0.1 – 0.2%)

Typically larger injection volumes (20 – 200 μ l)

Primary use is to measure molecular weight distribution

Resolution provided by **non-interactive** mechanism

Larger columns (300 x 7.5 mm industry standard)


Use of multiple columns in series (2 – 4)

Conventional GPC


Generating molecular weights

- Calibrate the column with a set of polymer standards
- Plot retention time (RT) versus peak log molecular weight ($\log M$)
- Calibration is used to generate molecular weights of unknowns
- **BUT** Molecular Weights are only **equivalent** to the Standards used
- So a Polystyrene calibration will give Polystyrene equivalent molecular weights for all samples analysed.

Chromatogram of narrow standards



GPC Calibration Curve



Molecular Weight Averages by GPC

$$\text{Number average } M_n = \frac{\sum N_i M_i}{\sum N_i}$$

M_n can be correlative with polymer colligative properties, e.g. freezing point depression

$$\text{Weight average } M_w = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$

M_w may be correlated with properties such as melt viscosity

$$\text{Z average } M_z = \frac{\sum N_i M_i^3}{\sum N_i M_i^2}$$

M_z may be correlated with properties such as toughness

$$\text{Polydispersity, } d = \frac{M_w}{M_n}$$

Polydispersity characterises the shape of the distribution

Introduction to GPC/SEC

GPC Separation Mechanism

To learn more about
Conventional GPC
please ask your local
Account Manager for the

Introduction to GPC Primer



An Introduction
to Gel Permeation
Chromatography and Size
Exclusion Chromatography


PRIMER

The Measure of Confidence

 Agilent Technologies

A Typical GPC/SEC System

Components of a GPC/SEC system



Columns

For polymer applications using organic solvents

PLgel series

- Covers all molecular weight ranges
- high pore volume
- high efficiency
- maximize resolution



Unequalled solvent compatibility

- easy transfer between polar and non polar eluents
- outstanding physical rigidity
- provides extended lifetimes that minimize downtime

Available in a range of column dimensions, mini-bore to prep

Columns

For polymer applications using aqueous solvents

PL aquagel-OH series

Available with mixed and individual pore sizes, and 5, 8 and 15 μm particle sizes, to cover a very wide range of molecular weights.

- chemically and physically stable
- “neutral” surface
- high performance analyses
- neutral, ionic & hydrophobic moieties

Available in a range of column dimensions, mini-bore to prep



Polymer Standards

Ideal reference materials

- generate accurate, reliable GPC/SEC column calibrations
- highest quality polymer standards
- extremely narrow polydispersity
- widest molecular weight range commercially available

Comprehensive range

- EasiVial
- EasiCal
- Traditional calibration kits
- Covers all molecular weight ranges
- Organic & aqueous GPC/SEC applications.



Scope of GPC

Application Range

GPC can be applied to an extremely broad range of polymers

- Wide range of solvents employed
 - Organic – Polar – Aqueous
- Wide range of operating temperatures
- Wide range of molecular weights

Dialkyl Phthalates

Columns: 2xPLgel 3 μ m 100Å, 300x7.5mm (PL1110-6320)


Eluent: THF

Flow Rate: 1.0ml/min

Detector: RI

KEY

1. Dioctyl phthalate
2. Dibutyl phthalate
3. Diethyl phthalate
4. Dimethyl phthalate
5. Toluene



Scope of GPC

Application Range


GPC can be applied to an extremely broad range of polymers

- Wide range of solvents employed
 - Organic – Polar – Aqueous
- Wide range of operating temperatures
- Wide range of molecular weights



Polyethylenes

Columns: 3xPLgel 10 μ m MIXED-B, 300x7.5mm (PL1110-6100)
Eluent: TCB
Flow Rate: 1.0ml/min
Temp: 160°C
Detector: RI



Agilent GPC/SEC Software

The fastest, easiest way to complete your GPC Workflow




Simplifies and facilitates your GPC/SEC process

- Everything required for GPC/SEC in a single package



Agilent GPC/SEC Software

Single solution for all your GPC/SEC Requirements

Control of Agilent GPC/SEC Systems



Collection from Agilent GPC/SEC Systems





Agilent GPC/SEC Software

Fast and Easy GPC/SEC

Baselines / integration regions added with single click

Narrow standard Column Calibrations easily created

Instant report or overlay data



Testing of material batches

Analysis of poly(styrene/butadiene) copolymers

Background:

A poly(styrene/butadiene) block copolymer mimics the properties of natural rubber.

Characteristics provided by:

- hard polystyrene chains surrounded by
- a network of rubbery polybutadiene

which provides strength and flexibility over a large temperature range.



Problem:

New batch of material is failing rheology testing

- Synthetic method to produce the copolymer has not changed
- End properties of the polymer significantly different to previous batches of the same material.

Solution:


Investigate molecular weight distribution using GPC

Testing of material batches

Analysis of poly(styrene/butadiene) copolymers

The GPC chromatogram shows evidence of low molecular weight material within the failed batch.

System : 1260 Infinity GPC System
Columns : 2 x PLgel 5 μ m MIXED-D
Temperature : 40 °C
Eluent : THF




This was attributed to residual homopolymer from the synthesis which was imparting different characteristics to that of the desired copolymer.


With the problem identified the cause was discovered and rectified.

Increasing the information from GPC/SEC

Most common detectors for GPC/SEC are *concentration* detectors:



RID



UV / DAD



Evaporative LSD


These provide information on the amount of polymer eluting from the column at any given time.

Increasing the information from GPC/SEC

There are two common *molecular weight sensitive* detectors.



Viscometer



Static Light Scattering Detector

These not only respond to the amount of polymer eluting but provide direct information on properties of the polymer relating to molecular weight.

Increasing the information from GPC/SEC

GPC/SEC provides critical information for the polymer chemist:

Distribution of chain lengths

(Relative molecular weights)

Further parameters can be determined by employing advanced detectors

The molecular **weight** (accurate or absolute)

The polymer's **size**

The polymer's **shape**

Can be used to investigate polymer **branching**

Polysaccharides

Investigation into structure

Complex polymers constructed from various sugar units.



Wide range of polysaccharides with large structural difference depending on source.

Pullulan, produced from starch by the action of a fungus, is composed of maltotriose units and has a linear structure.

Dextran, manufactured from sucrose by bacterial action has a highly branched structure.

Investigating their structure is of great interest for determining their properties and therefore end use application.

Light Scattering & Viscometry are employed to highlight molecular weight and structural differences



Summary

GPC/SEC is an indispensable analytical tool for the analysis of

- Polymers
- Polymer Additives

Applicable in a wide range of solvents allowing almost all commercial polymers to be analysed

The critical and often unique information attained is the population of polymer chain lengths in any given sample

The GPC/SEC experiment can be expanded to further investigate polymer properties – **more on this later**

Thanks for Listening

Any Questions



COMPREHENSIVELY BETTER POLYMER ANALYSIS

From research to quality control, each of your challenges in polymer characterization requires a unique approach. No matter which state-of-the-art solution you prefer, Agilent Technologies delivers proven answers. Experience the versatility, performance, and confidence that come with 35 years of GPC/SEC engineering excellence from Agilent Technologies—the Polymer Analysis People.



REPRODUCIBLE
Agilent 1260 Infinity GPC



ACCURATE
Agilent 1280 Infinity MDS GPC



DEDICATED
Agilent PL-GPC 50



POWERFUL
Agilent PL-GPC 220



EXTENSIVE
Agilent GPC/SEC Columns
and Calibrants

New
VERSATILE
Agilent GPC/SEC
Software

Learn more at
www.agilent.com/chem/gpc-solutions



© Agilent Technologies, Inc. 2012

The Measure of Confidence



Agilent Technologies

GPC/SEC Resources

Brochure	Agilent GPC/SEC solutions - comprehensively better polymer analysis	5990-8844EN
Application Compendia	Engineering polymers	5990-6970EN
	Polyolefin analysis	5990-6971EN
	Analysis of elastomers by GPC/SEC	5990-6866EN
	Biodegradable polymers	5990-6920EN
	Low molecular weight resins and prepolymers	5990-6845EN
	Excipient analysis	5990-7771EN
	Analysis of food additives by GPC/SEC	5990-8634EN
Primers	Introduction to GPC/SEC	5990-6969EN
	A guide to multi-detector GPC	5990-7196EN
Product guides	1260 Infinity GPC/SEC System	5990-9920EN
	1260 Infinity Multi-GPC/SEC System	5990-9921EN
	Agilent PL-GPC 50	5990-9937EN
	Agilent PL-GPC 220	5990-9926EN
	Agilent GPC/SEC Software	5991-0478EN
	Aqueous and polar GPC/SEC columns	5990-7995EN
	Organic GPC/SEC columns	5990-7994EN
Selection guide	GPC/SEC standards	5990-7996EN
	Quick guide for selecting columns and standards for GPC/SEC	5990-6868EN



Download at <http://www.agilent.com/chem/gpc-solutions>