

# Tracer Measurement by LC-MS

- Chemical Tracers such as Fluorobenzoic Acids (FBA) are used in reservoir fluids produced from the oil well being tested.
- Measurement of FBAs as tracers allows informed decision-making as to continue clean-up or begin production.
- FBAs are compatible in both aqueous fluid as salts and in organic base fluids as acids.
- Techniques such as GC, GC-MS, LC-UV, FTIR have been employed for the measurement of these tracers.
- Researchers continue to look for optimal methodologies with reduced sample preparation time, analysis time, multiplexing ability in newer technologies.
- LC-MS/MS using a QQQ system presents the ability to improvise the analytical method over existing technologies

# LC-MS/MS method for analysis of FBAs

- **Many of the FBAs are positional isomers hence they are isobaric. MS Resolving Power on its own is insufficient.**
- Challenge is to provide adequate resolution in the LC dimension for at least the critical FBAs (10-14 analytes).
- Sensitivity required in the low ppb to upper ppt for all analytes.
- Sample matrix is high in salt content.

# List of FBA Tracers


## MS/MS- Dynamic MRM Conditions

• Time	Name	Transition	CE	Fragmentor
• 1.78	2,6 DFBA	157-113	3	60
• 3.05	2 FBA	139-95	6	60
• 3.41	2,5 DFBA	157-113	6	60
• 3.77	2,4 DFBA	157-113	6	68
• 3.93	2 CBA	155-111	2	68
• 4.14	3 FBA & 4 FBA	139-95	10	88
• 4.41	2,4,5 TFBA	175-131	6	70
• 4.82	2,3,4 TFBA	175-131	6	58
• 5.02	3,4 FBA & 3,5 FBA	157-113	10	72
• 5.35	2 TFMBA	189-145	10	64
• 5.77	2,3,4,5 TFBA	193-149	6	64
• 5.96	3 CBA & 4 CBA	155-111	6	70
• 7.22	3 TFMBA	189-145	10	88
• 7.73	4 TFMBA	189-145	10	88



# LC-MS/MS of FBA Tracers

## 10 $\mu$ L injection of a 1 ppb standard mixture



# Results

- Chromatographic Separation of most critical tracers achieved using as UHPLC column. Requires Fast Acquisition rates in the MS – well suited to QQQ.
- Sensitivity limits of 10 ppt was achieved on tracers by using MRM analysis conditions on the QQQ MS.
- The samples were prepared by simple dilution and the LC-flow was diverted to waste for up to 1.5 min to eliminate the salt plug from the samples.
- Since samples are highly salty the source is cleaned more often than other methods – Cleaning the source on an Agilent LC-MS is practical and user-friendly.

# Conclusion

- A simple, sensitive, fast LC-MS/MS method with very little sample preparation has been developed for the analysis of tracers such as FBAs,
- The current analysis time is only **7 min** long in comparison to GC-MS or LC-UV methods that are ~ 20-30 min long/sample.
- Sample preparation for the LC-MS method involves simple dilution as opposed to GC-MS methods that require multi-step SPE and derivatization.