

What you don't see...
...CAN hurt you




Sample Prep for Today's
Analytical World

Tina Chambers
Christophe Deckers
Application Engineers, Sample Prep



Agilent Technologies



Traditional Reasons for Sample Prep



- Removal of interferences which would otherwise affect detection of analyte
- Concentration of an analyte to detectable levels
- Solvent switching into an analytically more compatible solvent



Today's World



**A
chromatographer's
dream: single peak
on flat baseline**




Agilent Technologies

Some samples are very conducive to reduced Sample Prep...


Single MS: SIM 283.8

100 fg HCB in Clean Matrix




SIM about equal to MS/MS in clean matrix


MS/MS: 283.8:213.9



300 fg HCB in Diesel




MS/MS 15x better than SIM in complex matrix




Agilent Technologies

...others not so much



Curtain plate after injection of 25 samples with extracts from raisins without cleanup

Common Contaminant #1: Particulates



Agilent Technologies

Why Filter?

Benefits of Filtration

Optimal Instrument Performance	Less System Downtime	Extend Column Lifetime	Greater Sample Integrity
			



Achieve lower detection limits

Agilent recommends filtering prior to chromatographic analysis to remove particulates from the sample



Why Filter the Sample?

Extreme Performance Requires Better Sample “Hygiene”




- Prevents blocking of capillaries, frits, and the column inlet (especially important for UHPLC)
- Results in less wear and tear on the critical moving parts of injection valves
- Results in less downtime of the instrument for repairs
- Produces improved analytical results by removing potentially interfering contamination
- “Functionalized” filtration can improve results over mechanical filtration by removing chemical interferences in addition to particulates



Agilent Technologies

Syringe Filter Benefits: Improved sub-2 micron LC Column Lifetime – With Human Plasma Extract



Agilent Captiva Filtration Products: Setting a New Standard in Filtration



Agilent Technologies

Captiva Premium Syringe Filters



- Certified to be free of UV detectable extractables on HPLC
 - PES & Glass Fiber also certified for LC/MS
- Color-coded boxes for easy identification
- Comprehensive portfolio to meet your application needs

Premium Syringe Filters						
Membrane	Diameter / Pore size					
	4 mm		15 mm		25 mm (* 28 mm)	
	0.2 µm	0.45 µm	0.2 µm	0.45 µm	0.2 µm	0.45 µm
PTFE	♦	♦	♦	♦	♦	♦
Nylon			♦	♦	♦	♦
PES	♦	♦	♦	♦	♦	♦
Regenerated Cellulose	♦	♦	♦	♦	♦	♦
Cellulose Acetate					♦	♦
Glass Microfiber			♦		♦*	
Depth filters: Glass / PTFE			♦	♦	♦	♦*
Depth filters: Glass / Nylon			♦	♦	♦	♦


Cellulose Acetate	Nylon	Glass Fiber/ PTFE	Glass Fiber	PES	PTFE	Regenerated Cellulose	Glass Fiber/ Nylon
-------------------	-------	-------------------	-------------	-----	------	-----------------------	--------------------



Cleanliness of Agilent Captiva PES Filters




ESI-Positive Mode



Captiva Premium Syringe Filters: Fastest Flow Rates



Agilent **RC**, **Nylon** and **PTFE** filters compared to competitors




Captiva Premium Syringe Filters: High Loading Capacity




- ✓ Due to nature of membrane, there is less clogging by particulates
- ✓ Assures more sample will be filtered

Capacity (volume) of 15 mm syringe filters over time (with Particulate-Laden Samples)



...Now back to our source.

Will filtration solve everything?




Curtain plate after injection of 25 samples with extracts from raisins without cleanup

Common Contaminant #2: Salts and Minerals

- Non-Volatile: will form residues in injection liner of ANY GC system
- Very water soluble: will carry through to detector in ANY LC system
- Ubiquitous in water (4g/L-30g/L) and body fluids (40-220 meq/L/day in urine)
- In LC/MS systems, ion suppression and precipitation in ion source/loss of sensitivity
- Easily eliminated by liquid-liquid extraction (conventional or solid-supported), but very polar compounds are difficult to extract



Salt Precipitation in Ion Source





Easily cleaned with MeOH, but downtime/run interruptions still exist



Agilent Technologies

Orthogonal ESI Ion Source Condition after 3000x Urine Dilute/shoot Injections



Nebulizer/Sprayer

Spray Shield/MS Inlet/Capillary



Agilent Technologies

6410 QQQ Sensitivity Results

Dilute/Shoot (1/10 dilution) versus SPE Sample Preparation

Opiates/Opioids

Compound	D/Shoot LLOQ (ng/ml)	SPE LLOQ (ng/ml)	ULOQ (ng/ml)
6-monoacetyl morphine	10	<1	1000
buprenorphine	10	1	1000
codeine	25	<1	1000
dihydrododeine	25	<1	1000
EDDP	10	<1	1000
fentanyl	1	<1	1000
heroin	10	<1	1000
hydrocodone	10	<1	1000
hydromorphone	5	<1	1000
meperidine	5	<1	1000
methadone	10	<1	1000
morphine	5	<1	1000
naloxone	5	<1	1000
naltrexone	10	<1	1000
N-desmethyltramadol	10	1	1000
norbuprenorphine	25	3	1000
norfentanyl	1	<1	1000
normeperidine	5	<1	1000
norpropoxyphene	5	<1	1000
o-desmethyltramadol	5	<1	1000
oxycodone	10	<1	1000
oxymorphone	5	<1	1000
propoxyphene	5	<1	1000
tapentadol	5	<1	1000
tramadol	1	<1	1000
trazodone	1	<1	1000

Sedatives/hypnotics

Compound	D/Shoot LLOQ (ng/ml)	SPE LLOQ (ng/ml)	ULOQ (ng/ml)
2-OH-ethylflurazepam	200	5	1000
7-aminoclonazepam	10	<1	1000
7-aminoflunitrazepam	5	<1	1000
alpha-OH-midazolam	10	<1	1000
alprazolam	10	<1	1000
a-OH-alprazolam	20	<1	1000
a-OH-triazolam	50	<1	1000
chlordiazepoxide	10	<1	1000
clonazepam	25 to 50	<1	1000
desalkylflurazepam	20	1	1000
diazepam	10	<1	1000
flunitrazepam	10	1	500
flurazepam	5	1	1000
lorazepam	50	20	1000
midazolam	10	<1	1000
nitrazepam	25	5	1000
nordiazepam	25	<1	1000
oxazepam	50	25	1000
temazepam	25	<1	1000
triazolam	5	<1	1000
zolpidem	5	<1	1000



Agilent Technologies

Sample Prep Solutions

- Liquid/Liquid or SSLE for non-polar compounds
- SPE for very complex matrices or large sample volumes
- QuEChERS for polar compounds (with or without dispersive SPE)

- Let's look at these in turn-



Sample Preparation – Supported LLE (SLE)

Hydromatrix™ - diatomaceous earth sorbent

- Composed of fossilized diatoms
- Purified at high temperatures
- High surface area for water adsorption
- Very polar surface



Chem Elut™ - pre-assembled cartridges with Hydromatrix

Combilut™- 96-well plate filled with Hydromatrix




Agilent Technologies

The SLE Process


Before Extraction



Apply Sample



Extract with Organic Solvent




Agilent Technologies

The Chem Elut Method

Aqueous
sample
being applied



Solid support
adsorbs water onto
high surface area
particles




Organic
extraction solvent



Agilent Technologies

The SPE mechanism




Agilent Technologies

The Bond Elut Method

Green = Blue and Yellow

Blue is more non polar than yellow

Blue is retained



Agilent Technologies

QuEChERS

First step extraction

Pictorial Representation of the QuEChERS Steps



- 1) Weigh sample, add water if needed, spike
- 2) Add 10ml ACN
- 3) Vortex
- 4) Add salt packet
- 5) Shake 1 minute
- 6) Centrifuge at 4,000 rpm for 5 minutes




Second Step – Dispersive SPE



7) Choose d-SPE kit based on matrix characteristics

8) Transfer 1-8ml aliquot, vortex 1 minute

9) Centrifuge



11) Analyze by GC/MS or LC/MS



Common Contaminant #3: Proteins and Lipids

- Low volatility in any GC systems, contaminate inlet liner, gold seal, front end of column
- Limited solubility of proteins in high organic reverse phase mobile phases in HPLC systems
- Both lipids and proteins coat column in any HPLC system: decrease/alteration of chromatographic separation
- Commonly encountered in human and animal blood and tissue samples, including mother's milk
- These will typically NOT show up as discreet contamination peaks



GC System Components Contamination with Biological Samples



Agilent Technologies

Sample Prep Solutions

- Protein Precipitation with or without Lipid Removal
- SPE
- QuEChERS for Fatty and Waxy Samples




Plasma contains high levels of proteins and lipids



Agilent Technologies

Manual Protein Precipitation



Faster or High Throughput with Filtration Plates

Captiva and Captiva ND are mechanical filters only


- Five times faster than manual (5 min vs. 26 min)
- No transfer steps
- Easy method transfer
- Filtration from 0.2um to 20um pore size



Agilent Technologies

Ion Suppression Regions

Protein precipitation sample
PCI with procainamide




Interference type	Salt/Polar ionics	Proteins/ Peptides	Lyso-phosphatidylcholines	Lipids and other hydrophobics
Typical Elution Conditions (C18 column)	At or near void with < 20% organic	10's of column volumes at 40% - 70% organic	10's of column volumes at 70% - 90% organic	10's to 100's of column volumes at > 90% organic
Short term effect (single injection)	Significant ion-suppression	Significant ion-suppression	Significant ion-suppression	Some ion suppression, however, usually retained on LC column)
Long term effect (multiple injections)	Unknown	Unknown	Decreased sensitivity, Increased variability	Decreased sensitivity, Increased variability
Likely long term causes	Ion source contamination	Ion source contamination	Ion source contamination, Some column build-up	Ion source contamination, Column build-up



Agilent Technologies

Captiva™ ND*Lipids* Lipid removal filtration plate



Features:

- Captiva particulate filter → removes protein interferences
- Proprietary Lipid Stripping Media → removes lipids
- Non-Drip Membrane → ease of use




Agilent Technologies

Lipid removal

Green = ppt only

Red = lipid-stripped ppt with CaptivaND^{Lipids}



Agilent Technologies

Sample Preparation Time Comparison

PPT (centrifugation) vs. Captiva ND Lipids



Captiva ND Lipids	Time (min)
Add 0.6 mL of MeOH and 0.2 mL of plasma sample to Captiva ND 96-well plate.	5
Mix each well with a pipette 5 times and apply vacuum for filtration.	
Directly transfer injection plate for analysis.	0
Total time required for sample preparation	5

This time comparison is based on the preparation of 96 samples.



Comparison to Protein Precipitation Antidepressants

Method Used:

Sample preparation with Captiva ND^{Lipids} filtration plates:

1. Add 600 µl methanol (with 0.1% formic acid) into the well of the Captiva ND^{Lipids} filtration plate.
2. Add 200 µl rat plasma (spiked, 1 – 128ng/mL) into each well, and process 6 replicates for each concentration, mix.
3. Process the plate by vacuum (15" Hg) for 4 min, and collect filtrate.

Sample preparation with standard protein precipitation:

1. Add 200 µl rat plasma (spiked, 1 – 128ng/mL) into 1.5 ml vial, 6 replicates for each concentration.
2. Add 600 µl methanol (0.1% formic acid) into the vial, vortex.
3. Centrifuge for 10 min at 14,000 rpm.
4. Take 600 µl aliquots of the supernatant from each vial, transfer to a 96-well sample plate.

Application note 01736



Agilent Technologies

Confidentiality Label
August 1, 2013

Absolute Recoveries; Sensitivity Enhancement Over Protein Precipitation

Recovery of antidepressant drugs from rat plasma


Compound	Extraction Method	Avg Rec%
Tranylcypromine Log P: 1.4	CaptivaND ^{Lipids}	67.1
	Protein Precipitation	29.1
Nomifensine Log P: 2.9	CaptivaND ^{Lipids}	77.4
	Protein Precipitation	38.1
Amoxapine Log P: 3.4	CaptivaND ^{Lipids}	77.8
	Protein Precipitation	38.8
Maprotiline Log P: 5.1	CaptivaND ^{Lipids}	82.8
	Protein Precipitation	71.1
Nefazodone Log P: 4.7	CaptivaND ^{Lipids}	75.0
	Protein Precipitation	63.1

Greatest enhancement is seen with lower logP compounds (<5) although enhanced sample cleanliness is observed with all



Protein and Lipid removal using the SPE approach

- Big proteins will be excluded
- Lipids adsorb and do not elute under normal elution conditions



Protein and Lipid Removal using QuEChERS

- Proteins precipitate during extraction step due to 50% ACN content
- dSPE kits for fatty/waxy samples available, contains C18 which adsorbs lipids
- Very fatty matrices (e.g. cooking oil) require addition of hexane in extraction step, see acrylamide application note for further details



Common Contaminant #4: Pigments and Dyes

- Various polarities and chemical structures, therefore various effects
- Anthocyanins (reds/blues) commonly found in berries, grapes, wine scavenged by PSA
- Chlorophyll (green) commonly found in foliage, vegetables, hay and carotenoids, scavenged by GCB
- SPE methods available for large number of chemical dyes
- dSPE (low capacity) or scavenging cartridges (high capacity)




Common Contaminant #4: Carbohydrates and polysaccharides

- Found in produce and beverages
- Very polar and often of high molecular weight (ex: starch)
- Viscous in organic solvents
- Best bet is with non-polar SPE, just like protein removal
- QUECHERS are a good option




Common Contaminant #6: Surfactants



- Can be polar, non-polar and ionic simultaneously.
- Often in high concentrations
- Can cause non specific binding on SPE cartridges
- Use ion-exchange SPE for ionic surfactants
- Use polar SPE sorbents like diol for non-ionic surfactants




What is SPE?




- Can be thought of as digital chromatography – compounds either bind or flow through
- Has a wide choice of sorbents with selectivities similar to sorbents used in HPLC, but it is not HPLC
- Many samples can be analyzed in parallel. Manifolds of 10 and 20 ports are readily available.
- Methods are simple to perform, but the process can be time-consuming, however **SPE** can be readily automated
- **SPE** methods give high selectivity, recovery and reproducibility






Conditioning



Clean-up



Before and after clean-up

SPE Manifolds



CaptiVac Vacuum Collar

96 Well

CaptiVac Collar: For use with Bond Elut 96 1mL and Captiva filtration plates



VacElut 12, 20 and SPS 24

Vacuum manifolds for SPE barrels.



Agilent Technologies

SPE Modes—"Digital Chromatography"

Analyte Adsorption (Bind-Elute)

Analyte(s) retained ($K_D \gg 1$)

Matrix unretained
and/or strongly retained ($K_D \sim 0$)
($K_D \gg 1$)

Preconcentration factor

Cleaner extracts

Load at 1-3 drops/sec (recovery \propto
1/flow)

Capacity issues may be more
important

Matrix Adsorption (Interference Removal)

Analyte(s) unretained ($K_D \sim 0$)

Matrix retained ($K_D \gg 1$)

No preconcentration advantage

Eluates may not be as clean

Sample loading often gravity fed

Used less often than analyte adsorption



Agilent Technologies

Solid Phase Extraction: Most efficient clean-up for complex matrices

Four Step of SPE:

- **Conditioning:** Preparation of the sorbent prior to sample addition
- **Retention:** Analytes of interest and other interferences adsorb onto the surface of the sorbent during sample addition (over 40 chemistries available!)
- **Washing:** Elimination of undesired interferences
- **Elution:** Selective desorption and collection of desired analytes from the sorbent



Retention Mechanisms

Polar

matrix is organic (ie organic phase from a liquid/liquid extraction)

analyte is water soluble

wash solvents are non-polar (hexane, methyl t-butyl ether etc)

elution solvents are polar (water, methanol, acetonitrile etc)

Non-polar

matrix is aqueous (foods, biological fluids)

analyte is organic soluble

wash solvents are aqueous

elution solvents are organic

Mixed mode

matrix is aqueous (foods, biological fluids)

analyte can be polar, hydrophilic, or hydrophobic

wash solvents are aqueous and organic

elution solvents are organic

sorbents are either mixed silica (such as C8/SCX) or polymer



Agilent Technologies

Retention Mechanisms

Cation exchange

matrix is aqueous (foods, biological fluids)

analyte is basic (cationic)

wash solvents are aqueous

elution solvents are :

- high ionic strength,
- pH is increased above the pKa of the target compound,
- competition with a cation (such as Na^+) with greater affinity for the sulfonic acid

Anion exchange

matrix is aqueous (foods, biological fluids)

analyte is acidic (anionic)



wash solvents are aqueous

elution solvents are :

- high ionic strength,
- pH is increased below the pKa of the target compound ,
- competition with an anion (such as SO_3^-) with greater affinity for the positively charge amine




Fundamental Steps for “Bind-Elute” SPE



Load



Wash



Elute



Agilent Technologies

Sample Pre-treatment or prewash



Agilent Technologies

Pre-extraction compound preparation


- Dilution
 - Reduce organic content before RP SPE
 - Reduce salt concentration before IEX SPE
 - Facilitate flow
 - Plasma is typically diluted 2:1 or 3:1
- Addition of internal standard
 - Deuterated or compound similar to analyte of interest
- Filtering, settling or centrifugation
 - Remove particulates to facilitate flow and reduce chance of clogging
- Buffering
 - Establish necessary pH
 - Neutralize analytes for RP SPE
 - Charge analytes for IEX SPE



Hydrolysis


For some clinical or drugs of abuse methods, many drugs, during metabolism, will take a glucuronide group to make the compound more polar allowing urinary excretion. Morphine is an example: morphine-3-glucuronide

Many GC applications will require chemical or enzymatic hydrolysis.



Large Sample Volumes

If sample has particulates (ground water, etc) sample needs to be filtered for cartridge applications. The smaller the bed mass diameter, the faster the SPE format will clog.



Protein Precipitation

- Used mostly for serum, plasma, and blood
- Precipitates cell wall and other proteins, which are then filtered out or centrifuged down
- Proteins are the single largest interference in biological samples
- Can also free protein-bound analytes
- Acetonitrile, methanol, and acetone are most commonly used.
- 3:1 to 10:1 ratio solvent:sample



Whole Blood

When ACN precipitation is not effective:

- Isotonic solutions with or without sonication
 - Dilution with water causes cells to burst from taking in water
 - Sonication facilitates the process
 - Pro: Does not add acids or salts
- See Certify Manual for other approaches
 - Acids
 - Zinc sulfates
 - others



Note – unlike many matrices, some compounds can be retained in the cell wall of the blood. Spiking whole blood samples, while still the best experimental approach doesn't simulate the difference between free and bound analytes.



Solid samples

- typically need to disrupt cell walls and/or separate bound analytes from the matrix
- Soil
 - Tumbling in a solvent is often used
 - Polar solvents (water, MeOH, acetone) used for polar and mid-polar compounds
 - Nonpolar solvents are used for more non-polar extractions
- Biological tissues and food
 - MSPD
 - Quechers
 - Homogenization with a solvent
 - freezing followed by milling (food)
 - esp when enzymes that can alter accurate analyte quantitation are present.



Non-aqueous samples

Oil samples are often diluted in hexane or other nonpolar solvents. Depending on the interferences either polar extraction techniques (diol, NH₂, etc) can be used




Agilent Sample Preparation Products



Agilent
Bond Elut Sample
Prep
Solutions
EXTRACTION

Bond Elut
Bond Elut Plexa
SPEC
Omix
SampliQ
Bond Elut QuEChERS



Chem Elut
Combilute
Hydromatrix

Agilent
Chem Elut Sample
Prep
Solutions
PURIFICATION



Agilent
Captiva Sample
Prep
Solutions
FILTRATION

Captiva ND Lipids
Captiva™



Technical Support



spp-support@agilent.com

1-800-227-9770, 3,3,3

www.agilent.com/chem/subscribe



Agilent Technologies



Thank You!



Agilent Technologies