

Rapid Low-level Identification and Quantitation of Host Cell Proteins




Dr. Shuai Wu
Application Scientist
Agilent Technologies



Dr. Linfeng Wu
Application Scientist
Agilent Technologies



Steve Madden
Software Product Manager
Agilent Technologies



HCP Introduction

Residual Host Cell Proteins (HCP) are process impurities remained in a purified drug product


- HCPs can influence product stability as well as other unwanted effects in final product
- FDA requires that HCP contaminants in the final product are measured and reported
- ELISA is still a widely accepted method for HCP quantification

Strengths

- Very sensitive
- High level of reproducibility
- High-throughput (plate format, automation)

Challenges

- Lack of specificity, no identification of individual HCPs
- Quantitation is based on a cohort of HCPs
- Lack of coverage for non-immunoreactive HCPs



Modified with courtesy of Wikipedia
https://en.wikipedia.org/wiki/ELISA#/media/File:ELISA_diagram.png

LC/MS as a Solution for Host Cell Protein Analysis




Advantages

- Doesn't require protein specific antibodies
- Improve early purification process
- Identify individual protein including immunogenic HCPs
- High analytical sensitivity (low ppm)
- Provides both qualitative and quantitative information


Challenge for LC/MS Analysis of HCPs

- Low abundant HCP peptides co-elute with very intense “product” mAb peptides
 - Need broad dynamic range and great separation
 - Reproducible chromatography needed for good quantitation




**Able to measure 5 orders
in-spectrum dynamic range**

6545XT AdvanceBio LC/Q-TOF




Intact Protein

Industry leading data quality and sensitivity at the intact level



Peptide Mapping and PTMs

Ready for routine confirmation or digging deeper with Iterative MS/MS




Glycan Profiling

Full released glycan workflow with BioConfirm B.09



Iterative MS/MS Acquisition

An Easy way to Dig Deeper



Automatically look for different peptides with each injection

Host Cell Protein Analysis

Dr. Linfeng Wu Dr. Shuai Wu



Experiment to Evaluate HCP Analysis

- ❑ Spike-in UPS2 standards in purified CHO-cultured mAb before digestion
 - UPS2 mix of 48 proteins spanning 6 orders of magnitude
 - standard protein levels from 0.0004 to 313 ppm
 - mAb without UPS2 was used as a negative control
- ❑ Automated sample preparation using AssayMAP Bravo
 - Denaturation, reduction, alkylation, digestion, desalting, fractionation
- ❑ 60min LC method on an AdvanceBio Peptide Plus column (2.1x150 mm)
- ❑ Data-dependent acquisition by 6545XT AdvanceBio LC/Q-TOF
 - Iterative MS/MS vs. Auto MS/MS

AssayMAP Bravo



6545XT
AdvanceBio
LC/Q-TOF



Iterative MS/MS Decision Engine Improves Protein Identification


Comparison of Unique Peptide Sequences between Iterative MS/MS and Auto MS/MS

Protein Accession	Protein Spiking Level (ppm)	3 Injections per Method	
		Iterative MS/MS	Auto MS/MS
mAb_HC	NA	419	382
mAb_LC	NA	201	186
ALBU_HUMAN_spike	313.0	46	43
CAH2_HUMAN_spike	137.3	19	15
CAH1_HUMAN_spike	135.6	9	9
LEP_HUMAN_spike	76.2	4	1
HBB_HUMAN_spike	74.8	12	10
HBA_HUMAN_spike	71.3	7	6
UBIQ_HUMAN_spike	50.0	6	6
CO5_HUMAN_spike	40.3	4	4
CATA_HUMAN_spike	28.1	2	2
SUMO1_HUMAN_spike	18.3	3	1
NQO1_HUMAN_spike	14.5	2	0
PRDX1_HUMAN_spike	10.4	3	0
PPIA_HUMAN_spike	9.5	4	4
MYG_HUMAN_spike	8.0	2	1

- Iterative MS/MS acquisition method identified more unique peptide sequences

Excellent Chromatography Reproducibility and Dynamic Range

Overlaid chromatograms of three LC-MS/MS runs



In-spectrum Dynamic Range
~ 4.3 orders of magnitude

Peptide	Mass Error (ppm)	Intensity	Intensity %RSD	Protein spiking Level (ppm)	Protein Name
ALELFR	-1.1	6.76E+03	10.3%	8	Myoglobin
TIAQDYGVLK	-1.8	1.51E+05	6.2%	10.4	Peroxiredoxin 1
SAVTALWGK	4.8	1.36E+05	6.0%	74.8	Hemoglobin subunit beta
EPQVYTLPPSR	1.0	1.38E+08	1.2%	NA	mAb


Report Directly from Protein Metrics Software -- Byologic

	[Sample] ←	mAb	mAb_spike ↓
Protein alias name			
mAb_HC		100	100
mAb_LC		41.2	31.6
ALBU_HUMAN_spike			0.0973
CAH2_HUMAN_spike			0.0586
CAH1_HUMAN_spike			0.0354
HBB_HUMAN_spike			0.0199
HBA_HUMAN_spike			0.0157
PPIA_HUMAN_spike			0.0146
LEP_HUMAN_spike			0.0133
UBIQ_HUMAN_spike			0.0073
PRDX1_HUMAN_spike			0.00319
CATA_HUMAN_spike			0.00266
NQO1_HUMAN_spike			0.00117
SUMO1_HUMAN_spike			0.000967
CO5_HUMAN_spike			0.000752
MYG_HUMAN_spike			0.000467


← 8 ppm spiking level

Identification of all the spiking proteins above 8 ppm coupling with Iterative MS/MS

Stacked Barchart of Normalized Protein per Sample



Label-free Quantification



Improve Identification by Fractionation using AssayMAP Bravo



Fractionation: Using AssayMAP

v1.0

A. Run Plate Stacking Utility


Number of Fractions

B. Application Settings

Number of Full Columns of Cartridges

Step	Conduct Step?	Volume (μL)	Flow Rate (μL/min)	Wash Cycles
Initial Syringe Wash	<input checked="" type="checkbox"/>			<input type="text" value="3"/>
Prime Cartridges	<input checked="" type="checkbox"/>	<input type="text" value="100"/>	<input type="text" value="300"/>	<input type="text" value="1"/>
Equilibrate Cartridges	<input checked="" type="checkbox"/>	<input type="text" value="50"/>	<input type="text" value="10"/>	<input type="text" value="1"/>
Load Sample	<input checked="" type="checkbox"/>	<input type="text" value="100"/>	<input type="text" value="5"/>	<input type="text" value="3"/>
Collect Flow Through	<input checked="" type="checkbox"/>			
Cartridge Cup Wash	<input checked="" type="checkbox"/>	<input type="text" value="50"/>		<input type="text" value="1"/>
Internal Cartridge Wash	<input checked="" type="checkbox"/>	<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="3"/>
Add Wash to Flow Through	<input checked="" type="checkbox"/>			
Predispense Elution Buffer	<input checked="" type="checkbox"/>	<input type="text" value="25"/>		
Elute Fraction 1		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Elute Fraction 2		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Elute Fraction 3		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Elute Fraction 4		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Elute Fraction 5		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Elute Fraction 6		<input type="text" value="25"/>	<input type="text" value="5"/>	<input type="text" value="1"/>
Final Syringe Wash	<input checked="" type="checkbox"/>			<input type="text" value="3"/>

C. Deck Layout




D. Labware Table

Deck Location	Labware Type
1	96AM Wash Station
2	96AM Cartridge & Tip Seating Station
3	96 AbGene 1127, 1 mL Deep Well, Square Well, Round Bottom
4	Stack of n*: 96 Greiner 650201, U-Bottom Standard PolyPro
5	96 Greiner 650201, U-Bottom Standard PolyPro
6	96 Greiner 650201, U-Bottom Standard PolyPro
7	Stack of n*: 96 Greiner 650201, U-Bottom Standard PolyPro
8	96 Greiner 650201, U-Bottom Standard PolyPro
9	96 Greiner 650201, U-Bottom Standard PolyPro

* The number of plates in a stack equals the Number of Fractions (0 to 6).

Improve Identification by Fractionation using AssayMAP Bravo



Improve Identification by Fractionation using AssayMAP Bravo




Protein Accession	Protein Spiking Level (ppm)	# Unique Peptide Sequences	
		1D LC-MS/MS	HPRP + LC-MS/MS
ALBU_HUMAN_spike	313.0	46	79
CAH2_HUMAN_spike	137.3	19	32
CAH1_HUMAN_spike	135.6	9	15
LEP_HUMAN_spike	76.2	4	6
HBB_HUMAN_spike	74.8	12	22
HBA_HUMAN_spike	71.3	7	14
UBIQ_HUMAN_spike	50.0	6	9
CO5_HUMAN_spike	40.3	4	6
CATA_HUMAN_spike	28.1	2	14
SUMO1_HUMAN_spike	18.3	3	11
NQO1_HUMAN_spike	14.5	2	8
PRDX1_HUMAN_spike	10.4	3	9
PPIA_HUMAN_spike	9.5	4	11
MYG_HUMAN_spike	8.0	2	2
CYB5_HUMAN_spike	7.6	0	2
EGR_HUMAN_spike	3.0	0	1
SYHC_HUMAN_spike	2.7	0	5
KCRM_HUMAN_spike	2.0	0	3


← 2 ppm

Identification of all the spiking proteins above 2 ppm with Iterative MS/MS acquisition

MS/MS Spectrum of a Peptide from a Protein Spiked at 2 ppm




Improved Identification of Endogenous HCPs by AssayMAP Bravo Fractionation



HCP Discovery

LC/Q-TOF
+ Protein
Metrics SW

- Ideal platform for HCP Discovery
- Excellent chromatographic reproducibility and dynamic range
- Easy processing of DDA data for simultaneous ID and semi-quant



Iterative
MS/MS

- Automated approach to easily dig deeper
- Identified all spiked proteins down to **8ppm**





AssayMAP
Bravo

- Personal sample prep workbench
- Fractionation allowed detection down to **2ppm**

From Discovery to Targeted Protein Quantification


Protein Discovery

6545XT
AdvanceBio
LC/Q-TOF




Targeted Quantification by MRM

6495B
Triple Quadrupole
LC/MS



Targeted Protein Quantification by Isotope Dilution Strategy and MRM



Skyline Automation

Skyline - HSA-Feb13-auto.sky

File Edit View Settings Tools Help

Targets

Replicates: opt

HSA

AAATECCQAADK

HPYFYAPPELLFFAK

LVNEVTEFAK

Tools

- SRM Collider
- QuaSAR
- Agilent Automation
- External Tools...
- Immediate Window
- Options...

opt Step-A_HSA-A1

Step -2

Step -1

AAATECCQAADK - 686.2870++

Step

Step 2

Intensity (10³)

Peak Area (10³)

CE Replicate

Opt CE

1600

1400

1200

1000

Skyline Automation - D:\chris\HSA-Feb13-auto.sky

Project setting

Template method D:\MassHunter\Methods\chris\A

Folder D:\MassHunter\Data\chris

Name demo

Timestamp

Action selections

- ☒ Step-A (Update Retention Times)
- ☒ Step-B (Optimize Collision Energy)
- ☒ Step-C (Export method, create worklist)
- ☒ Execute workli

Step-A Step-B Step-C

(Step-A)

Export method name: Demo

☒ Single metho

☐ One method per prot

☐ Multiple metho ☐ Ignore proteins

Max transitions per sample injection: 200

Optimizing: None

Method type: Standard

Dwell time (ms): 5

Create Project

Submit to Study Manager

Close

Project is created at: D:\MassHunter\Data\chris\demo.s. (3/4/2013 7:57:09 AM)

Total samples: 12

Update Retention Time

Optimize Collision Energy

Export Dynamic MRM Method

Step-A > (Step-B) > (Step-C)

Sample prefix Standard Start position P1-A1 All samples in same position

Edit	Sample Name	Sample Position	Method	Data File
1	Standard1	P1-A1	Demo.m	Step-A_Demo1.d

Step-B

Sample prefix Standard Start position P1-A1 All samples in same position


Edit	Sample Name	Sample Position	Method	Data File
1	Standard1	P1-A1	Demo_0001.m	Step-B_Demo_00011.d

Step-C



Sample prefix Samples Start position P1-A1 All samples in same position

Edit	Sample Name	Sample Position	Method	Data File
1	Samples1	P1-A1	Demo.m	Step-C_Demo1.d
2	Samples2	P1-B1	Demo.m	WorklistData2.d
3	Samples3	P1-B2	Demo.m	WorklistData3.d
4	Sample4	P1-B3	Demo.m	WorklistData4.d
5	Sample5	P1-B4	Demo.m	WorklistData5.d

MRM Assay to Evaluate Targeted HCP Quantification





Three Peptides Selected to Evaluate Targeted HCP Quantitation



Standard Curve by Heavy SIL Peptide Standard

Peptide sequence: LLLEYLEEK





13 min LC-MS analysis that includes sample on-line desalting

Standard curve ranges from 6.25 amol/ μ g to 125 fmol/ μ g (0.24 ppm ~ 4765 ppm)

Standard Curve by Heavy SIL Peptide Standard

Peptide sequence: VFDVIIR





13 min LC-MS analysis that includes sample on-line desalting

Standard curve ranges from 6.25 amol/ μ g to 125 fmol/ μ g (0.35 ppm ~ 7002.5 ppm)

Standard Curve by Heavy SIL Peptide Standard


Peptide sequence: DPGVLDR




13 min LC-MS analysis that includes sample on-line desalting

Standard curve ranges from 12.5 amol/μg to 125 fmol/μg (0.13 ppm ~ 1336.3 ppm)

Absolute Quantification Using Heavy SIL Peptides as Internal Standards



Absolute Quantification of Targeted Proteins



Targeted Protein	SUMO1_HUMAN	SYHC_HUMAN	Protein S100-A11 (G3HUU6)
Protein MW	38,815 Da	58,233 Da	11,241Da
Peptide sequence	LLLEYLEEEK	VFDVIIR	DPGVLDR
Spiked protein level	18.3 ppm	2.7 ppm	NA
Measured protein level	10.1 ppm	1.2 ppm	1.6 ppm

Summary of Workflows

HCP ID

Detection down
to 2 ppm

Sample Prep
**AssayMAP
Bravo**

Detection
**6545XT
AdvanceBio
LC/Q-TOF**

Iterative MS/MS

Processing
Protein Metrics

Targeted Quantitation

Sub-ppm monitoring

Sample Prep
**AssayMAP
Bravo**

Detection
**6495B Triple
Quad**

Skyline integration for
method development

Processing
**MassHunter
Quant**

Thank you!

