

PLOT TWIST: Using PLOT Columns for the Analysis of Gasses and Very Volatile Compounds

Porous Layer Open Tubular
columns

Mark Sinnott
Application Engineer
November 2024



Table of Boiling Point Fractions

Carbon #	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	C ₉	C ₁₀	C ₁₁	C ₁₂	C ₁₃	C ₁₄	C ₁₅	C ₁₆
Bpt of n-Paraffin at 760 mm Hg																
Centigrade	-161	-89	-42	-0.5	+36	69	98	126	151	174	196	216	235	253	270	287
Fahrenheit	-259	-127	-44	+31	97	156	209	258	303	345	384	421	421	488	519	548

Refinery off gas



Naphtha solvents



Cracking naphtha



Reformate



Gasoline



Kerosine, diesel



Jet fuel



Gas oil



Lube oil




Pitch asphalt



Wax




Residues



Wall Coated Open Tubular (WCOT) Columns

$$\text{WCOT} = \underline{W}_{\text{all}} \underline{C}_{\text{coated}} \underline{O}_{\text{open}} \underline{T}_{\text{ubular}}$$




$K_c \Rightarrow$ large

$K_c \Rightarrow$ small



Midrange Hydrocarbon Separation (C4 – C12)




Column: DB-DHA
100 m x 0.25 mm id, 1.0 μ m DB-1
5 meters of DB-5, 0.25 mm id, 1.0 μ m
Carrier: Helium at 24 cm/s
Oven: **5 °C** for 10 min
5 to 50 °C at 5 °C/min
50 °C for 52 min
50 to 200 °C at 1.4 °C/min
Injector: Split 1:100. 250 °C, 0.3 μ L unleaded gasoline
Detector: FID, 250 °C
Nitrogen make up gas at 30 mL/min

Ethylene Analysis on Nonpolar WCOT

ETHYLENE


1. Methane
2. Ethylene
3. Ethane
4. Propylene

Column: DB-Petro 100
100 m x 0.25 mm x 0.5 μ m
Carrier: He at 28.5 cm/s
Oven: **-40 °C (20 min)** $\rightarrow$ 120 °C (5 °C/min)
Injector: Split 1:10, 100 μ L, 200 °C
Detector: FID, 250 °C



- ☐ Column Nth = 350,000
- ☐ Oven start at -40 °C
- ☐ Poor separation of C2s
- ☐ Poor separation of C3s

Ethylene Analysis on Alumina PLOT



- Column Nth $\approx 75,000$
- Oven initial at 35 °C
 - No Cryo!
- Excellent separation of C2s/C3s/C4s

Column: **GS-Alumina**
50 m x 0.53 mm I.D.

J&W P/N: **115-3552**

Carrier: Helium at 11 mL/min, measured at 35°C

Oven: 35°C for 2 min
35-190°C at 6°/min
190°C for 3 min

Injector: Split 1:30, 200°C
0.2 mL of trace hydrocarbons in ethylene

Detector: FID, 200°C
Nitrogen makeup gas at 20 mL/min


1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Cyclopropane
6. Propylene
7. Isobutane
8. Acetylene
9. *n*-Butane
10. Propadiene
11. *trans*-2-Butene
12. 1-Butene
13. Isobutylene
14. *cis*-2-Butene
15. Isopentane
16. *n*-Pentane
17. Propyne
18. 1,3-Butadiene
19. 1-Pentene
20. *n*-Hexane

What Is a PLOT column?

- PLOT = **P**orous **L**ayer **O**pen **T**ubular
- Porous layer (surface) on the inner wall capillary
- Deposition of porous particles on the wall from suspension
- Porous layer serves as a stationary phase
- Gas-solid adsorption chromatography
Differences in analyte distribution between carrier gas and surface of the adsorbent (stationary phase)
- Gas-liquid chromatography – WCOT
Differences in solubility in liquid phase (stationary phase)

Advantages of GSC:

- Improved selectivity (α)
- Retention (k) of highly volatile solutes



Selectivity Interactions in PLOT Phases

Shape/size




Surfaces

Zeolites

Bonded carbon
Molecular sieves



Porous polymers



Al_2O_3

Bonded silica



*Molecular sieve by Bearfruitidea from the Noun Project

Creative Commons: Daniel Kulinski

Porous Layer Open Tubular (PLOT) Columns


Advantages

- Retention for high vapor pressure solutes
 - No cryogenics needed
 - Useful for multicolumn and valve switching, as well as heart cutting techniques
- Variety adsorbents
 - Porous polymers, alumina, molecular sieve, silica
 - Selectivity for isomeric compounds (HC isomers)
 - Gas and solvent separations
- Efficiency of a capillary
 - Preferred over packed columns
 - Packed columns have advantages in specific applications




PLOT – Some Limitations and Considerations

- Lower sample capacity than thick film (1 to 5 μm) liquid phase columns
The 0.53 mm PLOT columns offer highest loadability
- More limited maximum temperatures than siloxanes – limiting analyte scope
- Highly active – some compounds do not elute or require extensive column conditioning ($\text{CO}_2/\text{H}_2\text{O}$)
- Molecular rearrangements above Maximum Allowable Operating Temperature (MAOT) – selectivity might change
Do not exceed MAOT
- Particles may dislodge under mechanical/flow stress – but solutions do exist



PLOT – Voids/gaps

- Different manufacturing process than for WCOT columns
- Normal to find some Voids/gaps are normal
 - Areas of uneven coating
- Visual examination is performed to ensure no excessive number of voids/gasp
 - Picture to the right is showing unacceptable gaps/voids
- Each column is individually tested to ensure proper chromatographic performance



This is considered excessive!

Agilent J&W PLOT Column Portfolio

Porous polymers

HP-PLOT Q and U
PoraBOND Q and U
PoraPLOT Q, U, and S
GS-Q
PoraPLOT Q-HT

Zeolites

HP-Molesieve
CP-Molsieve 5A

Oxygenates

CP-Lowox
GS-OxyPLOT

Alumina


Al₂O₃/KCl, HP-PLOT Al₂O₃/KCl
Al₂O₃/Na₂SO₄, HP-PLOT Al₂O₃ S
GS-Alumina
Select MAPD

Porous silica

GS-GasPro - No Particles
SilicaPLOT - No Particles

Graphitised carbon

CarboPLOT P7
CarboBOND
GS-CarbonPLOT - No Particles





Compounds and PLOT Column Selection

- Permanent gases - Noble gases, O₂, N₂, H₂, CO, CO₂, CH₄	Molsieve	Select Permanent Gases
- Light hydrocarbons - C1 – C8, saturated, unsaturated, - Aromatics, natural gas, C2, C3, C4 streams	Alumina	Porous carbon
	Silica	Porous polymers
- Volatile sulfur compounds - H₂S, COS, mercaptans, sulfides, disulfides	Select Low Sulfur	Porous polymers
	Silica	
- Oxygenates - Alcohols, ketones, ethers	Lowox	OxyPLOT
- Solvents - Oxygenates, aromatics, alkanes, - Chlorinated hydrocarbons	Porous polymers	
Chlorinated and fluorinated hydrocarbons	Silica	Porous polymers

Porous Polymer PLOT Columns

- Analyses of polar and nonpolar volatile compounds
Oxygenates, gases, halogenated compounds, hydrocarbons C1 to C6, ketones, solvents
- Most versatile of adsorbent materials
- Suitable for aqueous injections and polars
(in split and direct mode through the inlet)
- Elution of water as a sharp peak, plus quantification of water
(on the polar U type)
- Recommended for column switching systems
- Styrene/divinyl benzene polymer = Nonpolar “Q type”
- Styrene-glycol methacrylate copolymer = Polar “U type”
- Divinylbenzene/vinylpyridine polymer = Midpolar “S type”


Styrene/divinylbenzene polymer



Porous Polymer PLOT columns



HP-PLOT Q, GS-Q, and PoraPLOT Q

- Bonded polystyrene-divinylbenzene
- Separation ethane/ethylene, propane/propylene
 - Percentage levels, not ppm impurities in C2, C3 matrix
- H₂S and COS
 - Refinery gas analysis (RGA)
 - 100 ppm - percentage, **but not low ppm due to strong surface interactions for H₂S**
- Precolumn for O₂, N₂, CH₄, CO separations on Molsieve column
- Poor C4 isomer separation
 - Alumina preferred
- Not sensitive to moisture, stable retention, no selectivity changes
- Good choice for solvents and CFCs



Refinery Fuel Gas, Porous Polymer Separation


HP- PLOT Q




To Molsieve column for further separation

Aqueous Injections on Porous Polymer Q-PLOT


Before



After 5 x 5 μ L water, **splitless** at 80 °C



Retention times
unchanged
for all compounds




Samples containing
water can be analyzed

Isothermal and short
cycle time


Agilent “BOND” Technology and Traditional Particle Traps

Particle PLOT Q




Particle size: 0.1 to 2 μm
Number: $> 10^{12}$

PoraBOND Q




Bonded polymer layer
Reduced need particle
trap



Solvent Analysis on PoraBOND Q

Conditions


Technique	: GC-wide-bore
Column	: Agilent PoraBOND Q, 0.53 mm x 25 m fused silica PLOT (df = 10 µm) (Part no. CP7354)
Temperature	: 100 °C (2 min) → 300 °C, 5 °C/min
Carrier Gas	: He, 25 kPa (0.25 bar, 3.5 psi)
Injector	: Split, T = 250 °C
Detector	: FID, T = 250 °C
Sample Size	: 0.5 µL
Concentration Range	: 0.1% per compound
Solvent Sample	: DMSO



Peak identification

1. methane
2. methanol
3. ethanol
4. acetonitrile
5. acetone
6. dichloromethane
7. 2-propanol (isopropanol)
8. dimethyl sulfide
9. diethyl ether
10. 1-propanol
11. pentane
12. 2-butanone
13. trichloromethane
14. tetrahydrofuran
15. ethyl acetate
16. 2-methoxyethanol
17. isobutanol
18. butanol
19. hexane
20. benzene
21. trichloroethylene
22. cyclohexane
23. 1,4-dioxane
24. pyridine
25. N,N-dimethylformamide
26. dimethyl sulfoxide



Integrated Particle Trap (PT)



Assurance and safety
for valves and detectors
in valve switching
and backflush modes

Proof of Agilent J&W PLOT PT Column Performance

- 150 °C → 250 °C at 20 °C/min
- Pressure 3x higher than optimal
- 15 cycles
- 10x carrier gas off/on



Agilent PLOT PT

- Zero particle shedding
- No unions and fittings – no blockage
- Multicolumn valve switching
- Particle traps are integrated on both ends – supports backflush and CFT technology
- Similar selectivity, efficiency, and peak shape performance to existing Agilent porous polymer PLOT columns
- Minimum method redevelopment required
- Available for the most popular porous polymer PLOT columns

PLOT-PT Chromatographic Performance

Column: PoraBOND Q PT, 30 m × 0.25 mm, 3 μm

Carrier : Helium, 43 cm/s at 90 °C

Oven : 90 °C to 140 °C at 10 °C/min


140 °C for 6 min

140 °C to 200 °C at 5 °C/min


200 °C for 10 min

Injection: Split, 250 °C, split ratio 1:160

Detector: MSD, 280 °C transfer line, full scan at m/z 30-350



Selectivity Comparison of Porous Polymer PLOT



GC Conditions	
Columns	Agilent J&W PoraPLOT Q, 25 m × 0.53 mm × 20 µm (p/n CP7554) Agilent J&W PoraPLOT S, 25 m × 0.53 mm × 20 µm (p/n CP7574) Agilent J&W PoraPLOT U, 25 m × 0.53 mm × 20 µm (p/n CP7584) Agilent J&W PoraPLOT Q, 10 m × 0.53 mm × 20 µm (p/n CP7553) Agilent J&W PoraBOND Q, 25m × 0.53 mm × 10 µm (p/n CP7354)
Carrier	Helium, constant flow, 5 mL/min
Oven	50 °C (1.0 min), Ramp 10 °C/min to 200 °C (2 min) *Ramp to 180 °C for the J&W PoraPLOT U
Inlet	SSL Inlet, split mode, 200 °C, split ratio 20:1
Inlet Liner	Ultra Inert, split, low pressure drop, glass wool (p/n 5190-2295)
GC/FID	Agilent 7890B GC equipped with FID
Sampler	Agilent 7693 automatic liquid sampler
FID Conditions	
Temperature	200 °C
Hydrogen	30 mL/min
Air	350 mL/min
Col + Make Up	25 mL/min

5994-3485EN

Alumina Adsorbent and PLOT Columns

Best selectivity for light hydrocarbon separations

Applications


- General C1 to C6 (C9) hydrocarbons
- Natural gas
- Ethylene streams, impurities
- Impurities in propylene
- Butylene streams, impurities
complex C4 composition
- Environmental hydrocarbon distributions
- Not suitable for polars
- Sample needs to be dry



Separation and conditions

- Degree of saturation
Elution: alkane, alkene, alkyne, (dialkenes)
- Selectivity: types of deactivation
KCl, Na₂SO₄ and proprietary
- Selectivity: column flow and temperature
- Moisture carrier gas (and sample)
Retention and selectivity impacted
Gas Clean moisture filters are recommended
- Not suitable for polar solutes

Selectivity Difference between KCl and Na₂SO₄



1. Methane
2. Ethane
3. Ethene
4. Propane
5. Propene
6. Cyclopropane
- 7. Ethyne
8. Isobutane
9. Propadiene
- 10. N-butane
11. T-2-butene
12. 1-butene
13. Isobutene
14. C-2-butene
- 15. 1,3-butadiene
- 16. Propyne

KCl for:

Impurities in ethylene

Impurities propylene

Impurities in 1,3 butadiene


Na₂SO₄ for C4 feeds

Alumina columns tend to exhibit overload for highly unsaturated analytes more quickly

Impurities in 1,3-Butadiene

50 m x 0.32 mm Al₂O₃ / KCl, 5 µm

100 °C → 200 °C, 6 °C/min



Absorption and Response of Hydrocarbons


The case for Al_2O_3 MAPD – Methyl Acetylene - PropaDiene



Be aware

Alumina columns exhibit partial irreversible absorption of highly unsaturated HC

–
No 100% level calibration



Standard Al_2O_3

Select Al_2O_3 MAPD

1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Propylene
6. Butane
7. Propadiene
8. Acetylene
9. tr-2-Butylene
10. Butylene
11. iso-Butylene
12. cis-2-Butylene
13. 1,3-Butadiene
14. Methyl acetylene


Improved Response on Select for MAPD



Select for MAPD
30 m x 0.53 mm

Peak Identification


1. Ethylene
2. Isobutane
3. n-butane
4. Propadiene
5. Cis-2-butene
6. 1,2-butadiene
7. 1,3-butadiene
8. Methyl acetylene
9. Ethyl acetylene
10. unknown



Courtesy: J. Luong, Dow Chemical Canada

Silica-based PLOT

Inorganic gases on J&W GC-GasPro
30 m x 0.32 mm




- GS-GasPro and CP-SilicaPLOT
- Light hydrocarbon separation: C1 to C4
- Extended hydrocarbon range compared to other PLOTs
 - Up to ~C12
- Inert enough for low ppm light sulfurs, H₂S, COS, mercaptans
- CFCs
- Not sensitive to moisture in carrier gas
- MSD compatible
- MAOT of 260/300 °C (GasPro)
- MAOT of 225/225 °C (SilicaPLOT)

GS-GasPro - Hydrocarbons

Column: GS-GasPro
30 m x 0.32 mm I.D.
J&W P/N: 113-4332
Carrier: Nitrogen at 29 cm/sec
Oven: 40°C for 2 min
40 - 200°C at 10°/min
Injector: 200 µL
Detector: FID


1. Methane
2. Ethane
3. Ethylene
4. Acetylene
5. Propane
6. Cyclopropane
7. Propylene
8. *cis*-2-Butene



GS-GasPro - Hydrocarbons



Column: GS-GasPro
60 m x 0.32 mm I.D.
P/N: 113-4362
Carrier: Helium at 40 cm/sec (3.3 mL/min),
measured at 80°C
Oven: 80°C for 0.5 min
80-175°C at 25°/min
175°C for 2 min
175-250°C at 25°/min
Injector: Split 1:17, 250°C
500 µL injection of 100 ppmV
SUMMA cannister mixture
Detector: FID, 275°C
Nitrogen makeup gas at 32 mL/min



27. 2-Methylpentane
28. 3-Methylpentane
29. Isoprene
30. *n*-Hexane
31. 4-Methyl-1-pentene
32. *trans*-2-Hexene
33. 2-Methyl-1-pentene
34. *cis*-2-Hexene
35. 2,4-Dimethylpentane
36. Methylcyclohexane
37. 2,3-Dimethylpentane
38. 2-Methylhexane
39. 3-Methylhexane
40. *n*-Heptane
41. Benzene
42. Isooctane (2,2,4-trimethylpentane)
43. 2,3,4-Trimethylpentane
44. 3-Methylheptane
45. 2-Methylheptane
46. *n*-Octane
47. Toluene
48. *n*-Nonane
49. Ethylbenzene
50. *m*-Xylene
51. *p*-Xylene
52. *o*-Xylene
53. Styrene
54. Isopropylbenzene (cumene)
55. *n*-Propylbenzene
56. 1,3,5-Trimethylbenzene
57. 1,2,4-Trimethylbenzene

GS-GasPro – Permanent Gases Cryogenic (-80 °C)

Column: GS-GasPro, 60 m x 0.32 mm I.D.

J&W P/N: 113-4362

Carrier: Helium, 20 psig (constant pressure)


Inlet: 200°C, 1:15 split

Sample: 100 µL, manual inject

Detector: TCD, 250°C, background compensated


Oven: -80°C for 3 min, -80° to 260°C at 25°/min, hold at 260°C

1. Hydrogen
2. Oxygen
3. Nitrogen
4. Carbon monoxide
5. Methane
6. Ethane
7. Carbon dioxide
8. Ethylene
9. Propane
10. Propylene
11. Isobutane
12. *n*-Butane
13. 1-Butene
14. *trans*-2-Butene
15. *cis*-2-Butene
16. Isopentane
17. *n*-Pentane



Molecular Sieves

- Noble gases, Kr, He, Ar, Xe
- N₂, O₂, H₂, CH₄, CO (not CO₂)
- HP-Molesieve
- CP-Molsieve 5 A
- Very sensitive to moisture/CO₂ in carrier gas/sample
 - Conditioning needed
 - Gas Clean moisture filters in carrier gas essential




Molecular Sieves



Permanent Gases

1. Neon
2. Argon
3. Oxygen
4. Nitrogen
5. Methane
6. Carbon monoxide



GS-CarbonPLOT – Permanent Gases

Column: **GS-CarbonPLOT**

30 m x 0.32 mm I.D., 3.0 μ m


P/N: **113-3133**

Carrier: Helium at 27 cm/sec (calculated off N₂ at 35°C)

Oven: 35°C isothermal

Injector: Split 1:30, 185°C, 250 μ L injection volume

Detector: TCD, 150°C



Oxygenates in Hydrocarbons

Oxygenates blended in gasoline


- Additives to boost octane content, prevent engine “knocking”
- MTBE, ETBE, ethanol
- **Percentage level oxygenates**
- **GC analyses on WAX or TCEP polar liquid phase columns**

Oxygenate in intermediates (monomers, naphthas)

- Lower effectiveness catalysts – lower yield
- Higher costs of catalysts
- More refinery downtime
- **Oxygenates at ppm level**
- **GC analyses on Lowox, OxyPLOT columns, Agilent exclusives**

Analyzing Oxygenates in Hydrocarbon Matrix


- Low ppm concentration level of oxygenates – low analyte capacity of stationary phase
- FID detection (MSD uncommon)
- Agilent columns with high selectivity hydrocarbons/oxygenates: CP-Lowox/OxyPLOT
- Very moisture sensitive 
- Use Gas Clean moisture filters in carrier gas




ASTM Trace Oxygenate Analysis Methods

Valve Configuration


Transfer of oxygenates
Valve off



Venting hydrocarbons
Valve on



Oxygenate Analysis




Select Low Sulfur – PLOT

- Super permeable and porous stationary phase
Proprietary material
- Unique selectivity characteristics for volatile sulfurs/hydrocarbons
No detector quenching for propane/propylene matrix
- Highly inert for volatile sulfurs, H₂S, methyl mercaptan
- One part number: 60 m x 0.32 mm, p/n CP8575
- Near 100% recovery for H₂S at 20 to 100 ppb levels
- Zero particle loss
- FPD, PFPD, SCD, and AED compatible

Sulfur Components in Propylene by PFPD

- H_2S and COS at 500 ppbv
- Same response
- No absorption H_2S



Technique: GC-PFPD
Column: Agilent J&W Select Low Sulfur, 60 m $\times$ 0.32 mm (p/n CP8575)
Oven: 65 °C isotherm
Carrier gas: Helium, constant flow, 2.0 mL/min
Injector: 200 °C, split 1:20
Detector: PFPD, 200 °C
Sample: Propylene matrix containing ~500 ppb H_2S , COS, and CH_3SH
Injection volume: 1 mL
Injection: Gas sampling valve


Considerations for PLOT Column Analysis

Columns

- Selectivity – consider elution pattern
- Capacity – overloading – 0.53 mm
- Inertness (low concentration sulfurs)
- Temperature limits


Contamination

- Efficiency loss, “ghost peaks”, increase in bleed
- Water, CO₂, high molecular weight hydrocarbons
- Carrier gas purifiers are highly recommended



Summary

- PLOT columns are an excellent option for achieving separation of solutes that are gasses at non-cryo/room temperature conditions.
- PLOT columns have a lower capacity than equivalent diameter WCOT columns.
- Due to their high degree of retention, be careful to not inject solutes that are too high boiling or incompatible with the phase.
- Be mindful of larger volumetric flow when interfacing to MSD (require TURBO pump or restrictor)
- Alumina columns are excellent for hydrocarbon separations, but not compatible with water or polar compounds.
- For water and polar compounds, porous polymer columns are ideal (PLOT-Q, PLOT-U).
- GS-GasPro/CP-SilicaPLOT are not compatible with alcohols (absorbed).
- Use a particle trap when interfacing with MSD or VALVES (“PT” columns are best).
- All columns are fully QC-tested to assure column-to-column reproducibility.



Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 option 3, option 3:

- Option 1 for GC or GC/MS columns and supplies
- Option 2 for LC or LC/MS columns and supplies
- Option 3 for sample preparation, filtration and QuEChERS
- Option 4 for spectroscopy supplies
- Option 5 for chemical standards (formerly ULTRA)



- gc-column-support@Agilent.com
- lc-column-support@agilent.com
- spp-support@agilent.com
- spectro-supplies-support@agilent.com
- chem-standards-support@agilent.com

CP-SilicaPLOT – Hydrocarbons

Column : Agilent CP-SilicaPLOT, 0.53 mm x 30 m, fused silica
PLOT CP-SilicaPLOT (df = 6 µm) (Part no. CP8570)

Temperature : 50 °C (5 min) → 225 °C, 5 °C/min

Carrier Gas : He, 20 kPa (0.2 bar, 2.9 psi)

Injector : Direct
T = 225 °C

Detector : FID
T = 250 °C

Sample Size : 2 µL

Concentration Range : ppm %

Peak identification

1. methane
2. ethane
3. ethylene
4. acetylene
5. propane
6. cyclopropane ●
7. propylene
8. isobutane
9. butane
10. 1-butene
11. propyne (methylacetylene)
12. 1,3-butadiene
13. isobutene
14. cis-2-butene
15. hexane

