Unlocking Pesticide Analysis Excellence: Proven Tactics for Optimal Results with GC/MS/MS Using Helium and Hydrogen

Becoming a Better Chromatographer Educational Webinar Series January 31, 2024

Anastasia A. Andrianova, Ph.D. GC/MS Applications Scientist





Pesticides in Food

- 5.6 billion pounds of pesticides are used worldwide annually
- 1 billion pounds in the US annually
- The global market for agricultural pesticides was estimated at \$106,530 million in the year 2022
- Projected to reach \$144,820 million by 2028
- Government agencies determine safe levels of pesticide in food



Rev Environ Health. 2009 Oct-Dec; 24(4): 303-309.




Pesticide Testing



- Consumer safety: ensure compliance within maximum residue limits (MRLs)
- Product confidence and quality control





In Today's Presentation...

I. Five key components of successful pesticide analysis with GC/TQ


II. Hydrogen as a carrier gas for GC/TQ analysis of pesticides






Five Key Components of Successful Pesticide Analysis


- 1. Effective sample extraction and matrix cleanup
- 2. Evaluation of the matrix in full scan
- 3. Midcolumn backflushing
- 4. A leak-free GC/TQ
- 5. Temperature-programmed MMI inlet with a 2 mm dimpled liner














1. Effective Sample Extraction and Matrix Cleanup



Adequate sample preparation (especially with the GC/TQ 7010C) Sometimes dilution may yield better sensitivity and enhance robustness



The Enhanced Matrix Removal (EMR) with Carbon S Passthrough Cleanup for Pigmented Food Matrices

- **Novel** Captiva EMR with Carbon S cartridges
 - Carbon S is an advanced hybrid carbon material
 - <u>Simplified</u> passthrough <u>cleanup</u> after QuEChERS extraction
 - Optimized formula based on different matrix complexities and pigment level
 - Improved method performance
- Compared to the traditional GCB included products
 - Efficient pigment removal with the compromise of sensitive analytes loss
 - Poor matrix cleanup for complex pigmented dry matrices









Captiva EMR Passthrough Selection Guide



Captiva EMR-HCF High Chlorophyll Fresh • Spinach, arugula, chard



Captiva EMR-GPD

General Pigmented Dry • Spices, seasoning, herbal medicine



Captiva EMR-GPF General Pigmented Fresh • Berries, peppers, broccoli



Captiva EMR-LPD Low Pigmented Dry • Nuts, tobacco, light pigmented spices



Simplified Passthrough Workflow

Baby spinach



Blackberry




Almond flour




Cinnamon powder





Crude ext. loading

Gravity elution or w/ low level external force



Captiva EMR cartridges

Collected eluent



Posttreatment before instrument analysis







Standard

QuEChERS extraction



Conclusions – Value of Captiva EMR with Carbon S

Data quality

Improved quantitation results quality

- · Improved sensitive pesticides recovery and reproducibility
- Reduced matrix interferences and matrix effect on targets, and thus improved method sensitivity and selectivity.
- Overall higher pass rate for multiclass multiresidue pesticides analysis, plus improved detection method robustness

Ease of use/	
throughput	Simplified preparation procedure

- One sample preparation method can serve both GC and LC applications
- Simplified passthrough cleanup provides easier and faster workflow saving 15 to 30% sample preparation processing time

Five new part numbers Simple and easy selection

- Three for fresh pigmented matrices, Captiva EMR-HCF 1 and 2, Captiva EMR-GPF
- Two for dry complex matrices, Captiva EMR-GPD, Captiva EMR-LPD








2. Evaluation of the Matrix in Full Scan




Adequate sample preparation (especially with the 7010C) Sometimes dilution may yield better sensitivity and enhance robustness



Evaluation of in-source loading of the matrix in full scan: Sometimes less is more







3. Midcolumn Backflushing



Adequate sample preparation (especially with the 7010C) Sometimes dilution may yield better sensitivity and enhance robustness




Evaluation of in-source loading of the matrix in full scan: Sometimes less is more



Backflushing



3. Midcolumn Backflushing



If no backflushing



3. Midcolumn Backflushing



Backflushing eliminates the need for bakeout:

- No deposition of high-boilers and GC column stationary phase into the EI source
- Increases GC column lifetime
- Significantly reduces the need for column trimming



Self Tightening Column Nut Installation – Inlet and Detectors <u>https://www.agilent.com/en/video/stcn-inlet-detector</u> Self Tightening Column Nut Installation – MS Interface <u>https://www.agilent.com/en/video/stcn-mass-spec</u> Gold-plated Flexible Metal Ferrules for GC – p/n G2855-28501



Keys to Unlocking Excellent Performance in Pesticides Analysis with GC/TQ

Evaluation of in-source loading of the matrix in full scan: Sometimes less is more



Adequate sample preparation (especially with the 7010C) Sometimes dilution may yield better sensitivity and

3



Keep the leaks out

Backflushing





Agilent

Keys to Unlocking Excellent Performance in Pesticides Analysis with GC/TQ



Evaluation of in-source loading of the matrix in full scan: Sometimes less is more



Adequate sample preparation (especially with the 7010C) Sometimes dilution may yield better sensitivity and enhance robustness







Backflushing

Keep the leaks out



Use of the temperature-programmed MMI with a 2 mm dimpled liner






https://www.agilent.com/cs/library/applications/application-sensitive-and-robust-detection-of-pesticides-in-dried-cannabis-8890-7010b-triple-quadrupole-gc-ms-5994-0568en-agilent.pdf



Results

App note 5994-4965EN





Food & Beverage Testing

<u>5994-4965EN</u>

Five Keys to Unlock Maximum Performance in the Analysis of Over 200 Pesticides in Challenging Food Matrices by GC/MS/MS



Authors

Abstract


Anastasia A. Andrianova and Limian Zhao Agilent Technologies, Inc.

This application note describes five best practices to enhance analytical performance in the analysis of over 200 pesticides in challenging matrices including spinach, walnut, and cayenne pepper. The novel Agilent Captiva EMR passthrough cleanup procedure following the Agilent QuEChERS extraction enabled a cleaner matrix background. The cleanup and extraction reduced matrix interferences with target analytes and extended the maintenance-free operation time of the instrument. Calibration performance was demonstrated over a wide dynamic range to over four orders of magnitude. It was shown that the Agilent 8890/7000E triple quadrupole GC/MS system achieved excellent linearity over a concentration range of 0.1 to 5,000 ppb. The Agilent 8890/7010C triple quadrupole GC/MS system demonstrated superior sensitivity yielding a higher signal-to-noise ratio at lower concentrations.




High Sensitivity for Deltamethrin and Pentachloronitrobenzene in Walnut

Deltamethrin: MRL in walnut 100 ppb



Pentachloronitrobenzene: No MRL in walnut But regulated in many vegetables and fruits (20 ppb – 1 ppm) Cannot be done with LC/MS




Calibration Performance Over a Wide Dynamic Range with the 7000E GC/TQ

Bifenthrin: encompassing varying MRLs in different commodities


Spinach








Calibration Over a Wide Dynamic Range with the 7000E and 7010C GC/TQ Pyriproxyfen and fludioxonil: widely different MRLs in spinach (3,000 and 10 ppb)





Calibration Over a Wide Dynamic Range with the 7000E and 7010C GC/TQ Pyriproxyfen and fludioxonil: widely different MRLs in spinach (3,000 and 10 ppb)



Fludioxonil: MRL in spinach 10 ppb



Calibration Performance for the 203 Pesticides with the 7000E and 7010C GC/TQ in Spinach, Walnut, and Cayenne Pepper

Number of compounds with R2>0.99 and their calibration ranges with the 7000E GC/TQ ranges with the 7010C GC/TQ 200 180 168 174 180 155 160 160 140 Number of target compounds Number of target compounds 136 140 130 120 120 100 90 00 80 80 64 60 60 40 33 40 17 17 17 20 20 13 Q 0 \cap Linear, 0.1- Linear, 0.1- Linear, 0.5-Quadratic, Linear, 0.1- Linear, 0.1- Linear, 0.1- Linear, 0.5-Quadratic, Quadratic, Quadratic, Linear, Linear, 5,000 ppb 1,000 ppb 5,000 ppb 0.1-5,000 0.5-5,000 1,000 ppb 500 ppb 250 ppb 1,000 ppb 0.1-1,000 other other other ppb ppb ranges ranges ranges ■ Spinach ■ Walnut ■ Cayenne Pepper Spinach Walnut Cavenne Pepper

Number of compounds with R2>0.99 and their calibration





Quadratic


other

ranges

5

ppb


Spinach QuEChERS Extract – for Longevity Study 7000E (pesticides spiked at 20 ppb), 10 minute analysis





Longevity: Area Normalized by the ISTD, 20 ppb in Spinach with 7000E (10 min) Normalized response of pesticides at 20 ppb spiked into QuEChERS EMR-HCF spinach

Area normalized by alpha-BHC-d6, %



• BHC-alpha (benzene hexachloride), 8% • BHC-beta, 9%

- Pirimiphos-methyl, 9%
- Bupirimate, 11%
- Metalaxyl, 10%


- - Bromophos-ethyl, 10%
 - Chlorthiophos , 11%

- Atrazine, 11%
- Prothiofos, 10%
- Fluquinconazole, 13%



Five Best Practices in Sample Preparation and GC/MS/MS Analysis

- Simplified and improved sample preparation achieved with the novel and improved Agilent Captiva EMR pass-through cleanup, following traditional Agilent QuEChERS extraction
- Evaluation of in-source loading of the matrix in full scan data acquisition mode
- Midcolumn backflushing
- Leak-free GC/triple quadrupole system enabled with Self Tightening collared column nuts and CFT gold-plated flexible metal ferrules
- Use of temperature-programmed multimode inlet with a 2 mm dimpled liner (no glass wool)
- Robustness over 700 injections



Application Note Food & Beverage Testing



7000E GC/TQ

Agilent


How About Hydrogen Carrier Gas?





Hydrogen As a Carrier Gas in GC/MS/MS


- Hydrogen is the preferred carrier gas for GC/MS
- Ongoing shortages and increased prices cause unpredictability







Agilent EI GC/MS Instrument Helium to Hydrogen Carrier Gas Conversion 5994-2312EN



- Contains detailed instructions for method conversion from He to H₂ carrier.
- **Essential read** before converting to H₂.



Hydrogen for Pesticide Residue Analysis: How to Convert GC/MS/MS Analysis from Helium to Hydrogen and Meet the MRLs

Webinar available ondemand: Link to watch

Separation Science

On-Demand Presentation

Hydrogen for Pesticide Residue Analysis: How to Convert GC/MS/MS Analysis from Helium to Hydrogen and Meet the MRLs

On-Demand


Key strategies for pesticide analysis utilizing hydrogen as the carrier gas, while maintaining sensitivity that allows for meeting the MRLs.



GC/MS Recipe for Success with Hydrogen

- H₂ carrier gas
- HydroInert Source optimized for use with hydrogen
- 20 × 20 m (0.18 mm × 0.18 µm) HP-5ms UI
- Starting column flows: 0.9 and 1.1 mL/min
- Retention time lock to chlorpyrifos methyl (9.143 min)
- <u>Same retention times</u> with H₂ as with the conventional He method (20 minutes)
- Ability to perform 10 minute analysis with excellent resolution





Hydrogen Carrier Gas for Analyzing Pesticides in Pigmented Foods with GC/MS/MS




Abstract

Anastasia A. Andrianova, Bruce D. Quimby, and Limian Zhao Agilent Technologies, Inc.

Authors


This application note describes the key strategies for pesticide analysis with gas chromatography/triple quadrupole mass spectrometry (GC/TQ) using hydrogen as the carrier gas while maintaining sensitivity to meet maximum residue limits (MRLs). The key aspects addressed in this work include the recommended column configuration, the optimized injection conditions, and the appropriate choice of the mass spectrometer (MS) electron ionization (EI) source hardware developed for use with hydrogen carrier gas. The 20 m × 20 m (0.18 mm × 0.18 µm) Agilent HP-Sms UI midcolumn backflush configuration allowed for maintaining the same retention times as with helium, leading to time savings associated with method translation. The resulting chromatographic resolution achieved under the optimal emedicine.




and the Agilent High Efficiency Source (HES) resulted in nearly identical spectra observed with hydrogen and helium, which allowed using the same multiple reaction monitoring (MRM) transitions and collision energies as with helium. The ability to use the same MRMs, collision energies, and retention times greatly simplified the transition from helium to hydrogen.



GC Configuration and Method Considerations



Cannot be used with H₂ carrier gas because of insufficient <u>inlet pressure</u>







GC Configuration and Method Considerations

20 x 20 m 0.18 mm x 0.18 μm

Compatible with H₂ carrier gas Analysis time **20 minutes**

Maintain the same RTs and elution order 30 m, 0.25 mm x 0.25 μ m \rightarrow 20 m, 0.18 mm x 0.18 μ m







GC Configuration and Method Considerations


20 x 20 m 0.18 mm x 0.18 μm


Compatible with H₂ carrier gas Analysis time **20 minutes**

10 m x 10 m 0.18 mm x 0.18 μm

Compatible with H₂ carrier gas Analysis time **10 minutes**

Maintain the same RTs and elution order 30 m, 0.25 mm x 0.25 $\mu m \rightarrow$ 20 m, 0.18 mm x 0.18 μm








What Happens to the Retention Times and Resolution?






Exact Retention Time Matching Achieved with Method Translation and Retention Time Locking






Superior Resolution with Hydrogen


EIC 163






Proof of Concept: <u>10 Minutes</u> with Hydrogen











10 Minutes with Hydrogen and Still Good Resolution



Electron Ionization MS Sources

Conventional El Extractor (Inert Plus) Source








40



Calibration Performance in Spinach QuEChERS Extract for 203 Pesticides





Summary of Analyzing Pesticides with H₂ Carrier Gas

- Over 90% of the evaluated pesticides could be detected at or below the default MRL of <u>10 ppb</u> in the vial
- However, just like with He, some pesticides may be more challenging to analyze at low concentrations
- Using **optimized GC conditions (injection, column set)** is essential when using H₂
- The optimized setup with H₂ improved chromatographic resolution and allowed for precise matching of retention times with He
- The same MRM transitions, with the same collision energies for the targets eluting at the same retention times as with He, could be used with H₂.

GC/TQ EI source considerations

The HydroInert EI source:


- Improved sensitivity compared to the standard EI with H₂
- Best spectral fidelity

The HES EI source:

• Best sensitivity observed with H₂



Agilent GC/MS Best Practices Infographic link



Dilute More, Inject Less The GC/MS High-Efficiency Source (HES) is designed for ultimate sensitivity where a little goes a long way.



interferences from

matrix background

horizontal line

GC/MS operation with excellent ultra-trace level sensitivity



Reduce

analysis time

Use Backflush

operation of the GC/MS system.



Extend column life



Significantly

reduce carryover


Save time cleaning, reconditioning, and recalibration

Achieve more consistent responses for target analytes

Implement JetClean

Keep the GC/MS source clean for increased sample throughput

Patented JetClean technology uses a gentle flow of hydrogen through the source.



Check for Leaks

hardware and software tools to aid in identification of leaks at their source.





Leak test software functionality for GC/TQ identifies the source of leaks in real time when using a leak testing gas.

🕂 Agilent



43

