

# Agilent Dako CoverStainer—the Efficient, Safety-Focused, and Accurate Choice for H&E



## Continuous delivery of stained slides in a lean workflow

Agilent Dako CoverStainer fully automates every step of the H&E process, from baking to drying, reducing manual processes and streamlining the workflow.

The instrument—combined with Dako Ready-to-Use (RTU) reagents, Dako validated, preoptimized protocols, and the Dako Reagent Management System (DakoRMS)—provides you with a solution that is efficient, safety-focused, and provides accurate staining results.

Dako CoverStainer enables you to achieve 240 stained slides per hour, with only 46 minutes until your first coverslipped slides are ready.

## The Dako H&E Solution includes:

### Validated Dako Ready-to-Use H&E reagents

- Three staining intensities
- Consistent for 3,000 slides or up to five days use (whichever comes first)

### Validated Dako H&E protocols

- Capacity for 19 customized protocols

### Dako Reagent Management System

- Constantly circulating reagents
- Optimizes consumption and prevents precipitation

### Unique rack design

- Ensures only slides are lowered into the reagent to minimize carryover
- Eliminates risk of cross contamination with slides side-by-side<sup>1</sup>

# Dako CoverStainer Specifications

**Table 1.** Instrument specifications.

Specification	Description
Instrument	<p>The Dako CoverStainer is a fully automatic working station that covers:</p> <ul style="list-style-type: none"> <li>– Baking</li> <li>– Dewaxing</li> <li>– Staining</li> <li>– Coverslipping</li> <li>– Drying</li> </ul>
Labeling	In Vitro Diagnostic Regulation (IVDR) (EU) 2017/746
Quality system	ISO13485
Instrument dimensions	<p>150 cm W x 67.2 cm D x 128.0 cm H</p> <p>105 kg (231.5 lbs)</p>
Electrical requirements/power	<p>100-120V~ 6A 50/60Hz</p> <p>Fuses: 2 x T6.3AL 250V</p> <p>Inlet voltage setting: 115V</p> <p>210-240V~ 3A 50/60Hz</p> <p>Fuses: 2 x T3.15AL 250V</p> <p>Inlet voltage setting: 230V</p>
Networking/LIS connection	<p>LAN: preferably with DHCP server</p> <p>Cable: minimum CAT-5 patch cable, maximum length: 3 m (10 ft)</p>
External ventilation	<p>110 mm (4") duct</p> <p>Airflow range: 80–140 m³/h required (47–82 ft³/min)</p>
Waste/drain	<p>Sink or drain required within 2 m (6.6 ft)</p> <p>Individual reagent waste drains back to original bottle</p> <p>One bulk water waste bottle which empties automatically (10 L)</p> <p>One backup bulk hazardous waste bottle, manually emptied, if required (10 L)</p>
Environmental conditions	<p>The work area for the Dako CoverStainer must be an environment with an ambient temperature between 18°C-26°C (64°F-79°F) and 25-60% relative humidity (RH). The instrument must not be exposed to direct sunlight.</p> <ul style="list-style-type: none"> <li>– Noise level &lt; 50 dB</li> <li>– Heat emission &lt; 220 W</li> </ul>
Water supply	<p>¾" threat tap</p> <p>Pressure range: 2-4 Bar</p> <p>Connection to tap water. Not all sources of water may be of sufficient quality for H&amp;E staining. It is recommended to use reagent-grade water corresponding to Clinical Laboratory Reagent Water (CLRW) standard as specified by CLSI, or water of similar quality.</p>

**Table 2.** Staining Module.

Specification	Description
Heating chamber capacity	12 racks at a time, maximum 120 slides
Heating chamber temperature	Minimum 60 °C and maximum 70 °C (140–158 °F)
Slide rack capacity	1 to 10 glass slides per slide rack
Loading of slide racks	Slide racks can be loaded and unloaded continuously
Capacity	Up to 240 slides per hour
Overnight run	Maximum 240 slides 120 slides coverslipped and dried 120 slides waiting in exit bath before coverslipping
Reagent stations/reagent bottles	18 dip tanks divided into two compartments which give 36 reagent stations  Dip tanks 1-2 and 14-16 (i.e. reagent stations 1-4 and 27–32) are reserved for xylene or xylene substitute due to special requirements to the solvent pumps  Dip tanks 17-18 (i.e. reagent stations 33–36) are reserved for running tap water  One reagent station contains approx. 300 mL 20 reagent bottle positions (1,000 mL)  All reagents circulate constantly during processing (DakoRMS) 20 membrane pumps (five pumps for xylene or xylene substitute, 15 pumps for other reagents) Each pump uses one reagent bottle and can go to one reagent station or two dip tank compartments
Reagent capacity	Dako CoverStainer reagents (Dako Hematoxylin, Dako Eosins, Dako Differentiation Solution and Dako Bluing Buffer) or user-defined reagents. Please note: We recommend reagents to be changed every five days or for every 3,000 slides (whichever comes first).
Protocols	20 protocols, of which 19 can be programmed by the user  It is possible to change number of steps, reagents, and time in the protocols  Each rack can be loaded with any of the 20 protocols as needed  Three validated, preoptimized protocols for both progressive and regressive staining
Alarms	Alarm notification when protocol is completed  Alarm notification, if error occurs (screen error message pops up)

**Table 3.** Coverslipping.

Specification	Description
Cover glass size	24 x 50 mm
Mounting medium	Validated and recommended by Agilent Technologies: <ul style="list-style-type: none"> <li>– CS703 Mounting Medium</li> <li>– CS705 Mounting Medium, Toluene-Free</li> </ul>
Programming of mounting parameters	Placement of both cover glass and amount of mounting medium applied to each slide can be programmed

## Dako CoverStainer

### The efficient, safety-focused, and accurate choice for H&E

- Walkaway solution designed for LEAN processes
- Processes up to 240 H&E slides per hour
- First batch of coverslipped slides in as little as 46 minutes
- Proven consistent staining quality—3,000 slides or five days' usage (whichever comes first)
- Three sets of validated Dako CoverStainer reagents
- Prefilled RTU reagent bottles minimize spills and reduce handling
- Cost efficient waste management

## Reference

1. Case Study on Dako CoverStainer for H&E - No Risk of Cross Contamination, **2019**, 38686)

[www.agilent.com](http://www.agilent.com)

D0111564\_1.00

This information is subject to change without notice.

© Agilent Technologies, Inc. 2016-2024  
Published in the USA, June 19, 2024  
38645 2024MAY23