



Agilent Technologies

Innovating the HP Way

1100 Series DAD Wavelength Calibration Test Procedure Using the Chemstation

DAD Wavelength Calibration Test

The detector uses the alpha (656.1 nm) and beta (486 nm) emission lines of the deuterium lamp for checking wavelength calibration. The sharp emission lines enable more accurate calibration check than is possible with holmium oxide. When calibration test is started, the 1 nm slit is moved into the light path automatically, and the gain is set to zero. On completion of the scan, the alpha- and beta-line deviations (in nm) are displayed. These values indicate how far the detector calibration deviates from the actual positions of the alpha and beta emission lines. To eliminate effects due to absorbing solvents, the test should be done with water in the flow cell (during the DAD Self Test, the flowcell should be removed). After calibration, the Holmium Oxide Test provides verification of wavelength accuracy at three additional wavelengths. The wavelength calibration test should be done after maintenance of the flow cell, lamp exchange, or after major repair. If test indicates the detector calibration is incorrect, Wavelength Recalibration should be done.

Caution: Please note that this test requires the spectral evaluation module G2180AA to be installed.

DAD Wavelength Calibration Test Results

Limits

Absorbance Maxima

486.0 nm

656.1 nm

Maximum deviation

Limits

485.5 - 486.5 nm

655.6 - 656.6 nm

-0.5nm to +0.5nm

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

Probable causes of test failure

Cause

Absorbing solvent or air bubble in flow cell.
with water.

Incorrect calibration

Dirty or contaminated flow cell.

Dirty or contaminated optical components
(achromat, windows).

Corrective action

Ensure the flow cell is filled

Recalibrate (see Wavelength
Calibration) and repeat the
test.

Run the Cell Test. If the test
fails, exchange the flow cell
windows.

Clean optical components
with alcohol and lint-free
cloth.